

VCHECK V200

PROGESTERONE TIMING CHART

results using different analyzers will not render useful comparisons. All dogs exhibit variability and deviation from the standard rules, in best practice, adjust the chart as needed for specific animals. Failure to achieve pregnancy may be due to

Test Frequency Recommendation

- If Value is between 1-3 ng/mL, test again in 3 days
- If value is between 3-8 ng/mL, test every 2 days
- If value is 7 or greater, test every day

Result Interpretation

Anestrus (out of heat) <1.0

2-3 ng/mL **LH Surge**

5-7 ng/ml **Estimated Ovulation**

7-10 ng/mL Egg Maturation Period

Allow two days post ovulation for maturation, a number in this range does not guarantee ovulation. Test again to confirm if in fertile

10-18 ng/ml Fertile Range

Natural or AI Inseminations

18-22 ng/ml Surgical-Al-Frozen

Reverse Progesterone

>3.0 ng/mL To Maintain Pregnancy

<3.0 ng/ml **Nearing Whelping**

<36 hours before whelping, observe temperature



Standard Timing By

Ovulation

Al With Fresh or Chilled Semen, and Natural Breed

2-4 Days After Ovulation

- *verify good sperm counts and forward motility
- *ovulation must be confirmed and result in fertile range required

Al or TCI With

Cryofrozen Semen

3-5 Days After Ovulation

- *verify health of sperm being used for breeding
- *instrument not recommend for surgical procedures

The ranges listed above has been compiled through a professional 720 test correlation with 3 other automated immunology analyzers. The ranges listed above has been compiled through a professional 720 test correlation with 3 other automated immunology analyzers. Peer group data has also been used. Numbers from other analyzers will always have different ranges since every analyzers. er has a different method of measurement. This chart may be adjusted slightly with the onset of new kit lot data provided. Intellectual property of MR Diagnostic Services 2022

MRDS Support@MRDiagnosticServices.com



Available Test Reagents

Cardiac

Detection and monitoring of heart disease in dogs and cats, including MMVD and DCM

Canine NT-proBNP

Feline NT-proBNP

Acute Inflammation

Earlier detection of acute inflammation in dogs and cats

Canine CRP

Feline SAA

Pancreatitis

Definite diagnosis of Pancreatitis in dogs and cats

 cPL (canine Pancreasspecific Lipase) fPL (feline Pancreasspecific Lipase)

Coagulation

Highly sensitive marker for thromboembolism

· Canine D-dimer

Hormones

A hormone marker for Cushing's and Addison's diseases

· Canine Cortisol

A hormone marker for diagnosis of thyroid gland diseases

• T4 (Thyroxine)

· Canine TSH

For determining optimal breeding dates

· Canine Progesterone

Specifications

Model Vcheck V200
Test capacity 1 test at a time

Power AC/DC adaptor 120V
Display 7" Color Touch Screen

Printer Built-in

Connectivity HL7 v2.6(PCD-01) / POCT1-A **Dimension** 7.8 " Wide x 9.4" Deep x 8" High

Weight 5.51bs

Simple Testing Procedure







 Dilute the sample

2. Add diluted sample into a test device

Check the test result

Product Name	Product Number	Product Type	Packing Unit
Vcheck Canine NT-pro BNP	VCF132DC	Device	5 Tests/Kit
Vcheck Feline NT-pro BNP	VCF130DC	Device	5 Tests/Kit
Vcheck Canine CRP 2.0	VCF109DD	Device	10 Tests/Kit
Vcheck Feline SAA 2.0	VCF128DD	Device	10 Tests/Kit
Vcheck cPL 2.0	VCF129DD	Device	10 Tests/Kit
Vcheck fPL 2.0	VCF127DD	Device	10 Tests/Kit
Vcheck D-dimer	VCF107DD	Device	10 Tests/Kit
Vcheck cCortisol	VCF105DD	Device	10 Tests/Kit
Vcheck T4	VCF106DD	Device	10 Tests/Kit
Vcheck CanineTSH	VCF118DC	Device	5 Tests/Kit
Vcheck cProgesterone	VCF122DD	Device	10 Tests/Kit



For More (Information: CanineP4.Com Info@CanineP4.com 844-673-7378





