

## **HEALVET 300**

### **PROGESTERONE TIMING CHART**


This Canine Progesterone Timing Chart is designed to work only with the Vet test kits. Comparing values to laboratory results using different analyzers will not render useful comparisons. All dogs exhibit variability and deviation from the standard rules, in best practice, adjust the chart as needed for specific animals. Failure to achieve pregnancy may be due to improper use of chart, incorrect testing techniques, or poor semen quality

Test Frequency Result Interpretation				
Recommendation	<1.0		Anestrus (out of heat)	ATT BOLLS
<ul> <li>If Value is between <u>1-3 ng/mL, test</u></li> </ul>	1.5-3 ng/	mL	LH Surge	
again in 3 days	4-6 ng/n	nl	Estimated Ovulation	
If value is between			Egg Maturation Period	Standard Timing By
<u>3-8 ng/mL, test</u> <u>every 2 days</u>	6-8 ng/mL		Allow two days post ovulation for maturation, a number in this range does not guarantee ovulation. Test again to confirm if in fertile range	Ovulation
• <u>If value is 8 or</u>				Al With Fresh or Chilled
greater, test every				Semen, and Natural Breed <u>2-4 Days After Ovulation</u>
day	8-15 ng/	ml	Fertile Range	*verify good sperm counts and
			Natural or AI Inseminations	forward motility *ovulation must be confirmed
			2-3 days past ovulation	and result in fertile range required
	15-20 ng/	/ml	Surgical-AI-Frozen	Al or TCl With
· · · ·			Late stage fertility, 4-5 days	Cryofrozen Semen
Reverse Progesterone *verify health of sperm being				
>5.0 ng/mL To Maintain Pregnancy			used for breeding	
		Noarir	ag W/bolping	*instrument not recommend for surgical procedures
<3.	·• ··6/ ····		ng Whelping	
This shart may be adjusted slightly with the	ansat of now kit lot data prov	vidad Intal	ectual property of MR Diagnostic Services 2023. Do not copy,	distribute or make alternations to this desument

This chart may be adjusted slightly with the onset of new kit lot data provided. Intellectual property of MR Diagnostic Services 2023. Do not copy, distribute, or make alterations to this document without expressed written consent from MR Diagnostic Services LLC

# MRDS Support@MRDiagnosticServices.com

#### **Ordinary Hormone Levels in Canines from Anestrus Through Pregnancy**



#### How do I track the ovulation date?

<u>Ovulation dates are best tracked by running another test again the next day</u> after her result places her in the "Estimated Ovulation" window. Bitches that ovulate above the recommendations, and ovulation stalling can hamper successful pregnancies if they occur unknowingly. For that reason, continuing testing well into the Green fertile range. Some bitches will <u>double</u>, or jump, 1-2 days after truly ovulating, and others may stall or progress more slowly. <u>You may see a large increase in a mere 24 hours, or in a week, depending on the female.</u>

If, after testing the next day, she is still in the ovulation range, she DID NOT ovulate the previous day and may be experiencing an ovulation stall, or split heat. Continue to monitor and test daily.

Having her result place her in the Estimated Ovulation range DOES NOT GUARANTEE she has actually ovulated as most bitches will ovulate based on a number of internal and environmental factors, and do not all ovulate at the exact same value. She may remain in the estimated ovulation for days in some cases.

#### **CONVERSION TO NMOL/L**

To Convert to nMol/L, simply multiply the result in ng/mL by 3.18.

IE, 10 ng/mL (times 3.18) = 31.80 nMol/L