



Quick Scan Mini™

PROGESTERONE TIMING CHART

This Canine Progesterone Timing Chart is designed to work only with the QuickScan Mini Vet test kits. Comparing values to laboratory results using different analyzers will not render useful comparisons. All dogs exhibit variability and deviation from the standard rules, in best practice, adjust the chart as needed for specific animals. Failure to achieve pregnancy may be due to improper use of chart, incorrect testing techniques, or poor semen quality

Result Interpretation

Test Frequency Recommendation	<1.0	Anestrus or Proestrus	Standard Timing By Ovulation
<ul style="list-style-type: none"> If result is between 0-1 ng/mL, test every 3 days 	1-3 ng/mL	LH Surge progesterone begins rising, onset of Estrus cycle	
<ul style="list-style-type: none"> If result is between 2-4 ng/mL, test every 1-2 days 	3-5 ng/ml	Estimated Ovulation Test again the following day(s) Confirm ovulation has occurred with a result at or beyond the Egg Maturation period prior to breeding	
<ul style="list-style-type: none"> If value is between 5-6 ng/mL, test again the next day 	5-6 ng/mL	Egg Maturation Period Egg cells require up to 48 hours to mature after fertilization. Once ovulation is confirmed, breed 2-4 days after the ovulation date. A result in this window confirms ovulation has occurred, the result will continue rising	
	6-9 ng/ml	Fertile Range Natural or AI Inseminations Breed 2-4 days after the confirmed ovulation date	Example A result of 3.4 ng/mL occurred on the 5th. On the 6th, the dog retested and resulted in 5.4 ng/mL. Ovulation has been confirmed, and breeding should be done between the <u>7th and 9th</u> (2-4 days later)

To Convert to nMol/L

Multiply ng/mL by 3.18




The ranges listed above has been compiled through a professional 720 test correlation with 3 other automated immunology analyzers. Peer group data has also been used. Numbers from other analyzers will always have different ranges since every analyzer has a different method of measurement. This chart may be adjusted slightly with the onset of new kit lot data provided. Intellectual property of MR Diagnostic Services 2022

Reverse Progesterone

>3.0 ng/mL To Maintain Pregnancy

<3.0 ng/ml Nearing Whelping

<24 hours before whelping, observe temperature drop and behavior



Individual Result Tracking Chart

[illegible]