

For Screening of Thyroid Diseases in Canines and Felines

An in vitro diagnostic test kit for the in-clinic, quantitative measurement of T4 (thyroxine) concentration in canine serum or feline serum, Vcheck's T4 test allows for the identification of this hormone biomarker. Canine hypothyroidism is a common disease related to thyroid function in dogs, and feline hyperthyroidism is the most common endocrine disease affecting geriatric cats. Vcheck's T4 test measures the concentration level of this hormone and can be used to screen for these diseases. The prognosis of hyperthyroidism and hypothyroidism is excellent when canines and felines are diagnosed at an early stage and appropriately treated and managed.

Clinical Applications

· Screening of hypothyroidism/hyperthyroidism and treatment monitoring

Specifications

Species	Canine, Feline	
Sample Type	Serum 50 µl	
Measurement	Quantitative	
Range	0.5 - 8 µg/dL	
Testing Time	20 minutes	
Storage Condition	2 - 8º C	

Simple Testing Procedure



to draw 50 ul

to the assay

diluent tube.



Mix & Incubate

Use a 50 µl pipette Use a tablet pipette to mix until the tablet is of serum and add completely dissolved (8 times or more). Incubate the mixture at room temperature for 10 min.

Lilv

Measure

Add 100 µl of incubated sample to the sample well of the test device and press [START].

Product Name	Product Number	Product Type	Packing Unit
Vcheck T4	VCF106DD	Device	10 Tests/Kit

BIONOTE

A Closer Look: T4


The Vcheck T4 test kit screens for T4, a major thyroid hormone that is important for normal regulation of metabolic rates and activity in various organs. When T4 concentration level is measured, it can detect hypothyroidism/hyperthyroidism in canines and felines.



Clinical Signs of Thyroid Diseases

Specific Clinical Application

The Vcheck T4 test kit is based on the competitive immunoassay method for the quantitative measurement of T4 concentration. T4 is one of the hormones secreted by the thyroid gland and is responsible for the regulation of metabolism. The Vcheck T4 can be used to assist in the monitoring of structural or functional abnormalities of the thyroid gland.



bionote.com customerservice@bionote.com 800-727-5169

from **BIONOTE**

All of Bionote's Vcheck biomarker tests are available for use on the Vcheck V200 and V2400 analyzers.

