



Eagle  
Dental

Feel the Difference



---


Endo Access Solutions

# Eagle Dental Endo Access Enables You to Perform with Confidence.

## ENDO ACCESS BURS

### Endo Access on Zirconia-Based Crown

One of the most daunting tasks faced by endodontists and dentists is a root canal through zirconia. Strauss has solved this problem with a unique line of Endo Access burs, the 801-018PF, 801-022PF, 807-012PF, and 807-016PF. These burs tunnel through the porcelain layer in seconds, penetrating through zirconia lithium disilicate based crowns with efficiency previously thought impossible.



To widen the access hole use the 856-018PF.

### Zirconia and Ceramic-Based Endo Access Burs Fast, safe, and without micro-fracturing.

For Endo Access it is recommended to move the handpiece from side to side with a light pressure while using coolant.

FG

Cat. No.	801-018PF	801-022PF	807-012PF	807-016PF	856-018PF
Head Size	018	022	012	016	018
Head Length			4	4	8

Pack of 6

### Endo Access to PFM Crown and Amalgam

- For PFM-based crowns use a Magic Touch round or inverted cone shape bur to cut through the ceramic until metal is showing. Then use one of the Barracuda Metal cutting burs to penetrate through the metal structure.
- For Endo Access through amalgam use one of the Barracuda metal cutting burs.

FG

Cat. No.	B1558	B331	B4
Head Size	012	010	014
Head Length	4.4	1.7	
FG	◆	◆	◆


Metal Cutting Burs

Pack of 25



Natural Teeth Endo Access


To extend your reach, use an FGOS or RAOS 2, 4, or 6 Surgical Length Carbide bur.



Plain Cut  
FG RA


Cat. No.	2	4	6
Head Size	010	014	018
Head Length			
OS	◆	◆	◆
RAOS	◆	◆	◆

Round Surgical Length



Pack of 10

Ideal for initial penetration to gain straight-line access to all canals and minimize chances of file separation.



FG


Cat. No.	801-012	801-014	801-018	801L-012	801L-014	801L-018	801XL-014	801XL-018
Head Size	012	014	018	012	014	018	014	018
Head Length								
Medium	▽	▽	▽	▽			▽	▽
Coarse	▽	▽	▽	▽	▽	▽	▽	▽
S. Coarse	▽	▽	▽	▽	▽	▽	▽	▽

Round      Surgical Length Round 24mm      Extra Long 28mm

Pack of 6

Round Safe-End

The Safe-End diamond bur is designed to widen and flare the pulp chamber safely without perforation, and to achieve an optimal straight-line access to the canal.

Safe End  
FG


Cat. No.	851-010	851-012	851-016	857-012	857-014	857-016
Head Size	010	012	016	012	014	016
Head Length	8.0	8.0	8.0	10.0	10.0	10.0
Fine	▽	▽			▽	
Medium	▽	▽	▽	▽	▽	▽

Round Safe-End Taper

Pack of 5

Endo Access

Special endodontic burs for reducing the pulp chamber roof and extending the access opening for successful root canal treatments. The long diamond-coated instrument neck facilitates producing straight access to the root canal and creating a funnel-shaped access to the pulp chamber. The special diamond grit guarantees rapid and gentle removal of the tooth structure.



FG


Cat. No.	802L-016	802L-019
Head Size	016	019
Head Length	10	10
Coarse	▽	▽

Endo Access

Pack of 6

Round End Tapered

The round end diamond bur is designed to widen the canal wall without the risk of gouging.



FG

Cat. No.	850-014	850-016	850L-012	850L-024	850XL-016
Head Size	014	016	012	24	016
Head Length	10	10	11	11	10
Medium	▽	▽	▽	▽	▽
Coarse	▽	▽	▽	▽	▽
S. Coarse	▽	▽	▽	▽	▽


Round End Tapered      XL 31mm

Pack of 6



Endo Access Z & Safe-End EZ


- The Endo Access Z burs have a **non-safe end**.
- The EZ is a specially designed carbide bur with a non-cutting safe-end tip.
- Safely widens pulp chamber without perforation and ledging.
- Creates funnel shape to allow easier access to root canals.
- Safe ended tip prevents penetration of the pulp chamber floor.



Micro Drills

Efficiently cutting special endodontic drill used for detecting all root canals in the case of multirouted teeth. With a high cutting efficiency and long neck, this drill ensures a clear view in deeper regions of the cavity for opening the access cavity of the root canal system. Use in a contra-angle with micromotor at 1,500 – 2,000 RPM, without cooling, dry and apply only light pressure during preparation.

Swiss-made by Edenta





**FG**

Cat. No.	40151Z	40199Z
Head Size	014	016
Head Length	11	11
Total Legth	23	23
FG	◆	◆


**Endo Access Z**



**FG**

Cat. No.	4219EZ	4152EZ
Head Size	014	016
Head Length	9	9
Total Legth	23	23
FG	◆	◆

**Endo Access Safe-End EZ**



**RA**

Cat. No.	H1SNL.206-006	H1SNL.206-008	H1SNL.206-010
Head Size	006	008	010
Total Length	34	34	34


**34mm Micro Drills**

Pack of 5

Gates Drills

Edenta Gates Drills are used to open the root canal and to further shape and taper the canal. Specially designed to bind while drilling in the root canal.






**RA**

Cat. No.	180.15.204-1	180.15.204-2	180.15.204-3	180.15.204-4
Head Size	050	070	090	110
Neck Length (L)	15	15	15	15

**28mm Gates Drills**



**RA**

Cat. No.	180.19.204-1	180.19.204-2	180.19.204-3	180.19.204-4
Head Size	050	070	090	110
Neck Length (L)	19	19	19	19


**32mm Gates Drills**

Pack of 6

Peeso Drills

Edenta Peeso Drills are specially designed to enlarge the cervical portion of the canal and to prepare the canal entrance. Useful in removing gutta percha in preparation for creating post space.






**RA**

Cat. No.	183.19.204-1	183.19.204-2	183.19.204-3	183.19.204-4
Head Size	070	090	110	130
Neck Length (L)	19	19	19	19

**32mm Peeso Drills**

Pack of 6





Eagle  
Dental

Feel the Difference