

Forza

Nitrile Chemo Exam Gloves



**WE MAKE
GLOVES™**



Specifications

NITRILE CHEMO EXAM GLOVES

Material	100% Butyronitrile
Brand	V12
Name	Forza
Use	Medical, Chemo rated
USA	ASTM D6978-5
EU	EN 374-5
Sizes	S, M, L, XL
SKU	S -12214201 / M -12214202 / L -12214203 / XL -12214204
Sterility	Non-sterile
Color	Violet
Thickness	Avg Finger 5mil, Avg Palm 4mil

Description

Our Violet color, V12 Forza Nitrile Chemo Exam Gloves conforms with ASTM D6978-5 and EN 374.

Product Packaging

NITRILE CHEMO EXAM GLOVES



BOX:

100 pcs/box

CASE:

10 boxes/case

Case size: 35.3x24x26cm

Weight: S 4.8kg / M 5.3Kg / L 5.6Kg / XL 5.8Kg

PALLET:

84 Cases/pallet

NW: S 403Kg / M 445Kg / L 470Kg / XL 487Kg



Testing. Development. Problem Solving.

November 12, 2019

TEST REPORT


PN 149986A

PHARMACEUTICAL SERVICES


Prepared For:

John Zhao
Intco Medical Industries Inc.
805 Barrington Avenue
Ontario, CA 91764

Prepared By:


Tiffany Heller
Manager, Pharmaceutical Services

Approved By:


Ana C Barbur, M.S.
Vice President, Analytical & Chemical Services

Rev 101218



An A2LA ISO 17025 Accredited Testing Laboratory — Certificate Numbers 255.01 & 255.02
ISO 9001:2015 Registered

ISO 9001:2015
Registered

Letters and reports are for the exclusive use of the clients to whom they are addressed and shall not be reproduced, except in full, without the written permission of Akron Rubber Development Laboratory, Inc. (ARDL). The information contained herein applies to the specific material, products or processes tested or evaluated. No warranty of any kind is herein construed or implied. The liability of ARDL, Inc. shall be limited to the amount of consideration paid for services. ARDL, Inc. is ISO 17025 accredited by A2LA for the test methods listed on the referenced certificates.

www.ardl.com | 2887 Gilchrist Rd. | Akron, Ohio 44305 | answers@ardl.com | Toll Free (800) 830-ARDL
Fax (330) 794-6610 | Worldwide (330) 794-6600



Testing. Development. Problem Solving.

November 12, 2019

John Zhao
Intco Medical Industries Inc.

Page 2 of 5
PN 149986A

SUBJECT: Permeation testing per ASTM D 6978 on sample submitted by the above company.

RECEIVED: One (1) light blue glove type identified as; Disposable Exam Nitrile Gloves, Blue Color, Powder-Free, Non-sterile, Lot# ISCT201907.

TEST CHEMICALS:

Table 1. List of the Testing Drugs and their Sources

TESTING CHEMOTHERAPY DRUGS	DRUG SOURCE
Bortezomib (Velcade), 1 mg/ml (1,000 ppm)	LC Labs; Lot# BBZ-116; Expiration 4/2025
Carboplatin, 10 mg/ml (10,000 ppm)	TEVA; Lot# 18F06MA; Expiration 06/2020
Cyclophosphamide, 20 mg/ml (20,000 ppm)	Sandoz; Lot# 18050125; Expiration 04/2020
Etoposide (Toposar), 20 mg/ml (20,000 ppm)	Teva; Lot# 31325485B; Expiration 07/2021
Fluorouracil (5 Flu), 50 mg/ml (50,000 ppm)	MedChem; Lot# 22009; Expiration 12/2019
Gemcitabine, 38 mg/ml (38,000 ppm)	LC Labs; Lot# GMC-105; Expiration 1/6/2025
Ifosfamide, 50 mg/ml (50,000 ppm)	Baxter Healthcare; Lot# 8C015G; Expiration 03/2021
Irinotecan, 20 mg/ml (20,000 ppm)	LC Labs; Lot# RCN-105; Expiration 03/2024
Mitoxantrone, 2 mg/ml (2,000 ppm)	Sigma Aldrich; Lot# MKBR2210V; Expiration 04/2021
Oxaliplatin, 2 mg/ml (2,000 ppm)	LC Labs; Lot# XAP-111; 12/2019
Paclitaxel, 6 mg/ml (6,000 ppm)	Actavis; Lot# 8AZ5292; Expiration 05/2021
Pemetrexed, 25 mg/ml (25,000 ppm)	LC Labs; Lot# PEM-106; Expiration 04/2024
Thiotepa, 10 mg/ml (10,000 ppm)	USP; Lot# R11380; Expiration 10/2020
Vincristine Sulfate, 1 mg/ml (1,000 ppm)	USP; Lot # R046R0; Expiration 10/2019

TESTING CONDITIONS:

Standard Test Method Used:	ASTM D 6978
Deviation from Standard Test Method:	Used 1" Permeation Cell
Analytical Method:	UV/VIS Spectrometry
Testing Temperature:	35.0°C ± 2.0
Collection System:	Closed Loop
Specimen Area Exposed:	5.067 cm ²
Selected Data Points:	25/test
Number of Specimens Tested:	3/test
Location Sampled From:	Cuff area

*ARDL is ISO 17025 accredited by A2LA for the test methods listed on the certificates referenced on page one. Unless specified, the current specification version is used.
NOTE: Non-ISO 17025 accredited test methods are designated with the ^ symbol to differentiate from ISO 17025 accredited methods in the body of the test report.*

www.ardl.com | 2887 Gilchrist Rd. | Akron, Ohio 44305 | answers@ardl.com | Toll Free (800) 830-ARDL
Fax (330) 794-6610 | Worldwide (330) 794-6600



November 12, 2019

John Zhao
Intco Medical Industries Inc.

Page 3 of 5
PN 149986A

COLLECTION MEDIA:

Table 2. Collection Media for Test Chemicals

TESTING CHEMOTHERAPY DRUGS	COLLECTION MEDIUM
Bortezomib (Velcade), 1 mg/ml (1,000 ppm)	Distilled Water
Carboplatin, 10 mg/ml (10,000 ppm)	Distilled Water
Cyclophosphamide, 20 mg/ml (20,000 ppm)	Distilled Water
Etoposide (Toposar), 20 mg/ml (20,000 ppm)	Distilled Water
Fluorouracil (5 Flu), 50 mg/ml (50,000 ppm)	9.20 pH Sodium Hydroxide Solution
Gemcitabine, 38 mg/ml (38,000 ppm)	Distilled Water
Ifosfamide, 50 mg/ml (50,000 ppm)	Distilled Water
Irinotecan, 20 mg/ml (20,000 ppm)	Distilled Water
Mitoxantrone, 2 mg/ml (2,000 ppm)	Distilled Water
Oxaliplatin, 2 mg/ml (2,000 ppm)	Distilled Water
Paclitaxel, 6 mg/ml (6,000 ppm)	30% Methanol Aqueous Solution
Pemetrexed, 25 mg/ml (25,000 ppm)	Distilled Water
Thiotepa, 10 mg/ml (10,000 ppm)	Distilled Water
Vincristine Sulfate, 1 mg/ml (1,000 ppm)	Distilled Water

DETECTION METHOD OF CHEMICAL PERMEATION:

UV/VIS ABSORPTION SPECTROMETRY:

Instrument: Perkin Elmer UV/VIS Spectrometer Lambda 25

UV/VIS Absorption Spectrometry was used to measure the absorbance of test chemicals, which permeated through the specimens into the collection medium. The collection medium was circulated in a closed loop at 11 ml/minute of flow rate through the testing period. Data collection was performed according to the programmed schedule by means of UV Winlab software from the Perkin Elmer Corporation. The list of the characteristic wavelengths is shown below.

Table 3. Characteristic Wavelengths used in UV/VIS Absorption Spectrometry

TESTING CHEMOTHERAPY DRUGS	WAVELENGTH (nm)
Bortezomib (Velcade), 1 mg/ml (1,000 ppm)	206
Carboplatin, 10 mg/ml (10,000 ppm)	192
Cyclophosphamide, 20 mg/ml (20,000 ppm)	200
Etoposide (Toposar), 20 mg/ml (20,000 ppm)	205
Fluorouracil (5 Flu), 50 mg/ml (50,000 ppm)	269
Gemcitabine, 38 mg/ml (38,000 ppm)	202
Ifosfamide, 50 mg/ml (50,000 ppm)	200
Irinotecan, 20 mg/ml (20,000 ppm)	200
Mitoxantrone, 2 mg/ml (2,000 ppm)	242
Oxaliplatin, 2 mg/ml (2,000 ppm)	199
Paclitaxel, 6 mg/ml (6,000 ppm)	231
Pemetrexed, 25 mg/ml (25,000 ppm)	196
Thiotepa, 10 mg/ml (10,000 ppm)	199
Vincristine Sulfate, 1 mg/ml (1,000 ppm)	220

*ARDL is ISO 17025 accredited by A2LA for the test methods listed on the certificates referenced on page one. Unless specified, the current specification version is used.
NOTE: Non-ISO 17025 accredited test methods are designated with the ^ symbol to differentiate from ISO 17025 accredited methods in the body of the test report.*

November 12, 2019

John Zhao
Intco Medical Industries Inc.

Page 4 of 5
PN 149986A

SAMPLE CHARACTERISTICS:

Table 4. Cuff Thickness characteristics for the tested specimens: Light blue glove identified as: Disposable Exam Nitrile Gloves, Blue Color, Powder-Free, Non-sterile, Lot# ISCT201907.

Testing Chemotherapy Drugs	Thickness (mm)			Average (mm)
	Sample 1	Sample 2	Sample 3	
Bortezomib (Velcade)	0.059	0.062	0.060	0.060
Carboplatin	0.058	0.057	0.061	0.059
Cyclophosphamide	0.070	0.073	0.060	0.068
Etoposide (Toposar)	0.058	0.067	0.057	0.061
Fluorouracil (5 Flu)	0.059	0.062	0.057	0.060
Gemcitabine	0.068	0.060	0.059	0.063
Ifosfamide	0.070	0.061	0.059	0.064
Irinotecan	0.054	0.060	0.056	0.057
Mitoxantrone	0.055	0.060	0.056	0.057
Oxaliplatin	0.056	0.069	0.060	0.062
Paclitaxel	0.056	0.059	0.062	0.059
Pemetrexed	0.069	0.056	0.059	0.061
Thiotepa	0.063	0.063	0.060	0.062
Vincristine Sulfate	0.066	0.060	0.057	0.061
Weight/Unit Area (g/m2)	64.0			

*ARDL is ISO 17025 accredited by A2LA for the test methods listed on the certificates referenced on page one. Unless specified, the current specification version is used.
NOTE: Non-ISO 17025 accredited test methods are designated with the ^ symbol to differentiate from ISO 17025 accredited methods in the body of the test report. *



November 12, 2019


John Zhao
Intco Medical Industries Inc.


Page 5 of 5
PN 149986A

RESULTS:

Table 5. Permeation Test Results on: Light blue glove identified as; Disposable Exam Nitrile Gloves, Blue Color, Powder-Free, Non-sterile, Lot# ISCT201907.

TEST CHEMOTHERAPY DRUG AND CONCENTRATION	MINIMUM BREAKTHROUGH DETECTION TIME (Sample 1,2,3) (Minutes)	STEADY STATE PERM. RATE (Sample 1,2,3) (µg/cm ² /minute)	OTHER OBSERVATIONS
Bortezomib (Velcade), 1 mg/ml (1,000 ppm)	>240	N/A	Slight swelling and no degradation
Carboplatin, 10 mg/ml (10,000 ppm)	>240	N/A	Slight swelling and no degradation
Cyclophosphamide, 20 mg/ml (20,000 ppm)	>240	N/A	Slight swelling and no degradation
Etoposide (Toposar), 20 mg/ml (20,000 ppm)	>240	N/A	Slight swelling and no degradation
Fluorouracil (5 Flu), 50 mg/ml (50,000 ppm)	>240	N/A	Slight swelling and no degradation
Gemcitabine, 38 mg/ml (38,000 ppm)	>240	N/A	Slight swelling and no degradation
Ifosfamide, 50 mg/ml (50,000 ppm)	>240	N/A	Slight swelling and no degradation
Irinotecan, 20 mg/ml (20,000 ppm)	>240	N/A	Slight swelling and no degradation
Mitoxantrone, 2 mg/ml (2,000 ppm)	>240	N/A	Slight swelling and no degradation
Oxaliplatin, 2 mg/ml (2,000 ppm)	>240	N/A	Slight swelling and no degradation
Paclitaxel, 6 mg/ml (6,000 ppm)	>240	N/A	Slight swelling and no degradation
Pemetrexed, 25 mg/ml (25,000 ppm)	>240	N/A	Slight swelling and no degradation
Thiotepa, 10 mg/ml (10,000 ppm)	125.8 (128.0,134.6,125.8)	0.6 (0.6,0.5,0.6)	Slight swelling and no degradation
Vincristine Sulfate, 1 mg/ml (1,000 ppm)	>240	N/A	Slight swelling and no degradation

Prepared By: 
Tiffany Heller
Manager, Pharmaceutical Services

Approved By: 
Ana C Barbur, M.S.
Vice President, Analytical & Chemical Services

*ARL is ISO 17025 accredited by A2LA for the test methods listed on the certificates referenced on page one. Unless specified, the current specification version is used.
NOTE: Non-ISO 17025 accredited test methods are designated with the ^ symbol to differentiate from ISO 17025 accredited methods in the body of the test report.*



Testing. Development. Problem Solving.

October 30, 2019

TEST REPORT

PN 149986B

PHARMACEUTICAL SERVICES

Prepared By:

Tiffany Heller
Manager, Pharmaceutical Services

Approved By:

Ana C Barbur, M.S.
Vice President, Analytical & Chemical Services

Rev 101218



An A2LA ISO 17025 Accredited Testing Laboratory — Certificate Numbers 255.01 & 255.02
ISO 9001:2015 Registered

ISO 9001:2015
Registered

Letters and reports are for the exclusive use of the clients to whom they are addressed and shall not be reproduced, except in full, without the written permission of Akron Rubber Development Laboratory, Inc. (ARDL). The information contained herein applies to the specific material, products or processes tested or evaluated. No warranty of any kind is herein construed or implied. The liability of ARDL, Inc. shall be limited to the amount of consideration paid for services. ARDL, Inc. is ISO 17025 accredited by A2LA for the test methods listed on the referenced certificates.

www.ardl.com | 2887 Gilchrist Rd. | Akron, Ohio 44305 | answers@ardl.com | Toll Free (800) 830-ARDL
Fax (330) 794-6610 | Worldwide (330) 794-6600

October 30, 2019

Page 2 of 5
PN 149986B

SUBJECT: Permeation testing per ASTM D 6978 on gloves submitted by the above company.

RECEIVED: One (1) glove type identified as; Disposable Exam Nitrile Gloves, Blue Color, Powder-Free, Non-sterile, Lot# ISCT201907.

TEST CHEMICALS:

Table 1.1 List of the Testing Drugs and their Sources

TESTING CHEMOTHERAPY DRUGS	DRUG SOURCE
Bleomycin Sulfate, 15 mg/ml (15,000 ppm)	TEVA; Lot# 31321906B; Expiration 09/2019 (Tested prior to exp.)
Busulfan, 6 mg/ml (6,000 ppm)	Sigma; Lot# BCBS8240V; Expiration 12/2021
Cetuximab, 2 mg/ml (2,000 ppm)	Lilly; Lot# C1700079; Expiration 04/2020
Cisplatin, 1 mg/ml (1,000 ppm)	Accord; Lot# PY02454; Expiration 04/2021
Cyclosporin A, 100 mg/ml (100,000 ppm)	USP; Lot# J0M382; Expiration 11/2019
Cytarabine, 100 mg/ml (100,000 ppm)	Sigma Aldrich; Lot# LRAB3688; Expiration 12/2021
Dacarbazine, 10 mg/ml (10,000 ppm)	Teva; Lot# 31325414B; Expiration 09/2021
Daunorubicin, 5 mg/ml (5,000 ppm)	Sigma Aldrich; Lot# 125M4750V; Expiration 03/2020
Docetaxel, 10 mg/ml (10,000 ppm)	Sigma Aldrich; Lot# LRAB0750; Expiration 12/2021
Doxorubicin HCl, 2 mg/ml (2,000 ppm)	WestWard; Lot# BJ0044; Expiration 10/2020
Epirubicin HCl (Elevance), 2 mg/ml (2,000 ppm)	Actavis; Lot 7U15152; Expiration 10/2020
Fludarabine, 25 mg/ml (25,000 ppm)	USP; Lot# H1K220; Expiration 12/2019
Idarubicin, 1 mg/ml (1,000 ppm)	Sigma Aldrich; Lot# R080E0; Expiration 12/2019
Melphalan, 5 mg/ml (5,000 ppm)	USP; Lot# R086P0; Expiration 04/2021
Methotrexate, 25 mg/ml (25,000 ppm)	Mylan; Lot# 7801774; Expiration 04/2020
Mitomycin C, 0.5 mg/ml (500 ppm)	Sigma Aldrich; Lot# MKCD6056; Expiration 03/2020
Paraplatin (Carboplatin), 10 mg/ml (10,000 ppm)	TEVA; Lot# 18F06MA; Expiration 06/2020
Retrovir, 10 mg/ml (10,000 ppm)	GlaxoSmithKline; Lot# C819185; Expiration 04/2020
Rituximab, 10 mg/ml (10,000 ppm)	Hetero Healthcare; Batch# RB1921B; Expiration 08/2021
Topotecan, 1 mg/ml (1,000 ppm)	USP; Lot# R093L0; Expiration 11/2020
Trisenox, 1 mg/ml (1,000 ppm)	Sigma Aldrich; Lot# 129K0039V; Expiration 12/2020

TESTING CONDITIONS:

Standard Test Method Used:	ASTM D 6978
Deviation from Standard Test Method:	Used 1" Permeation Cell
Analytical Method:	UV/VIS Spectrometry
Testing Temperature:	35.0°C ± 2.0
Collection System:	Closed Loop
Specimen Area Exposed:	5.067 cm ²
Selected Data Points:	25/test
Number of Specimens Tested:	3/test
Location Sampled From:	Cuff Area

*ARDL is ISO 17025 accredited by A2LA for the test methods listed on the certificates referenced on page one. Unless specified, the current specification version is used.
NOTE: Non-ISO 17025 accredited test methods are designated with the ^ symbol to differentiate from ISO 17025 accredited methods in the body of the test report.*

October 30, 2019

Page 3 of 5
PN 149986B

COLLECTION MEDIA:

Table 2 Collection Media for Test Chemotherapy Drugs

TEST DRUG AND CONCENTRATION	COLLECTION MEDIUM
Bleomycin Sulfate, 15 mg/ml (15,000 ppm)	Distilled Water
Busulfan, 6 mg/ml (6,000 ppm)	Distilled Water
Cetuximab, 2 mg/ml (2,000 ppm)	Distilled Water
Cisplatin, 1 mg/ml (1,000 ppm)	Distilled Water
Cyclosporin A, 100 mg/ml (100,000 ppm)	Distilled Water
Cytarabine, 100 mg/ml (100,000 ppm)	Distilled Water
Dacarbazine, 10 mg/ml (10,000 ppm)	Distilled Water
Daunorubicin, 5 mg/ml (5,000 ppm)	Distilled Water
Docetaxel, 10 mg/ml (10,000 ppm)	Distilled Water
Doxorubicin HCl, 2 mg/ml (2,000 ppm)	Distilled Water
Epirubicin HCl (Ellence), 2 mg/ml (2,000 ppm)	Distilled Water
Fludarabine, 25 mg/ml (25,000 ppm)	Distilled Water
Idarubicin, 1 mg/ml (1,000 ppm)	Distilled Water
Melphalan, 5 mg/ml (5,000 ppm)	Distilled Water
Methotrexate, 25 mg/ml (25,000 ppm)	Distilled Water
Mitomycin C, 0.5 mg/ml (500 ppm)	Distilled Water
Paraplatin (Carboplatin), 10 mg/ml (10,000 ppm)	Distilled Water
Retrovir, 10 mg/ml (10,000 ppm)	Distilled Water
Rituximab, 10 mg/ml (10,000 ppm)	Distilled Water
Topotecan, 1 mg/ml (1,000 ppm)	Distilled Water
Trisenox, 1 mg/ml (1,000 ppm)	Distilled Water

DETECTION METHOD OF CHEMICAL PERMEATION:

UV/VIS ABSORPTION SPECTROMETRY:

Instrument: Perkin Elmer UV/VIS Spectrometer Lambda 25

UV/VIS Absorption Spectrometry was used to measure the absorbance of test chemicals, which permeated through the specimens into the collection medium. The collection medium was circulated in a closed loop at 11 ml/minute of flow rate through the testing period. Data collection was performed according to the programmed schedule by means of UV Winlab software from the Perkin Elmer Corporation. The list of the characteristic wavelengths is shown below.

Table 3 Characteristic Wavelengths used in UV/VIS Absorption Spectrometry

TESTING CHEMOTHERAPY DRUGS	WAVELENGTH (nm)
Bleomycin Sulfate, 15 mg/ml (15,000 ppm)	290
Busulfan, 6 mg/ml (6,000 ppm)	197
Cetuximab, 2 mg/ml (2,000 ppm)	199
Cisplatin, 1 mg/ml (1,000 ppm)	199
Cyclosporin A, 100 mg/ml (100,000 ppm)	199
Cytarabine, 100 mg/ml (100,000 ppm)	272
Dacarbazine, 10 mg/ml (10,000 ppm)	320
Daunorubicin, 5 mg/ml (5,000 ppm)	269
Docetaxel, 10 mg/ml (10,000 ppm)	231
Doxorubicin HCl, 2 mg/ml (2,000 ppm)	232
Epirubicin HCl (Ellence), 2 mg/ml (2,000 ppm)	233 & 253
Fludarabine, 25 mg/ml (25,000 ppm)	261
Idarubicin, 1 mg/ml (1,000 ppm)	257
Melphalan, 5 mg/ml (5,000 ppm)	260
Methotrexate, 25 mg/ml (25,000 ppm)	303
Mitomycin C, 0.5 mg/ml (500 ppm)	217
Paraplatin, 10 mg/ml (10,000 ppm)	192
Retrovir, 10 mg/ml (10,000 ppm)	266
Rituximab, 10 mg/ml (10,000 ppm)	192
Topotecan, 1 mg/ml (1,000 ppm)	254
Trisenox, 1 mg/ml (1,000 ppm)	197

*ARDL is ISO 17025 accredited by A2LA for the test methods listed on the certificates referenced on page one. Unless specified, the current specification version is used.
NOTE: Non-ISO 17025 accredited test methods are designated with the ^ symbol to differentiate from ISO 17025 accredited methods in the body of the test report.*



October 30, 2019

Page 4 of 5
PN 149986B

SAMPLE CHARACTERISTICS:

Table 4. Cuff thickness for glove identified as: Disposable Exam Nitrile Gloves, Blue Color, Powder-Free, Non-sterile, Lot# ISCT201907.

Testing Chemotherapy Drugs	Thickness (mm)			Average (mm)
	Sample 1	Sample 2	Sample 3	
Bleomycin Sulfate	0.068	0.062	0.066	0.065
Busulfan	0.066	0.058	0.069	0.064
Cetuximab	0.060	0.069	0.065	0.065
Cisplatin	0.063	0.073	0.062	0.066
Cyclosporin A	0.065	0.064	0.066	0.065
Cytarabine	0.065	0.065	0.063	0.064
Dacarbazine	0.067	0.064	0.068	0.066
Daunorubicin	0.067	0.061	0.068	0.065
Docetaxel	0.065	0.063	0.068	0.065
Doxorubicin HCl	0.061	0.066	0.070	0.066
Epirubicin HCl (Ellence)	0.068	0.065	0.075	0.069
Fludarabine	0.072	0.074	0.063	0.070
Idarubicin	0.068	0.065	0.073	0.068
Melphalan	0.056	0.061	0.065	0.060
Methotrexate	0.064	0.062	0.081	0.069
Mitomycin C	0.060	0.071	0.068	0.066
Paraplatin	0.056	0.068	0.064	0.063
Retrovir	0.060	0.064	0.064	0.063
Rituximab	0.073	0.063	0.075	0.070
Topotecan	0.073	0.066	0.063	0.067
Trisenox	0.074	0.068	0.063	0.069
Weight/Unit Area (g/m²)	64.4			

*ARL is ISO 17025 accredited by A2LA for the test methods listed on the certificates referenced on page one. Unless specified, the current specification version is used.
NOTE: Non-ISO 17025 accredited test methods are designated with the ^ symbol to differentiate from ISO 17025 accredited methods in the body of the test report.*


October 30, 2019


Page 5 of 5
PN 149986B

RESULTS:

Table 5 Permeation Test Results on Testing of: Disposable Exam Nitrile Gloves, Blue Color, Powder-Free, Non-sterile, Lot# ISCT201907.

TEST CHEMOTHERAPY DRUG AND CONCENTRATION	MINIMUM BREAKTHROUGH DETECTION TIME (Specimen 1/2/3) (Minutes)	STEADY STATE PERM. RATE (Specimen 1/2/3) ($\mu\text{g}/\text{cm}^2/\text{minute}$)	OTHER OBSERVATIONS
Bleomycin Sulfate, 15 mg/ml (15,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Busulfan, 6 mg/ml (6,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Cetuximab, 2 mg/ml (2,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Cisplatin, 1 mg/ml (1,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Cyclosporin A, 100 mg/ml (100,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Cytarabine, 100 mg/ml (100,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Dacarbazine, 10 mg/ml (10,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Daunorubicin, 5 mg/ml (5,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Docetaxel, 10 mg/ml (10,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Doxorubicin HCl, 2 mg/ml (2,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Epirubicin HCl (Ellence), 2 mg/ml (2,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Fludarabine, 25 mg/ml (25,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Idarubicin, 1 mg/ml (1,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Melphalan, 5 mg/ml (5,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Methotrexate, 25 mg/ml (25,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Mitomycin C, 0.5 mg/ml (500 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Paraplatin, 10 mg/ml (10,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Retrovir, 10 mg/ml (10,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Rituximab, 10 mg/ml (10,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Topotecan, 1 mg/ml (1,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation
Trisenox, 1 mg/ml (1,000 ppm)	No breakthrough up to 240 min.	N/A	Slight swelling and no degradation

Prepared By: 
Tiffany Heller
Manager, Pharmaceutical Services

Approved By: 
Ana C Barbur, M.S.
Vice President, Analytical & Chemical Services

*ARDI is ISO 17025 accredited by A2LA for the test methods listed on the certificates referenced on page one. Unless specified, the current specification version is used.



Certifications

FDA U.S. FOOD & DRUG ADMINISTRATION Follow FDA | En Español **SEARCH**

[Home](#) [Food](#) [Drugs](#) [Medical Devices](#) [Radiation-Emitting Products](#) [Vaccines, Blood & Biologics](#) [Animal & Veterinary](#) [Cosmetics](#) [Tobacco Products](#)

Establishment Registration & Device Listing

[FDA Home](#) [Medical Devices](#) [Databases](#)

New Search Back To Search Results

Establishment:
V12 Health LLC
Business Trade Name:
V12 Health LLC

11400 NW 34 TH STREET
Doral , FL 33178
Registration Number: 3016842185
FEI Number*: 3016842185
Status: Active
Initial Distributor/Importer: Yes
*Note Firm May Have Additional Establishment Types.
Please Review Listings For Further Information.
Date Of Registration Status: 2021

Owner/Operator:
V12 Health LLC
11400 NW 34 TH STREET
Doral, FL US 33178
Owner/Operator Number: [10073104](#)

Official Correspondent:
Tarek Kirschen
11400 NW 34 TH STREET
Doral, FL 33178
Phone: 1-305-8909900

* Firm Establishment Identifier (FEI) should be used for identification of entities within the imports message set

Establishment Registration & Device Listing

[FDA Home](#) [Medical Devices](#) [Databases](#) Results per Page [New Search](#)

1 to 3 of 3 Results for **Owner Operator**
Number : [10073104](#)

Establishment Name	Registration Number	Current Registration Yr
V12 Health LLC FLUSA	3016842185	2021
<ul style="list-style-type: none"> Medical Gloves With Chemotherapy Labeling Claims - Test For Use With Chemotherapy Drugs - FORZA Specification Developer Polymer Patient Examination Glove - PISTA Specification Developer Latex Patient Examination Glove - Nero Specification Developer 		
V12 Health LLC TURKEY	3018501299	2021
<ul style="list-style-type: none"> Non-Surgical Isolation Gown - V12 Isolation Gown Manufacturer Gown, Examination - V12 Surgical Gowns Manufacturer Cap, Surgical - V12 Bouffant Cap Manufacturer Cover, Shoe, Operating-Room - V12; V12 Shoe Cover Manufacturer 		
V12 Health LLC CHINA	No number listed	2021
<ul style="list-style-type: none"> Latex Patient Examination Glove - Nero Manufacturer Medical Gloves With Chemotherapy Labeling Claims - Test For Use With Chemotherapy Drugs - FORZA Manufacturer Polymer Patient Examination Glove - PISTA Manufacturer Patient Examination Glove - CORSA Manufacturer 		

Can't find what you're looking for? [Try a new search](#)



WE MAKE GL♥VES

V12 Health, LLC

33 SW 2nd Ave, PH2
Miami, FL 33130
U.S.A.

Federal Tax ID # 85-0688159 | FDA registration # 3016842185

"We Make Gloves" is a trademark of V12 Health, LLC

"Everyone Is Essential" is a trademark of V12 Health, LLC

+1 (305) 507 5007 | v12health.com | info@v12health.com