4G health watch

EN

1. Package Contents

Watch	Charging Cable	Quick Start Manual
		

2. How to Operate

2.1 Powering on Your Watch

Long press the power key to turn on your watch . Please make sure to fully charge the watch first if it Doesn' t turn on.

2.2 Charging your watch

Properly attach the magnetic charging base to the back of the watch,align the metal contacts on the back of the watch with the gold pin on the charging base while placing, then plug the USB end into a USB port on computer, charging dock or power bank for charging.

Input current: <0.6A

Input voltage: DC 5V

Charging time: about 2 hours

LED indicator

Charging: Purple indicator light flashes continuously.


Fully charged: Purple indicator light stays on constantly.

Notes:

- ① It is recommended to use the universal charger with certification mark on the market.
- ② Please keep the charging port dry before charging.
- ③ When the watch voice broadcasts "Low battery please charge," it indicates that the watch battery is critically low. To prevent it from shutting down, please recharge the watch promptly.

2.3 How to insert or remove the SIM card?

- (1) Turn off the watch, then use your fingernail to rotate and open the SIM card cover.
- (2) When inserting the SIM card, please follow the markings on the SIM card (chip facing up, notch facing inward) and gently push the SIM card into the slot.
- (3) Rotate the SIM card cover back into place and ensure it is securely closed.


Notes:

- ① The slot is designed to prevent reverse insertion to a certain extent, but it does not prevent oblique insertion, so please make sure to align the direction when inserting the card.
- ② The watch adopts a pop-up card slot. When removing the card, simply push the SIM card lightly into the slot, then release it, and it will pop out. Do not force it out.

2.4 Wearing Your Watch

- (1) For optimized tracking of the measured values, we recommend wearing the watch with one finger's width below your wrist bone.
- (2) Please make sure the watch is fairly snug against your skin and does not slide up or down your wrist during exercise.


2.5 How to Operate

	Long Press	Short Press	Power Key 
Power Key	● Power on/off ● Voice prompt	/	
	Emergency SOS	Voice prompt	
	/	● End the call ● Voice prompt ● Health measurement	
	Make a call to a contact	● Answer the call ● Voice prompt	

2.6 Replacing the Strap

Please choose the strap with a width of 27mm if you want to replace it.

- (1) Remove the strap from the watch by sliding the snap lock on the strap.
- (2) Align the new strap with the watch and buckle the strap in.

2.7 APP

(1) Downloading

IOS users can download the app from the App Store, Android users can download the 'CareFit' app from the Google Play Store, or scan the QR code below for download.

(2) Pairing

- 1) Open APP, register an account, and truthfully fill in your personal information (height, weight, date of birth);
- 2) Click "Device" > "4G Watch"
- 3) Scan the QR code attached to the packaging box to bind the watch.

Notes:

- ① To ensure successful data transfer to the App, please insert a functional SIM card into the smartwatch first (a SIM card is required for time synchronization).
- ② The default language of the watch is English, and you can switch languages through the app.

(3) Firmware Upgrade

When there is an new firmware version, there will be a notification in the app. Go to the "Device" interface of the app and select firmware upgrade.

3. Main Feature

3.1 Health Measurement

To wear the watch correctly, briefly press the watch's voice button and follow the voice prompts to enter health measurement mode (the watch will sequentially measure heart rate, blood oxygen, blood pressure, stress level, and body temperature). During the measurement process, the watch's indicator light will remain blue. After the measurement is complete, follow the voice prompts to sync the data with the app and view the measurement results.

Notes:

For the first blood pressure measurement, please calibrate the watch with a home blood pressure monitor.

3.2 Heart Rate Monitoring

- (1) Built-in PPG heart rate sensor coupled with HR algorithm, the watch can accurately monitor your heart rate during wearing.
- (2) Heart Rate Alert: Supports setting heart rate alert range through the app. Once set successfully, the watch will emit a voice prompt when the detected heart rate value is outside the set range.

Notes:

The light signal transmission might be blocked if your skin is too dark or with excess hair, or improper wearing also may lead to measurement failure.

3.3 Heart Health and Scientific Sleep

- (1) Open the app, tap "Device" > "24h Heart Health", and toggle the feature on or off as needed.
- (2) When this function is turned on, the watch will automatically monitor the heart health at night. If you wear the watch for more than 6 hours, you can view a scatter chart of your heart health after synchronizing your data with the watch on the App. The measurement results are for reference only, and can not be used as a basis for clinical diagnosis.
- (3) When this function is turned on, if you wear the watch to sleep between 8:00pm and 9:00 am the next morning, you can check REM sleep and other sleep data after synchronizing your data with the watch on the APP, so as to know about your sleep status.

Notes:

- ① When this function is turned on, the heart rate indicator light will remain on during sleep, which may shorten the device's battery life.

- ② The sleep monitor will stop recording within 5-10 minutes after you wake up to move.
- ③ The watch will not record daytime sleep data.

3.4 Voice Communication

The watch supports making and receiving calls independently with a SIM card. Press the call button on the watch, and follow the voice prompts to select the contact you want to call.

Long-pressing the call button directly will default to dialing the first contact number.

Note:

The contacts of the watch address book need to be set in the APP. Open the App> "Device"> "Contacts"> "+">, add contacts. (A maximum of 8 contacts can be added.)

4. More Information

4.1 Water Resistance Instruction

The water resistance rating: **IP67**

The water resistance performance of the device is not permanently valid, it may diminish as time passes. The device can be used during washing hands, rain or swimming in the shallow water, but not support hot water shower, diving, sauna, etc. Do not let the product come into contact with sea water.

4.2 Device Care

- (1) Do not press the key on the device while it is underwater.
- (2) Do not drop or impact the watch.
- (3) Do not use any sharp objects to clean the device.
- (4) Avoid exposing the device to extreme temperatures that are either too cold or too hot for an extended period, as this may cause permanent damage.
- (5) Prevent the product from coming into contact with seawater.

- Gently wipe the device using a cotton cloth with a mild neutral detergent.
- After each training session, please rinse the watch with clean water.

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