MSLA-Modeling – Gray (MMD-Roo1GY)

Check

- Vat and LCD screen are clean
- Shake well before pouring, 60 sec
- Vat & platform are tighten properly
- Do Not overfill vat

Setting recommendation Chart MSLA Dental & MSLA Gray

Recommended Setting Parameters


				Exposure
		Normal	Bottom	Time (sec)
	Layer	Exposure	(Burn-In)	for Bottom
Printer	Thickness	Time (sec)	Layers	(Burn-In)
Photon	50 μm	10	8	45
Photon	100 μm	15	10	50
Photon S	50 μm	5	8	40
Photon S	100 μm	8	10	45
Shuffle	50 μm	10	8	40
Shuffle	100 μm	15	10	45
Shuffle XL	50 μm	8	8	40
Shuffle XL	100 μm	12	10	45
Sonic	50 μm	2.5	8	30
Sonic	100 μm	4	10	30
Sonic Mini	50 μm	4	6	30
Sonic Mini	100 μm	6	6	30

Each printer is unique, light intensity varies from printer to printer eventhough they are the same model and make. Please adjust normol exposure time from our recommended setting ± 3 (± 0.5 for sonic) second for your printer and print geometry.


For size calibration, download the stl file: https://www.thingiverse.com/thing:1586206

Setting Example:

AnyCubic Photon:



For Phrozen Shuffle:





Washing

Submerge print in clean IPA for up to 5 minutes, gently shake-off / blow-off excess IPA, set print in a shaded airy place to dry before curing. (extended period in IPA causes deformation)

Dried print might be a little tacky to touch

Curing

Post-curing required for optimal material performance: Light wavelength: 405 nm Total light energy per unit area required on the print: 162,000 mJ/cm²

Example:

Light power on the print = 40mW/cm^2 , (higher power per unit area shortens curing time) 162,000 mJ / 30 mW / 60 sec = 90 minutes of curing time.

Controlled curing chamber temperature between: 30~40°C / 86~104°F

Storage

Don't keep resin in the vat for more than **2 days**. Filter out all debris in the case of fail print.

Note: Slight pigment settlement is normal. Gentally mix remaining resin in resin vat with soft wiper for color evenness.