


# WRIST ELECTRONIC BLOOD PRESSURE MONITOR **procare**

Instruction Manual  
MODEL: 240389



CE 0413  
FDA Approved

Distributed by: Arise Medical LLC  
Address: 654 Westminster Road,  
Wilkes-Barre, PA 18702, USA  
Customer Service: Tel: +1-866-277-7168  
Website: www.arisemedical.com

## TABLE OF CONTENTS

|  |    |
|--|----|
| INTRODUCTION.....                              | 2  |
| NOTES ON SAFETY.....                           | 2  |
| ABOUT BLOOD PRESSURE.....                      | 6  |
| PRECAUTIONS BEFORE USE.....                    | 8  |
| FEATURES OF THE PRODUCT.....                   | 9  |
| PARTS IDENTIFICATION.....                      | 10 |
| INSERT OR REPLACE BATTERIES.....               | 11 |
| TIME AND OF SYSTEM SETUP.....                  | 12 |
| UNIT CONVERSION mmHg/kPa DISPLAY.....          | 13 |
| WHO BLOOD PRESSURE CLASSIFICATION DISPLAY..... | 14 |
| ATTACHING THE ARM CUFF.....                    | 15 |
| HOW TO MEASURE BLOOD PRESSURE.....             | 16 |
| CARE AND MAINTENANCE.....                      | 17 |
| SPECIFICATIONS.....                            | 19 |
| TROUBLESHOOTING.....                           | 20 |

## INTRODUCTION

The Monitor uses the oscillometric method of blood pressure measurement. Measurement Automatic Electronic Blood Pressure Monitor is intended for use by medical professionals or at home to monitor and display diastolic, systolic blood pressure and pulse rate, with an air wrist cuff buckled around one's wrist according to the instructions in the "ATTACHING THE WRIST CUFF". The expected life of the product is 5 years.

The product complies with the electromagnetic compatibility requirement of EN60601-1-2 and safety standards of IEC60601-1 and performance of IEC 80601-2-30 as specified in EEC directive 93/42/EEC.

## NOTES ON SAFETY

- \* The warning signs and sample icons shown here are listed for your safe and correct use of the unit, so as to prevent injuries or damages to the device.
- \* The icons and meanings are as follow.

### Examples of signs

- The icon indicates prohibitions (what you should not do). Matters involving actual prohibitions are indicated by text or pictures in or near . The left icon refers to "general prohibition".

- The icon indicates something that is compulsory (what must always be observed). Matters involving actual compulsory actions are indicated by text or pictures in or near . The left icon refers to "general compulsion".
- The icon indicates something that is disassembled or "Don't disassemble". Matters involving actual compulsory actions are indicated by text or pictures in or near . The left icon refers to "general prohibition".
- Type BF Applied part. The following symbol indicates that the device is MR-unsafe.
- Please refer to the instructions for use.
- Indicates a medical device that needs to be protected from moisture.
- Marking of electrical and electronic equipment in accordance with Article 11(2) of Directive 2002/96/EC (WEEE).



- Patient must follow doctor's instruction and should not perform self-judgment and self-treatment by the measuring result. Self-diagnosis of measured results and treatment are dangerous. The device should not be used to judge illness, first aid and continuously monitor measuring.
- This device can not be used for Patient transport and surgical care. It can be used in household or fixed places only. Please press "on/off" button to stop work when you feel uncomfortable with the wrist, or if the air is inflating abnormally without stop.
- Do not let a child below 12 years old and the people who can't express one's intention. When it is used by the people of 12~18 years old, it should accompany by the Adult. May cause accident or trouble.
- Do not use the unit for purpose other than measuring blood pressure. May cause accident or trouble.
- Please do not use mobile phone around the device. Please do not use the device around the magnetic field.
- The device is prohibited from being used during movement.
- Do not use the equipment in outdoor or shower rooms.
- Do not disassemble, repair, or remodel the main unit or the wrist cuff of the blood pressure monitor. Will cause the unit to function erroneously.

## Requests from Manufacturer

- Make sure there is no connection tubing kinking before start measuring to avoid any injury to patient.
- For any patient, do not measure more than 3 times continuously, it should be at least above 5 minutes of interval rest between any two measurements, otherwise will cause extravasated blood.
- Do not measure your blood pressure over 6 times each day.
- Do not apply the cuff over a wound as this can cause further injury.
- Do not measure on the wrist which is on the side of a mastectomy, otherwise it could cause injury.
- Observe the air pressure value from the LCD display.
- When measuring, it could not exceed 280 mmHg, otherwise Please press "on/off" button to stop
- Do not use force to bend the wrist cuff or the air tube.
- Do not knock or drop the main unit.
- Always use the specified accessories in the manual. The use of other parts not approved by the manufacturer may cause faults or injuries.
- For service information, parts list etc., please contact the dealer.


- The PATIENT is an intended OPERATOR.
- Not servicing and maintenance while the ME EQUIPMENT is in use.
- The user can maintain the product, the maintenance method is described in the maintenance instructions of manual.
- Stop using the equipment immediately, if it is in contact with water.

## ABOUT BLOOD PRESSURE


1. What is blood pressure?  
Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts. Diastolic pressure occurs when the heart expands.  
Blood pressure is measured in millimeters of mercury (mmHg). One's natural blood pressure is represented by the fundamental pressure, which is measured first thing in the morning while one is still at rest and before eating.

2. What is hypertension and how is it controlled?  
Hypertension, an abnormally high arterial blood pressure, if left unattended, can cause many health problems including stroke and heart attack. Hypertension can be controlled by altering lifestyle, avoiding stress and with medication under a doctor's supervision.  
To prevent hypertension or keep it under control:

- Do not smoke
  - Exercise regularly
  - Reduce salt and fat intake
  - Have regular physical checkups
  - Maintain proper weight
3. Why measure blood pressure at home?  
Blood pressure measured at a clinic or doctor's office may cause apprehension and produce an elevated reading, 25 to 30 mmHg higher than that measured at home. Home measurement reduces the effects of outside influences on blood pressure readings, supplements the doctor's readings and provides a more accurate, complete blood pressure history.
4. WHO blood pressure classification  
Standards for assessment of high blood pressure, without regard to age, have been established by the World Health Organization (WHO), and is shown in chart below.
5. Blood pressure variations  
An individual's blood pressure varies greatly on a daily and seasonal basis. It may vary by 30 to 50 mmHg due to various conditions during the day. In




hypertensive individuals, variations are even more pronounced.  
Normally, the blood pressure rises while at work or play and falls to its lowest levels during sleep. So, do not be overly concerned by the results of one measurement.  
Take measurements at the same time every day using the procedure described in this manual, and know your normal blood pressure. Many readings give a more comprehensive blood pressure history.  
Be sure to note date and time when recording your blood pressure. Consult your doctor to interpret your blood pressure data.



- For people with irregular or unstable peripheral circulation problems due to diabetes, liver disease, hardening of the arteries, etc., there may be fluctuations in blood pressure values measured at the upper arm versus at the wrist.
- Measurements may be impaired if this device is used near televisions, microwave ovens, X-ray, mobile phone equipment or other devices with strong electrical fields. To prevent such interference, use the monitor at a sufficient distance from such devices or turn them off.
- Before using, should wash your hands.
- Do not measure on the arm which simultaneously used monitoring ME Equipment, otherwise it could cause loss of function.
- Consult your doctor if the unexpected readings are obtained, also please refer to "Trouble shooting" of the manual.
- The reading is probably a little lower than measured in the hospital due to the steady mood at home.
- Cuff pressure range 0-299mmHg


- WHO blood pressure classification display.
- Easy to use, Press a button to automatically measure, record the measurement values and measurement time.
- Automatically turns off (within 1 minute) to save power.

## PARTS IDENTIFICATION




- ### FEATURES OF THE PRODUCT
- Memory can store 90 measurements.
  - Large and clear LCD display.

- Press "SET" key to turn on.
- Press and hold "SET" key until the year number displays and flashes on LCD to enter setting mode.




- Press "MEM" key to adjust the year, then press "SET" key again to save your setting and enter the month setting mode.
- Press "MEM" key to adjust the month. Following the same steps to adjust date/hour/minute.




## UNIT CONVERSION mmHg/kPa DISPLAY

The goods have mm Hg(mmHg), kPa (kPa) two kinds of blood pressure display units(mmHg factory to express).  
Press "ON / OFF" button for 10 seconds to display unit switching interface, then press "MEM" key to select mmHg / kPa, press "ON / OFF" button to exit.


1. Fasten the wrist cuff
- 1) Wrap the wrist cuff around your wrist about (1-2)cm above your hand as shown in the figure at the right.
- 2) Fasten the wrist cuff tightly by using the Velcro Strip. For proper measurements, fasten the wrist cuff tightly and measure on a bare wrist.



2. How to take proper measurements
- For best accuracy in blood pressure measurement:  
Sit comfortably at a table. Rest your wrist on the table.  
Relax for about 5 to 10 minutes before measurement.  
Raise your hand so that the wrist cuff is at the same level as your heart.  
Remain still and keep quiet during measurement.  
Do not measure right after physical exercise or a bath.  
Measure your blood pressure at about the same time every day.



1. Fasten the wrist cuff according to the instructions in "ATTACHING THE WRIST CUFF".
2. Press the "ON/OFF" button. All icons appear two seconds on DISPLAY, then switch to measurement, and display "0" for last measurement record.



3. Start measurement, the cuff in the strap will automatically inflate. The mark will flash on LCD. When complete, the results will be displayed.



## READ MEMOR Y

Press "MEM" button, a memory reading out the latest measurements, "MEM" for the buttons(UP), "SET" button for the memory (DOWN)  
Power Measurement closure or after the end of the state, can press the "MEM" button read out the latest measurement of memory.

## DELETE MEMOR Y


The state read out the memory press the (memory) button five seconds, the LCD display "0" has been to delete all memory.



## CARE AND MAINTENANCE

To keep your digital blood pressure monitor in the best condition and protect the unit from damage, follow the directions listed below.  
Keep the monitor in the storage case when not in use.  
Do not fold the arm cuff too tightly.  
The fabric fastener could touch the inner surface of the arm cuff and damage it.

Clean the monitor and cuff with a soft dry cloth.  
Do not use any abrasive or volatile cleaners.



- CAUTION
- \* Do not submerge the device or any of the components in water.
  - \* Do not subject the monitor to extreme hot or cold temperatures, humidity or direct sunlight.
  - \* Store the device and the components in a clean, safe location.
  - \* Do not subject the monitor to strong shocks, such as dropping the unit on the floor.
  - \* Remove the batteries if the unit will not be used for three months or longer. Always replace all the batteries with new ones at the same time.
  - \* This product is designed for use over an extended period of time; however, it is generally recommended that it be inspected and calibrated every two years to ensure proper function and performance.
  - (\* Replacement is done by EU representatives)

## SPECIFICATIONS


|                      |   |
|----------------------|---|
| Measuring Method     | Oscillometric Measurement   |
| Indication           | Digital LCD display   |
| Measuring Range:     | Pressure: (30~280)mmHg<br>Pulse: (40~199)Beat/min                               |
| Accuracy:            | Static Pressure: ±3mmHg Pulse: ±5%  |
| Memory:              | 90 Memories   |
| Power supply:        | 2x1.5V Batteries(LR03 or AAA)<br>use alkaline battery, measure above 200 times. |
| Operating condition: | +5°C~+40°C, 15%RH~93%RH<br>Atmospheric pressure: 70kPa~106kPa                   |
| Storage condition:   | -20°C~+55°C, 0%RH~93%RH<br>Atmospheric pressure: 50kPa~106kPa                   |
| Dimensions:          | Approx: 62(W)X78(H)X31(D)mm   |
| Weight:              | Approx: 130g, excluding batteries   |
| Classification       | Type BF   |
| Wrist circumference  | (13~19.5)cm   |

\* Specifications may be changed without notice in the event of improvement being made.


## HOW TO TAKE PROPER MEASUREMENTS

1. Fasten the wrist cuff according to the instructions in "ATTACHING THE WRIST CUFF".
2. Press the "ON/OFF" button. All icons appear two seconds on DISPLAY, then switch to measurement, and display "0" for last measurement record.

3. Start measurement, the cuff in the strap will automatically inflate. The mark will flash on LCD. When complete, the results will be displayed.



## WHO BLOOD PRESSURE CLASSIFICATION DISPLAY



## TROUBLESHOOTING

If you have trouble in using the unit please check the following points first.

| ERROR DISPLAY  | POSSIBLE CAUSE   | HOW TO CORRECT  |
|--|--|---|
| Nothing is displayed                                 | No battery installation  | Insert batteries  |
| When you push the POWER button or Battery icon flash | Battery worn out<br>The polarities of batteries placed wrongly | Replace new batteries<br>Insert battery in the correct polarities |

## Guidance and manufacturer's declaration - electromagnetic immunity

The Model 240389 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 240389 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.

| Immunity test                                | IEC 60601 test level              | Compliance level                  | Electromagnetic environment-guidance   |
|--|-----------------------------------|-----------------------------------|--|
| Electrostatic discharge (ESD) IEC 61000-4-2  | ±8 kV contact<br>±2 kV, ±4 kV air | ±8 kV contact<br>±2 kV, ±4 kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Power frequency magnetic field IEC 61000-4-8 | 30 A/m, 50/60Hz                   | 30 A/m, 50/60Hz                   | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.      |

NOTE U<sub>i</sub> is the a.c. mains voltage prior to application of the test level

## Guidance and manufacturer's declaration - electromagnetic emissions

The Model 240389 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 240389 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.

| Immunity test              | IEC 60601 test level   | Compliance level | Electromagnetic environment-guidance   |
|----------------------------|--|------------------|--|
| Conducted RF IEC 61000-4-6 | 3 Vrms/50 kHz to 80 MHz<br>6 Vrms 80 kHz to 80 MHz outside ISM bands | N/A              | Portable and mobile RF communications equipment should be used no closer to any part of the Model 240389 Series Electronic Blood Pressure Monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance |

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a) The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are: 13,567 MHz to 13,567 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are: 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b) The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2,7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.

## Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model 240389 Series Electronic Blood Pressure Monitor is used exceeds the applicable RF compliance level above, the Model 240389 Series Electronic Blood Pressure Monitor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Model 240389 Series Electronic Blood Pressure Monitor.

d) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the Model 240389 Series Electronic Blood Pressure Monitor

| Rated maximum output of transmitter | 150 kHz to 80 MHz              | 80 MHz to 800 MHz              | 800 MHz to 2,7 GHz           |
|-------------------------------------|--------------------------------|--------------------------------|------------------------------|
| W                                   | $d = \frac{3.5}{U_i} \sqrt{P}$ | $d = \frac{3.5}{U_i} \sqrt{P}$ | $d = \frac{2}{U_i} \sqrt{P}$ |
| 0.1                                 | 0.12                           | 0.12                           | 0.23                         |
| 0.1                                 | 0.38                           | 0.38                           | 0.73                         |
| 1                                   | 1.2                            | 1.2                            | 2.3                          |
| 10                                  | 3.8                            | 3.8                            | 7.3                          |
| 100                                 | 12                             | 12                             | 23                           |

For transmitters separated at a maximum power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.