


MOIST/DRY NON-AUTO SHUT OFF HEATING PAD

80-800-000-L

BV
MEDICAL



Controller
(On/Low, Med, High, Off)

HOW TO USE YOUR HEATING PAD

1. The sponge may be lightly moistened or used dry.
2. Insert sponge and heating pad into cloth cover.
3. Insert plug into 110 - 120 Volt electrical outlet.
4. Press the button on the controller for desired heat setting of Low, Medium, or High.
5. Turn off by pressing button until setting indicator is off.

DO NOT SIT OR LIE ON THE HEATING PAD

WARNING – DO NOT USE OVER INSENSITIVE SKIN OR IF YOU HAVE POOR CIRCULATION. PERSONS WITH DIABETES SHOULD CONSULT THEIR PHYSICIAN BEFORE USE. PEOPLE WHO ARE INCAPACITATED AND CHILDREN SHOULD ONLY USE THE HEATING PAD UNDER STRICT SUPERVISION AND SHOULD NOT BE LEFT UNATTENDED. NEVER USE THE HEATING PAD AS A BED OR FOOT WARMER.

SERVICE CHECKS

If your heating pad fails to heat at any time, please check the following:

1. See that plug fits properly and firmly into the electric outlet.
2. If heating pad still does not work, unplug it and wait 30 seconds, then try plugging it into the electric outlet.
3. Make sure the electrical cord is not broken and fits firmly into the heating pack & plug.

MOIST HEAT THERAPY

A frequently prescribed therapy in medicine today is the application of moist heat. Blood vessels react to moist heat by dilating, which in turn, increases blood flow to the area of treatment.

Back pain caused by muscle strain leads to muscle spasm. Moist heat tends to be useful in treating back pain. Rest, as well as the use of the heating pad will gradually relax the muscle(s) in spasm and increase flexibility, thereby reducing pain and discomfort.

Other common conditions that moist heat may be useful for are: Arthritic and Rheumatic conditions.

WARNING: This heating pad does not have an automatic shut off feature. Do not leave unattended while in use.

LIMITED WARRANTY: This heating pad is warranted against defective material or workmanship for 3 years from date of purchase from the manufacturer. It will be repaired or replaced as long as it is used in accordance with safety instructions and at the sole discretion of BV Medical. This warranty does not cover misuse or abuse of the product.

To obtain maximum results while using the heating pad, treatment should not last longer than 15 minutes, unless prescribed by your physician. Excessive treatment may cause adverse effects. In the case of any unusual symptoms, discontinue use and call your physician immediately.

INSTRUCTIONS FOR CARE OF CLOTH COVER

Remove cloth cover over the heating pad. To keep the cover fresh and clean, wash it with mild soap in lukewarm water and drip dry to keep it fluffy. DO NOT MACHINE WASH OR DRY.

If cover becomes discolored, scorched or worn, please replace it with a new one. The sponge inside is hand wash only in lukewarm water and drip dry.

It is to be expected that the cover of the heating pad unit may eventually become discolored. Excessive scorching may cause the cover to deteriorate.

IMPORTANT SAFETY INSTRUCTIONS DO NOT DESTROY

1. All instructions must be read.
2. Use while sleeping is prohibited.
3. Not for use on small children or infants.
4. Check to see that skin is not red or inflamed during use.
5. DO NOT, use this pad on a person who cannot feel the heat, such as persons who are paralyzed or unconscious or who have poor circulation or insensitivity in the area being treated.
6. The cloth cover must be on the unit before use.
7. Do not use pins or metallic fasteners on this pad anywhere over your body.
8. Do not sit or lay on top of the heating pad.
9. Do not fold or crush the heating pad while in use.

SAVE THESE INSTRUCTIONS

Manufactured by:
BV Medical
28W206 Commercial Ave., Unit B
Lake Barrington, IL 60010
info@bvmedical.com | www.bvmedical.com
(888) 822-8293

Made in China