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MOUTHGUARD FITTING INSTRUCTIONS

- 1. Gather the necessary materials: your mouthguard, scissors, a small pot, a slotted spoon, a small bowl, and a towel. Fill the pot with water and bring to a boil. Fill the small bowl with ice water
- 2. Place the mouthguard in your mouth up against your top row of teeth, and make sure it fits comfortably. It may be necessary to trim the back of your mouthguard at the ends with scissors if it pokes into the back of your jaw. (note: mouthguards are mostly intended to protect your front teeth, and many athletes prefer shorter fits that don't entirely cover their molars)
- 3. Using a slotted spoon, dip your mouthguard into the pot of boiling water. Keep it fully submerged for 30-60 seconds
- 4. Remove mouthguard from boiling water and pat dry with a towel. Mouthguard will be warm to the touch, but shouldn't be hot enough to burn you
- 5. Place your mouthguard up against your top teeth and push with your thumbs. Then bite on the mouthguard to form smaller indentations with your bottom teeth. Suck the guard tightly onto your top teeth while applying pressure with your tongue for about 20 seconds. Be sure not to chew the mouthguard or move it around
- 6. After about 20 seconds, remove the mouthguard and place it in the ice water bowl
- 7. Once your mouthguard is cool, it's good to go. If the fit isn't comfortable, you can repeat these steps to re-fit the mouthguard

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