

### ★ Storage

Store at -20°C.

#### Contents

- Product Manual
- **amfiRivert** Reverse Transcriptase, 1 ul/rxn
- **amfiRivert** Reverse Transcriptase 5X Reaction Buffer

ALL PRODUCTS SOLD BY GenDEPOT ARE INTENDED FOR RESEARCH USE ONLY UNLESS OTHERWISE INDICATED. THIS PRODUCT IS NOT INTENDED FOR DIAGNOSTIC OR DRUG PURPOSE

### ★ Storage

Ship with ice pack and dry ice.

### ★ Introduction

**amfiRivert** Reverse Transcriptase is provided in quantities sufficient for 100 first-strand cDNA synthesis reactions of 20ul each. **amfiRivert** Reverse Transcriptase enables robust, full-length cDNA synthesis for the reproducible analysis of rare or long messages. The **amfiRivert** Reverse Transcriptase can be used to reverse transcribe total RNA, poly(A)+ mRNA or synthetic transcript RNA templates.

### ★ amfiRivert Reverse Transcriptase 5X Reaction Buffer

250mM Tris-HCl (pH 8.3 at 25°C), 375mM KCl and 50mM DTT, 17.5mM MgCl<sub>2</sub>

### ★ Quality Control Assays

#### First-Strand cDNA Synthesis

Reverse Transcription is performed on a mixture of RNA templates ranging in length from 0.5 to 10kb primed with oligo(dT). The ranging is compared to a reaction using a reference lot of the enzyme. Prominent bands of 8 different cDNA products are observed by gel electrophoresis and autoradiography.

#### Amplification

When 0.25 zeptomoles (approximately 100 copies) of 1.2kb Kanamycin Positive Control RNA is reverse transcribed at 42°C and amplified (40 cycles), the result is a clear, discrete 323bp DNA product as visualized on a agarose gel by ethidium bromide staining.

### ★ Protocol

#### Prepare RNA Target and Primer

1. Use sterile, nuclear-free, thin-walled tubes, rechilled on ice.
2. For each 20ul reverse transcription (RT) reaction, combine.

Components	Volume
RNA template	up to 1ug
Primer	20pmol or 0.5ug
Nuclease-free water to a final volume of	5ul

3. Incubate at 70°C for 5 minutes (option).
4. Quick-chill at 4°C for 5 minutes and hold on ice.

#### Prepare Reverse Transcription Mix

1. For each 20ul RT reaction, combine

Components	Final Conc	Volume
Water, DEPC Treated		Xul
<b>amfiRivert</b> 2X Reaction Buffer	1X	4ul
dNTP mix (10mM each)	0.5mM	1ul
RNase Inhibitor Plus	1ul	20ul
Vortex the mixture <b>amfiRivert</b> Reverse Transcriptase		1ul
Final Volume RT Mix per 20ul reaction		15ul

2. Vortex gently to mix.
3. Dispense 15ul aliquots into reaction tubes.

#### Note:

- **amfiRivert** Reverse Transcriptase 5X Reaction Buffer contains a fixed MgCl<sub>2</sub> concentration of 3.5mM (1X). However, higher concentrations may be achieved by add

ing additional MgCl<sub>2</sub>.

- We recommend keeping the RT reaction mix chilled on ice prior to incubation.

#### Reverse Transcription

1. Anneal at 25°C for 5 minutes.

#### Note:

- The 5 minutes, 25°C annealing step is suggested for using oligo dT and Random Hexamer.

- This step can be skipped for using gene specific primer.

- If you need to improve specificity, minimize nonspecificity, try annealing at a more elevated temperature.

2. Extend the first strand for 60 minutes at 42°C. The extension temperature may be optimized between 37°C-55°C.

3. Heat-inactivate the **amfiRivert** Reverse Transcriptase by incubating at 70°C for 15 minutes.

4. Analyze cDNA, proceed with PCR or store frozen.

#### Note:

- Annealing conditions may require optimization step. The extension

### ★ Related Products

Description	Cat No
dNTP mixture, 10mM each	D0610
Oligo(dT) <sub>18</sub>	O1024
RNase Inhibitor Plus	R2808
RNazor, RNase Decontamination Reagent, Spray bottle	R7000
DEPC Treated Water	W0805