

★ Storage

Store at 2-8 °C.

📁 Contents

- Product manual
- Gelatin Solution 0.1%, Cell culture tested

ALL PRODUCTS SOLD BY GenDEPOT ARE INTENDED FOR RESEARCH USE ONLY UNLESS OTHERWISE INDICATED. THIS PRODUCT IS NOT INTENDED FOR DIAGNOSTIC OR DRUG PURPOSE

★ Introduction

0.1% Gelatin Solution is used to coat cell culture flask. The use of coated flasks improves cell attachment for certain types of primary cells as certain continuous cell lines. This sterile-filtered solution contain 0.1% gelatin in water.

★ Handling / Directions For Use

NOTE: The gelatin solution may be cloudy upon receipt due to exposure to low temperatures during shipping. If this occurs, place bottle in a 37°C water bath for 1 hour (maximum of 2 hours). Ensure the temperature of the water bath is ≥ 35°C and < 40°C. The gelatin will dissolve, resulting in a clear solution. Product performance will not be affected.

TO COAT CULTURE VESSELS:

1. Aseptically dispense enough gelatin into the dish or flask to completely cover the bottom of the vessel. Recommended volumes are:
 - 3 mL per 25 cm² flask or 60 mm tissue culture dish
 - 7 - 8 mL per 100 mm tissue culture dish
2. Leave gelatin solution in contact with the culture ware for at least 20 minutes at room temperature (15 - 25°C).
3. Aspirate the excess gelatin solution and leave culture ware sitting open in the biosafety cabinet until dry.

NOTE: Gelatin-coated dishes or flasks can be stored at room temperature (15 - 25°C) for up to 2 weeks, or at 2 - 8°C until required.

★ Related Product

Product Name	Cat No
Ceraplex Defined Serum Replacement(10X)	C3317
Ceraplex, Premium Cell Culture Supplement(10X)	C3318