

☆ Storage

Store at -15 ~ 25 °C. Avoid exposure to frequent temperature changes. See the expiration date on the Product Information Label.

☆ Contents

- Product Manual
- DNase I (Pharmaceutical Grade)

ALL PRODUCTS SOLD BY GenDEPOT ARE INTENDED FOR RESEARCH USE ONLY UNLESS OTHERWISE INDICATED. THIS PRODUCT IS NOT INTENDED FOR DIAGNOSTIC OR DRUG PURPOSE

☆ Shipping Condition

Frozen on dry ice.

☆ Introduction

DNase I is an endonuclease derived from bovine pancreas that will degrade double-stranded DNA in the presence of divalent cations, producing 3'-OH oligonucleotides. Mg^{2+} in the reaction solution causes the enzyme to produce nicks in double-stranded DNA, while in the presence of Mn^{2+} , DNase I cleaves both strands of the DNA. DNase I is useful in nick translation for introducing single-stranded nicks that serve as primer sites for initiation of DNA synthesis and for cloning random DNA fragments by cleaving double-stranded DNA. DNase I is a chromatographically pure preparation. It is offered as lyophilized powder. For greatest stability, it is important that DNase be dissolved at a concentration of at least 1mg/ml in 50% Glycerol with 20mM Tris-HCl, pH 7.5 and 1mM $MgCl_2$. This solution can be stored at -20°C for at least a year.

☆ Recommended 1X Reaction Buffer

40mM Tris-HCl (pH 8.0), 10mM $MgSO_3$, 10mM $CaCl_2$.

☆ Heat Inactivation

10 minutes at 65°C in the presence of Stop Solution.

☆ Inhibitors

EGTA; EDTA; Salt concentration >100mM will reduce DNase activity.

☆ Requirement

Ca^{2+} and Mg^{2+} or Mn^{2+}

☆ Source

Bovine pancreas.

☆ Stop Solution

20mM EDTA (pH 8.0)

☆ Unit Definition

One Kunitz unit causes as increase in absorbance at 260nm of 0.001 per minute per ml at 25°C pH 5.0 when acting on highly polymerized DNA in the presence of ionized magnesium and calcium.

☆ Procedure

1. Add to an RNase-free PCR tube:

1ug of RNA sample

49ul of 1X Reaction Buffer:

1ul of DNase I, 1 unit/ml*

* Refer to the Certification of Analysis for the lot specific activity.

To dissolve DNase I at a concentration of at least 1unit/ml, storage buffer is recommended: 50% Glycerol, 20mM Tris-HCl, pH 7.5, and 1 mM $MgCl_2$.

This solution can be stored at -20°C for at least one year.

2. Incubate for 10 ~ 15 minutes at 37°C.

3. To stop the reaction, add 1ul of Stop Solution to bind calcium and magnesium ions and to inactivate the DNase I.

Note: The Stop Solution (20mM EDTA) must be added before heating to prevent metal (Mg/Ca) ion catalyzed hydrolysis of the RNA.

Heat at 70°C for 10 minutes to denature both the DNase I and the RNA. This product should not be used for digestions longer than 15 minutes or for digestions at temperatures higher than 37°C, or the residual contaminating RNase activity will begin to degrade the RNA.

☆ Usage Note

1. This DNase solution does not contain an RNase inhibitor.

2. Under different buffer conditions the amount of DNase required to completely digest a given amount of DNA may need to be empirically determined. For example, salt concentrations >100mM will reduce DNase activity.