

★ Storage

Store at 2-8°C.

★ Contents

- Product Manual
- West-Ez Blocking Buffer, 3% BSA, 100ml contains 3% (w/v) Bovine Serum Albumin in 25mM Tris, 150mM NaCl, pH 7.4 with Kathon Antimicrobial Agent

ALL PRODUCTS SOLD BY GenDEPOT ARE INTENDED FOR RESEARCH USE ONLY UNLESS OTHERWISE INDICATED. THIS PRODUCT IS NOT INTENDED FOR DIAGNOSTIC OR DRUG PURPOSE

★ Shipping Condition

Ship with ice pack.

★ Introduction

West-Ez Blocking Buffer, 3% BSA is a preformulated IgG-Free, Protease-Free BSA solution for blocking excess binding sites in ELISA, Western blotting, Immunohistochemistry and other Immunochemical applications.

★ Procedure

1. Add the blocking solution to the ELISA well, blotting membrane or immunohistochemical slide.
2. Incubate 30 minutes to 2 hours at room temperature or 37°C
3. Continue with the standard protocol for the specific application.

★ Related Products

Product Name	Cat No
Albumin, (IgG, Fatty Acid and Protease Free)	A0100
AffiSelect Cox IV Loading Control Antibody	A0017
AffiSelect GAPDH Loading Control Antibody	A0039
AffiSelect Beta-Actin Loading Control Antibody	A0042
AffiSelect a-Tubulin Loading Control Antibody	A0050
NP-40 Lysis Buffer (2X)	N1200
Affiselect Total Protein Extraction Solution	A0710
RIPA Cell Lysis Buffer (1X) with EDTA	R4100
Xert Protease Inhibitor Cocktail Solution (100X)	P3100
Xpert Prestained Protein Marker (6.5-240 kDa)	P8502
Xpert 2 Prestained Protein Marker (10-240 kDa)	P8503
TBS Buffer, 20X, pH 7.4	T8054
West-Ez Stripping buffer	S2100
Goat anti-Rabbit IgG(H+L)-HRP, 1mg/ml	SA002
Rabbit anti-Goat IgG(H+L)-HRP, 1mg/ml	SA007
Tween 20, Molecular Biology Grade	T9100
West-Q Pico Dura ECL Solution	W3653
West-Q Femto ECL Solution	W3700
West-Ez Blocker, 5% Non-Fat Kit	W3710