PT137 Platinum Thoracolumbar Support - Instructions for Use

INDICATIONS: Back sprains and strains, osteoarthritis, mild osteoporosis, posture control

Take a snug hip measurement and choose the appropriate size from the sizing chart below.

1. Unroll the garment and release the hook and pile closures.

2. Center the garment on the back, ensuring the bottom edge sits at the distal margin of the sacrum.

3. Attach the front sections together, ensuring top and bottom are level with one another.


4. Pull the top two side elastic straps and attach to the front. Do the same with the bottom straps.

5. Criss-cross the shoulder straps in the back and bring them to the front to attach.

6. When replacing the plastic stays with contoured stainless steel stays (model 02031-17), remove the garment, shape the steel stays to required patient contours and insert them into the garment casings.

Reapply on patient and make adjustments as needed. Teach the patient on how to remove and reapply it.

Heat Moldable Rigid Back Panel available separately. Model PT742.

Sizing guide and selection:

Measure HIP circumference

Model No.	Size	HIP circumference
PT137-02	S	75-85 cm / 29.5-33.5 in
PT137-03	M	85-95 cm / 33.5-37.5 in
PT137-04	L	95-105.5 cm / 37.5-41.5 in
PT137-05	XL	105.5-115.5 cm / 41.5-45.5 in
PT137-06	XXL	115.5-126 cm / 45.5-49.5 in
PT742	One size	Heat moldable rigid insert

Care instructions: Remove the stays. Close all hook and pile closures on themselves before washing. Hand wash or machine wash on delicate cycle with mild detergent in cold water. Lay flat to dry. Do not dry clean, bleach or iron; do not use fabric softener.

Caution: If a rash or redness occurs, please discontinue use and consult your practitioner. For single use only.

Warranty: 90 days. For warranty claim, contact the practitioner from whom you made the purchase.

NOTE: Spinal braces should be worn over a close fitted cotton shirt or undershirt.

QUESTIONS

Contact Customer Service at:

USA

Tel: +1 800 492 1088

Fax: +1 800 245 3765

Email: info-usa@trulife.com

Canada

Tel: +1 800 267 2812

Fax: +1 613 392 4139

Email: infocanada@trulife.com

Trulife has appointed Medical Device Safety Service (MDSS) of Hannover, Germany to act as our EU authorized representative. They may be contacted at:

MDSS GmbH

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30175 Hannover

Germany

Phone (+49)-511-6262 8630

FAX (+49) -511-6262 8633

LIMITED WARRANTY

Trulife warrants that this product will be free from defects in material and workmanship from the date of installation.

This warranty will not apply if the product has been damaged by misuse, abuse, neglect, improper care, failure to follow instructions, abnormal wear and tear, or in the event that this product has been modified/repaired by persons unauthorized by Trulife.

If a defect in material or workmanship is found during the warranty period, Trulife will, at Trulife's option, either repair or replace the product. If it is not possible to repair or replace the product, Trulife will be limited to refunding the purchase price.

Trulife will not be liable under any legal theory for any direct, indirect, special, incidental or consequential damages arising from the use of or inability to use this product.

The application guidelines for this Trulife product are for the use of and by certified, qualified practitioner only. Patients are not to attempt to apply or adjust the item unless expressly instructed to do so by the practitioner responsible for the prescription and/or initial fitting of the device. All patient questions should be referred to the practitioner and not to the manufacturer. The manufacturer warrants only that the enclosed product has been inspected for quality and can be effective for certain indications, but final decisions and ongoing monitoring must be made by the medical professional(s) prescribing and/or fitting the device to determine its effectiveness for an individual patient. Patient compliance is an integral part of the entire protocol and must be adhered to in order to avoid potential problems and to maximize the effectiveness of the prescribed product.

As a condition of the sale of any Trulife product, this product is restricted to a "Single Patient Use Only" by the originally fitted patient in order to protect the care provider and the patient against potentially adverse consequences of infectious disease transmission, material instability in adapting to the configuration of the original user and/or decrease in effectivity. Any express or implied warranties are voided if the product is reused or fitted to another patient. Additionally, a license of right to use under any relevant patents pertaining to the product is terminated with the cessation of use by the original patient. As with all Trulife products, this product must be prescribed and applied by a qualified practitioner to determine it meets the needs of the particular patient and accomplishes the desired results.

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