C €⁰¹²³

Dental Electric Motor Instruction

Please carefully read this Manual before.

Thank you for purchasing MT2 Dental Electric Motor manufactured by Guilin Woodpecker Medical Instrument Co., Ltd. To ensure the correct use of device, it is recommended to carefully read the content on Installation, Operation, Maintenance, etc. in the Manual. For easier check, it is recommended to place the Manual in a position easy to access.



ZMN-SM-641 V1.4- 20230221

Guilin Woodpecker Medical Instrument Co., Ltd.

Content

Forward	1
1 Introduction	1
2 Basic technical parameters	3
3 Product performance structure and composition	4
4 Main unit interface	7
5 Function and operation	9
6 Safety precautions	-12
7 Clean, disinfection, and sterilization	-13
8 Troubleshooting	-17
9 Storage and transport	-17
10 After-sales service	-18
11 Environment protection	-18
12 Symbol instruction	-19
13 EMC-Declaration of comformity	-19

Forward

Guilin Woodpecker Medical Instrument Co., Ltd is a professional manufacturer researching, developing, and producing dental products. Woodpecker owns a sound quality control system and four brands, Woodpecker, DTE, DBA and RTA. Its main products include Ultrasonic Scaler, Curing light, Apex locator, Ultrasurgery, Endo Motor, Dental electric Motor, etc.

1 Introduction

Dental electric motor for driving dental handpieces for dental surgery, Mainly used for dental aesthetics restoration, crown breaking, open marrow preparation, deburring, polishing and other aspects of power. The device must be operated in hospital and clinic by the qualified dentists.

1.1 Precautions before operation

Dangers:

To prevent electric shock, do not use a wet hand to pull the power cord, and please prevent the water from entering the control circuit.

Keep away from explosives and flammable materials. Do not use the dental electric motor for patients who are anesthetized with nitrous oxide.

Warnings:

This dental electric motor may malfunction when used in an environment where electromagnetic interference occurs. This dental electric motor cannot be installed near the device that releases the magnetic wave. When using an ultrasonic vibrating device or an electrode knife in the vicinity, turn off the switch on the dental electric motor control panel.

MT2 requires special precautions for EMC and needs to be installed and put into use according to the EMC environment.

Device with electromagnetic transmitting will affect the normal operation of MT2. Please do not run both devices at the same time.

• Do not use it in operating rooms that contain a mixture of potentially flammable gases.

To avoid possible injury or damage to the dental electric motor, make sure that the motor handpiece (hereinafter referred to as the motor) is completely stopped when changing the contra-angle. (And the contraangle tool should be replaced by the pedal controller.)

• A severe impact, such as a drop from high position, can result in damage to the dental electric motor.Do not try to disassemble the controlling penal or motor.

• After use, please immediately clean, lubricate and disinfect the dental handpiece (hereinafter referred to as the handpiece).

• Do not lubricate the motor. The lubricant can cause overheating and damage the motor.

• Do not use a solution with dissolving ability to clean the control panel.

• Do not remove the motor cable from the main unit.

• After each operation, turn off the power supply.

• APPLIED PART: The contra-angle

• Any serious incident that has occurred in relation to the device should be reported to manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The ratio displayed on the interface of the main unit is the speed ratio of the contra-angle handpiece, please select the contra-angle handpiece with the corresponding speed ratio.

1.2 Intended use

Electrical drive including media supply for dental handpieces in the field of preventative dentistry, conservative dentistry for tasks such as the preparation of cavities and prosthodontics for tasks such as the preparation of crown.

1.3 Model

MT2

1.4 Contraindications

None known.

1.5 Safety requirements

Guilin Woodpecker Medical Instrument Co., Ltd. will not be liable for any direct or indirect damages and losses under the following conditions:

• The device is used to the unmentioned usages or the usages outside the scope of application.

• The operator did not use the device with the method in accordance

with the procedures and requirements stipulated in the Instruction Manual.

• The wiring system of the room where the device is used does not meet the requirements of appropriate standard and other proper requirements.

• Assembling, operating, and repairing the equipment without the authorization of the Woodpecker.

• The environmental conditions in which the equipment is located or stored do not meet the requirements mentioned in the section on technical requirements.

• After each use, pull the contra-angle out from the motor in parallel, place it in the corresponding position, and powers off the main unit.

• During the process of use, the operator is in direct contact with the contra-angle. The contra-angle is made of copper.

2 Basic technical parameters

2.1 Specification of main unit

Model: MT2 Software version: MT2-V1.0.0 Size: 165.5mm×129.7mm×77.6mm

2.2 Specification of power

MODEL No.: UES90-300300SPA1 Power supply input: 100-240V~ 50/60Hz 1.5A Power supply output: DC 30V 3.0A

2.3 Motor specification

Model: E-MT

Rotation speed: 2000-40000 rpm; Voltage input: DC 24V Size: $\Phi 22 \times 76.7$ mm Tail cord length: 1.8m Spray water source: water pressure (2 bar ~ 5 bar), water flow > 50ml / min Spray gas source: air pressure (2.5 bar ~ 5 bar), air flow > 1.5L / min

2.4 Use environment

2.4.1 Ambient temperature: $+5^{\circ}C \sim +40^{\circ}C$

- 2.4.2 Relative humidity: $30\% \sim 75\%$
- 2.4.3 Atmospheric pressure: $70kPa \sim 106kPa$

2.5 Device safety classification

2.5.1 Type of protection against electric shock: Class II equipment

2.5.2 Degree of protection against electric shock: B type applied part

2.5.3 Degree of protection against harmful ingress of water: Ordinary equipment (IPX0). Not waterproof.

2.5.4 Classified by operation mode: Intermittent operating device

2.5.5 Degree of safety application in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide: Equipment cannot be used in the presence of a flammable anesthetic mixture with air, oxygen, or nitrous oxide.

2.5.6 The equipment cannot be used in the nuclear magnetic resonance environment.

3 Product performance structure and composition

3.1 Safety requirements during installation

Danger: Equipment is installed on the premise that the installation must comply with appropriate standards and associated electric safety requirements.

Danger: Never install the unit in an explosive atmosphere and do not operate in areas with flammable gases (anaesthetic mixtures, oxygen, etc.).

Danger: The installation site should be protected from shocks and splashes of water or other liquids.

DANGER: Do not install the unit near or above the heat source. It must be installed in a place where the surrounding air is sufficiently circulated. There is enough space around it, especially for the exhaust fan and the back position.

WARNING: Do not place the parts directly under sunlight or ultraviolet light.


WARNING: This equipment is movable, but please handles with care.

WARNING: Make sure the connection parts are dry before connecting the wires to the unit. If necessary, blow to dry it with an air gun.

WARNING: Detachable power cords should not be installed in locations that are difficult to disconnect.


WARNING: When using manual mode, the contra-angle handpiece is less than 3 minutes per contact.

3.2 See the packing list for the machine configuration.



It mainly consists of main unit, motor, power adapter, etc.








Figure 4

3.5.1 Connect/Disconnect motor and its tail cord

A. Connect, align and insert the motor and motor tail connections, twist the nut of the motor tail cord.

B. Remove the motor tail cord from the motor, unscrew and disconnect the motor wire nut, and gently exit the motor tail connector;

3.5.2 Connect/Disconnect motor and contra-angle

A. Insert the motor connecting shaft into the contra-angle, and then turn the contra-angle till you heard a "click" sound, to ensure accurate positioning;

B. When removing the contra-angle, pull the contra-angle out of the motor in parallel.

3.6 Installation steps

3.6.1 Open the package, check whether the items of the equipment are complete according to the packing list, and place the main unit on the stable surface.

3.6.2 Connect the four-hole handpiece tube to the control box through the four-hole handpiece connector and tighten it.

3.6.3 Connect the power adapter to the DC adapter connector and connect the socket to the DC adapter with the power cord.

4 Main unit interface

4.1 Main interface



Figure 5

4.2 Setting interface

	Setting	
	DEMO MODE DI/bar FACT SET	Lang
		L
	Figure 6	
4.3 Manual mode in	nterface Demo Mode	
•	Demo mode	
	200 000 RPM	+
	@16:1 @ 1:1 @ 1:5	
	Start FWD Light	
	Figure 7	

4.4 Interface of calibration mode



Figure 9

Enter the setting interface, click the calibration mode touch button, enter the air pressure calibration interface, click "Start", and the interface will pop up as shown in Figure 9. Fully press the foot pedal to until it displays 100%, and then release the pedal. Till then the calibration is successful.

4.5 Restore the factory setting



Figure 10

Enter the setting interface, and click the restore factory settings touch button to enter the restore factory settings interface for confirmation. In conforming interface, click OK, the interface shown in Figure 10 will pop up. Select "OK" to restore the original factory settings parameters or select "Cancel" to quit restoring factory settings.

4.6 Language selecting interface

语言设置 Langua	age se	tting
	请选择 se select	语言 language
	中文	Chinese
	英文	English
确定		取 消 Cancel

Figure 11

Enter the setting interface, click the language selection touch button so that the interface shown in Figure 11 will pop up, and select the desired language. When you click "OK", the corresponding language selected. When you select "Cancel", the original language setting would be maintained.

5 Function and operation

5.1 Install the product correctly according to the product installation steps, and the operator should face the screen.

5.2 Turn on the power switch on the main unit, the screen will display and enter the main control interface (Figure 1).

5.3 The electric motor operation is controlled by the foot pedal of the dental chair.

5.4 Make sure that the pedal control calibration is performed before using the unit for the first time.


5.5 Icon description

Icon	Name	Function	
P1 P2 P3	Mode	Select the preset fixed speed. (P1/P2/P3)	
@16:1 (@1:1) (@1:5)	Speed ratio	Select the speed of contra- angle (16:1/1:1/1:5)	
+	Speed adjustment	Increase speed	

	Speed adjustment	Decrease speed	
	Store	Store the set parameters	
FWD REV	Forward/Reverse rotation	Control the forward and reverse rotation of motor	
SPEED 5000 RPM	Speed	Display the set operating speed	
Ť	Setting	Enter the setting interface	
DEMO MODE	Manual mode	Enter manual mode adjusting interface	
psi/bar	Calibration	Enter calibration interface	
FACT SET	Restore factory setting	Restore the system to factory setting.	
Lang	Language selection	Enter language setting interface	
	Exit	Exit the submenu setting mode	
Start	Start	Start to activate the motor.	
Stop	Stop	Stop the motor	
-¥- Light	LED	Turn on/off the motor LED.	

200 000 RPM	Speed adjustment	Adjust the speed formerly set in the manual mode. Touch the progress bark to increase or decrease the speed.
-------------	------------------	--

5.6 Basic function adjustment on main unit controlling interface



rotation

5.6.1 Select speed ratio through touching key "1". The color of the key will change when touching the key.

If using a 16:1 contra-agnle, touch the 16:1 speed ratio button. If using a 1:1 contra-agnle, touch the 1:1 speed ratio button.

If using a 1:5 contra-agnle, touch the 1:5 speed ratio button.

5.6.2 Control the ON and OFF of LED through key "2".

5.6.3 5Adjust the speed of the motor by touching key "3". Increase or decrease speed with key "+" and "-".


5.6.4 Switch to the forward rotation or reverse rotation by touching key "4".

5.6.5 Select the working speed of corresponding mode by touching key "5".

5.6.6 Store the set speed ratio, speed, mode, forward rotation/reverse rotation, and the ON/OFF of LED by touching key "6".

5.6.7 Enter setting interface by touching key "7"..

5.7 Basic funtion adjustment in setting interface



5.7.1 In setting interface, touch key "8" Manual Mode key to enter manual mode interface. Under the manual mode state, the rotation of motor can be directly controlled without foot pedal. Touch Speed Ratio key to select corresponding speed ratio. Touch Speed key to adjust the speed. Touch Start key, Stop key, Forward/Reverse Rotation key, and LED key to control the output of motor. Touch the Exit key to exit the manual mode interface.

5.7.2 In setting interface, touch key "9" Calibration key, touch Start key, and the interface shown in Figure 8, 9would pop up. Fully step on the foot pedal until it displays 100%, and release the foot pedal to finish the calibration.

5.7.3 In setting interface, touch key "10" Restore Factory Setting key, click OK to enter interface shown in Figure 10, and decide whether to restore the factory setting by clicking "OK" or "Cancel" key.

5.7.4 In setting interface, touch key "11" to enter interface shown in Figure 7 for language selection.

5.7.5 In the setting interface, touch key "12" to exit submenu setting mode.

6 Safety precautions

Cautions:

6.1 For repairs and purchase of spare parts, please contact our authorized supplier.

6.2 The accuracy of the speed monitoring depends on the highprecision performance of the handpiece installed on the micro motor. If the handpiece of other manufacturers is used, the actual speed value may not be displayed correctly. To ensure the actual matching display speed, please use the matching handpiece.

6.3 Read this operating manual before use and fully understand the functions of each part.

6.4 Check the operating status of the dental electric motor before use to confirm that there is no abnormality.

6.5 Test the dental electric motor before use to ensure accurate operation.

6.6 If the dental electric motor is permanently malfunctioning (excessive vibration, noise and heat generation, etc.), please immediately close it and return it to the authorized dealer.

6.7 Clean the control panel with a damp cloth and turn off the power before cleaning.

6.8 The temperature of the motor handpiece surface of applied part may reach 43 °C, During use, attention should be paid to intermittent rest to avoid scald.

7 Clean, disinfection, and sterilization

7.1 The cleaning and sterilization of motor is as follows

Warning	Warning: 1. Before sterilization, the hose connector should be removed. 2. Do not spray any cleaning solution or lubricant into the motor. 3. Except the motor , the other parts such as main unit, power adapter and tail cord are not allowed to be sterilized. 4. Do not clean the motor in a washer-disinfectant unit. 5. The products may not be exposed to temperature above 137°C. 6. The cleaning and sterilization must be conducted before use.
Advice	The products have been designed for a large number of sterilization cycles. The materials used in manufacture were selected accordingly. However with every renewed preparation for use, thermal and chemical stresses will result in ageing of the products. The allowed maximum times of sterilization for motor is 250 times.
Processing principles	It is only possible to carry out effective sterilization after the completion of effective cleaning. Please ensure that, as part of your responsibility for the sterility of products during use, only sufficiently validated equipment and product-specific procedures are used for cleaning and sterilization, and that the validated parameters are adhered to during every cycle. Please also observe the applicable legal requirements in your country as well as the hygiene regulations of the hospital or clinic, especially with regard to the additional requirements for the inactivation of prions.
Preparation at the Point of Use	After use, the Contra-angel handpiece motor should be disassembled from the motor and then the motor should be disassembled from the motor pigtail. Immediately after disassembly, use cold water (<40°C) to remove any heavy dirt from the unit. Store the units in a humid environment to prevent dirt from remaining.
Transporta- tion	Safe storage and transportation to the reprocessing area to avoid any damage and contamination to the environment.

Preparation for Decontami- nation	The products must be reprocessed in a disassembled state.
Pre-clean	 Tools: tray, syringe, hose, clean and dry soft cloth Install the bearing protection sleeve. Install the bearing protection sleeve. ② Use a syringe to draw 500ml of purified water to flush the waterway and atomization airway to remove contaminants in the inner cavity. Do not clean the cooling air circuit, otherwise it will damage the motor. ³ Wipe the water stains on the interface with clean cloth to prevent water from entering the inner cavity and damaging the bearing. ³ Install the interfaceprotection sleeve. ⁵ Flush with running water for 2 min to remove surface contaminants. Be careful not to flush the port directly. ⁶ Flush with 500ml purified water to remove surface contaminants.

	Notice: Please do not use an automatic Washer-disinfector to clean the motor,
	otherwise the water and cleaning agents will cause damage to the motor.
	Notice: Never reprocess this medical device in an ultrasonic device . Otherwise
	cause malfunction and material damage.
	Do not allow any liquid enter internally into the motor either by the nose or hose
	connector.
	Step1: Clean the inner cavity
	(1) Remove the interface protection sleeve. (2) Use a suring to draw 500 m of determent to
	2) Use a syringe to draw 500ml of detergent to rinse the waterway and atomization airway to
	remove contaminants in the inner cavity.
	3 Wipe the water stains on the interface with clean cloth to prevent water from
	entering the inner cavity and damaging the bearing.
	(4)Install the interface protection sleeve.
	Step2: Wipe the surface
	① Wet the clean soft cloth in the detergent and wipe the surface of the test sample
	thoroughly for 5 times. Change a clean soft cloth after each wipe. If there are still
	visible contaminants, wipe repeatedly until there are no visible contaminants. 2 Wet the clean soft cloth in the detergent and wipe the surface of the test sample
	thoroughly for 5 min. If there are still visible contaminants, wipe repeatedly until there
	are no visible contaminants.
Cleaning	③ cleaning condition:
Citaning	Cleaning agent:RUHOF ENDOZIMER AW PLUS WITH APA
	Dilution ratio:1:270 Temperature:<60 °C
	Step3: Scrub the surface
	Scrub the test sample thoroughly with an instrument brush stained with detergent for 3
	min.
	Step4: Rinse
	(1) Remove the interface protection sleeve.
	2) Use a syringe to draw 500ml of purified
	water to rinse the waterway and atomization $e^{\circ} \circ \circ \circ$
	airway to remove contaminants in the inner cavity.
	③ Wipe the water stains on the interface with clean cloth to prevent water from
	entering the inner cavity and damaging the bearing.
	④ Install the interface protection sleeve.
	Step5: Rinse the surface
	① Rinse the sample surface with 500ml purified water to remove the residual
	detergent. Note: Do not rinse the port directly to avoid damaging the internal bearing.
	The boling international ordering.
Disinfection	N/A, The motor must eventually be sterilized, so the sterilization process is not
	applicable.

Drying	Drying should be conducted after cleaning. Wipe off the residual water stains on the surface of the test sample with dry absorbent soft cloth.
Maintenance	Visual inspection: Before packaging and auto sterilization, make sure that the motor has been maintained according to manufacturer's instruction. Visual check the integrity of the dental handpiece and the cleanliness of the device: Check the motor: 1. If there is still visible stain on the products after cleaning, the entire cleaning process must be repeated. 2. If the appearance of the dental handpiece is obviously damaged, smashed, detached, corroded or bent, it must be scrapped and not allowed to continue to be used.
Packaging	Pack the instruments in an appropriate packaging material for sterilization. The packaging material and system refer to EN ISO11607.
Sterilization	Sterilization of instruments by applying a fractionated pre-vacuum steam sterilization process (according to EN 285/EN 13060/EN ISO 17665) under consideration of the respective country requirements. Minimum requirements: 4 min at 134 °C (in EU: 5 min at 134 °C) Maximum sterilization temperature: 137 °C Put the motor in a high pressure steam sterilization bag, and seal it.And then sterilize the motor under the temperature of 134 °C (273 °F) and the pressure of 2.0 bar ~ 2.3 bar (0.20MPa~ 0.23MPa)) for not less than 4 minutes (in EU: 5 min at 134 °C) Allow a maximum sterilization time of 20 minutes at 134°C. Drying is needed after sterilization.
Precautions	 Only products that have been effectively cleaned or cleaned and disinfected are allowed to be sterilized; Before using the sterilizer for sterilization, read the Instruction Manual provided by the equipment manufacturer and follow the instructions. Do not use hot air sterilization and radiation sterilization as this may result in damage to the product; Please use the recommended sterilization procedures for sterilization. It is not recommended to sterilize with other sterilization procedures such as ethylene oxide, formaldehyde and low temperature plasma sterilization. The manufacturer assumes no responsibility for the procedures that have not been recommended. If you use the sterilization procedures that have not been recommended, please adhere to related effective standards and verify the suitability and effectiveness. *Fractional pre-vacuum procedure = steam sterilization with repetitive pre-vacuum. The procedure used here is to perform steam sterilization through three pre-vacuums.
Storage	Storage of sterilized instruments in a dry, clean and dust free environment at modest temperatures, refer to label and instructions for use.
Reprocessing validation study information	It is the duty of the user to ensure that the e reprocessing processes including resources, materials and personnel are capable to reach the required results. State of the art and often national law requiring these processes and included resources to be validated and maintained properly.

7.2 Clean, disinfection, and sterilization

7.2.1 Before each use, wipe the surface of the machine and the tail cord of the motor handpiece with a soft cloth or paper towel soaked in 75% medical alcohol. Repeat the wipe for at least 3 times.

7.2.2 After each use, wipe the surface of the device and the tail cord of

the motor handpiece with a soft cloth soaked in clean water (distilled or deionized water) or a clean disposable wipe. Repeat the wipe for at least 3 times.

8 Troubleshooting

Fault	Cause	Solution
Error 01	furning on the device so that the	Release the foot pedal and reboot.
Error 02		Check whether the power adapter is right.
Error 03	Motor abnormality	Check whether the motor is well connected or replace the motor.

If the problem still cannot be solved, please contact our local distributor or our company.

Note: The user must use the original accessories. Please contact our local dealer or the company for purchase. It is forbidden to use related accessories of other brands, so as to avoid damage to the electric motor or other dangers.

Warning: The equipment shall not be changed or modified without the express consent or authorization of Guilin City Woodpecker Medical Instrument Co., Ltd.

9 Storage and transport

9.1 The device should be handled carefully and lightly. Be sure that it is far from the vibration, and installed or kept in a cool, dry, and ventilated place.

9.2 Do not store the machine together with articles that is poisonous, combustible, caustic, or explosive.

9.3 This machine should be stored in a room where the relative humidity is $10\% \sim 93\%$, atmospheric pressure is 70kPa ~ 106 kPa, and the temperature is -20° C $\sim +55^{\circ}$ C.

9.4 Excessive impact and shake should be prevented during transport. Lay it carefully and lightly.

9.5 Do not put it together with dangerous goods during transport.

9.6 Avoid being exposed to sun, rain, and snow during transport.

10 After-sales service

Since the date of sales, for the device that cannot work normally for quality problem, with Warranty Card, our company is responsible for the repair. Please refer to the Warranty Card for details of warranty.

	Toxic or harmful substances or elements					
Part	(Pb)	(Hg)	(Cd)	(Cr6+)	(PBB)	(PBDE)
Main unit	0	0	0	0	0	0
Motor handpiece	0	0	0	0	0	0
Power adapter	0	0	0	0	0	0
Dental contra-angle	0	0	0	0	0	0
Mechanical elements, including bolts, nuts, washers, etc.	0	0	0	0	0	0

11 Environment protection

•: Indicates that the content of the toxic substance in all homogeneous materials of the part is below the limit requirement stipulated in SJ/T-11363-2006 Limit Requirements for Toxic and Hazardous Substances in Electronic Information Products.

 \times : indicates that the content of the toxic substance in at least one of the homogeneous materials of the part exceeds the limit requirement specified in SJ/T-11363-2006.

(This product meets EU RoHS environmental protection requirements; there is currently no mature technology in the world to replace or reduce the content of lead in electronic ceramics, optical glass, steel and copper alloy.)

According to the Administrative Measures on the Restriction of the Use of Hazardous Substances in Electric and Electronic Products and the Regulations on the Management of the Recycling of Waste Electric and Electronic Products and related standards, please observe the safety and precautions of the products, and after use, please recycle or dispose this product after according to the methods in local laws and regulations.

12 Symbol instruction

C € 0123	CE marked product	X	Appliance compliance WEEE directive
M	Date of manufacture		Manufacturer
Ŕ	Type B applied part		ClassIIequipment
	Used indoor only	8	Follow instructions for use
10%	Humidity limitation	106kPa	Temperature limitation
134℃ \$\$\$\$	Can be autoclaved	IPX0	Ordinary equipment
-20°C-++55°C	Atmospheric pressure for storage	ECREP	Authorised Representative in the EUROPEAN COMMUNITY

13 EMC-Declaration of comformity

The device has been tested and homologated in accordance with EN 60601-1-2 for EMC. This does not guarantee in any way that this device will not be effected by electromagnetic interference Avoid using the device in high electromagnetic environment.

Serial number	Cable name	cable length (m)	Cable type
1	Power cord (input)	1.2m	Unshielded parallel line
2	Power cord (output)	1.2m	Unshielded parallel line
3	Handle tail	1.8m	Unshielded parallel line

Technical Description Concerning Electromagnetic Emission

Table 1: Declaration - electromagnetic emissions

Guidance and manufacturer's declaration - electromagnetic emissions

The model MT2 is intended for use in the electromagnetic environment specified below. The customer or the user of the model MT2 should assure that it is used in such an environment.

RF emissions CISPR 11	Group 1	The model MT2 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF emissions CISPR11	Class B	The model MT2 is suitable for used in	
Harmonic emissions IEC 61000-3-2	Class A	all establishments, including domestic establishments and those directly connected to the public law valtage payor supply	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	to the public low-voltage power supply network that supplies buildings used for domestic purposes.	

Technical Description Concerning Electromagnetic Immunity

Table 2: Guidance & Declaration - electromagnetic immunity

Guidance & Declaration — electromagnetic immunity	
The model MT2 is intended for use in the electromagnetic environment specified below. The customer or the user of the model MT2 should assure that It is used in such an environment.	at

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2, ±4, ±8, ±15kV air	±8kV contact ±2, ±4, ±8, ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electric fast	±2kV for power		Mains power quality
transient/burst	supply lines	±2kV for power	should be that of a
IEC 61000-4-4	±1kV for Input/	supply lines	typical commercial or
	output lines		hospital environment.
	$\pm 0.5, \pm 1$ kV line	$\pm 0.5, \pm 1$ kV line to	Mains power quality
Surge	to line	line	should be that of a
1EC 61000-4-5	$\pm 0.5, \pm 1, \pm 2kV$	$\pm 0.5, \pm 1, \pm 2kV$	typical commercial or
	line to earth	line to earth	hospital environment.

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U _T (>95% dip in U _T .) for 0.5 cycle <5 % U _T (>95% dip in U _T .) for 1 cycle 70% U _T (30% dip in U _T) for 25 cycles <5% U _T	<5 % U _T (>95% dip in U _T .) for 0.5 cycle <5 % U _T (>95% dip in U _T .) for 1 cycle 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (>95 % dip in U _T)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the models MT2 requires continued operation during power mains interruptions, it is recommended that the models MT2 be powered from an uninterruptible		
	(>95 % dip in U _T) for 250 cycles	for 250 cycles	battery.		
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.		
NOTE U_T is the a.c. mains voltage prior to application of the test level.					

Table 3: Guidance & Declaration - electromagnetic immunity concerning **Conducted RF & Radiated RF**

Guidance & Declaration - Electromagnetic immunity					
The model MT2 is intended for use in the electromagnetic environment					
specified below. The customer or the user of the models MT2 should assure that					
it is used in such an environment.					
T	nmunity test IEC 60601 Compliance Electromagnetic environment - test level level guidance				
Immunity test	test level	level	guidance		

Immunity test		Compliance	Electromagnetic environment -
	test level	level	guidance

Conducted RF IEC 61000-4-6 Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 6 Vrms ISM frequency band 3 V/m 80 MHz to 2.7 GHz	3V 6V 3V/m	Portable and mobile RF communications equipment should be used no closer to any part of the models MT2, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.2 \times P^{1/2}$ $d=2 \times P^{1/2}$ 80 MHz to 800 MHz $d=2.3 \times P^{1/2}$ 800 MHz to 2.7 GHz where P is the maximum output power rating of the transmitter In watts (W) according to the transmitter manufacturer and d Is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,a should be less than the compliance level in each frequency range. b Interference may occur In the vicinity of equipment marked with the following symbol: $(((\cdot))))$	
NOTE I At 80 MHz end 800 MHz. the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

a Field strengths from fixed transmitters, such as base stations for radio (cellular/ cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the model MT2 is used exceeds the applicable RF compliance level above, the model MT2 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the model MT2. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

 Table 4: Recommended separation distances between portable and mobile

 RF communications equipment and the model MT2

Recommended separation distances between			
portable and mobile RF communications equipment and the model MT2			
The model MT2 is intended for use in electromagnetic environment in			
which radiated RF disturbances is controlled. The customer or the user of the			
model MT2 can help prevent electromagnetic interference by maintaining a			
minimum distance between portable and mobile RF communications equipment			
(transmitters) and the model MT2 as recommended below, according to the			
maximum output power of the communications equipment.			

Rated maximum	Separation distance according to frequency of transmitter m			
output power of transmitter W	150kHz to 80MHz d=1.2×P ^{1/2}	80MHz to 800MHz d=1.2×P ^{1/2}	800MHz to 2,7GHz d=2.3×P ^{1/2}	
0,01	0.12	0.12	0.23	
0,1	0.38	0.38	0.73	
1	1.2	1.2	2.3	
10	3.8	3.8	7.3	
100	12	12	23	

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) accordable to the transmitter manufacturer.

NOTE I At 80 MHz and 800 MHz. the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Scan and Login website for more information





Guilin Woodpecker Medical Instrument Co., Ltd. Information Industrial Park, Guilin National High-Tech Zone, Guilin, Guangxi, 541004 P. R. China

Tel: Europe Sales Dept.: +86-773-5873196 North/South America & Oceania Sales Dep.:+86-773-5873198 Asia & Africa Sales Dep.:+86-773-5855350 E-mail: woodpecker@glwoodpecker.com, sales@glwoodpecker.com Website: http://www.glwoodpecker.com Website: http://www.glwoodpecker.com MedNet EC-Rep GmbH Borkstrasse 10 · 48163 Muenster · Germany

Dental Electric Motor Warranty Card

Name of Customer		
Address Details		
Postal Code		
Tel		
Model		(I) For Distributor
Main Unit No.		
Handpiece No.		
Contra-angle No.		
Purchase Date		
Contact Person		
Date	Maintenance Record	Repairer
L		



Guilin Woodpecker Medical Instrument Co., Ltd. Information Industrial Park, Guilin National High-Tech Zone, Guilin, Guangxi, 541004 P. R. China

Europe Sales Dept. Tel: +86-773-5873196 North America, South America & Oceania Sales Dept. Tel: +86-773-5873198

+30-//3-36/3198 Asia & Africa Sales Dept. Tel: +86-773-5855350 Fax: +86-773-5822450 E-mail: woodpecker@glwoodpecker.com, sales@glwoodpecker.com Website: http://www.glwoodpecker.com

Distributor:	
	Seal

Dental Electric Motor Warranty Card

Name of Customer		(II) Return to Manufacturer
Address Details		
Postal Code		
Tel		
Model		
Main Unit No.		
Handpiece No.		
Contra-angle No.		
Purchase Date		
Contact Person		
Date	Maintenance Record	Repairer



Guilin Woodpecker Medical Instrument Co., Ltd. Information Industrial Park, Guilin National High-Tech Zone, Guilin, Guangxi, 541004 P. R. China

Europe Sales Dept. Tel: +86-773-5873196 North America, South America & Oceania Sales Dept. Tel: +86-773-5873198

+30-//3-36/3198 Asia & Africa Sales Dept. Tel: +86-773-5855350 Fax: +86-773-5822450 E-mail: woodpecker@glwoodpecker.com, sales@glwoodpecker.com Website: http://www.glwoodpecker.com

Distributor:

Seal

Warranty Instruction

I Period validity:

Within one year from the date of sale, the main unit, handpiece and contra-angle can be repaired for free by providing warranty card.

II Range of warranty:

Within the warranty period of validity, we are responsible for any troubles caused by quality problems or products technique and structure.

III The following are beyond our warranty:1. The damage caused by disobeying the operation instruction or lack of the needed condition.

2. The damage caused by unsuitable operation or disassembly without authorization.

3. The damage caused by unadvisable transportation or preservation.

4. There isn't the seal of distributor or the warranty card isn't filled in completed.

Warranty Instruction

I Period validity:

Within one year from the date of sale, the main unit, handpiece and contra-angle can be repaired for free by providing warranty card.

II Range of warranty:

Within the warranty period of validity, we are responsible for any troubles caused by quality problems or products technique and structure.

III The following are beyond our warranty:1. The damage caused by disobeying the operation instruction or lack of the needed condition.

2. The damage caused by unsuitable operation or disassembly without authorization.

3. The damage caused by unadvisable transportation or preservation.

4. There isn't the seal of distributor or the warranty card isn't filled in completed.