



DISPOSABLES

**KN95 MASK**



[WWW.SPAZIOMEDICAL.COM](http://WWW.SPAZIOMEDICAL.COM)

## Description

Non-Medical mask ideal for blocking bacteria, droplets and dust.


Mask body, nose clip and a mask band made of non-woven ultrasonic compound.

Material:

- 40g non-woven fabric surface layer
- 40g layer of cotton
- 40g layer of meltdown cloth
- 20g non-woven fabric inner layer

Intended for single use only

**Pollen prevention**  
**Allergy free**  
**Dustproof**  
**Sanitary material grade**



## Specifications

Name + Code	KN95 Mask	SPA/KN95.MASK
Material	Non - woven ultrasonic compound	
Packaging:	Box of 10	
Objects to be blocked	Bacteria, dust , droplets	
Dimensions	W.10.5 x H10.5 Earloop. 22cm	

## Instructions



Put the mask on your face covering your chin with the nosepiece on the top of your nose.



Pull the straps over your ears and adjust the mask until you feel comfortable.



Mold the nosepiece around the bridge of your nose to fit the mask to your face.



Adjust the mask into its final position.

1. Wash hands before putting on your mask, also wash hands after taking your mask off.
2. The facial mask must be placed so that metal nose piece faces outwards and at the top of the mask.
3. Put your mask on your face from the chin first, then pull the straps over your ears to adjust the mask properly.
4. Make sure there is no hair, facial hair, jewellery or clothing between your face and the mask. Adjust the the face mask until you feel comfortable.
5. Use both hands to mold the metal nosepiece around the bridge of your nose.
6. Perform a seal test check by placing both hands over your mask and breathe out, if there is any leakage around the nose, go back to step 5. If air leakage occurs around the edges of the face mask re-adjust the straps.