osseotouch
TOTAL CONTROL

magnetic MALETT

Medical device for surgery and dentistry - Instructions for use

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TABLE OF CONTENTS

1.0. Important information for the user	4
2.0. Device Illustrations	5
2.1. Control elements, indicators, connections and functions	7
3.0. Symbols - Labelling	8
4.0. Device description	10
4.1. Device general description	10
4.2. Description of the coupling parts	11
5.0. Safety Guidelines	12
5.1. Warning and Caution	12
5.2. Intended use	14
5.3. User's qualifications	14
5.4. Use conditions	14
5.5. Safety guidelines for the installation environment	14
5.6. Operating safety guidelines	15
6.0. Transport and storage	15
6.1. Transport and storage conditions	15
7.0. Installation and use	16
7.1. Unpacking	16
7.2. Installation and connection	16
7.3. Activation	16
7.4. Surgery instrument insertion on the handpiece OSTEOTOME function	17
7.5. Release of surgery instrument	17
7.6. How to activate the handpiece	18
7.7. Calibration of the force	18
8.0. Maintenance	19
8.1. Cleaning and maintenance	19
8.2. Disinfection and sterilization of the parts	19
8.3. Maintenance	21
8.4. Repairs	21
8.5. Fuses replacement	21
9.0. Disposal of waste	22
9.1. Responsibility	22
9.2. Warranty	22
9.3. Technical Literature	22
10.0. Technical description	23
10.1. Failures and malfunctions	23
11.0. Technical data	25
12.0. Electromagnetic compliance declaration	26
13.0. Crown remover handpiece	29
13.1 Important information	29
13.2 Illustrations	30
13.2.1 Handpiece parts	30
13.3 Intended use crown remover handpiece	30
13.3.1 Standard supply	31
13.3.2 Safety rules	31
13.4 Instructions for use	32
13.4.1 Hook insertion on crown remover handpiece	32
13.4.2 Hook release from crown remover handpiece	32
13.5 Sterilization instructions	33
13.6 Maintenance - Repairs	33

1.0. IMPORTANT INFORMATION FOR THE USER

This operating manual was made to help you to properly install and connect the Magnetic Mallet device. All the useful details for a proper use of this device are contained in this manual. You should read it very carefully and store it in its slipcase in a dry and clean place in order to gather any useful information in the future.


Attention: Please be aware you are referring to the last update of the instructions for use as the manufacturer reserves the right to make improvements at any time without prior warning. The images and drawings in this instructions for use are provided for illustrative purposes only.

This instruction for use complies with the EU 2017/745 of the European Commission of april 5, 2017 for medical devices, modifying the Regulations (EC) 178/2002 and 1223/2009 and repelling the directives 90/385/EEC and 93/42/EEC


2.0. DEVICE ILLUSTRATIONS


Magnetic Mallet⁽¹⁾

1) Magnetic Mallet can be also supplied
with the Crown Remover Handpiece (on request)


Front view

Rear view

Osteotome Handpiece

2.1. Control elements, indicators, connections and functions

CONTROL UNIT

- A. Socket
- B. Mains plug
- C. ON/OFF switch
- D. Control unit
- E. Manual control push button
- F. Knob for the calibration of the force
- G. Display showing the force level
- H. Power cable
- I. Fuse holder
- L. Footswitch
- M. Handpiece connector
- N. Footswitch connector

OSTEOTOME HANDPIECE

- O. Osteotome function handpiece
- P. Chuck
- Q. Ring Nut
- R. Instrument

3.0. EXPLANATION OF THE USED SYMBOLS - LABELLING

Symbols placed on the magnetic mallet chassis

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LC - 0000

Symbols labelling

	Information on the manufacturer
	Device serial number
	AC (Alternating current)
	Coupled part "BF type" (according to IEC 606011 rules)
	Prior to setting the device going, carefully read the instructions for use
	Caution hazardous voltage
	Store in a dry and clean place
	Caution! The improper use of this device can cause injury
	You should not use the device if the packaging is damaged
	Storage temperature
	CE mark including the identification number of the Notified Body
	Disposal of special waste (electric and electronic devices)
	Device Class II - IEC 4175172
	Medical Device
	Steam Sterilization
	UDI number
	Lot number

4.0. DEVICE DESCRIPTION

4.1. Device general description

Magnetic Mallet is an electromedical device for surgery and dentistry practice with a double function according to the coupled part (handpiece O):

- **OSTEOTOME** version to obtain the bone plastic deformation thanks to the available forces and their application timing.
- **CROWN REMOVER** version to obtain the crown and bridge removal from the abutment or the implant.

This device is basically made up by a control unit, a footswitch and two interchangeable handpieces which can be used depending on the needed practice.

The control unit case is made by self-extinguishing plastic. On the control unit you can find:

- The ON/OFF switch (C)
- Two different connectors. One for the footswitch and the other one for the handpiece. The two connectors are identifiable by the following graphic symbols:
 - ▶ Handpiece metal connector (M)
 - ▶ Footswitch plastic connector (N)
- Manual control pushbutton (E)
- Plastic knob to calibrate the level of the force (F)
- Display showing the selected force (G)

The control unit, the control/operational electronics together with the electronic/electric safeties are housed into the control unit case. The handpiece is powered by the control unit with a 58V c.c. maximum voltage. The handpiece is connected to the control unit by a cordset which can be manually unplugged. The handpiece can be activated by the pressure of both footswitch (L) and manual control push-button (E).

Both for the **OSTEOTOME** and for the **CROWN REMOVER** function the force applied to the handpiece through the pressure of the footswitch or of the push button is not continuous. Each pressure of the switch conveys a single force pulse

Turning the knob (F) you can select the desired level of the force among the four possible options (4); at the same time the selected force will be showed on the liquid crystal Display (G) and by the number selected on the knob. **Magnetic Mallet** has been designed to give 1 pulse per second. You should not exceed the maximum threshold of 200 consecutive pulses without having a 40 minutes break.

4.2. Description of the coupling parts

Osteotome Function Handpiece with ring nut

For the OSTEOTOME function you must use the Handpiece with the ring nut (O) equipped with the appropriate surgery instruments.

INTENDED USE

Medical device to perform a wide array of dentistry surgeries such as bone compaction, teeth and roots extraction and die cast/bridges removal (with the crown remover handpiece). The pulse conveyed to the surgery instrument allows you to achieve the plastic deformation of the spongy bone in order to create the site for the implant.

Through the handpiece the surgery instrument can be energized with 4 different forces so, that the surgeon can replace the energy carried out by the human force in the manual practice (state of art).

You can select the desired level of the force amongn the 4 possible options (1-4) by turning the knob (F) placed on the control unit.

In the premarketing tests we carried out, the most used forces have been the 2nd and the 3rd one; nevertheless it's important to remember that the whole spectrum of the forces has a very close progressive range. That's why they can all be considered suitable for surgery. It's entirely up to the surgeon to determine the more suitable one to achieve the desired result.

! In order to familiarize yourself with the device, in the first stage you should always start from the lighter force (No.1). It's also advisable to keep the surgery instrument between your fingers selecting the different levels of the forces. In this way you shall be able to assess the corresponding pulse.

Surgery instrument action is led both by the pressure exerted by the surgeon on the handpiece and by the axial forwards of the energy. In any case the maximum forwards of the surgery instrument is of about 1 mm.

Attention: The surgeon must always determine the pressure to be exerted on the bone in order to achieve the desired result.

5.0. SAFETY GUIDELINES

5.1. Warning and Caution

Please carefully read this Manual for use and rigorously adhere to the instructions. The captions "Caution", "Warning" and "Attention" have a precise meaning. You should carefully read the related instructions in order to grant a safe and effective working of this device.

Caution: indicates a danger/risk for the patient or the surgeon. The failure in comply with the indication can carry injuries to the patient or to the surgeon.

Warning: indicates the maintenance to comply with in order to avoid any possible damaging in the appliance.

Attention: indicates special information related to the use of the Magnetic Mallet device or other important information.

Caution: Prior to activating the device carefully read the manual for use. Pay close attention to the sections related to the device working and to safety guidelines in order to avoid any risk for the patient, your staff and yourself.

Caution: Install the Magnetic Mallet device keeping it out of the reach of the patient.

Caution: Always check the cables. If the cable is somehow damaged you should immediately replace it calling the authorized service provider.

Caution: Prior to activating the device you should always have a functional operational test (FOT) pressing the footswitch or the manual control push-button.

Caution: Prior to activating the device make sure you have connected the proper handpiece according to the intended practice/treatment.

Caution: It's strictly forbidden to open the control unit – **Electric shock risk!** Any opening of the control unit by non authorized persons nullifies the warranty.

Caution: You should move away the device from the patient every time there is stalling in the device (for example: power failure during the use).

Caution: You should unplug the device from the main power every time you goes on working on the control unit (for example : disinfection/cleaning, unplugging of the connector).

Caution: It's forbidden to put fingers or any tool inside the footswitch or the handpiece connectors.

Caution: Handle with care the packaging and avoid all sort of strikes or damages.

Caution: 1 pulse/sec for a maximum of 200 consecutive pulses then a 40 minutes break.

Warning: Keep dry! Make sure no liquid gets into the control unit and into the hand-piece. Do not lean liquid holders on it.

Warning: You should install the device in dry environments only and keep it dry. Take care no liquid gets into the control unit since this can lead to malfunction.

Warning: The device must be only supplied with the voltage which is stated on the rating plate placed on the rear panel of the control unit.

Warning: Use fuses having a value corresponding to the one stated on the line filter rating plate.

Warning: Handle and use with care the power cable and the handpiece/footswitch cordsets. In order to avoid cable stress use them in wide leeway situations.

Attention: The warranty coverage is nullified by improper usage of the device.

Attention: According to an official test the Magnetic Mallet meets the requirements of Regulation (UE) 2017/745 and 2014/30/EU concerning electromagnetic compatibility and doesn't emit any interfering radiation.

Attention: Do not use on patients:

- Less than 14 years old
- With health problems
- With systemic chronic dependency
- With problems in maxillary sinus
- With handicaps
- With coagulation disorders
- With periodontal infection in alveolar bone
- Smokers
- Alcohol or drugs addicted

The use of this device is unadvisable on patients with pace-maker without the prior opinion of his attending physician.

Attention: In the first stage you should start by a low level force and increase it step by step, as necessary.

Attention: It's up to the surgeon to determine the more suitable pressure to be applied on the bone in order to achieve the desired result.

Attention: It's up to the surgeon to determine the more suitable force to achieve the desired result.

ATTENTION: You should be knowledgeable about the device and its working before using it on a patient.

5.2. Intended use

The **Magnetic Mallet** is a medical device to be used in dental surgery. According to the coupled handpiece it can have a different use:

- OSTEOTOME function - HANDPIECE with ring nut
- CROWN REMOVER function - HANDPIECE with joint

Intended use: Osteotome function

Medical device to be used to perform dentistry surgeries. The pulse action plastically deforms the bone in order to implant partial denture.

The **Magnetic Mallet** in **OSTEOTOME function must be used with surgery instruments supplied by Meta Ergonomica only.**

*For the Crown-Remover function you should refer to the dedicated Instruction for use you can find in Section 13 of this manual

For your own safety you should not use the device in other areas of applicability.

Attention: For your own safety you should not carry out alterations on the device.

5.3. User's qualifications

The Magnetic Mallet should be used in medical environments only and performed by a dentistry surgeon cognizant on the device use.

5.4. Use conditions:

Temperature -5°/+45°

5.5. Safety guidelines for the installation environment

The device should be used in dentist's surgery or medical environments only.

Warning: The device should be installed in dry environments only and kept dry.

Warning: The device must be only supplied with the voltage which is stated on the rating plate placed on the rear panel of the control unit.

5.6. Operating safety guidelines

Prior to your first treatment on a patient you should be cognizant on the working of this medical device.

-  **Warning:** 1 pulse/sec for a maximum of 200 consecutive pulses then a 40 minutes break.
-  **Warning:** Prior to activating the device you should always carry a functional operational test (F.O.T.).
-  **Warning:** Instruments are not provided in a sterile package and must, therefore, be sterilised before their first use (see “Disinfection and sterilization of the parts” section). Carefully read the specific Instruments Instruction for use.
-  **Warning:** Move away the device from the patient every time there is a stalling in the device working (for example: power failure).
-  **Warning:** Prior to activate the device make sure you have plugged into the device the proper handpiece.
-  **Warning:** Prior to activate the device make sure you corerctly plugged the handpiece into the corresponding plug (M) and the footswitch in the corresponding plug (N).
-  **Warning:** Make sure no liquids penetrate into the handpiece.

6.0. TRANSPORT AND STORAGE

6.1. Transport and storage conditions

-  **Caution:** Handle with care. Take care the package is not damaged.

Take care the external labelling is undamaged.

On the delivery please inspect the device into the packaging. Verify the device is undamaged and the wholeness of the supply (see “Unpacking” section).

7.0. INSTALLATION AND INSTRUCTIONS FOR USE

7.1. Unpacking

Avoid any damaging to the labelling in pulling the device and the accessories out from the package. Verify no article is missing or damaged.

Avoid any damaging to the labelling.

Should any article be missing or be damaged please immediately contact the manufacturer or the supplier.

7.2. Installation and connection

 **Warning:** Install the device in dry environments only. Take care no liquid gets into the control unit.

 **Caution:** The device should be installed and kept away from the patient.

1. The power source and the Magnetic Mallet should have the same nominal voltage. The nominal voltage is indicated on the rating panel placed on the rear panel of the control unit.
2. Connect the power cable (H) to the device through the socket (A)
3. Plug in the power cable (H) of the device
4. Jack in the proper handpiece (O) into the connector on the power supply (D)
5. Jack in the footswitch (L), into the connector on the power supply (D)


 **Caution:** Prior to the activation, take sure you've plugged in the proper handpiece according to the required use/treatment.

7.3. Activation

1. Switch on the device pushing the ON/OFF push-button (C) placed on the rear panel of the control unit (D)
2. The display (G), placed on the front panel of the control unit, lights up.


 **Attention:** When the display (G) lights up it will show a number (from 1 to 4) according to the positioning of the knob (F) on the front panel.

7.4 Surgery instrument insertion into the handpiece

1. Insert the instrument (R) into the handpiece
2. Insert the ring nut (Q) on the instrument
3. Screw the ring nut (Q)
4. Ensure the surgery instruments is tightly fastened in order to avoid any possible ejection of the instrument during surgery. This could cause severe risk for the operator, the assistant and the patient.

7.5. Release of surgery instrument

1. Unscrew the ring nut (Q)
2. Pull out the instrument (R)

7.6. How to activate the handpiece

1. To activate the handpiece, press the footswitch (L) or the manual control (E) which is placed on the front panel of the control unit (D).

Each pressure of the footswitch (L) or of the manual control (E) conveys a single pulse.

Attention: In the first stage you should always start by the lighter force (No.1).

Caution: 1 pulse/sec for a maximum of 200 consecutive pulses then a 40 minutes break

7.7. Calibration of the forces

1. Turn the knob (F) on the control unit front panel in order to calibrate the force level (from 1 to 4).
2. The display (G) shows the calibrated force level according the knob (F) positioning.

Attention: The most suitable level of the force is determined by the surgeon according to the practice/treatment.

The surgery instrument action is mainly due to the pressure exerted by the surgeon on the handpiece and by the axial forwards movement of the energy.

Attention: It's up to the surgeon to determine the pressure to be exerted on the bone in order to achieve the desired result.

8.0. MAINTENANCE

8.1. Cleaning and maintenance

 **Caution:** You should unplug the device from the main power every time you go working on the control unit (for example: cleaning, connectors unplugging).

 **Caution:** The silicone cable with the connector and the handpiece should be sterilized by autoclave before any usage.

 **Caution:** You should never put your fingers or any tool inside the connector of the handpiece or into the footswitch one.

 **Warning:** Make sure no liquids enter into the control unit. You should not hold liquids on it. The control unit case should be routinely disinfected with cleaning products which aren't harmful for plastics.
Make sure no liquid enter into the handpiece. Do not soak into liquid disinfectants.

8.2. Disinfection and sterilization of the parts

 **Caution:** You should unplug the device from the main power every time you go working on the control unit (for example: cleaning, connectors unplugging).

DISINFECTION AND STERILIZATION BY AUTOCLAVE OF THE HANDPIECE WITH ITS CORDSET AND CONNECTOR, THE RING NUT AND THE SURGERY INSTRUMENT.

The handpiece (O) with its cordset and metal connector (M), the ring nut (Q) and the surgery instrument (R) **must be submitted to autoclave sterilization before their first use and before all treatment.**

After having uncoupled the parts (see the section 7.5) you should carry out the cleaning process. Brush instruments off under flowing tap water, then proceed with the sterilization by autoclave.

For instruments cleaning refer to the Section 4.0 of the specific instruction for use.

AUTOCLAVE DIRECTIONS FOR USE:

Achievement of the dew point through the proper combination of temperature and pressure values which should be the following ones:

- 135°C for 2 bar pressure

ITEM	WASHING/DISINFECTION	STERILIZATION BY AUTOCLAVE
Handpiece with cordset and connector	Avoid any entering of liquids into the handpiece. Do not use liquid disinfectants for the handpiece	Items should be sterilized by autoclave before their first use as they're non-sterile supplied. Items should be sterilized by autoclave before all treatment
Chuck, joint	Prior to sterilization by autoclave you should disinfect the chuck and the joint. Soak a cotton bud with an hospital grade disinfectant and wipe any visible debris from both inside and outside the piece. To get a proper disinfection of the parts you should use hospital grade disinfectants. You should carefully read the instructions for use on the disinfectant packaging in order to determine the suitability of the treatment	Items should be sterilized autoclave before their first use as they're non-sterile supplied. Items should be sterilized by autoclave before all treatment
Instruments, crown remover hook, inserter, cutter	Carry out the cleaning process as per instruction contained in the specific booklet. Brush off under flowing tap water then proceed with the sterilization by autoclave. Avoid any use of IUPAC (benzalklonium chloride) to clean instruments.	Items should be sterilized by autoclave before their first use as they're non-sterile supplied. Items should be sterilized by autoclave before all treatment.

Caution:

Prior to disinfection and sterilization, carefully read the instructions contained in the Instruments Instruction for Use - Section 4.0

8.3. Maintenance

Caution: Do not open the device! Electric shock danger! Any intervention on the device by unauthorized people nullifies the warranty.

NO LUBRICATION AND/OR MAINTENANCE IS NEEDED ON THIS DEVICE **NO INTERVENTION ON THIS DEVICE IS ALLOWED.**

In the event of any failure you should immediately contact the manufacturer/distributor (see section "Failures and malfunctions"). The only maintenance you can carry out is the fuses replacement (see section "Fuses replacement").

8.4. Repairs

Caution: Do not open the device! Electric shock risk! Any opening of this device by unqualified people nullifies the warranty.

Warning: FOR YOUR OWN SAFETY YOU SHOULD NOT MAKE ALTERATIONS TO THIS DEVICE.

Caution: The returned material should be previously cleaned and disinfected. I The handpiece and instruments should be sterilized too and packed into sterile bags. The material not suitably prepared will be rejected.

IN THE EVENT OF FAILURES OR MALFUNCTIONS CALLING FOR REPAIRS YOU SHOULD PROMPTLY ADVICE THE MANUFACTURER/DISTRIBUTOR SO THAT THE DEVICE WILL BE CALLED IN.

The distributor/manufacturer will carry out the failures test and the repair works.

8.5. Fuses replacement

Warning: You should replace fuses with some new one having the amperage value indicated in the rating plate on the fuse holder.

1. Take off the fuse holder (I) from the line filter using a little screwdriver.
2. In the fuse holder you will find two 1,6 A fuses.
3. Verify if one or both fuses are burnt out.
4. Pull out the burnt out fuse/fuses.
5. Pull into the holder the new fuse/fuses. Make sure they have the same amperage value indicated in the rating plate placed on the fuse holder.
6. Verify the correct working of the device.

9.0. DISPOSAL OF WASTE

This device meets the requirements of **Directive 2012/19/UE** for electric and electronic devices disposal - **RAEE**. At the end of the period of use the device and the accessories, must be sent for recycling of the materials or for disposal in a manner which poses no threat to humans or the environment. The manufacturer is responsible for the compliance with National requirements. The manufacturer is responsible for the compliance with National requirements. For instruments refer to the Section 6.0 of the Instruments instruction for use.

9.1. Responsibility

The manufacturer is liable for safety, reliability and performances of this device only if:

- The installing has been performed in order to adhere rigorously to the instructions.
- All necessary alterations or repairs have been carried out by authorized repair service only.
- The device has been used in accordance with the instructions for use and its intended use.

9.2. Warranty

The manufacturer undertakes to provide the final customer of this device with a warranty of satisfactory functions freedom from faults in both and manufacturing process for the duration of 24 months from the delivery date. In case of justifiable complaints the manufacturer will provide repairs and/or spare parts free of charge. Nevertheless, the manufacturer will charge the final customer with shipment costs and it is not accountable for risks arising from the shipment itself. For other instances the manufacturer will refer to the warranty indicated in the trade general conditions.

Any opening, repair or alteration carried out by unauthorized persons relieves the manufacturer of all responsibility concerning the safe working of the device and nullifies the warranty.

The life cycle of the device is of 8 years. After such a time the manufacturer will not ensure the possibility to refurbish the device.

9.3. Technical Literature

The manufacturer will furnish on request circuit diagram, the components list, all descriptions and information useful to the technical assistance in order to carry the authorized repairs.

10.0. TECHNICAL DESCRIPTION

10.1. Failures and malfunctions

 **Caution:** Every time you go on working on the control unit you should unplug the device from the main power.

Failure / Malfunction

Break down without sonic alarm.

Possible cause

- Failure of electrical supply.
- Line filter fuse burnt out.
- Internal circuit failure.

Possible solution

- Verify the electrical supply presence.
- Fuse/fuses replacement (see section "Maintenance"). If the malfunction is still present after the fuse/fuses and power cable replacement contact the authorized service provider.
- Contact the authorized service provider.

Failure / Malfunction

The display shows the capital letter "E", the sonic alarm is ringing and you notice a mild heating of the handpiece.

Possible cause

- Internal circuit failure.

Possible solution

- Switch off the device.
- Contact the authorized service provider.

Failure / Malfunction

When you press the footswitch no pulse is conveyed to the instrument.

Possible cause

- The footswitch is damaged.
- Possible wrong connection of the footswitch in the corresponding plug.
- Internal circuit failure.

Possible solution

- Verify you correctly connected the footswitch in the plug.
- Verify if the pulse is conveyed to the instrument by pressing the manual control push-button. If so contact the authorized service provider and ask for a new footswitch.
- Even if the footswitch doesn't work you can follow through the treatment using the manual control pushbutton.

Failure / Malfunction

If you press the manual control push-button no pulse is conveyed to the.

Possible cause

- The manual control push-button is damaged.
- Internal circuit failure

Possible solution

- Verify if the pulse is conveyed by pressing the footswitch. If so contact the authorized service asking for the manual control switch repair/replacement.
- Contact the authorized service.

Failure / Malfunction

The display is switched off or it shows some non descript number.

Possible cause

- The display is damaged.
- Internal circuit failure.

Possible solution

- Take note that the level of the force is all the same indicated by the knob positioning (1-2-3-4). In that event you can follow through the treatment before calling the authorized service for display replacement or repairs.

Failure / Malfunction

Mild heating of the handpiece

Possible cause

- Improper use. You have exceeded the maximum threshold of 200 consecutive pulses without having had the necessary break.

Possible solution

- Switch off the device. Wait 40 minutes before switching on the device again.

Failure / Malfunction

The O-ring type 108 of the handpiece is no more working as antiunscrewing for the ring nut

Possible cause

- O-ring consumption or breaking

Possible solution

- Replace the O-ring with the spare one supplied with the handpiece (see plastic container).
- If the problem endures, contact the authorized service provider.

11.0. TECHNICAL DESCRIPTION

Supply voltage: 230 V
Frequency: 50 Hz
Current input (stand by): 0,070 A
Current input (operative condition): 0,38 A
Fuse: F - 5x20 - 250V - 1.6A

**Medical device class IIa in accordance with MDD 93/42/CEE
and MDR (UE) 2017/745**

Medical device Class II – BF type in accordance with IEC 60601-1

Dimensions 226.5 mm x 223 mm x 109 mm
(l x h x d)
Control unit weight Kg. 2,600
Handpiece weight Kg. 0,390

Transport/storage conditions

Storage temperature range -20°C / +40°C
Humidity 10% . 90%
(relative humidity, without condensation)
Air pressure: +700hPa / +1060 hPa

Conditions for use

Temperature range for use -5°C / + 45°C
Humidity 15% . 85%
(relative humidity, without condensation)
Air pressure: +700hPa / +1060 hPa

12.0. ELECTROMAGNETIC COMPLIANCE DECLARATION

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2, 4, 8, 15kV air	Compliant	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast Transient/Burst IEC 61000-4-2	±2kV for input power ports ±1kV for I/O ports	Compliant	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	Input power ports: 0.5 and 1.0kV (line to line) 0.5, 1.0 and 2.0kV (line to earth) Signal I/O: 2kV (line to earth)	Compliant	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, Short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Dips: >95% UT for 0.5 and 1 cycle 30% UT for 25 cycles (50Hz) 30% UT for 30 cycles (60 Hz) Interruption: >95% UT for 250 cycles (50Hz) >95% UT for 300 cycles (60Hz)	Compliant	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m	Compliant	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment guidance
Conducted RF IEC 61000-4-6	3 Vrms 150KHz to 80MHz	3 Vrms	<p>Portable and mobile RF communication equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance.</p> $d=1,167 * \text{sqrt}(P)$ <p>$d=1,167 * \text{sqrt}(P)$ 80 MHz to 800 MHz $d=2,2333 * \text{sqrt}(P)$ 800 MHz to 2,5 GHz</p> <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey³, should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>
Radiated RF IEC 61000-4-3	3V/m 80 MHz to 2,5 GHz	3 V/m	

Note 1: at 80 MHz and 800 MHz, the higher frequency range applies

Note 2: these guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

The equipment is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled.

The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum power of communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 KHz to 80 MHz $d=1,17\sqrt{P}$ m	80 MHz to 800 MHz $d=1,17\sqrt{P}$ m	800 MHz to 2,5 GHz $d=2,33\sqrt{P}$ m
0,01	0,117	0,117	0,233
0,1	0,370	0,370	0,740
1	1,17	1,17	2,33
10	3,70	3,70	7,40
100	11,7	11,7	23,3

For transmitters rated at maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

13.0. CROWN REMOVER HANDPIECE

13.1 Essential information

The removal of the permanent cements for crown and bridges retention entails that it's always mandatory for the clinician to previously assess the inherent difficulties of the operation.

The mechanical action needed to smash the cement should be always evaluated according to the alveolar root firmness and the tilt angle degree of the abutment. If the angle goes from 3° to 6° the operation will be easier, while if the abutment shaping didn't take care about the proper tilt angle degree, the operation will be surely harder and hazardous.

This is the reason why the clinician assessment is always though, subtle and mandatory to determine if the crown should be removed by a mechanical pulse or not and that it should be always considered the possibility to encounter the abutment fracture or removal.


A further element to consider is if the screw, in the meantime, caused peri-implantitis, if any.

The Magnetic Mallet conveys a mechanical, adjustable pulse with 4 different force levels which is kept changeless in the timing of application. This pulse, directly conveyed to the prosthetic structures, allows the smashing of the retention cement and the subsequent loosening of the crown.

The conclusion is that, in case the surgeon assesses the feasability of the surgery, the use of Magnetic Mallet with the Crown Remover Handpiece will allow the surgeon to use a device more effective than the manual technique and the automatic devices of the state of the art.

13.2 Handpiece illustrations

13.2.1 Handpiece parts

R. Nipple

U. Crown Remover hook

V. Silicone cap for the Crown Remover handpiece

LC. Handpiece

13.3 Intended use: crown remover handpiece

Medical device for dentistry to be used in crowns and bridges removal from abutments or implants.

The action of the adjustable force allows to break the retentions cement used to fix the crown or the bridge, getting its consequent loosening and removal.

The Magnetic Mallet in the CROWN REMOVER function should be equipped with hooks produced by ASA DENTAL only.

PART NUMBER:

ASA DENTAL REF. 2306 – 2307 – 2308

13.3.1 Standard supply

- 1 Handpiece with junction - Crown remover function
- 3 Crown remover hooks
- 2 Silicone caps
- 1 Replacement junction

13.3.2 Safety guidelines

For your own safety you should not use the device in other areas of applicability.

Attention: For your own safety you should not carry out alterations on the device.

Attention: The improper use of this device can cause injuries

The hook action on retention cement is mainly determined by its impact on the prosthesis and incidentally by the extraction strength exerted by the surgeon.


In order to familiarize yourself with the device, in the first stage you should always start from the lighter force (nr.1).

It's also advisable to keep the crown remover hook between your fingers selecting the different levels of the forces. In this way you shall be able to assess the corresponding pulse.

13.4 Instruction for use

13.4.1 Insertion of the hook in the handpiece

1. Screw the nipple (R) on the rod thread (Z) of the handpiece (LC)
2. Push the silicone protective cap (V) on the nipple (R) all the way to the bottom such it adheres
3. Choose the proper hook (U) and couple it to the handpiece inserting it in the nipple (R) you've assembled to the handpiece (LC)
4. Screw the crown remover hook (U) tightening it.
5. Ensure the crown remover hook is tightly fastened then press the footswitch (L) or the manual control (E)

13.4.2 Release of the crown remover hook

1. Unscrew the crown remover hook (U).
2. Pull out the crown remover hook (U).
3. Pull out the silicone protective cap (V).
4. Unscrew the nipple (R).

The hook action is mainly determined by its impact on the prosthesis and incidentally by the extraction strength exerted by the surgeon.

DISCLAIMER: It's up to the surgeon to determine the more suitable force to be applied by the hook to the prosthesis in order to achieve the desired result. In order to avoid any breakage, the surgeon should always axially pull, in comparison to the abutment position.

PRECAUTIONS:

1. The removal of prosthesis always presents a risk of damage to the prosthesis itself as well as to the patients' hard, soft and anchorage tissues.
2. Take appropriate measures to prevent swelling and/or aspiration of loosen parts and substances during treatment.
3. Crown and bridges that are cemented into place with certain type of cement, shouldn't be mechanically treated for the increased risk of breakage.
4. The practice performed with the Magnetic Mallet and its crown remover hand-piece is very gentle. Nevertheless, the stability of the root should be carefully verified in case of advanced periodontitis in order to prevent the possible extraction of the root.

13.5 Instructions for the sterilization

Sterilization by autoclave of the handpiece (LC) with its cordset and connector (M), nipple (R), silicone cap (V) and crown remover hook (U).

Attention:

The handpiece (LC) with its cordet and connector (M), nipple (R), silicone cap (V), and the crown remover hook (U) should be sterilized by autoclave **before its first use** as they're **non-sterile supplied**.

All the parts must be sterilized **before and after each treatment**.

1. After having uncoupled the parts, (see the section 13.4.2) you should carry out the cleaning process.
2. Brush hooks and nipple off under flowing tap water then proceed with the sterilization by autoclave.
3. Soak a cotton bud with a hospital degree disinfectant and wipe any visible debris from both inside and outside the nipple.
4. Clean the hooks with hospital degree disinfectant then put all the parts with the handpiece into the autoclave for sterilization.

AUTOCLAVE DIRECTIONS FOR USE:

Achievement of the dewpoint through the proper combination of temperature and pressure values which should be the following ones:

135°C – 2 Bars

13.6. Maintenance and repairs

In the event of failures or malfunctions calling for repairs, you should promptly inform the manufacturer and/or the distributor so that the device can be called in.

The handpiece doesn't require lubrication and/or maintenance.

No technical intervention or alteration should be carried out by non-authorized personnel.

Any opening of the product nullifies the warranty.

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Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2, 4, 8, 15kV air	Compliant	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast Transient/Burst IEC 61000-4-4	±2kV for input power ports ±1kV for I/O ports	Compliant	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	Input power ports: 0.5 and 1.0kV (line to line) 0.5, 1.0 and 2.0kV (line to earth) Signal I/O: 2kV (line to earth)	Compliant	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips, Short interruptions and voltage variations on power supply input lines IEC 61000-4-11	Dips: >95% UT for 0.5 and 1 cycle 30% UT for 25 cycles (50Hz) 30% UT for 30 cycles (60 Hz) Interruption: >95% UT for 250 cycles (50Hz) >95% UT for 300 cycles (60Hz)	Compliant	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m	Compliant	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

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MANUFACTURER: Meta Ergonomica Srl