MICROSCOPE KIT

INSTRUCTION MANUAL

MODEL #44121

ENGLISH

INTRODUCTION

Congratulations on your Celestron microscope purchase. Your new microscope is a precision optical instrument, made of high quality materials to ensure durability and long life. It is designed to give you a lifetime of enjoyment with a minimal amount of maintenance.

Before attempting to use your microscope, please read through the instructions to familiarize yourself with the functions and operations to maximize your experience. Refer to the microscope diagrams to locate the parts discussed in this manual.

The microscope provides magnification from 40x to 600x. It is ideally suited for examining specimen slides of yeasts and molds, cultures, plant and animal parts, fibers, bacteria, etc.

You can also examine small and thin objects at low powers such as coins, rocks, insects, various materials, etc. You can also explore the exciting microscopic world by creating your own specimen slides. An excellent book to get you started, which describes various projects and activities along with techniques, is "The World of the Microscope" by Usborne.

The final section of this manual provides simple care and maintenance tips to ensure that your microscope provides you with years of quality performance.

Note: This product is designed and intended for use by those 14 years of age and older.

Standard Accessories Included with your Microscope

FIGURE 1

SPECIFICATIONS — MODEL # 44121

Stage	Plain Stage with metal clips — 74 mm x 70 mm (2.9" x 2.8")
Zoom Eyepiece	Glass optics — power continuous from 10x to 20x
Focuser	Coarse focus — dual knobs
Objectives	All glass optics — see magnification chart for powers
Illuminator — Top	Pen light style
Illuminator — Bottom	Uses 2AA batteries (user supplied)
Nosepiece	Triple with click stop
Weight/Dimensions (with batteries)	17 oz (482 g) — 5.25" x 3.13" x 9.75" (133 mm x 79 mm x 248 mm)

MAGNIFICATION TABLE

Use the following table to determine the magnification for the different eyepiece/objective lens combination of your microscope.

Objective Lens	4x	15x	30x
10x on Zoom Eyepiece	40x	150x	300x
20x on Zoom Eyepiece	80x	300x	600x

SETTING UP YOUR MICROSCOPE

1. Take the Styrofoam container out of the carton.
2. Remove the tape from the Styrofoam container holding the various parts in place.
3. Carefully remove the microscope and other parts from the container and set them on a table, desk, or other flat surface.
4. Remove the plastic bag covering the microscope.
5. Remove the plastic cap from the zoom eyepiece Fig. 1 (1).
6. Install two AA batteries (user supplied) in the base Fig. 1 (11) of the microscope for the bottom illuminator. See Figure 3A, showing the battery compartment closed at the back of the microscope. Figure 3B shows the batteries being installed in the compartment

(see the inside of the battery compartment door to locate where the positive (+) and negative (-) ends of the batteries go). Pull out on the door to open it while holding the base firmly and push firmly to close it after the batteries are installed.

FIGURE 3A

FIGURE 3B

You are now ready to use your microscope for looking at specimen slides or small objects through the zoom eyepiece!

MICROSCOPE OPERATION

VIEWING A SPECIMEN

Carefully place a specimen slide under the stage clips Fig. 1 (8) and center the specimen directly over the hole in the center of the stage Fig. 1 (7). Figure 4A shows the stage area with the hole in the center and Figure 4B shows a specimen slide centered over the hole in the stage. It will take some experimenting to place slides or objects in the center of the stage as the image you see is upside

down and reversed, but after some usage you will have an easy time centering. Read the sections below on Focusing, Changing Power, and Illumination before proceeding.

You are now ready to focus and view the specimen, but first you must take some precautions so you do not damage a specimen slide or valuable object. When using the higher powers while you are focusing, make sure that the objective lens does not hit the slide or object being viewed.


FIGURE 4A

FIGURE 4B

FIGURE 5

FOCUSING & CHANGING POWER (MAGNIFICATION)

Now that the specimen slide (or object) is placed directly under the objective lens, use the focus knob Fig. 1 (9) to focus on the specimen. **Note: For very small objects, you should set them on the clear slide with a recessed hold in the center.**

1. Always start with the lowest power (4x objective lens) and have the zoom eyepiece at the 10x position, all the way counterclockwise (when you are facing the microscope from the front) until it stops, so that the total power is 40x – this is the most useful power for solid objects.
2. For slightly higher power, you can rotate the knurled ring on the zoom eyepiece (see Figure 5) clockwise to obtain powers of 40x to 80x (or anywhere in between) as you continue rotating to the 20x eyepiece position. Note that you will have to refocus whenever you rotate the eyepiece to obtain a sharp focus.
3. For much higher powers, you will have to rotate the nosepiece Fig. 1 (5) to change the objective lens to 15x (provides total power of 150x to 300x depending on what position you have the zoom eyepiece in or 30x provides total power of 300x to 600x). You rotate the nosepiece by holding the microscope above the nosepiece with one hand and rotating the nosepiece with the other hand until it clicks at the position. Be cautious not to let the objective lens touch the specimen slide or object when changing to higher powers. You should turn the focus knob first to lower the stage to a low position. **Note: The power range of the objective lens you are using is shown on the nosepiece after it clicks into position.**
4. At the highest powers, your views will be greatly magnified, but somewhat darker. The most enjoyable views can be at the lower powers, which have a wider field of view and brighter illumination.

FIGURE 6A


FIGURE 6B

ILLUMINATION

To get the sharpest and best views, the illumination (lighting) will have to be adjusted.

1. The top illuminator Fig. 1 (3) is used only for solid objects (not specimen slides) so that light shines down onto the object. Push the button on the top of the illuminator to turn it on/off. After some usage, you can determine the best way of adjusting the light to provide the most pleasant views.
2. The bottom illuminator Fig. 1 (10) is used for specimen slides. The light shines up through the hole in the stage through the slide. You can see a close up of the illuminator in Figure 6A below. The illuminator is turned on by rotating it so the light goes up through the hole. You turn off the illuminator by rotating it down so that the mirror is on the upper part (the mirror is not useful with this microscope since you have the much better electric illuminator).
3. The light from the bottom illuminator can be increased or decreased by rotating it with very slight movements.
4. The bottom illumination may be too bright with some specimen slides. Included with your microscope is a light diffuser, which reduces the brightness and glare somewhat and can make the views sharper with a higher contrast level. In Figure 6A, the diffuser is the small black piece. The diffuser fits over the bulb area by press fitting it on. Figure 6B shows what the illuminator looks like with the light diffuser in place. It will take experimenting to get this right.

CARE, MAINTENANCE, AND WARRANTY

Your Celestron microscope is a precision optical instrument and should be treated with care at all times. Follow these care and maintenance suggestions and your microscope will need very little maintenance throughout its lifetime.

- When you are done using your microscope, remove any specimens left on the stage.
- Turn off the top and bottom illuminators when you are done using the microscope.
- If you will not be using your microscope for a long period of time (30 days or longer), remove the batteries.
- Always place the dust cap over the eyepiece when not in use or when being stored.
- Store the microscope in a clean, dry place.
- Be very careful if using your microscope in direct sun light to prevent damage to the microscope or your eyes.
- When moving your microscope, carry it by the "arm" with one hand.
- Clean the outside surfaces with a damp cloth.
- Never clean optical surfaces with cloth or paper towels as they can scratch optical surfaces easily.
- Remove dust from optical surfaces with a camel's hair brush or air blower.
- To remove fingerprints from optical surfaces, use a lens cleaning agent and lens tissue available at most photo outlets. When cleaning, do not rub in circles as this may cause streaks or scratches to occur.
- Never disassemble or clean internal optical surfaces. This should only be done by qualified technicians at the factory or other authorized repair facilities.
- When handling specimen slides, use care as the edges can be sharp.

YOUR MICROSCOPE HAS A TWO YEAR LIMITED WARRANTY.

FOR ASSISTANCE WITH THIS PRODUCT, PLEASE VISIT CELESTRON'S ONLINE TECHNICAL SUPPORT CENTER AT

[HTTPS://WWW.CELESTRON.COM/PAGES/TECHNICAL-SUPPORT](https://www.celestron.com/pages/technical-support). HERE, YOU MAY SEARCH THROUGH A COMPREHENSIVE DATABASE OF FREQUENTLY ASKED QUESTIONS OR SUBMIT A REQUEST FOR ASSISTANCE.

FCC NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

2835 Columbia Street • Torrance, CA 90503 U.S.A.
[celestron.com/pages/technical-support](https://www.celestron.com/pages/technical-support)

Telephone: 1(800) 421-9649 • ©2019 Celestron
All rights reserved. • 02-19

Product design and specifications are subject to change without prior notification.

This product is designed and intended for use by those 14 years of age and older.

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