Harvest Poly Splint™

Acrylic milling blanks for removable & temporary bite appliances.
Fräsrohlinge zur Herstellung von herausnehmbaren Aufbisschienen
Acrílico son piezas en bruto para fresado de férulas dentales removibles.

Harvest Dental Products, LLC
905 Columbia St.
Brea, CA 92821 USA
Tel: 1+714-674-7400
Fax: 1+714-674-7402
Em: hello@harvestdental.com
Web: www.harvestdental.com

INSTRUCTIONS FOR USE – ENGLISH
INTENDED USE Harvest Poly Splint™ Acrylic milling blanks for removable & temporary bite appliances.
TECHNICAL DATA E-Module: 1.900 MPa Bending Strength: >90 MPa
COMPOSITION _ Polymethyl methacrylate (PMMA)
SHADE Clear (transparent)
HAZARD WARNINGS Product quality and durability can be compromised if handled improperly. Dust from grinding can irritate eyes, skin or mucous membranes, take protective measures.
STORAGE No special measured required.
WORKING INSTRUCTIONS Mill resin discs while taking into account parameters for milling acrylic polymers. Check milled framework for obvious chips or cracks, use only fault free framework for bite splint appliance. Ensure that object is dust free.
FINISHING For trimming and separating frameworks, use a cross-cut tungsten carbide or separating disc suitable for acrylics.
RESTORING Bite-Splint appliances can be restored and extended at any time using commercially available MMA-based, cold-curing resins.
GUARANTEE Application recommendations, no matter whether verbally, in writing, or by means of practical instructions are based on our experience and tests and are thus only a recommendation. Our products are subject to continuous further development. We reserve the right to alter the products in terms of design or composition.
SPECIAL INFORMATION Enter the LOT No. with each procedure that requires identification of the material.

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