ENGLISH

DRUG TEST

Self-Test for the detection of 6 drugs in urine:
Morphine, Amphetamine, Cocaine, Marijuana, Methamphetamine and Methadone

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Test	Cut-off (ng/mL)	Time permanence in urine after exposure
Amphetamine (AMP1000)	1000	4 to 24 hours
Cocaine (COC 300)	300	24 to 48 hours
Kokain (COC 300)	300	3 to 10 days
Marijuana	50	3 to 10 days
Methadone (MTD 300)	300	3 to 5 days
Methamphetamine (MET 500)	500	2 to 10 days
Morphine (MOP 300)	300	3 to 5 days

PRINCIPLE OF THE TEST

The DRUG TEST Panel is a rapid urine screening test that can be performed without the use of an instrument. The test utilizes monoclonal antibodies to specifically detect elevated levels of specific drugs in urine. The speed and sensitivity of the immunochromatographic techniques have made themselves the most relevant methods for drugs of abuse screening. The DRUG TEST Panel will generate a colored line in the specific test region because of the dispstick because of the region's drug-positive urine specimen will saturate all the binding sites of the antibody. Therefore, the colored line will not form in the test region's drug competition, while a drug-negative urine specimen will generate a line in the test region because of the absence of drug competition. To serve as procedural control, a colored line will always appear at the control region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

KIT CONTENT

Materials provided

1 Test cassette

1 Instructions for Use

Materials needed but not provided

Urine collection cup

Timer

- 1) The DRUG TEST Panel provides only a qualitative, preliminary analytical result.
- 2) There is a possibility that technical or procedural errors, as well as interfering substances in the urine specimen may cause a positive result.
- 3) A positive result does not indicate level or intoxication or concentration in urine.
- 4) Negative result may not necessarily indicate drug-free urine. Negative results can be obtained when drug is present but below the cut-off level of the test.
- 5) This test does not distinguish between drugs of abuse and certain medications.

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PRECAUTIONS

1. The assay is designed for use with human urine only. Fresh urine specimens should be used, and the urine does not require any special handling or pretreatment.
2. Carefully read the instructions before performing the test. The test is only interpretable if the instructions are carefully respected.
3. Keep out of the reach of children.
4. The test Panel should remain in the sealed pouch until use.
5. Do not use after the expiry date printed on the label and on the protective pouch.
6. Make sure you correctly follow the "time" instructions when carrying out the test and observing the results.
7. After use, all components can be disposed according to local regulations; all specimen should be considered potentially hazardous and handled in the same manner as an infectious agent.
8. The test kits should be kept away from direct sunlight, moisture and heat.
9. In case of difficulties in color identification (such as Daltonism), ask for help in test reading.
10. External use only. DO NOT SWALLOW.

MULTI-DROEGE RAPID TEST

DEUTSCH

A rapid self-test for the simultaneous, qualitative detection of multiple drugs and drug metabolites in human urine. The DRUG TEST is a rapid chromatographic immunoassay for the qualitative detection of multiple drugs and drug metabolites in urine. Their permanence in urine and their detectable levels are expressed in the following table:

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PRECAUTIONS

1. The test is designed for use with human urine only. Fresh urine specimens should be used, and the urine does not require any special handling or pretreatment.
2. Carefully read the instructions before performing the test. The test is only interpretable if the instructions are carefully respected.
3. Keep out of the reach of children.
4. The test Panel should remain in the sealed pouch until use.
5. Do not use after the expiry date printed on the label and on the protective pouch.
6. Make sure you correctly follow the "time" instructions when carrying out the test and observing the results.
7. After use, all components should be disposed according to local regulations; all specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
8. The test kits should be kept away from direct sunlight, moisture and heat.
9. In case of difficulties in color identification (such as Daltonism), ask for help in test reading.
10. External use only. DO NOT SWALLOW.

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
SPECIMEN HANDLING AND STORAGE

The urine specimen should be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitates should be centrifuged, filtered, or allowed to settle to obtain a clear specimen for testing...

TEST PROCEDURE

Allow the test and the urine specimen, and/or controls to reach room temperature (15–30 °C) before testing.

- 1) Bring the pouch to room temperature before opening it. Remove the test panel from the sealed pouch and use it within one hour. 2) Remove the cap. 3) With the arrow pointing toward the urine specimen, immerse the test panel vertically in the urine specimen for at least 10 to 15 seconds...

RESULTS INTERPRETATION

RESULT SHOULD BE READ AT 5 MINUTES. DO NOT INTERPRET THE RESULT AFTER 10 MINUTES. The shade of the colored lines(s) in the Test region (T) may vary.

The result should be considered negative whenever there is even a faint line. Please refer to the illustration

NEGATIVE

A colored line appears in the Control region (C) and no line appears in the Test region (T). This negative result means that the concentrations in the urine sample are below the designated cut-off levels for a particular drug tested.

POSITIVE

A colored line appears in the Control region (C) and NO line appears in the Test region (T). The positive result means that the drug concentration in the urine sample is greater than the designated cut-off for a specific drug.

NON VALID

No line appears in the Control region (C). Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for Control line failure. Read the directions again and repeat the test with a new test card. If the result is still invalid, contact your manufacturer.

Amphetamine (AMP) Amphetamine is a potent central nervous system (CNS) stimulant that is used in the treatment of attention deficit hyperactivity disorder (ADHD), narcolepsy, and obesity but it is also available on the illicit market. Amphetamines are chemically related to the human body's natural catecholamines: epinephrine and norepinephrine.

Cocaine (COC) Cocaine is a potent central nervous system stimulant and a local anesthetic. Initially, it brings about extreme energy and restlessness while gradually resulting in tremors, over-sensitivity and spasms.

Marijuana (THC) THC (delta-9-tetrahydrocannabinol) is the primary active ingredient in cannabis (marijuana). When smoked or orally administered, THC produces euphoric effects. Users have impaired short-term memory and slowed learning.

Methadone (MTD) Methadone is a narcotic analgesic prescribed for the management of moderate to severe pain and for the treatment of opiate dependence (Heroin, Vicodin, Percocet, Morphine). The pharmacology of oral methadone is very different from IV methadone.

Methamphetamine (MET) Methamphetamine is an addictive stimulant drug that strongly activates certain systems in the brain. Methamphetamine is closely related chemically to amphetamine, but the central nervous system effects of methamphetamine are greater.

Morphine (MOP) Opiate refers to any drug that is derived from the opium poppy, including the natural products, morphine and codeine, and the semi-synthetic drugs such as heroin. Opioid is more general, referring to any drug that acts on the opioid receptor.

The effects of Methamphetamine generally last 2-4 hours and the drug have a half-life of 9-24 hours in the body. Methamphetamine is excreted in the urine primarily as Amphetamine, and oxidized and deaminated derivatives.

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