
# DURONIC

Creative Consumer Electronics

## BPM120 - Blood Pressure Monitor

[duristic.com](http://duristic.com)

[duristic](https://www.youtube.com/duristic)

[@duristicmedia](https://www.instagram.com/duristicmedia)

[Duristic](https://www.facebook.com/Duristic)

[@duristicmedia](https://www.twitter.com/duristicmedia)

## Contents

|   |            |
|---|------------|
| <b>English Manual</b>                     | <b>5</b>   |
| <b>Manuel d'utilisation en français</b>   | <b>31</b>  |
| <b>Gebrauchsanleitung auf Deutsch</b>     | <b>53</b>  |
| <b>Manual de Instrucciones en Español</b> | <b>75</b>  |
| <b>Manuale d'istruzioni In Italiano</b>   | <b>97</b>  |
| <b>Instrukcja obsługi. Polski</b>         | <b>119</b> |

**English Manual**

## Contents

|   |           |
|---|-----------|
| <b>INTRODUCTION</b> .....   | <b>7</b>  |
| Safety information  |           |
| Product symbols   |           |
| Cautions  |           |
| LCD display   |           |
| Monitor components  |           |
| Box contents  |           |
| <b>BEFORE YOU START</b> .....   | <b>11</b> |
| The choice of power supply  |           |
| Installing and replacing the batteries                                    |           |
| Setting date and time   |           |
| Select the User   |           |
| <b>TAKING MEASUREMENTS</b> .....  | <b>15</b> |
| Tie the cuff  |           |
| Start a measurement   |           |
| <b>MANAGING THE RECORDS</b> .....   | <b>17</b> |
| Recall the records  |           |
| Delete the records  |           |
| <b>INFORMATION FOR USER</b> .....   | <b>19</b> |
| Tips for accurate measurements  |           |
| Maintenance   |           |
| <b>ABOUT BLOOD PRESSURE</b> .....   | <b>21</b> |
| What are systolic pressure and diastolic pressure?                        |           |
| What is the standard blood pressure classification?                       |           |
| Irregular Heartbeat Detector  |           |
| Why does my blood pressure fluctuate throughout the day?                  |           |
| Why do I get a different blood pressure at home compared to the hospital? |           |
| Is the result the same if measuring on the right arm?                     |           |
| <b>TROUBLESHOOTING</b> .....  | <b>24</b> |
| <b>SPECIFICATIONS</b> .....   | <b>25</b> |
| <b>CONTACT INFORMATION</b> .....  | <b>26</b> |
| <b>COMPLIED EUROPEAN STANDARDS LIST</b> .....                             | <b>26</b> |
| <b>AUTHORISED COMPONENTS</b> .....  | <b>27</b> |
| <b>EMC GUIDANCE</b> .....   | <b>27</b> |
| <b>WARRANTY</b> .....   | <b>28</b> |

Thank you for selecting Duronic BPM120 arm blood pressure monitor. The monitor features blood pressure measurement, pulse rate measurement and auto-saves the results recorded. This model provides you with at least two years of reliable service. Readings taken by the blood pressure monitor are equivalent to those obtained by a trained observer using the cuff and stethoscope auscultation method. This instruction manual contains important safety and care information, and provides step-by-step instruction for using the product. Please read the manual thoroughly before using the product.

**BPM120 Features:**

- 60x41mm Digital LCD display
- Holds a maximum of 60 records
- Uses inflating measuring technology

**Safety Information**

**Product Symbols**

The symbols shown below may appear in the user manual, labelling or other components. They are the requirement of standard and using.

|   |   |   |   |
|---|---|---|---|
|    | THE OPERATION GUIDE<br>MUST BE READ   |    | Type BF Applied Parts   |
| CE 0123   | CE Mark: Conforms to essential requirements of the Medical Device Directive 93/42/EEC.  |    | DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.  |
|  | Manufacturer  |  | Direct Current  |
| SN  | Serial Number   |  | Authorized Representative in the European Community   |
|  | Manufacture Date  |  | RECYCLE: Electrical waste products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice. |
|  | The Green Dot is the license symbol of a European network of industry-funded systems for recycling the packaging materials of consumer goods. |   |   |

## Cautions

This device is intended for adult use only. This device is intended for non-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on any extremities other than the arm or for any functions other than obtaining a blood pressure measurement. Do not confuse self-monitoring with self-diagnosis. This machine allows you to monitor your blood pressure. Do not begin or end medical treatment without asking a doctor for treatment advice.

If you are taking medication, consult your doctor to determine the most appropriate time of day to measure your blood pressure. Never change a prescribed medication without consulting your doctor. When the device is used to measure patients, who have common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, the best result may occur with deviation. Please consult your doctor about the result.

If the cuff pressure exceeds 40 kPa (300 mmHg), the unit will automatically deflate. Should the cuff not deflate when pressure exceeds 40 kPa (300mmHg), detach the cuff from the arm and press the START /STOP button to stop inflation.

The equipment is not AP/APG equipment and is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide. The operator shall not touch output of batteries and the patient simultaneously. To avoid measurement errors, please avoid exposing the machine to strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.

The user must check that the equipment functions safely and see that it is in proper working condition before being used. This device is contraindicated for any female who may be suspected of, or is pregnant. Besides providing inaccurate readings, the effects of this device on the fetus are unknown.

Manufacturer will make available on request circuit diagrams, component parts list etc.

This unit is not suitable for continuous monitoring during medical emergencies or operations. Otherwise, the patient's arm and fingers will become numb, swollen and even purple due to a lack of blood.

Please use the device under the environment which was provided in the user manual, otherwise, the performance and lifetime of the device will be impacted and reduced.

During use, the patient will be in contact with the cuff. The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. It will not cause any potential sensitization or irritation reaction.


Please use the accessories and detachable parts specified/authorised by Duronic. Otherwise, it may cause damage to the unit or danger to the user/patients.

The device should not need to be calibrated within the two years of reliable service. When disposing of the accessories, detachable parts, or any of the medical equipment, do so according to the local refuse/ recycling guidelines.

If you have any problems with this device, such as setting up, maintaining or using it, please contact the Customer Care Team at Duronic. Do not open or repair the device by yourself.


When cleaning, use the soft cloth to clean the whole unit. Don't use any abrasive or volatile cleaners.

### LCD Display


| SYMBOL               | DESCRIPTION                    | EXPLANATION  |
|----------------------|--------------------------------|--|
| <b>SYS</b>           | Systolic blood pressure        | High blood pressure  |
| <b>DIA</b>           | Diastolic blood pressure       | Low blood pressure   |
| <b>PUL</b><br>/min   | Pulse display                  | Pulse in beats per minute  |
| AVG                  | Average value                  | The average value of blood pressure  |
| 1/3                  | Memory                         | Indicates memory mode and memory group   |
|                      | Motion indicator               | Motion may result in an inaccurate measurement                                 |
| <b>mmHg</b>          | mmHg                           | Measurement Unit of the blood pressure (1mmHg=0.133kPa)                        |
| Lo +                 | Low battery                    | Batteries are low and need to be replaced                                      |
|                      | Irregular heartbeat            | Blood pressure monitor is detecting an irregular heartbeat during measurement. |
|                      | Blood pressure level indicator | Indicate the blood pressure level  |
| D M Y<br>00/00/00 PM | Current Time                   | Year/Month/Day, Hour/Minute  |
|                      | Heartbeat                      | Blood pressure monitor is detecting a heartbeat during measurement.            |
|                      | User 1                         | Start measurement for User 1   |
|                      | User 2                         | Start measurement for User 2   |

## Monitor Components

## Box Contents


1. Blood Pressure Monitor

2. Cuff (22~42cm)  
(Type BF applied part)

3. 4x AAA Batteries


4. User Manual

## Choice of Power Supply

1. Battery powered mode: 6VDC 4\*AAA batteries
2. AC adaptor powered mode: 5V  $\overline{\text{---}}$  1A (not included)  
Corresponding AC adapters with USB cables are widely available.

If in doubt, contact Duronic for advice.


**Caution:** In order to get the best effect and protect your monitor, please use the right battery and special power adapter. Do not use different types of batteries together.

## Installing and Replacing the Batteries

1. Slide off the battery cover.
2. Install the batteries by matching the correct polarity, as shown.
3. Replace the cover.

Replace the batteries whenever the below happens:

- The  symbol shows
- The display dims
- The display does not light up

**Caution:** Remove batteries if the device is not likely to be used for some time. Old batteries are harmful to the environment - do not dispose of with other daily rubbish. Remove the old batteries from the device and follow your local recycling guidelines. Do not dispose of batteries in fire - batteries may explode or leak.

## Setting Date and Time

It is important to set the clock before using your blood pressure monitor, so that a time stamp can be assigned to each record that is stored in the memory.


(year: 2018—2058 time:24H/12 H)

1. When the monitor is off, press the “SET” button and hold for 3 seconds to enter the mode for year setting.

2. Press the “MEM” repeatedly to change the [YEAR].


3. When you get the correct year, press “SET” to save and it will turn to next step automatically.

4. Next, repeat steps 2 and 3 to set the [MONTH] and [DAY].


5. Repeat step 2 and 3 to set the [TIME FORMAT] between 12h and 24h.

6. Repeat steps 2 and 3 to also set the [HOUR] and [MINUTE].


6. After the [MINUTE] is set, the screen will show “DONE”, it will show you the date/time you have set and then it will turn off automatically.


## Select the User


1. When the monitor is off, press and hold the “MEM” button to enter the user setting mode. The user ID will blink showing the User 1 symbol.

2. Then press “MEM” button again to select between User 1 and User 2.

3. After selecting the suitable user ID, press the “SET” button to confirm and save. It will display the selected user with the word “DONE”. It will then turn off automatically.


## Tie the Cuff

Remove all jewellery such as watches, bracelets from your left arm. If your doctor has diagnosed you with poor circulation in your left arm, use your right arm instead.


1. Roll or push up your sleeve to expose the skin. Make sure your sleeve is not too tight.

2. Wrap the cuff around your upper arm, then position the tube off-centre toward the inner side of arm, in-line with the little finger.

3. The cuff should be snug but not too tight. You should be able to insert one finger between the cuff and your arm.

4. Sit comfortably with your left arm resting on a flat surface. Place your elbow on a table so that the cuff is at the same level as your heart. Turn your palm upwards. Sit upright and take 5-6 deep breaths.


### ***Tips for Accurate Measurements***


- Rest for 5 minutes before measuring.
- Wait at least 3 minutes between measurements. This allows your blood circulation to recover.
- For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a doctor.
- Sit comfortably. Do not cross your legs and keep your feet flat on the ground. Keep your back against the backrest of the chair.

## Start a Measurement


1. When the monitor is off, press the “START/STOP” button to turn on the monitor. It will start and finish the whole measurement automatically.


LCD Display when switched on.

Adjusts the zero automatically.

Inflating and measuring automatically.


Displays and saves the results automatically.

2. If you need to turn off the monitor during the measurement process, press the “START/STOP” button to turn it off, otherwise it will turn off automatically within 1 minute when it completes the measurement.





## Recall the Records

1. When the monitor is off, press the “MEM” button to show the most recent record.

2. Press the “MEM” or “SET” buttons to move through the historic records to find the record you want.

The record will be shown in this order: first showing the record number, followed by the date and then the time.


**Caution:** The most recent record (1) is shown first. Each new measurement is assigned to the first (1) record. All other records are pushed back one digit (e.g., 2 becomes 3, and so on), and the last record (60) is dropped from the list.


## Delete the Records

If you did not get the correct measurement, you can delete either single records individually, or delete the whole set of records altogether.


### A: To delete one individual record:

1. Press the "MEM" button once to enter Memory Recall Mode. Press the "SET" button to find the record you wish to delete.


2. Press and hold down the "SET" button until the screen above shows, with the words "del" and "one".


3. Press and hold down the "START/ STOP" button until the screen above shows, with the words "del" and "done".


### B: To delete all of the records:


1. Press the "MEM" button once to enter Memory Recall Mode. The above screen will show.

2. Press the "MEM" and hold down until the above screen shows, with the words "del" and "all".

3. Press and hold down the "START/ STOP" button until the screen above shows, with the words "del" and "done".


4. When there is no record, the screen display above will show on the monitor.

Note: To exit out the delete mode press "Start/Stop" button shortly.

## Tips for Accurate Measurements

Measurements may be inaccurate if taken during the following circumstances:

Within 2 hours of eating or drinking.


Immediately after drinking tea or coffee, or after smoking.

Within 20 minutes of taking a bath.

When talking or moving your fingers.


In a very cold environment.

When you want to discharge urine.

## Maintenance


In order to get the best performance, please follow the below instructions.


Store in a dry place and avoid direct sunlight.


Avoid the monitor touching water. If the machine gets wet, dry immediately afterwards.

Avoid intense shaking and collision.

Avoid dusty environments and unstable temperatures.

Use a wet (well wrung) cloth to remove any dirt.


Do not clean the cuff with water and never immerse the cuff in water or any other liquid.


## What are systolic pressure and diastolic pressure?

When heart ventricles contract and pump blood out of the heart, blood pressure reaches its maximum value, the highest pressure in the cycle, which is known as systolic pressure. When the ventricles relax between heartbeats, the lowest blood pressure is known as diastolic pressure.

## What is the standard blood pressure classification?

The blood pressure classification mode by World Health Organization (WHO) and International Society of Hypertension (ISH) in 1999 is illustrated below:

| Level \ Blood Pressure (mm Hg) | Optimal | Normal  | High-normal | Mild    | Moderate | Severe |
|--------------------------------|---------|---------|-------------|---------|----------|--------|
| SYS                            | <120    | 120-129 | 130-139     | 140-159 | 160-179  | ≥180   |
| DIA                            | <80     | 80-84   | 85-89       | 90-99   | 100-109  | ≥110   |

**Caution: Only a doctor can tell you your normal blood pressure range and the point at which you are at any risk. Contact your doctor for advice if the measurements taken with these monitor falls outside the range.**

## Irregular Heartbeat Detector

An irregular heartbeat is detected when a heartbeat rhythm varies while the monitor is measuring the systolic pressure and diastolic pressure. During each measurement, the monitor will keep a record of all the pulse intervals and calculate the average value of them.

If there are two or more pulse intervals, the difference between each interval and the average is more than the average value of +25% then the irregular heartbeat symbol will appear on the display with the measurement result.

If there are four or more pulse intervals with the difference between each interval and the average is more than the average value of +15%, then the irregular heartbeat symbol will appear on the display with the measurement result.

**Caution:** The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

## Why does my blood pressure vary throughout the day?

Individual blood pressure varies naturally multiple times during the day. It is also affected by the way you tie your cuff around your arm and the position of your arm during measurement. For this reason, we advise you take measurement at the same time of day and in the same position.

Please note: If the user takes medication the pressure will vary more.

*If you wish to measure a second time, wait at least 3 minutes before taking another measurement.*

## **Why do I get a different blood pressure at home compared to the hospital?**

Blood pressure is different even throughout the day due to weather, emotions, exercise etc.

Also, there is the “white coat” effect, which means blood pressure usually increases in clinical settings due to anxiety or nerves.

What you need to pay attention to when you measure your blood pressure at home:


- If the cuff is tied properly.
- If the cuff is too tight or too loose.
- If the cuff is tied on the upper arm.
- If you feel anxious.

Advice: Take 2-3 deep breaths before beginning a measurement as this will help get a more accurate result. If you are particularly anxious or upset, relax yourself for 4-5 minutes until you calm down.

## **Is the result the same if measuring on the right arm?**

It is fine to use the cuff on either arm, but some people will experience a difference in measurement results.

For this reason, we recommend using the same arm every time.

## Troubleshooting

This section includes a list of error messages and frequently asked questions for problems you may encounter with your blood pressure monitor. If the monitor is not operating as you think it should, check here before arranging for servicing.

| Problem                | Symptom  | Check This  | Remedy  |
|------------------------|--|---|---|
| <b>No power</b>        | Display is dim or will not light up.   | Batteries are exhausted.  | Replace with new batteries.   |
|                        |  | Batteries are inserted incorrectly.   | Insert the batteries correctly.   |
|                        |  | AC adapter is inserted incorrectly.   | Insert the AC adapter tightly.  |
| <b>Low batteries</b>   | <br>Shows on the display. | Batteries are low.  | Replace with new batteries.   |
| <b>Error message</b>   | 'E01' shows  | The cuff is not secure, either too tight or too loose.                              | Refasten the cuff and then measure again.   |
|                        | 'E02' shows  | The monitor has detected motion, talking, or the pulse is too poor while measuring. | Relax for a moment and then measure again.  |
|                        | 'E03' shows  | The measurement process does not detect the pulse signal.                           | Loosen the clothing on the arm and then measure again.  |
|                        | 'E04' shows  | The treatment of the measurement failed.  | Relax for a moment and try again.   |
|                        | 'EExx' shows   | A calibration error occurred.   | Retake the measurement. If the problem persists, contact Duronic for further assistance. Refer to the warranty for contact information. |
| <b>Warning message</b> | 'out' shows  | Out of range measurement.   | Relax for a moment. Refasten the cuff and then measure again. If the problem persists, contact your doctor.                             |

## Specifications

|   |  |
|---|--|
| Power supply                              | Battery powered mode:<br>6VDC 4*AAA batteries<br>AC adaptor powered mode: 5V  1A<br>Can be supplied by AC adaptor model<br>(Not Included) |
| Display mode                              | Digital LCD V.A.60*41mm  |
| Measurement mode                          | Oscillographic testing mode  |
| Measurement range                         | Rated cuff Pressure 0mmHg~299mmHg(0kPa<br>39.9kPa) Measurement pressure: SYS:<br>60mmHg~230mmHg (8.0kPa~30.7kPa) DIA:<br>40mmHg~130mmHg (5.3kPa~17.3kPa)<br>Pulse value: (40~199) beat/minute                              |
| Accuracy                                  | Pressure:<br>5°C~40°C within $\pm 0.4$ kPa(3mmHg)<br>pulse value: $\pm 5\%$  |
| Normal working condition                  | Temperature: +5°C to +40°C. Relative<br>humidity range of 15% to 90%. Atmospheric<br>pressure: 700 hPa to 1060hPa  |
| Storage & transportation<br>condition     | Temperature: -20°C to 60°C<br>Relative humidity: $\leq 93\%$<br>Atmospheric pressure: 50kPa to 106kPa  |
| Measurement perimeter of<br>the upper arm | About 22cm~42cm  |
| Weight                                    | Approx.162g(Excluding the dry cells)   |
| External dimensions                       | Approx.93x130x32,5mm   |
| Attachment                                | 4*AAA batteries,user manual  |
| Mode of operation                         | Continuous operation   |
| Degree of protection                      | Type BF applied part   |
| Protection against ingress of<br>water    | IP21   |
| Software Version                          | A01  |

## Contact Information

### Distributed by:

Company: SHINEMART LTD

Address: Unit 1, Spilsby Road, Harold Hill, Romford, Essex, RM3 8SB, UK.

### Manufactured by:

Company: Guangdong Transtek Medical Electronics Co., Ltd.

Address: Zone A, No.105 ,Dongli Road, Torch Development District, Zhongshan, 528437, Guangdong, China.

### Authorized European Representative:

Company: MDSS - Medical Device Safety Service GmbH Address: Schiffgraben 41, 30175 Hannover, Germany.


## Complied European Standards List

|  |  |
|--|--|
| <b>Risk Management</b>                 | EN ISO 14971:2012 / ISO 14971:2007 Medical devices - Application of risk management to medical devices   |
| <b>Labelling</b>                       | EN ISO 15223-1:2016 / ISO 15223-1:2016 Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. Part 1 : General requirements   |
| <b>User Manual</b>                     | EN 1041:2008 +A1:2013 Information supplied by the manufacturer of medical devices  |
| <b>General Requirements for Safety</b> | EN 60601-1:2006+A1:2013/ IEC 60601-1:2005+A1:2012 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance<br>EN 60601-1-11:2015/ IEC 60601-1-11:2015 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment                      |
| <b>Electromagnetic Compatibility</b>   | EN 60601-1-2:2015/ IEC 60601-1-2:2014 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests  |
| <b>Performance Requirements</b>        | EN ISO 81060-1:2012 Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type<br>EN 1060-3:1997+A2:2009 Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems<br>IEC 80601-2-30:2009+A1:2013 Medical electrical equipment- Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers |
| <b>Clinical Investigation</b>          | EN 1060-4:2004 Non-invasive sphygmomanometers - Part 4: Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers<br>ISO 81060-2:2013 Non-invasive sphygmomanometers - Part 2: Clinical validation of automated measurement type  |
| <b>Usability</b>                       | EN 60601-1-6:2010+A1:2015/IEC 60601-1-6:2010+A1:2013 Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability<br>IEC 62366-1:2015 Medical devices - Part 1: Application of usability engineering to medical devices   |
| <b>Software Life-cycle Processes</b>   | EN 60601-1-6:2010+A1:2015/IEC 60601-1-6:2010+A1:2013 Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability<br>IEC 62366-1:2015 Medical devices - Part 1: Application of usability engineering to medical devices   |
| <b>Bio-compatibility</b>               | ISO 10993-1:2009 Biological evaluation of medical devices- Part 1: Evaluation and testing within a risk management process<br>ISO 10993-5:2009 Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity<br>ISO 10993-10:2010 Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization   |

## Authorised Components

Please use the Duronic authorised power adapter.

(Not Included)

Adapter  
 Input: 100~240V, 50~60Hz,400mA  
 Output: 5V  $\square$  1000mA  
 Type: BLJ06L0501000U-V  
 BLJ06L0501000U-B

## EMC Guidance

This medical equipment monitor is suitable for home healthcare environments.

Warning: Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with any other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antennas) should be used no closer than 30cm (12 inches) to any part of the equipment TMB-1872-A, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

**Table 1 Guidance and MANUFACTURER's declaration – ELECTROMAGNETIC EMISSIONS- for all ME EQUIPMENT and ME SYSTEMS**

| Guidance and manufacturer's declaration - electromagnetic emissions |             |
|---|-------------|
| Emissions test  | Compliance  |
| RF emissions<br>CISPR 11  | Group 1     |
| RF emissions<br>CISPR 11  | Class [ B ] |
| Harmonic emissions<br>IEC 61000-3-2                                 | Class A     |
| Voltage fluctuations/<br>flicker emissions<br>IEC 61000-3-3         | Comply      |

**Table 2 Guidance and MANUFACTURER's declaration – electromagnetic IMMUNITY – for all ME EQUIPMENT and ME SYSTEMS**

| <b>Guidance and manufacturer's declaration – electromagnetic Immunity</b>                              |  |   |
|--|--|---|
| <b>Immunity Test</b>   | <b>IEC 60601-1-2 Test level</b>  | <b>Compliance level</b>   |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2   | ±8 kV contact<br>±2 kV, ±4kV, ±8 kV, ±15 kV air  | ±8 kV contact<br>±2 kV, ±4kV, ±8 kV,<br>±15 kV air  |
| Electrical fast transient/burst<br>IEC 61000-4-4   | ±2 kV for power supply lines<br>±1 kV signal input/output<br>100 kHz repetition frequency  | ±2 kV for power supply lines<br>±1 kV signal input/output<br>100 kHz repetition frequency   |
| Surge<br>IEC61000-4-5  | ±0.5 kV, ±1 kV differential mode<br>±0.5 kV, ±1 kV,±2 kV common mode   | ±0.5 kV, ±1 kV differential mode<br>±0.5 kV, ±1 kV,±2 kV common mode  |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11 | 0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°.0 % UT;<br>1 cycle and 70 % UT; 25/30 cycles;<br>Single phase: at 0°.0 % UT; 250/300 cycle | 0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°.<br>0 % UT; 1 cycle and 70 % UT;<br>25/30 cycles; Single phase: at 0°.<br>0 % UT; 250/300 cycle |
| Power frequency magnetic field<br>IEC 61000-4-8  | 30 A/m<br>50Hz/60Hz  | 30 A/m<br>50Hz/60Hz   |
| Conducted RF<br>IEC61000-4-6   | 3 V<br>0,15 MHz – 80 MHz<br>6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz z<br>80 % AM at 1 kHz   | 3 V<br>0,15 MHz – 80 MHz<br>6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz z<br>80 % AM at 1 kHz  |
| Radiated RF<br>IEC61000-4-3  | 10 V/m<br>80 MHz – 2,7 GHz<br>80 % AM at 1 kHz   | 10 V/m<br>80 MHz – 2,7 GHz<br>80 % AM at 1 kHz  |
| NOTE U <sub>T</sub> is the a.c. mains voltage prior to application of the test level.                  |  |   |

**Table 3 Guidance and MANUFACTURER’s declaration – electromagnetic IMMUNITY – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING**

| Guidance and manufacturer’s declaration - electromagnetic Immunity   |                      |            |  |                                  |                |              |                           |
|--|----------------------|------------|--|----------------------------------|----------------|--------------|---------------------------|
| Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment) | Test Frequency (MHz) | Band (MHz) | Service  | Modulation                       | Modulation (W) | Distance (m) | IMMUNITY TEST LEVEL (V/m) |
|  | 385                  | 380-390    | TETRA 400  | Pulse modulation b) 18Hz         | 1.8            | 0.32         | 7                         |
|  | 450                  | 430-470    | GMRS 460, FRS 460  | FM c) ± 5kHz deviation 1kHz sine | 20             | .3           | 28                        |
|  | 710                  | 704-787    | LTE Band 13, 17  | Pulse modulation b) 217Hz        | 0.20           | .3           | 9                         |
|  | 745                  |            |  |                                  |                |              |                           |
|  | 780                  |            |  |                                  |                |              |                           |
|  | 810                  | 800-960    | GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5         | Pulse modulation b) 18Hz         | 2              | 0.32         | 8                         |
|  | 870                  |            |  |                                  |                |              |                           |
|  | 930                  |            |  |                                  |                |              |                           |
|  | 1720                 | 1700-1990  | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4,25; UMTS | Pulse modulation b) 217Hz        | 2              | 0.32         | 8                         |
|  | 1845                 |            |  |                                  |                |              |                           |
|  | 1970                 |            |  |                                  |                |              |                           |
|  | 2450                 | 2400-2570  | Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7           | Pulse modulation 217 Hz          | 2              | 0.32         | 8                         |
|  | 5240                 | 5100-5800  | WLAN 802.11 a/n  | Pulse modulation 217 Hz          | 0.2            | 0.39         |                           |
|  | 5500                 |            |  |                                  |                |              |                           |
|  | 5785                 |            |  |                                  |                |              |                           |

# Warranty

1 Year Guarantee from Shine-Mart Ltd, Trademark owners of Duronic Products.

NOTE: THESE STATEMENTS DO NOT AFFECT YOUR STATUTORY RIGHTS AS A CONSUMER

This product has been manufactured under the strictest or quality control procedures, and using the highest quality of materials, to ensure excellent performance and reliability. It will give very good and long lasting service, provided it is properly used and maintained. The product is guaranteed for 1 years from the date of original purchase. If any defect arises due to a faulty materials or poor workmanship, the faulty product must be returned to the original place of purchase. Refund or replacement is at the discretion of that company.

Duronic Products are offered with a 1 year guarantee under the following conditions:

1. The product must be returned to the retailer with original proof of purchase.
2. The product must be installed and used in accordance with the instructions contained in this operating manual.
3. It must be used only for domestic purposes.
4. It does not cover wear and tear, malicious damage, misuse, neglect, inexpert repairs or consumable parts.
5. Shine-Mart Ltd has no responsibilities for accidental or consequential loss or damage.
6. Shine-Mart Ltd are not liable to carry out any type of servicing work, under the guarantee
7. Valid in the UK only

CE 0123

Information on waste Disposal for Consumers of  
Electrical & Electronic Equipment.

This mark on a product and/or accompanying documents indicates that when it is to be disposed of, it must be treated as Waste Electrical & Electronic Equipment, (WEEE).

Any WEEE marked products must not be mixed with general household waste, but kept separate for the treatment, recovery and recycling of the materials used. For proper treatment, recovery and recycling; please take all WEEE marked waste to your Local Authority Civic waste site, where it will be accepted free of charge.

If customers dispose of Waste Electrical & Electronic Equipment correctly, they will be helping to save valuable resources and preventing any potential negative effects upon human health and the environment, of any hazardous materials that the waste may contain.

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1 Rok Gwarancji firmy Shine-Mart LTD, marki handlowej Duronic Products.
5. Shine-Mart Ltd. nie ponosi odpowiedzialności za przypadkowe albo celowe zgubienie produktu lub jego uszkodzenie.
6. Firma Shine-Mart Ltd nie jest zobowiązana do przeprowadzania jakichkolwiek prac serwisowych w ramach gwarancji.