Blood Pressure Monitor BPM490

PLEASE READ ALL THE INSTRUCTIONS BEFORE USE

duronic.com

[duronic](https://www.youtube.com/duronic)

[@duronicmedia](https://www.instagram.com/duronicmedia)

[Duronic](https://www.facebook.com/Duronic)

[@duronicmedia](https://www.twitter.com/duronicmedia)

Contents

English Manual	5
Manuel d'utilisation en français	33
Gebrauchsanleitung auf Deutsch	61
Manual de Instrucciones en Español	89
Manuale D'istruzioni In Italiano	117

English Manual

Contents

INTRODUCTION

Safety Information
LCD Display Signal
Monitor Components

BEFORE YOU START

The Choice of Power Source
Setting Date, Time and Measurement Unit
Pair-up the Blood Pressure Monitor with Your Device

MEASUREMENT

Tie the cuff
Start a measurement

THE OPERATION OF RECALLING RECORDS

Recall the Records
Delete the Records
Data Transmission

INFORMATION FOR USER

Tips for operation
Maintenance

ABOUT BLOOD PRESSURE

What are systolic pressure and diastolic pressure?
What is the standard blood pressure classification?
Why does my blood pressure fluctuate throughout the day?
Why do I get a different blood pressure at home compared to the hospital?
Is the result the same if measuring both arms?

TROUBLESHOOTING

SPECIFICATIONS

AUTHORIZED COMPONENT

CONTACT INFORMATION


COMPILED EUROPEAN STANDARDS LIST

EMC GUIDANCE

Thank you for selecting our arm type blood pressure monitor. The monitor features blood pressure measurement, pulse rate measurement and auto-save the result. The design provides you with two years of reliable service. Readings taken by the blood pressure monitor are equivalent to those obtained by a trained observer using the cuff and stethoscope auscultation method. This manual contains important safety and care information, and provides step by step instruction for using the product. Read the manual thoroughly before using the product.

Safety Information

The below signs might be in the user manual, labeling or other component. they are the requirement of standard and using.

	THE OPERATION GUIDE MUST BE READ		Type BF Applied Parts
	CE Mark: conforms to essential requirements of the Medical Device Directive 93/42/EEC.		DISPOSAL: Do not dispose this product as unsorted municipal waste. Collection of such waste separately for special treatment is necessary.
	Manufacturer		Direct Current
SN	Specifies Serial Number		Authorized Representative in the European Community
	Manufacture Date		Symbol for "Including RF transmitter"

CAUTION

This device is intended for adult use only.

This device is intended for non-invasive measuring and monitoring of arterial blood pressure. It is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.

Do not confuse self-monitoring with self-diagnosis. This unit allows you to monitor your blood pressure. Do not begin or end medical treatment without asking a physician for treatment advice. If you are taking medication, consult your physician to determine the most appropriate time to measure your blood pressure. Never change a prescribed medication without consulting your Physician.

When the device was used to measure patients who have common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation, the best result may occur with deviation. Please consult your physician about the result.

If the cuff pressure exceeds 40 kPa (300 mmHg), the unit will automatically deflate. Should the cuff not deflate when pressures exceeds 40 kPa (300 mmHg), detach the cuff from the arm and press the corresponding user button to stop inflation.

The equipment is not AP/APG equipment and not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.

The operator shall not touch output of batteries and the patient simultaneously.

To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.

The user must check that the equipment functions safely and see that it is in proper working condition before being used.

This device is contraindicated for any female who may be suspected of, or is pregnant. Besides providing inaccurate readings, the affects of this device on the fetus are unknown.

Manufacturer will make available on request circuit diagrams, component parts list etc.

This unit is not suitable for continuous monitoring during medical emergencies or operations.

Otherwise, the patient's arm and fingers will become anaesthetic, swollen and even purple due to a lack of blood.

Please use the device under the environment which was provided in the user manual. Otherwise, the performance and lifetime of the device will be impacted and reduced.

During use, the patient will be in contact with the cuff. The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. It will not cause any potential sensitization or irritation reaction.

Please use ACCESSORIES and detachable partes specified/ authorised by MANUFACTURE.

Otherwise, it may cause damage to the unit or danger to the user/patients.

The device doesn't need to be calibrated within two years of reliable service.


Please dispose of ACCESSORIES, detachable parts, and the ME EQUIPMENT according to the local guidelines.

If you have any problems with this device, such as setting up, maintaining or using, please contact the SERVICE PERSONNEL of SHINEMART. Don't open or repair the device by yourself.

Please report to SHINEMART if any unexpected operation or events occur.


Please use the soft cloth to clean the whole unit. Don't use any abrasive or volatile cleaners.

LCD display signal

SYMBOL	DESCRIPTION	EXPLANATION
SYS	Systolic Blood Pressure	High blood pressure
DIA	Diastolic Blood Pressure	Low blood pressure
Pul/min	Pulse	beat/minute
	User 1	Start measurement for user 1 and transmit the measuring result automatically.
	User 2	Start measurement for user 2 and transmit the measuring result automatically.
	Data Pending to Transmit	Measurement data stored in the equipment.
	Data Transmitting	Data Transmission succeeds.
	Shocking reminder	Shocking will result in inaccurate
	Low Battery	Low battery and please replace the batteries.
mmHg	Unit	Measurement unit of blood pressure
kPa	Unit	Measurement unit of blood pressure
	Memory	Display the serial number of the measurement
	Current Time	Month:Day (Hour:Minute)
	Deflating	Exhaust the air in the cuff
	Irregular Heartbeat Detector	Please refer to Page 15 for more information

Monitor Components

List


1. Blood Pressure Monitor

2. Cuff (Type BF applied part) (22cm~32cm)

3. 4*AA Batteries


4. User Manual

The Choice of Battery

1. Battery powered mode:
6VDC 4×AA alkaline batteries
2. AC adaptor powered mode: 6V==1A

(Can be supplied by AC adaptor model UE08WCP-060100SPA !)(Not Included) Please unplug the adaptor to depart from the using utility power.

Note:

The adaptor socket is located on the right side of the monitor. Make sure the socket is not obstructed.

Caution: In order to achieve the best performance and protect your monitor, please use the authorized / specified battery and power adaptor.

The Installment And Replacement Of Battery

1. Open the battery door.
2. Insert the batteries according to the polarity indications.
3. Close the battery door.


Replace the batteries under following circumstances:

displays on the LCD.

The LCD display dims.

When powering on the monitor, the LCD does not light up.


Caution

Remove batteries if the device is not likely to be used for some time. Worn batteries are harmful to the environment. Do not dispose with daily garbage. Remove the old batteries from the device following your local recycling guidelines. Do not dispose of batteries in fire. Batteries may explode or leak.


Setting the date, time and measurement unit

Please proceed to time setting before your initial use so as to ensure each piece of record are labeled with a time stamp. (Default Value: Date: 01-Jan-2012; Time: 00:00; Measurement Unit: mm Hg)


1. When the monitor is OFF, press and hold "User 1" button to enter Time Setting Mode.

2. As pictured in the right, the blinking numeral "2012" representing [YEAR]. Press "Query" button to change the numeral. Each press will increase the numeral by one in a cycling manner.




3. Press "User 1" button again to confirm [YEAR]. Then the numeral representing [MONTH] blinks.


4. Repeat step 2 and 3 to confirm [MONTH] and [DAY].

5. Repeat step 2 and 3 to confirm [HOUR] and [MINUTE].

Hour Format: 24-Hour

6. Repeat step 2 and 3 to confirm the measurement unit.


7. After confirming the measurement unit, the LCD will display "dONE" and the monitor will shut off automatically.


Pair-up the Blood Pressure Monitor with Your Device

1. Download the free app "TransHealth" on the App Store or Google Play Store and install it on your phone or tablet.
2. Turn on Bluetooth and the app. Make sure both are ON when pair-up is proceeding.


When the monitor is OFF, press and hold the  button (User 2) to start pair-up. The symbol  and the symbol  will be shown on the LCD alternatively, indicating pair-up is proceeding.

If SUCCESSFUL, symbol  will be shown on the LCD.


If UNSUCCESSFUL, symbol  be shown on the LCD.

3. The monitor will shut off automatically after Pair-up process is complete.


Tie The Cuff

1. Tie the cuff on your upper arm, then position the tube off-center toward the inner side of arm in line with the little finger.

2. The cuff should be snug but not too tight. You should be able to insert one finger between the cuff and your arm.

3. Sit comfortably with your left arm resting on a flat surface.

Rest for 5 minutes before measuring.

Wait at least 3 minutes between measurements. This allows your blood circulation to recover.

For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.


Start The Measurement

Press the  (User 1) to turn on the monitor and it will finish the whole measurement automatically, and then save the measure data for User 1. The same with  (User 2).


In this manual, take the measurement of user 1  as an example.

1. Press the  to turn on the monitor.

LCD Display


Adjust to zero automatically.

Inflating and measuring automatically.

Display and save the results.
The data transmission will proceed automatically.

2. Press the  to power off , otherwise it will turn off automatically within one minute after last operation.

TIPS

A. When the whole measurement is complete, press the other user button,the monitor will start measuring again.


B. Maximum 60 records are both for user 1 and user 2.

Recalling the records

1. When the monitor is OFF, press “Q” button to access the memory.


2. The LCD will display the latest measuring result of the user ID which completes the last measurement.

* The record number, measuring date and measuring time will be displayed alternatively.


3. Press “Q” button to rotate the history records.

4. When in the memory mode, press the  (User 1) to recall the measurement history of User 1, or press the  (User 2) button to recall the measurement history.

5. When no history stored for the specific user in the monitor, press “Q” button and the LCD will display as pictured to the right.

CAUTION

The most recent record (1) is shown first. Each new measurement is assigned to the first (1) record. All other records are pushed back one digit (e.g., 2 becomes 3, and so on), and the last record (60) is dropped from the list.



Delete the Records

1. When under the query mode, press and hold “ Q ” button for 3 seconds to clear the memory.

2. When the LCD displays “dEL ALL”, press “ Q ” button to confirm.

3. The LCD will display “dEL dOnE” and then shut off automatically.

4. If you wish to stop clearing the memory, you may press the other button, rather than “ Q ” button to turn off the monitor, or wait until the monitor shuts off automatically.

Tips For Measurement

Measurements may be inaccurate if taken in the following circumstances.

Wait at least 1 hour after dinner or drinking

Immediate measurement after tea, coffee, smoking

Wait at least 20 minutes after taking a bath

When talking or moving your fingers


In a very cold environment

When you need the toilet

Maintenance

In order to get the best performance, please follow the below instructions.

Put in a dry place and
avoid the sunshine

Avoid immersing it in the water.
Clean it with a dry cloth in case.

Avoid intense shaking and collisions

Avoid dusty environment and
unstable temperature surrounding.


Use the slightly damp cloth
to remove the dirt.

Avoid washing the cuff.

What are systolic pressure and diastolic pressure?


When ventricles contract and pump blood out of the heart, blood pressure reaches its maximum value, the highest pressure in the cycle is known as systolic pressure. When the heart relaxes between heartbeats, the lowest blood pressure is diastolic pressure.

What is the standard blood pressure classification

The blood pressure classification published by World Health Organization (WHO) and International Society of Hypertension (ISH) in 1999 is as follows:

Level		Normal			Mild	Moderate	Severe
Blood Pressure (mmHg)							
SYS		< 120	120~129	130~139	140~159	160~179	≥ 180
DIA		< 80	80~84	85~89	90~99	100~109	≥ 110

Only a physician can tell your normal BP range. Please contact a physician if your measuring result falls out of the range. Please note that only a physician can tell whether your blood pressure value has reached a dangerous point.

Irregular Heartbeat Detector

An irregular heartbeat is detected when a heartbeat rhythm varies while the unit is measuring the systolic and diastolic blood pressure. During each measurement, this equipment records the heartbeat intervals and works out the standard deviation. If the calculated value is larger than or equal to 15, the irregular heartbeat symbol appears on the symbol when the measurement results are displayed.

Caution:

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. Please note that the device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

Why does my blood pressure fluctuate throughout the day?

1. Individual blood pressure varies multiple times everyday. It is also affected by the way you tie your cuff and your measurement position, so please take the measurement under the same conditions.
2. If the person takes medicine, the pressure will vary more.
3. Wait at least 3 minutes for another measurement.

Why do I get a different blood pressure at home compared to the hospital?

The blood pressure is different even throughout the day due to weather, emotion, exercise etc. Also, there is the "white coat" effect, which means blood pressure usually increases in clinical settings.

What you need to pay attention to when you measure your blood pressure at home:

- If the cuff is tied properly.
- If the cuff is too tight or too loose.
- If the cuff is tied on the upper arm.
- If you feel anxious.
- Taking 2-3 deep breaths before beginning will be better for measuring.
- Advice: Relax yourself for 4-5 minutes until you calm down.

Is the result the same if measuring on the right arm?

It is ok for both arms, but there will be some different results for different people. We suggest you measure the same arm every time.

This section includes a list of error messages and frequently asked questions for problems you may encounter with your blood pressure monitor. If the products not operating as you think it should, check here before arranging for servicing.


Problem	Symptom	Check This	Remedy
No power	Display is dim or will not light up.	Batteries are exhausted.	Replace with new batteries.
		Batteries are inserted incorrectly.	Insert the batteries correctly.
Low batteries	 Show on the display	Batteries are low.	Replace with new batteries
Error message	E1 shows	Communication error	Check if both the APP and the Bluetooth are on, operate and send the data again.
	E3 shows	The cuff is not secure.	Re-adjust the cuff and relax for a moment and then measure again.
	E10 or E11 shows	The monitor detected motion while measuring.	Movement can affect the measurement. Relax for a moment and then measure again.
	E20 shows	The measurement process does not detect the pulse signal.	Loosen the clothing on the arm and then measure again.
	E21 shows	The treatment of the measurement failed.	Relax for a moment and then measure again.
	EExx,shows on the display.	A calibration error occurred.	Retake the measurement. If the problem persists, contact the retailer or our customer service department for further assistance. Refer to the warranty for contact information and return instructions.

Power supply	Battery powered mode: 6VDC 4*AA batteries AC adapter powered mode: 6V==1A (Can be supplied by AC adaptor model UE08WCP-060100SPA only!) (Not Included)
Display mode	Digital LCD V.A.78mm×92mm
Measurement mode	Oscillographic testing mode
Measurement range	Rated cuff pressure: 0kPa - 40kPa (0mmHg~300mmHg) Measurement pressure: 5.3kPa-30.7kPa (40mmHg-230mmHg) pulse value: (40-199) beat/minute
Accuracy	Pressure: 5°C-40°C within ±0.4 kPa (3 mmHg) Pulse Value: ±5%
Normal working condition	Temperature:5 °C-40°C Relative Humidity: ≤85%RH Atmospheric Pressure: 86kPa-106 kPa
Storage & transportation condition	Temperature:-20 °C-60°C Relative Humidity: 10%RH-93%RH Atmospheric Pressure: 50kPa-106 kPa
Measurement perimeter of the upper arm	About 22cm~42cm
Weight	Approx.385g(Excluding the dry cells)
External dimensions	Approx.120*160*69mm
Attachment	4*AA batteries,user manual
Mode of operation	Continuous operation
Degree of protection	Type BF applied part
Protection against ingress of water	IP22
Software Version	V01

WARNING: No modification of this equipment is allowed.

Authorized Component

1. Please use the
Duronic authorized
adaptor.
(Not Included)

Adaptor
Type: UE08WCP-060100SPA Input:
100-240V, 50-60Hz,400mA
Output: 6V==1A
(Expected Service Life: 50,000 Hours)

Contact Information

For more information about our products, please visit
UNIT 1 SPILSBY ROAD, HAROLD HILL, ROMFORD
Manufactured by: Guangdong Transtek Medical Electronics Co., Ltd.
Company: Guangdong Transtek Medical Electronics Co., Ltd.
Address: Zone A, No.105 ,Dongli Road, Torch Development District,
Zhongshan,528437,Guangdong,China
Authorized European Representative:
Company: MDSS - Medical Device Safety Service GmbH
Address: Schiffgraben 41, 30175 Hannover, Germany

Complied European Standards List

Risk Management	ISO/EN 14971:2012
Labeling	ISO/EN 15223-1:2012
User Manual	EN 1041:2008
General Requirements for Safety	EN 60601-1:2006 IEC/EN 60601-1-11:2010
Electromagnetic Compatibility	IEC/EN 60601-1-2:2007
Performance requirements	EN 1060-1:1995+A2:2009 EN 1060-3:1997+A2:2009
Clinical investigation	EN 1060-4:2004
Usability	IEC/EN 60601-1-6: 2010 IEC/EN 62366: 2007
Software life-cycle processes	IEC/EN 62304:2006+AC: 2008

EMC Guidance

**Table 1 – Guidance and MANUFACTURER'S declaration
ELECTROMAGNETIC EMISSIONS – for all ME EQUIPMENT and ME
SYSTEMS**

Guidance and manufacturer's declaration – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 2	The device must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Table 2 – Guidance and MANUFACTURER'S declaration – electromagnetic IMMUNITY – for all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	<5 % U_T (>95 % dip in U_T) for 0,5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Table 4 – Guidance and MANUFACTURER'S declaration – electromagnetic IMMUNITY – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
IMMUNITY test	IEC 60601 TEST LEVEL	Compliance level	Electromagnetic environment – guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>3 Vrms 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2,5 GHz</p>	<p>3 Vrms</p> <p>3 V/m</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.167 \sqrt{P}$ <p>$d = 1.167 \sqrt{P}$ 80 MHz to 800 MHz</p> <p>$d = 2.333 \sqrt{P}$ 800 MHz to 2,5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.</p> <p>NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p>a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.</p> <p>b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

Table 6 – Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM – for ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the device			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.167 \sqrt{P}$	80 MHz to 800 MHz $d = 1.167 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2.333 \sqrt{P}$
0,01	0.117	0.117	0.233
0,1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.690	3.690	7.378
100	11.67	11.67	23.33
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Warranty

1 Year Guarantee from Shine-Mart Ltd, Trademark owners of Duronic Products.

NOTE: THESE STATEMENTS DO NOT AFFECT YOUR STATUTORY RIGHTS AS A CONSUMER

This product has been manufactured under the strictest or quality control procedures, and using the highest quality of materials, to ensure excellent performance and reliability. It will give very good and long lasting service, provided it is properly used and maintained. The product is guaranteed for 1 years from the date of original purchase. If any defect arises due to a faulty materials or poor workmanship, the faulty product must be returned to the original place of purchase. Refund or replacement is at the discretion of that company. Duronic Products are offered with a 1 year guarantee under the following conditions:

1. The product must be returned to the retailer with original proof of purchase.
2. The product must be installed and used in accordance with the instructions contained in this operating manual.
3. It must be used only for domestic purposes.
4. It does not cover wear and tear, malicious damage, misuse, neglect, inexpert repairs or consumable parts.
5. Shine-Mart Ltd has no responsibilities for accidental or consequential loss or damage.
6. Shine-Mart Ltd are not liable to carry out any type of servicing work, under the guarantee
7. Valid in the UK only

Imported by Shine-Mart Ltd, RM3 8SB

Information on waste Disposal for Consumers of Electrical & Electronic Equipment.

This mark on a product and/or accompanying documents indicates that when it is to be disposed of, it must be treated as Waste Electrical & Electronic Equipment, (WEEE).

Any WEEE marked products must not be mixed with general household waste, but kept separate for the treatment, recovery and recycling of the materials used.

For proper treatment, recovery and recycling; please take all WEEE marked waste to your Local Authority Civic waste site, where it will be accepted free of charge.

If customers dispose of Waste Electrical & Electronic Equipment correctly, they will be helping to save valuable resources and preventing any potential negative effects upon human health and the environment, of any hazardous materials that the waste may contain.

Follow Us

duronic.com

[duronic](https://www.youtube.com/duronic)

[@duronicmedia](https://www.instagram.com/duronicmedia)

[Duronic](https://www.facebook.com/duronic)

[@duronicmedia](https://twitter.com/duronicmedia)

Distributed by: Model: LS802-B Ref. No.: BPM490
 SHINEMART LTD Unit 1, Spilsby Road, Harold Hill,
 Romford, Essex, RM3 8SB UK

 **Guangdong Transtek Medical Electronics Co., Ltd.**
 Zone A, No.105 ,Dongli Road, Torch Development District,
 Zhongshan,528437,Guangdong,China

EC	REP
----	-----

 MDSS - Medical Device Safety Service GmbH
 Schiffgraben 41,30175 Hannover, Germany

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5. Shine-Mart Ltd non è responsabile per Danni accidentali o consequenziali.
6. Shine-Mart Ltd non è tenuta a eseguire qualsiasi tipo di riparazione, ai sensi della garanzia.
importata da Shine-Mart Ltd, RM3 8SB