



1. Before using the orthosis, clip or file the claw of the spur, so that it does not irritate the leg. If necessary, a hole can be made in the orthosis to relieve the spur.




2. The orthosis is easiest to apply while the animal lies down and the carpal joint is in a natural / stretched position.



3. Place the limb on the open orthosis with the wrap just above the back of the main footpad*.



4. Tighten so that the orthosis is shaped around the leg. Fasten the orthosis with the elastic straps.



5. If the overlap is too wide (5a) or the wrap is too long (5b), cut them with scissors. If you are planning to use only the orthosis without the stabilizing nylon straps, attach the elastic straps tightly over the outside of the wrap. If the straps are too long (5c), cut them with scissors and burn the edges gently with a lighter to prevent fraying.




6. Apply the orthosis. Note! The Paw must have free movement and the orthosis must not cause pressure or chafing on top of the paw.

For additional stability and support Use the included stabilizing nylon straps.

- Apply the stabilizing straps after applying the orthosis and before the animal is putting weight on the leg.
- Set the desired joint angle / load position before applying the stabilization straps.
- The stabilizing bands are to be applied to the outside of the orthosis where the animal needs additional stability (eg carpal hyperextension).
- Cut the stabilizing straps to the desired length.
- Tighten the outer elastic straps over both the orthosis and the stabilizing straps.

Examples of stabilization with nylon straps.

NOTE! Fasten nylon straps when the orthosis is already applied on the animal!



Important information

- Examine the leg 2-4 times / day for swelling of the toes, as well as heat gain and chafing.
- Contact your vet if exuding wounds, odours or other problems occur. It is always easier to treat complications that are noticed early, so let them know if there are skin irritations, preferably before it has become a wound or drainage of fluid.
- The orthosis always affects the gait. Does it seem uncomfortable for the patient? If the animal licks, bites or is not using the leg, please contact your vet. Be aware of changing behavior, for example, if the mood is not as it used to, or if the animal does not want to eat. Contact your veterinary clinic if you have questions!

Care instructions

- Remove fur and debris from the Velcro for best adhesion.
- Do not wash the straps.
- Clean the orthosis with mild detergent by hand, maximum 30 °C.

Contact us for questions regarding the orthosis
and for any problems that may occur.

Tel. +46 (0)31-21 41 28

Email: info@scandiop.se

www.scandiop.se



Scandiop

© 2016 Scandiop