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COVID-19 Antigen Rapid Test Clinical Sensitivity and Specificity Study Report

1. Objective

The CLUNGENE[®] COVID-19 Antigen Rapid Test (hereinafter referred to as the CLUNGENE Device) manufactured by Hangzhou Clongene Biotech Co., Ltd. is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasopharyngeal swab and oropharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider.

This study is intended to evaluate the clinical performance, between the CLUNGENE Device and the comparator RT-PCR assay.

2. Method

A study of 770 direct nasopharyngeal swabs was performed. The specimens were prospectively collected from patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care at 5 locations and tested at a single central laboratory.

Two nasopharyngeal swabs were collected from individual symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19. At all locations, one nasopharyngeal swab was tested directly with the COVID-19 Antigen Rapid Test according to product instructions for use, and the other swab was eluted in 3 mL viral transport media and tested with RT-PCR assay for detection of SARS-CoV-2. Swabs were randomly assigned to testing with the COVID-19 Antigen Rapid Test or RT-PCR assay and were tested by operators who were blinded to the RT-PCR test result.

The positive percent agreement (PPA) was calculated as $100\% \times (\text{True Positive} / [\text{True Positive} + \text{False Negative}])$. The negative percent agreement (NPA) was calculated as $100\% \times (\text{True Negative} / [\text{True Negative} + \text{False Positive}])$. Accuracy was calculated as $100\% \times ([\text{True Positive} + \text{True Negative}] / \text{Total sample Qty})$. The 95% (two-sided) confidence interval (CI) was calculated using the Wilson Score Method.

3. Comparator method

Real-Time Fluorescent RT-PCR Kit for Detecting SARS-CoV-2, manufactured by BGI Genomics Co. Ltd., is a real-time reverse transcription polymerase chain reaction (rRT-PCR) test. This product has got CE, NMPA certification and FDA Emergency Use Authorized. A specimen is positive for SARS-CoV-2 if the Ct value of ORF1ab gene is not higher than 37 and the Ct value of human housekeeping gene β -Actin is not higher than 35.

4. Enrollment criteria (inclusion/exclusion criteria)

4.1 Inclusion criteria

- Patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care.
- Symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19.

4.2 Exclusion criteria

- Unable to obtain samples of information needed for the experiment
- Samples that have been contaminated or contaminated during sample storage
- Samples with inappropriate storage conditions

5. Result

The results are summarized in the following table.

The RT-PCR cycle threshold (Ct) is the relevant signal value. Lower Ct value indicate higher viral load. The

sensitivity was calculated for the different Ct value range (Ct value \leq 33 and Ct value \leq 37).

COVID-19Antigen		RT-PCR (Ct value \leq 33)		Total
		Positive	Negative	
CLUNGENE[®]	Positive	145	2	147
	Negative	3	593	596
Total		148	595	743

PPA (Ct \leq 33): 98.0% (145/148), (95%CI: 94.2% ~99.3%)

NPA: 99.7% (593/595), (95%CI: 98.8% ~99.9%)

Accuracy: 99.3% ((145+593)/743), (95%CI: 98.4% ~99.7%)

COVID-19Antigen		RT-PCR (Ct value \leq 37)		Total
		Positive	Negative	
CLUNGENE[®]	Positive	161	2	163
	Negative	14	593	607
Total		175	595	770

PPA (Ct \leq 37): 92.0% (161/175), (95%CI: 87.0% ~95.2%)

NPA: 99.7% (593/595), (95%CI: 98.8% ~99.9%)

Accuracy: 97.9% ((161+593)/770), (95%CI: 96.6% ~98.7%)

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

6. Conclusion

Taken together, the CLUNGENE[®] COVID-19 Antigen Rapid Test had a positive percent agreement (sensitivity) of 98.0% (95%CI: 94.2% ~99.3%) with specimens of a Ct count \leq 33, 92.0% (95%CI: 87.0% ~95.2%) with specimens of a Ct count \leq 37, negative percent agreement (specificity) of 99.7% (95%CI: 98.8% ~99.9%), and accuracy of 99.3% (95%CI: 98.4% ~99.7%) with specimens of a Ct count \leq 33, 97.9% (95%CI: 96.6% ~98.7%) with specimens of a Ct count \leq 37.

COVID-19 Antigen Rapid Test Cassette (Nasal Swab) Clinical Sensitivity and Specificity Study Report

1. Objective

The CLUNGENE[®] COVID-19 Antigen Rapid Test Cassette (Nasal Swab) (hereinafter referred to as the CLUNGENE Device) manufactured by Hangzhou Clongene Biotech Co., Ltd. is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasal swab from individuals who are suspected of COVID-19 by their healthcare provider.

This study is intended to evaluate the clinical performance, between the CLUNGENE Device and the comparator RT-PCR assay.

2. Method

A study of 617 direct nasal swabs was performed. The specimens were prospectively collected from patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care at 5 locations and tested at a single central laboratory.

Two nasal swabs were collected from individual symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19. At all locations, one nasal swab was tested directly with the COVID-19 Antigen Rapid Test Cassette (Nasal Swab) according to product instructions for use, and the other swab was eluted in 3 mL viral transport media and tested with RT-PCR assay for detection of SARS-CoV-2. Swabs were randomly assigned to testing with the COVID-19 Antigen Rapid Test Cassette (Nasal Swab) or RT-PCR assay and were tested by operators who were blinded to the RT-PCR test result.

The positive percent agreement (PPA) was calculated as $100\% \times (\text{True Positive} / [\text{True Positive} + \text{False Negative}])$. The negative percent agreement (NPA) was calculated as $100\% \times (\text{True Negative} / [\text{True Negative} + \text{False Positive}])$. Accuracy was calculated as $100\% \times ([\text{True Positive} + \text{True Negative}] / \text{Total sample Qty})$. The 95% (two-sided) confidence interval (CI) was calculated using the Wilson Score Method.

3. Comparator method

Real-Time Fluorescent RT-PCR Kit for Detecting SARS-CoV-2, manufactured by BGI Genomics Co. Ltd., is a real-time reverse transcription polymerase chain reaction (rRT-PCR) test. This product has got CE, NMPA certification and FDA Emergency Use Authorized. A specimen is positive for SARS-CoV-2 if the Ct value of ORF1ab gene is not higher than 37 and the Ct value of human housekeeping gene β -Actin is not higher than 35.

4. Enrollment criteria (inclusion/exclusion criteria)

4.1 Inclusion criteria

- Patients in community meeting Department of Public Health definition of a suspected COVID-19 case and being tested for SARS-CoV-2 part of routine medical care.
- Symptomatic (within 7 days of onset) or asymptomatic patients who were suspected of COVID-19.

4.2 Exclusion criteria

- Unable to obtain samples of information needed for the experiment
- Samples that have been contaminated or contaminated during sample storage
- Samples with inappropriate storage conditions

5. Result

The results are summarized in the following table.

The RT-PCR cycle threshold (Ct) is the relevant signal value. Lower Ct value indicate higher viral load. The sensitivity was calculated for the different Ct value range (Ct value ≤ 30 and Ct value ≤ 37).

COVID-19 Antigen		RT-PCR (Ct value \leq 30)		Total
		Positive	Negative	
CLUNGENE[®]	Positive	117	3	120
	Negative	3	462	465
Total		120	465	585

PPA (Ct \leq 30):97.5% (117/120), (95%CI: 92.9%~99.2%)

NPA: 99.4% (462/465), (95%CI: 98.1%~99.8%)

COVID-19 Antigen		RT-PCR (Ct value \leq 37)		Total
		Positive	Negative	
CLUNGENE[®]	Positive	139	3	142
	Negative	13	462	475
Total		152	465	617

PPA (Ct \leq 37):91.4% (139/152), (95%CI: 85.9%~94.9%)

NPA: 99.4% (462/465), (95%CI: 98.1%~99.8%)

PPA - Positive Percent Agreement (Sensitivity)

NPA - Negative Percent Agreement (Specificity)

6. Conclusion

Taken together, the CLUNGENE[®] COVID-19 Antigen Rapid Test Cassette (Nasal Swab) had a positive percent agreement (sensitivity) of 97.5% (95%CI: 92.9%~99.2%) with specimens of a Ct count \leq 30, 91.4% (95%CI: 85.9%~94.9%) with specimens of a Ct count \leq 37, negative percent agreement (specificity) of 99.4% (95%CI: 98.1%~99.8%).

EC DECLARATION OF CONFORMITY

Name and address of the manufacturer: **Hangzhou Clongene Biotech Co., Ltd.
No.1 Yichuang Road, Yuhang Sub-district
Yuhang District
311121 Hangzhou
China**

We declare under our sole responsibility that

the medical device: **COVID-19 Antigen Rapid Test**

of class: **Other**
according to article 9 of directive 98/79/EC

meets the provisions of the directive 98/79/EC and its transpositions in national laws which apply to it. The declaration is valid in connection with the "final inspection report" of the device.

Conformity assessment procedure: **Directive 98/79/EC Annex III**

Applicable standards: **EN ISO 13485:2016 EN ISO 15223-1:2016
EN ISO 23640:2015 EN13612:2002/AC:2002
EN 13975:2003 EN ISO 14971:2012
EN ISO 18113-1:2011 EN ISO 18113-2:2011
EN 62366-1:2015**

Name and address of the authorized representative: **Shanghai International Holding Corporation GmbH (Europe)
Eiffestrasse 80
20537 Hamburg
Germany**

Hangzhou, July.15.2020



Place, date

Shujian Zheng, Legal representative

Name and function


COVID-19 Antigen Rapid Test Cassette

Hangzhou Clongene Biotech Co., Ltd.
en.clongene.com

COVID-19 & SARS-CoV-2

COVID-19 is an acute respiratory infectious disease caused by novel coronavirus (SARS-CoV-2), and people are generally susceptible. Based on the current epidemiological investigation, the incubation period is 1 to 14 days, mostly 3 to 7 days. Novel coronavirus includes four typical structural proteins: Spike Protein, Envelope Protein, Membrane Protein and Nucleocapsid Protein.

Nucleocapsid (N) protein is the most abundant protein with highly conserved in SARS-CoV-2. N protein is used as the core raw material of rapid diagnostic reagent for immunology in the market.

Clongene has developed the COVID-19 Antigen Rapid Test Cassette. The COVID-19 Antigen Rapid Test is a lateral flow immunoassay intended for the qualitative detection SARS-CoV-2 nucleocapsid antigens in nasal swab, nasopharyngeal swab and oropharyngeal swab from individuals who are suspected of COVID-19 by their healthcare provider.

Kit Contents

Contents				
Extraction reagents	Work Station	Test Cassette	Sterilized Swab	Extraction Tube & Dropper Tip
				

Product Features

 CE certification

 Easy to collect samples


 No equipment required

 Instant result at 15 minutes

 Results are clearly visible

 Suitable for large-scale rapid screening

Principle


The COVID-19 Antigen Rapid Test is a lateral flow immunoassay based on the principle of the double-antibody sandwich technique. If the specimen contains SARS-CoV-2 antigen, a colored test line (T) would be visible in the result window. Absence of the T line suggests a negative result. The control line (C) is used for procedural control, and should always appear if the test procedure is performed properly.


Specimens

The detect specimens include nasal swab, nasopharyngeal swab and oropharyngeal swab.


Nasal Swab

Nasopharyngeal Swab

Oropharyngeal Swab


Inadequate specimen collection, improper specimen handling and/or transport may yield false results; therefore, training in specimen collection is highly recommended due to the importance of specimen quality to obtain accurate test results.

Sample Collection


Nasal Swab

Gently rotating the swab, insert swab about 2.5 cm (1 inch) into nostril until resistance is met at turbinates. Rotate the swab several times and repeat in other nostril using the same swab.


Nasopharyngeal Swab

Insert the swab through the nostril parallel to the palate (not upwards) until resistance is encountered or the distance is equivalent to that from the ear to the nostril of the patient, indicating contact with the nasopharynx. Gently rub and roll the swab.

Oropharyngeal Swab

Insert swab into the posterior pharynx and tonsillar areas. Rub swab over both tonsillar pillars and posterior oropharynx and avoid touching the tongue, teeth, and gums.

Test Procedure


1

Add all of the extraction reagent into an extraction tube.


2

Insert the swab specimen into the extraction tube. Roll the swab at least 5 times and leave the swab in the extraction tube for one minute.


3

Remove the swab while squeezing the sides of the tube to extract the liquid from the swab.


4

Cover the extraction tube with a dropper tip tightly.


5

Transfer 3 drops (approximately 100µL) to the specimen well of the test cassette.


6

Interpret the test results at 15 minutes. Do not read results after 20 minutes.

Interpretation of Results

Performance Characteristics

Clinical Performance

■ For nasopharyngeal swab:

770 nasopharyngeal swabs were collected from individual symptomatic patients. The swabs were detected by COVID-19 Antigen Rapid Test of Clongene and the RT-PCR. Summary data was showed as below:

COVID-19 Antigen		RT-PCR (Ct value ≤ 33)		Total
		Positive	Negative	
CLUNGENE [®]	Positive	145	2	147
	Negative	3	593	596
Total		148	595	743

PPA (Ct ≤ 33): 98.0% (145/148), (95%CI: 94.2% ~ 99.3%)

NPA: 99.7% (593/595), (95%CI: 98.8% ~ 99.9%)

COVID-19 Antigen		RT-PCR (Ct value \leq 37)		Total
		Positive	Negative	
CLUNGENE [®]	Positive	161	2	163
	Negative	14	593	607
Total		175	595	770

PPA (Ct \leq 37): 92.0% (161/175), (95%CI: 87.0% ~ 95.2%)

NPA: 99.7% (593/595), (95%CI: 98.8% ~ 99.9%)

- PPA - Positive Percent Agreement (Sensitivity)
- NPA - Negative Percent Agreement (Specificity)

■ For nasal swab:

617 nasal swabs were collected from individual symptomatic patients. The swabs were detected by COVID-19 Antigen Rapid Test of Clongene and the RT-PCR. Summary data as below:

COVID-19 Antigen		RT-PCR (Ct value \leq 33)		Total
		Positive	Negative	
CLUNGENE [®]	Positive	132	3	135
	Negative	4	462	466
Total		136	465	601

PPA (Ct \leq 33): 97.1% (132/136), (95%CI: 92.7% ~ 98.9%)

NPA: 99.4% (462/465), (95%CI: 98.1% ~ 99.8%)

COVID-19 Antigen		RT-PCR (Ct value ≤ 37)		Total
		Positive	Negative	
CLUNGENE®	Positive	139	3	142
	Negative	13	462	475
Total		152	465	617

PPA (Ct ≤ 37): 91.4% (139/152), (95%CI: 85.9% ~ 94.9%)

NPA: 99.4% (462/465), (95%CI: 98.1% ~ 99.8%)

- PPA - Positive Percent Agreement (Sensitivity)
- NPA - Negative Percent Agreement (Specificity)

Limit of Detection (Analytical Sensitivity)

The study used cultured SARS-CoV-2 virus, which is heat inactivated and spiked into nasopharyngeal swab specimen. The Limit of Detection (LoD) is 5.7×10^2 TCID₅₀/mL.

Cross Reactivity (Analytical Specificity)

We have evaluated 32 commensal and pathogenic microorganisms that may be present in the nasal cavity and no cross-reactivity was observed.

Interference

17 potentially interfering substances with different concentration were evaluated and found no affect to the test performance.

High-dose Hook Effect

The COVID-19 Antigen Rapid Test was tested up to $1.0 \times 10^{5.67}$ TCID₅₀/mL of inactivated SARS-CoV-2 and no high-dose hook effect was observed.

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Common List of COVID-19 rapid antigen tests recognized in EU

EUROPEAN COMMISSION
DIRECTORATE-GENERAL FOR HEALTH AND FOOD SAFETY

Public health, country knowledge, crisis management
Health Security and Vaccination

III. Rapid antigen tests of which the test results are mutually recognised

- Hangzhou Clongene Biotech, Clongene COVID-19 Antigen Rapid Test Kit

EC DECLARATION OF CONFORMITY

Name and address of the manufacturer: **Hangzhou Clongene Biotech Co., Ltd.**
No.1 Yichuang Road, Yuhang Sub-district
Yuhang District
311121 Hangzhou
China

We declare under our sole responsibility that

the medical device: **COVID-19 Antigen Rapid Test Cassette (Nasal Swab)**

of class: **Other**
according to article 9 of directive 98/79/EC

meets the provisions of the directive 98/79/EC and its transpositions in national laws which apply to it. The declaration is valid in connection with the "final inspection report" of the device.

Conformity assessment procedure: **Directive 98/79/EC Annex III, excluding Section 6**

Applicable standards: **EN ISO 13485:2016** **EN ISO 15223-1:2016**
EN ISO 18113-1:2011 **EN ISO 18113-2:2011**
EN ISO 23640:2015 **EN ISO 14971:2019**
EN 13641:2002 **EN13612:2002/AC:2002**
EN 13975:2003 **EN 62366-1:2015**
ISO 17511:2020

Name and address of the authorized representative: **Shanghai International Holding Corporation GmbH (Europe)**
Eiffestrasse 80
20537 Hamburg
Germany

Hangzhou, February 04, 2021.

Place, date

Shujian Zheng, Legal representative

Name and function

Certificate

The Certification Body of
TÜV Rheinland LGA Products GmbH

hereby certifies that the organization
**Hangzhou Clongene Biotech
Co., Ltd.**
**No. 1 Yichuang Road, Yuhang Sub-district
Yuhang District
311121 Hangzhou
P.R. China**

has established and applies a quality management system for medical devices
for the following scope:

**Design and Development, Manufacture and Distribution of
In-vitro Diagnostic Rapid Test of Fertility, Drug of Abuse,
Infectious Diseases, Tumour Markers and Cardiac Markers**

Proof has been furnished that the requirements specified in

EN ISO 13485:2016

are fulfilled. The quality management system is subject to yearly surveillance.

Effective Date: 2020-11-25
Certificate Registration No.: SX 60152722 0001
An audit was performed. Report No.: 15073650 008
This Certificate is valid until: 2023-11-18

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Hangzhou Clongene Biotech Co., Ltd.
No.1 Yichuang Road, Yuhang Sub-district, Yuhang District,
311121 Hangzhou, China

EC REP
Shanghai International Holding Corp. GmbH
(Europe) Eiffestrasse 80, D-20537 Hamburg,
Germany

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July 6, 2021

To whom it may concern,

We, as the manufacturer of the CLUNGENE COVID-19 Antigen Rapid Test, assure that the following major SARS-CoV-2 variants which seem to spread more easily and quickly than other variants and are circulating globally, can be detected by CLUNGENE COVID-19 Antigen Rapid Test:

- ✓ **B.1.1.7 (Alpha)**: This variant was initially detected in the UK;
- ✓ **B.1.351 (Beta)**: This variant was initially detected in South Africa in December 2020;
- ✓ **P.1 (Gamma)**: This variant was initially identified in travelers from Brazil, who were tested during routine screening at an airport in Japan, in early January;
- ✓ **B.1.617.2 (Delta)**: This variant was initially detected in India.

The recombinant N proteins of these variants can be detected by the CLUNGENE COVID-19 Antigen Rapid Test, at concentrations same as that of the ordinary recombinant N protein without variation (amino acid: Met1-Ala419, accession # QHO62884.1).

Yours sincerely,

Shujian Zheng,

Managing Director of Hangzhou Clongene Biotech Co., Ltd.