



IBI SCIENTIFIC

9861 Kapp Court Peosta, IA 52068
Tele: (563) 690-0484 Toll Free: 1-800-253-4942
Fax: (563) 690-0490 E-Mail: info@ibisci.com
Web: www.ibisci.com

ISO 9001:2008 Certified

ACCLG0002 – Replacement Belly Dancer Legs

Important! – For proper Belly Dancer operation, when a unit requires new legs it is important to **replace all four legs**. Replacing only one or two legs will result in poor performance and premature failure of the other Belly Dancer legs.

- 1.) Prior to performing any type of service to the Belly Dancer unit, ensure the unit is unplugged from any/all electrical sources.
- 2.) Now, using a side-cutters or sharp knife cut the eight (8) tie straps shown at locations “A” and “B”. Use caution when removing straps! Remove the cut straps from unit.
- 3.) Next, supporting the platform remove the tubing from the upper and lower standoffs in all four (4) positions. Discard all four (4) used legs.
- 4.) Remove four (4) new legs from packaging and replace legs onto standoffs in the same manner the old legs were removed. Ensure the legs are fully seated onto standoffs.
- 5.) Now remove the eight (8) tie straps from packaging and wrap around the legs where they meet the standoffs, see positions “A” and “B”. Tighten straps around legs to secure properly.
- 6.) Cut excess tie strap end away from legs and discard.
- 7.) Reconnect unit’s power cord to input entry module and standard power source. Turn control dial on front of unit and ensure unit will now function properly.

