Applicator Tampon User Manual

Hey, Ü! Our bodies & our planet
deserve organic cotton.

Read me carefully

Our tampons are made of 100% certified organic cotton and plastic-free cardboard applicator. so no nasty synthetics inside this box or our bodies!

What you should know about toxic shock syndrome

TSS is a rare but serious illness that can occur among men, women and children. It is caused by toxins produced by the bacteria *Staphylococcus aureus*. This bacteria is commonly found on the skin, and inside the nose and vagina.

Around half the cases of TSS are by menstruating teens and women under thirty years of age. TSS is commonly associated with tampon use but can also occur when tampons are not used.

Symptoms

If during your period you get any of the following symptoms immediately remove your tampon and contact your doctor. The initial symptoms are very similar to the flu. They usually appear very quickly and are often severe.

Symptoms include: High fever (above 39 degree/102 f), vomiting, diarrhea, a sunburn like rash, dizziness and / or fainting.

Tell you doctor that have been using tampons and that you are concerned about TSS. TSS can have fatal consequences, but can be easily treated at an early stage. Please consult a doctor as soon as you can.

Seek medical advice before resuming the use of tampons if you have had TSS warning signs in the past. You can reduce the risk of TSS by using tampons with the lowest absorbency necessary and by changing your tampon regularly (every 4 to max. 8 hours).

Top tips

Can tampons be worn at night? They can be worn at night but make sure you change yours right before you go to bed and first thing when you wake up.

How many tampons can I use at one time? Only ever use one tampon at a time. Make sure to remove yours at the end of your period and don't use them when you aren't on your period.

Should I be able to feel my tampon? If your tampon is inserted correctly, you shouldn't be able to feel it. If you can, you may not have pushed it in far enough or you may be using an absorbency that is too high for your needs causing it to feel dry.

How do I know when to change my tampon? Don't wear your tampon for longer than 8 hours. If you start to feel your tampon, this is a sign it's time to change it as it will be full.

How do I know which absorbency to choose? It's best to use a tampon with the lowest absorbency that meets your needs. On lighter days, perhaps towards the end of your period, choose a lighter absorbency. On heavy flow days, you can use a higher absorbency tampon.

Use the guide below to select which
absorbency is right for you.

Regular  6-9g

Super  9-12g

How to use our tampons (written with love by LUÜNA women)

Preparation

- Wash your hands.
- Check your tampon paper wrapper. Is it completely sealed? If not, don't use it.
- Open the paper wrapper at the bottom of the tampon (the Ü on the wrapper will be facing up when it's held the correct way up).
- Remove the whole wrapper, holding the bottom of the cardboard applicator.

Inserting the tampon

- Find a comfortable position. You can sit on the toilet with your legs open or stand with one leg raised (best to rest your foot on something like the toilet - unless you have exceptional balance).
- Relax. This is super important. If you don't relax, your muscles will be tense which will make it difficult to insert your tampon.
- Gently hold the skin around your vaginal opening.
- Insert the tip of the cardboard applicator into your vaginal opening until the outer tube is almost completely inserted.
- Then, with your finger, push the inner tube until the tampon is released from the applicator completely.
- You'll know your tampon is in the right place because you won't be able to feel it anymore. If you do, push it in a little further.
- Leave the string hanging outside of your body.
- Wash your hands.

Removing the tampon

- You should change your tampon regularly and not leave it in for more than 8 hours. We recommend changing every 4-6 hours.
- If you start to feel your tampon, this is a sign it's time to change it.
- First, wash your hands.
- Firmly tug the string and pull the tampon out.
- If this feels uncomfortable and dry, it might mean you are using a tampon with an absorbency that is too high for you. Try switching to a lower absorbency.
- Throw your tampon away (not down the toilet!)
- Wash your hands.

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