



Thank you for purchasing the EZ Lite Cruiser Adjustable Leg Rests

Follow these instructions step-by-step for installation.

For any questions, please e-mail:

support@ezlitecruiser.com

Or call us at: (888) 544-6054

Operators are available anytime between 9 AM and 5 PM Pacific Time, Monday through Friday.

Watch a VIDEO of the Installation:

<http://www.ezlitecruiser.com/ez10>


Here's are the parts included with your order:



Please identify these components, according to the numbers below, and lay them out for easy access during installation. Please also see the photos to the right for a detailed identification of the sub-parts.

If anything is missing, please let us know.

1. Two (2) Frame Clamps
2. Two (2) Mounting Rods
3. Two (2) Leg Rests
4. One (1) Hex Wrench



MOUNTING PINS
TOOTH

MOUNTING HOLES
LATCH

Detailed
Identification of
Components

Step 1: Attach Clamps & the Mounting Rods

A. Loosen the screws of the clamps enough to fit the larger side of the clamp around the frame. After attaching the clamp to the frame, moderately tighten the clamp so it rests on the frame, but don't tighten it completely yet.

B. Insert the Mounting Rod into the smaller side of the clamp. When facing the chair, the rounded end ("Tooth") of the Mounting Rod should be facing inward, and parallel to the seat, with the Mounting Pins facing to the outside of the chair.

****Note**, there is a Left and Right side Mounting Pin and Leg Rest, so don't get them mixed up.



Repeat Step 1 process "A" & "B" for both the Left and Right side.



Then, once both of the Clamps and the Mounting Rods are firmly secured to the frame of the chair, Attach each Leg Rest by fitting the Mounting Holes (shown above, right) to the Mounting Pins (shown above, left) on the Mounting Rod.



Step 2: Swing Leg Rests Inward and Latch them to the Mounting Rods.

Swing the Leg Rest Inward slowly.

If you find that the Mounting Rod is spinning as you swing the leg rests inward, then either:

- A) The Clamps are not secured tight.
- B) The wrong end of the Clamp is attached to the frame.
- C) The Clamp is attached to the Rod in the wrong place (too high.)

After the Latch meets the "Tooth" of the Mounting Rod, gently push the latch so it locks on to the "Tooth" and secures the Leg Rests to the Mounting Rods. Do not force the Latch to try to lock by continuing to swing the leg rests inward! This will cause the Mounting Rod to spin out of alignment.

