CoaguChek PT Controls

REF		SYSTEM
06679684190	4 bottles each of lyophilized control plasma Level 1 & 2 <p>8 diluent-filled bulb pipettes (red) <p>1 quality control code chip</p></p>	CoaguChek® Pro II

English

Intended use

CoaguChek PT Controls are used for system checks and quality control testing of prothrombin time with the CoaguChek Pro II meter and CoaguChek PT Test strips. The CoaguChek PT Controls are intended for near-patient testing. Not for self-testing.

Reagents - working solutions

Each bottle of control contains lyophilized anticoagulated rabbit plasma. Each diluent-filled bulb pipette contains calcium chloride and preservatives in water.

Target values and ranges

The code chip provided in the pack contains all relevant lot-specific information to perform a quality control test.

The CoaguChek Pro II meter displays the control range, the reading obtained and an evaluation of the result (inside or outside of the control range). The result is automatically saved as a control reading in the memory of the meter.

The control intervals and limits should be adapted to each laboratory’s individual requirements. Values obtained should fall within the defined limits. Each laboratory should establish corrective measures to be taken if values fall outside the defined limits.

Follow the applicable government regulations and local guidelines for quality control.

Precautions and warnings

For in vitro diagnostic use for health care professionals. Exercise the normal precautions required for handling all laboratory reagents.

Infectious o microbial waste:

Warning: handle waste as potentially biohazardous material. Dispose of waste according to accepted laboratory instructions and procedures.

Environmental hazards:

Apply all relevant local disposal regulations to determine the safe disposal.

Safety data sheet available for professional user on request.

Handling

Use scissors to cut the tip off the bulb pipette very close to the welding. Make sure there is no liquid contained in the cut-off section. Also make sure that no liquid is left in the bulb itself.

- Dissolve the content of one bottle carefully by adding the complete content of 1 bulb pipette (red).
- Mix carefully by gently swirling the bottle. Avoid foam formation.
- Allow to stand closed for 5 minutes to reconstitute.
- Use up the reconstituted control material within 30 minutes.

Storage and stability

Store at 2-8 °C.

The lyophilized control plasma is stable up to the stated expiration date. The reconstituted control material is stable for 30 minutes.

Materials provided

- 4 bottles of lyophilized control plasma Level 1
- 4 bottles of lyophilized control plasma Level 2
- 8 diluent-filled bulb pipettes (red)
- 1 code chip

Materials required (but not provided)

- REF 07210841xxx, CoaguChek Pro II meter (with W-LAN) or
- REF 07237944xxx, CoaguChek Pro II meter (without W-LAN)
- REF 06688721xxx, CoaguChek PT Test

Assay

- Prepare the meter and a test strip according to the instructions given in the operator’s manual and in the method sheet for the test strips.
- Insert the quality control code chip into the meter.
- Use the bulb pipette to apply one hanging drop of the control solution to the sample application area of the test strip.
- Save the remaining control solution until the test result is obtained for repeating the quality control check if necessary.
- Use the control solution within 30 minutes after reconstitution.

System checks

The CoaguChek Pro II meter has a number of inbuilt system checks. Please refer to the operator’s manual for details. The CoaguChek PT Test strip has an inbuilt quality control function.

Performance data

Reproducibility of CoaguChek PT Test was determined using CoaguChek PT Controls at 4 external study sites. Controls were measured on 21 days in 2 runs per day and with 2 or 3 lots per study site. Reproducibility figures were calculated by ANOVA (Analysis of Variance).

Reproducibility				
PT Control Level	Mean (INR)	SD (INR)	CV (%)	
Level 1	1.28	0.04	3.2	
Level 2	2.94	0.09	3.1	

For further information, please refer to the appropriate Meter Manual and the Method Sheets of all necessary components.

A point (period/stop) is always used in this Method Sheet as the decimal separator to mark the border between the integral and the fractional parts of a decimal numeral. Separators for thousands are not used.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or the patient is established.

The Summary of Safety & Performance Report can be found here: https://ec.europa.eu/tools/eudamed

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© 2022, Roche Diagnostics
COAGUCHEK is a trademark of Roche.
All other product names and trademarks are the property of their respective owners.
Roche Diagnostics GmbH, Sandhofer Strasse 116, D-68305 Mannheim