

# SONICTEC

healthcare made easier



## SARS-Cov-2 & Influenza A/B Combo Rapid Test Cassette (swab)

An Antigen rapid test for the detection of SARS-Cov-2 and influenza A/B virus in nasal swab. For self-testing use.

Read the instructions carefully before taking the test.

REF: K751416D

English

Australia Sponsor: Sonictec Pty Ltd

Customer Support Number: 02 8328 1008

Address: 17 Chisholm St, Wollli Creek NSW 2205

Hours: 9am-7pm AEST / 9am-8pm AEDT, 7 days per week

Email: [info@sonictec.com.au](mailto:info@sonictec.com.au) Website: [www.sonictec.com.au/sonictec3in1](http://www.sonictec.com.au/sonictec3in1)




Scan the QR code for information on how to use the SARS-Cov-2 & Influenza A/B Combo Rapid Test Cassette (swab).


### MATERIALS PROVIDED




Test device



Package insert




Extraction Tube with buffer



Sterilized Swab (The information about the swab can be found on its packaging.)



Biohazard Specimen Bag



Tube stand (The 25 tests/kit package contains the tube stand, the 1 test/kit, 2 tests/kit and 5 tests/kit package use the test box itself as tube stand.)

Materials required but not provided

Timer

### TEST PROCEDURE

<p><b>Step 1</b> Wash or clean your hands and make sure they are dry before starting the test.</p>	<p><b>Step 2</b> Read the instruction for use carefully.</p>	<p><b>Step 3</b> Take out one extraction tube, pull off the sealed aluminum foil on the extraction tube. Place extraction tube into tube stand or box tube stand.</p>	<p><b>Step 4</b> Unpack the swab. <b>Caution: Do not touch swab tip when handling the swab.</b></p>	
<p><b>Step 5</b> Tilt your head back slightly. Insert the swab about 2 cm - at least with the entire soft swab tip - into the left nostril. Gently rotate the swab at least five times against the nasal wall.</p>		<p><b>Step 6</b> Insert the same swab about 2 cm - at least with the entire soft swab tip - into the right nostril. Again, gently rotate the swab at least five times against the nasal wall. Remove the swab from the second nostril. <b>Caution: If the swab stick breaks during specimen collection, please use a new swab.</b></p>		
<p><b>Step 7</b> Dip the soft swab tip into the liquid. Rotate the swab for at least 15 seconds while pressing the head against the inside of the tube to dissolve the specimen in the liquid.</p>	<p><b>Step 8</b> Remove the swab from the extraction tube by squeezing the sides of the tube together and pulling the swab out to ensure most of the buffer remains in the tube. Discard swab in biohazard specimen bag.</p>	<p><b>Step 9</b> Screw on and tighten the nozzle onto the extraction tube.</p>	<p><b>Step 10</b> Shake the extraction tube vigorously to mix the specimen and the sample extraction buffer.</p>	
<p><b>Step 11</b> Open the foil pouch and take out the test cassette. Place the checked test cassette on a flat, clean surface. <b>CAUTION: Perform the test within 60 minutes after the foil pouch is opened.</b></p>	<p><b>Step 12</b> Add 3 drops of the solution from the specimen collection tube to each sample well of the test cassette.</p>	<p><b>Step 13</b> Set timer for 15 minutes. <b>CAUTION: Do not read the result beforehand, even if a line has already appeared at control region C. The test results will be invalid if read after 20 minutes.</b></p>	<p><b>Step 14</b> Please dispose of the test materials in a biohazard specimen bag with the household refuse. If there are local regulations, please follow them.</p>	<p><b>Step 15</b> Wash hands thoroughly after test completion.</p>

### INTERPRETATION OF RESULTS

Positive	Negative	Invalid
<p><b>SARS-Cov-2:</b> Two red lines appear. One red line appears in the control region (C), and one red line in the test region (T).</p> <p><b>Influenza A and Influenza B:</b> Three distinct colored lines appear. One colored line should be in the control region (C) and 2 colored lines should be in the Influenza A region (A) and Influenza B region (B).</p> <p><b>Influenza B:</b> Two distinct colored lines appear. One colored line should be in the control region (C) and another colored line should be in the Influenza B region (B).</p> <p><b>Influenza A:</b> Two distinct colored lines appear. One colored line should be in the control region (C) and another colored line should be in the Influenza A region (A).</p> <p><b>SARS-Cov-2, Influenza A and Influenza B:</b> Five distinct colored lines appear. One red line appears in the SARS-Cov-2 control region (C), and one red line in the test region (T). One colored line should be in the Influenza A region (A) and Influenza B region (B). The shade may vary, but if even a faint line appears, it should be considered positive. If you have a SARS-CoV-2 POSITIVE result, follow the guidance from your local State or Territory Health Department for guidance on confirmation testing if necessary, and if unwell seek medical assistance. If you have a Influenza POSITIVE result, individuals with a positive result or who are unwell are advised to consult a medical practitioner</p>	<p>Only one red line appears in the control region (C), and no line in the test region (T/A/B). The negative result indicates that there are no Novel coronavirus particles and influenza A/B in the sample, or the number of viral particles is below the detectable range. However, a negative result does not rule out COVID-19 and influenza A/B. <b>If you have symptoms like fever, cough and/or shortness of breath. Please retest in 1-3 days. You must continue following the applicable hygiene and distancing rules even with a negative result. If symptoms persist or if unwell please consult a medical practitioner for follow-up clinical care.</b></p>	<p>No red line appears in the control region (C). The test is invalid even if there is a line on test region (T/A/B). Insufficient sample volume or incorrect procedural techniques are the most likely reasons for control line failure. <b>Review the test procedure and repeat the test using a new test device. If the problem persists, discontinue using the test kit immediately and contact your local distributor.</b></p>