AquaVial™ PRO500 is designed as a water quality screening tool for water treatment professionals, for monitoring microbial levels, early detection of biofilm formation, and as a Legionella prevention screening tool.

AquaVial™ PRO500 kit is designed to detect any waterborne bacteria, or fungi, and a number of other waterborne microorganisms. To date, the kit has been tested and proven to detect *Mycobacteria spp.*, *E. coli* and other coliform bacteria such as *Klebsiella spp.*, *Legionella spp.*, *Bacillus spp.*, *Pseudomonas spp.*, *Staphylococcus spp.*, *Salmonella spp.*, and fungi such as *Aspergillus spp.*, *Candida spp.*

AquaVial™ PRO500 is not a substitution for *E. coli* or coliforms testing in areas where the government requires drinking water testing to be done at a certified lab.

AquaVial™ PRO500 detection limit is 500 colony forming units per milliliter (CFU/ml).

- Test kit contents:**
- 20 ml Disposable syringe
  - AquaVial™ PRO500 Disposable filter cartridge
  - Test vial (pink, black cap)
  - Saline solution vial (clear, white cap)

- Before you test:**
- ☑ Wear gloves to eliminate the chance of contaminating samples.
  - ☑ Check colour of test vial (pink, black cap) and shake. If not pink, do not use.

- Disposal/Recycling:**
- Filter cartridge and syringe are not reusable.
  - Dispose all components in the recycle bin or according to local regulations.

## Instructions


Follow exactly as written, any deviation may cause false positives or negatives!

### 1 Collect and filter sample


(i) Collect **40 ml** of sample water from the water source to be tested in a clean container.

40 ml sample water →

(ii) Unwrap the test kit and filter cartridge. Screw/attach the filter cartridge onto the syringe tip.


(iii) Position the filter cartridge tip into the water sample previously collected. Pull back the syringe plunger and draw in **20 ml** of sample water into the syringe.

(iv) Take the tip of the filter cartridge completely out of the sample water and continue to pull back the plunger completely to fill any remaining space inside the syringe with air. This step will ensure accurate results.


(v) Unscrew the filter cartridge completely from the syringe. Slowly press the plunger down to empty **all** the sample water and air from the syringe. Dispose this water as you will no longer need it for testing.

### 2 Wash sample

(iii) Position the tip of the filter cartridge into the saline solution vial. Pull back the syringe plunger and draw in **all** of the saline solution into the syringe.

NOTE: Washing with saline solution will not affect the bacterial count, but ensure your sample is chemical free and will yield accurate results.


### 3 Test sample

(i) Remove the cap of the test vial, and place it face up so that the inside is not in contact with any surface. Do not touch the inside of the vial as this could result in a contaminated sample.

### 3 Test sample

(ii) Slowly press the plunger down to add **all** the saline solution from the syringe into the test vial. Once the syringe is empty, slightly pull back the plunger to pull in **only a few drops of solution**, and then push down on the plunger to release the drops. Repeat to ensure liquids are thoroughly mixed. Do not repeat more than twice. Do not draw in all the test solution. You should not see solution inside the syringe during this process.

(iii) Recap the test vial. Do not reopen test vial after sealing. Use a Sharpie marker to label the test vial appropriately, indicating the water source. This is important when testing various sources simultaneously.

### 4 Observe colour and record results

Let the test vial sit at room temperature (20°C–30°C), and perform a quick check after 15 minutes. If no color change is observed after 15 minutes, leave the test vial for an additional 15 minutes. Immediately compare the colour within the test vial with the **Results Interpretation Chart** below and record test results.

#### Results Interpretation Chart | Bacteria level

Clean water ✓	Contaminated water ✗	Highly contaminated water ✗
< 500 CFU/ml	500 – 2,000 CFU/ml	> 2,000 CFU/ml
Pink	Purple	Clear

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