**KIT CONTENTS:**

- SCE-002-1: 2 Safe Check vials + Testing instructions
- SCE-024-1: 24 Safe Check vials + Testing instructions

**IMPORTANT SAFETY NOTES:**

- The media in the vials is non-toxic and poses no significant hazard when product is used as directed.
- DO NOT INGEST.
- For best results, test before the expiration date.

**DISPOSAL:** Vials with no colour change can be thrown in the trash. **TO RECYCLE:** Empty contents into the sink and recycle the cap and vial. Vials with a colour change can be thrown in the trash. Do not open vials. **TO RECYCLE:** Autoclave vials first, then empty contents into the sink and recycle the vial and cap.

**DIRECTIONS FOR USE:**

Safe Check is designed specifically for in-office dental waterline testing. Detect microbial levels that exceed 500 CFU/mL within 24 hours of incubation at room temperature. Safe Check is compatible with most dental waterline cleaners.

1. Twist the cap counter-clockwise to open the vial. Place the cap open side up on a clean surface. Do not touch the inside of the cap or vial to avoid contamination.
2. Fill the vial with sample water collected from the waterline. Do not exceed the 5 mL fill line.
3. Do not leave the vial open. Twist the cap clockwise to close the vial.
4. Shake the vial until the dry media deposited on the bottom of the vial has dissolved. The water sample will turn blue.

\* *TIP: Use a permanent marker to label the vial appropriately indicating the water source and date and time. This is important when testing multiple water sources simultaneously.*

**WARNINGS:**

- ▶ If the water sample does not turn blue after shaking, do not use it to test.
- ▶ If the blue colour changes within the first 5–10 minutes, this indicates chemical interference from a waterline cleaner, and the test cannot provide accurate results. Purge the system for at least 30 minutes to flush out the chemicals. If you have a cartridge/straw, remove it prior to purging. Test again the following day using a new vial.

5. Incubate the vial in an upright position for 24 hours at room temperature (20–30°C / 68–90°F).


*NOTE: Read results between 24–28 hours, afterwards the detection limit will decrease with time, reflecting a limit less than the recommended 500 CFU/mL.*

\* *TIP: Take a photo at the beginning and the end of the incubation time. This provides a visual record and allows for easier BEFORE and AFTER colour comparison.*

6. After the 24 hour incubation, check for changes in the vial colour. Compare the color of the test vial against the colour diagram below. If the test vial changes colour to purple, pink or turns clear, the microbial level exceeds 500 CFU/mL. This test is intended for rapid screening only.

7. Record results. Scan the QR code or visit [www.duwltesting.com](http://www.duwltesting.com) to access The Hub — your free online tool for scheduling, testing and recording results.

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