


Products for Dentistry



OSUNG Catalogue 2022•2023



DIAMOND BUR

Page. 224~285

We will always try to make products that are loved by dentists.

Products for Dentistry

OSUNG Catalogue 2022•2023

Products for Dentistry



OSUNG Catalogue 2022•2023



OSUNG MND CO.,LTD.

NEW Product



Contraster P.022



NiTi-periotome P.085



Palatal Wedge P.118



Micro Saw Shield P.141



Narrow Bone Removing Bur P.145



Ridge Split Bur P.145



Lateral Approach Bur P.145



Calibration Instrument P.181



Zirconia Removing Bur P.296



Orthodontic Tweezer P.331



Crimpable Hook 331



Orthodontic Cassette 349

Products for Dentistry

OSUNG Catalogue 2022/2023



Designed Instruments by a clinician.

Stellar is a new brand name for the worldwide market.
The core value of the stellar instrument is user-friendly design. It is selected,
modified, and evaluated by the dentists' group.

Technological Innovation by OSUNG

**We study every technology related to dental instruments.
And we aim for the top.**

We have a clear goal.
It is becoming to the top in quality, service and technology.

And we keep doing our endeavors for
mechanization of manual labor process
enhancement of manufacturing system
development of user-friendly design
fusion with new concept & idea
application of the latest technology

And we will grow up as
a specialist
a representative manufacturer
a technological leading company
an intrinsic value firm
in the dental industry.

Contents

09	Diagnostic
37	Periodontal
73	Surgery
129	Implant
169	Restorative
201	Endodontic
221	Prosthodontic
307	Orthodontic
343	Instrument management
361	Laboratory
369	Index

DIAGNOSTIC

Diagnostic	Explorer	012
	Mirror Handle	016
	Mirror	018
	Oversized Mirror	018
	Plastic Mouth Mirror	018
	Spoon Mirror	019
	Para Mirror	019
	Photo Mirror(Glass)	020
	Photo Mirror(Matal)	021
	Handle Photo Mirror(Matal)	021
	Contraster	022
	FF-Photo Mirror	023
	Probe	026
	Implant Probe	027
	EX-Probe	030
	Locking Plier	031
	Tweezer	031
	Retractor	032
	Mouth Prop	036

PERIODONTAL

Ultrasonic Scaling	Ultrasonic Scaler Tip	040
Root Planing & Curettage	Sickle Scaler	041
	Micro Sickle Scaler	044
	Gracey Curette	046
	Rigid Gracey Curette	054
	Mini Five Curette	056
	Mini Five Rigid Curette	058
	After Five Rigid Curette	058
	Universal Curette	060
	Special Curette	062
	Implant Curette	063
	Chisel Scaler	063
	Hoe Scaler	064
	File Scaler	065
	Periodontal File Scaler	065
Option	Sharpening Stone	066
	Perio Scaling Kit	066
Periodontal Treatment	Manual	067

SURGERY

Preparation for Surgery	Surgical Suction Tip	076
	Extension Hose	076
	Extension Hose Adapter	076
	Frazier Suction Tip	076
	Suction Tip	077
	Surgical Drape	078
	Wrapping Cloth	078
	Towel Clamp	078
	Anesthesia Syringe	079
	Scalpel Handle	080
	Composite Scalpel Handle	081
Oral Surgery	Periosteal Elevator	082
	NiTi-Periotome	085
	Periotome	086
	Root Picker	087
	Luxating Elevator	088
	Elevator	090
	Extraction Forceps(Adult)	092
	Extraction Forceps(Pedo)	096
	Surgical Curette	098
	Bone Rongeur	102
	Nipper	102
	Bone File	103
	Mallet	103
	Hemostat	104
	Needle Holder	105
	Anatomic Dressing Forceps	106
	Tissue Plier	107
	Scissors	107
Periodontal Surgery	Periodontal Knife	110
	Periodontal Chisel	110
	Periodontal Surgical Curette	111
Maxillofacial Surgery	Periosteal Elevators for Maxillofacial Surgery	112
	V-Notch Periosteal Elevator	112
	Channel Retractor	113
	Spatula Periosteal Chisel	113
	Retractor	114
	Tunneling Instrument	116
	Palatal Wedge	118
Simple Extraction	Manual	119
Excision of Torus	Manual	123

IMPLANT

Fixture Implantation	Lindemann Drill	132
	Implant Depth Gauge	134
	Caliper	134
Bone Graft	Bone Spreader	135
	Bone Expander Hand Kit	136
	Bone Expander Engine Kit	137
	Micro Saw Shield	140
	Micro Saw	142
	Trephine Bur	144
	Narrow Bone Removing Bur	145
	Ridge Split Bur	145
	Lateral Approach Bur	145
	Surgi-Drill Stand	145
	Convex Osteotome	146
	Concave Osteotome	147
	Bone Scraper	148
	Block Bone Clamp	148
	Bone Collect Chisel	149
	Bone Collector	149
	Hexa Wrench	149
	Bone Mill	150
	Bone Crusher	150
	Bone Crusher Mallet	150
	Bone Syringe	151
	Bone Well	152
	Bone Carrier	152
	Bone Packer	152
	Membrane Forceps	153
	Sinus Rongeur	153
Sinus Lift	Crestal Approach Kit	154
	Lateral Approach Kit	156
	Sinus Lift	158
	Bone Screw	160
	Bone Tack	161
	Bone Tack Offset Holder	161
Implant 2nd Surgery	Hand Tissue Punch	162
	Tissue Punch	163
Implant Crown Setting & Maintenance	Screw Removal Kit	164
	Implant Curette	166
Option	PRF & GRF Box	167
	New Product	168

RESTORATIVE

Cavity Preparation	Excavator	172
	Gingival Retractor	175
	Margin Trimmer	176
Amalgam	Amalgam Carrier	177
	Amalgam Well	177
	Amalgam Plugger	178
	Carver	179
	Amalgam Burnisher	180
Composite Resin	Measuring Instrument	181
	Placement	181
	Composite Instrument	182
	Composite Instrument Kit	190
Amalgam Filling	Manual	191
Resin Filling	Manual	197

ENDODONTIC

Cavity Preparation	Intraligamentary Syringe	204
	Endodontic Explorer	204
	Broach Holder	204
	Endodontic Excavator	205
Endodontic	Spreader	206
	Endo Locking Plier	206
	Endo Ruler	206
	Endo Box	206
	Root Canal Plugger	207
Rubber Dam Instrument	Rubber Dam Set	208
	Rubber Dam Punch	209
	Rubber Dam Plier	209
	Rubber Dam Frame	209
	Rubber Dam Clamp Stand	210
	Rubber Dam Clamp	211
	OrthoMTA Carrier	212
	OrthoMTA Syringer	212
	OrthoMTA Plugger	212
Root Canal Treatment	Manual	213

PROSTHODONTIC

Dental Diamond Bur	Dental Diamond Bur	224
	Bur Kit	253
	My Bur Kit Case	279
	Bur Block	286
Impression	Spatula	287
	Paper Holder	287
	GingiCord Packer	288
	Gingimaster Injector	289
	Impression Tray	290
	Agar Syringe	295
Crown Removing & Setting	Zirconia Removing Bur	295
	Crown Remover	296
	Crown Forceps	298
	Crown Gripper	298
Articulator	Occlusal Plane Plate	299
	Willis Gauge	299
	Occlusal Rim Plate	300
Prosthodontic Treatment	Manual	301

ORTHODONTIC

Orthodontic Diagnosis	Photo Mirror(Glass)	310
	Photo Mirror(Matal)	311
	Handle Photo Mirror(Matal)	311
	FF-Photo Mirror	312
	Orthodontic Strip	314
	Strip Holder	314
	Orthodontic Arch Wire	316
	Coil Spring	317
	Orthodontic Wire	317
	Splint PET	319
	Bracket Positioning Gauge	320
	Bracket Positioning Height Gauge	320
	Band Preparation Instrument	321
	Ligature Tucker Instrument	321
Orthodontic Instrument	Hook-Crimping Plier	322
	Bracket Remover	322
	Wire Bending Plier	323
	Band Remover	324
	Tying and Holding Plier	327
	Wire Cutting Instrument	328
	Aligner Plier	329
	Orthodontic Tweezer	331
	Crimpable Hook	331
	Orthodontic Instrument Cassette	332
	Bos Sunny Orthodontic Plier Kit	334
	Bos Sunny Surgical Instrument Kit	335
	Metal Strip Holder	336
	Band Cutting Scissors	336
	Fixator	336
Orthodontic Treatment	Manual	337

INSTRUMENT MANAGEMENT

Instrument Sterilization	Instrument Cassette	346
Instrument Storage	Instrument Tray	351
	Instrument Color-Coding Item	352
	Chairside Management	353
Instrument Sterilization	Endo Ruler	354
	Endo Box	354
	Surgical Drape	355
	Wrapping Cloth	355
	Instrument Pouch	355
	Scaler Tip Stand	356
	Scaler Tip Torque Wrench	356
	Sharpening Stone	356
	Bur Block	357
	My Bur Kit Case	357
	Surgi-Drill Stand	358
Unit Chair Accessory	Cotton Pellet Dispenser	359
	Cotton Pellet Push Device	359

LABORATORY

LAB Products	Casting Machine	364
	P.K.Thomas	366
	Waxing & Carving Instrument	367
	Spatula	368

INDEX 369

The structure and feature of the catalog

Focused on better use & search availability and user convenience.

The instruments are classified by clinical field and the basic preparation is shown as an example.

- 1 **Way to find a certain instrument**
Instruments are classified by clinical field. Refer to the chapter title.
- 2 **Feature of instrument**
The main features of each instrument are described with pictures.
- 3 **How to order**
Refer product code on your order always.
- 4 **Detailed use of instrument**
Video clip provided for detailed instruction.
- 5 **New or recommended product**
New product or recommended product are indicated with an icon.
- 6 **Magnified picture**
Enlarged picture of working part provided for detailed look.
- 7 **Technical information**
The information of basic technology used for instrument manufacturing is shown.

Periodontal Sickle Scalers

- Has cutting edge at both side as remover of supra-gingival calculus
- Tip end is pointed.
- There are curved & straight types.

Curved Sickle Scaler

- both cutting edges are focused to one point according to shape of rounded curved blade.

Straight Sickle Scaler = Jacquette Scaler


- both cutting edges are focused to one point according to shape of straight blade.
- Jacquette scaler

Sickle Scaler, Silicone Handle

Autoclavable

NEW PRODUCT
2LSH5-33

To remove calculus of interproximal & cervical in anterior



Video Clip

NEW PRODUCT
2LSH6-H7

Anterior, Premolar
to remove calculus of interproximal

NEW PRODUCT
2LSH6-33

Anterior
One side is standard size and the other side is small size


Science & Technology

The selection is important

For case and heat treatment is very important.
Make a good product if there is no analytical metal crystal structure and heat treatment?

Surprisingly, even the famous instrument manufacturers in developed countries are not able to secure their analytical skills. However, we have world-class technology and know-how in analytical engineering for metal as a result of many efforts for a long time.

Figure. SAM image for checking the crystal grain size, solid solubility of carbide and etc...



Products for Dentistry

OSUNG Catalogue 2020/2021

Diagnostic


Products for Dentistry

OSUNG Catalogue 2022°2023



DIAGNOSTIC

Explorer	012
Mirror	016
Photo Mirror	020
Probe	024
EX-Probe	028
Tweezer	029
Retractor / Lip Wider	030
Mouth Prop	032
Suction Tip	032




Explorers

Detail of Explorer



It is to detect dental caries or calculus with a sharp point at the end

Type of Explorers



Straight Type
Subgingival calculus and caries




Shepherd Hook Type
Subgingival calculus and caries




Curved Type
Deep pocket and furcation



Urban-Type
Root surface on anterior detection and facial & lingual on posterior



11/12-Type
Calculus on anterior and posterior. Inspect root surface.



Endo Type
Root canal entrance



Pigtail Type & Cowhorn Type
Detection for periodontal pocket not deeper than a third dental root cervical or gingival sulcus calculus.



Explorers

The explorer provides the tactile information to the clinician's fingers and is used to locate calculus deposits, tooth surface irregularities, defective margin on restorations, decalcified areas and carious lesions.

Explorer_Silicone Handle



BEST

2EXD5H

- Shepherd Hook Type + Urban Type
- Rigid type.
- For subgingival calculus and caries




2EXD5W

- Shepherd Hook Type + Urban Type



2EXD17-23

- Shepherd Hook Type + Urban Type



2EXD5-8

Explorer, EXD5-8



2EXDG16

- Used to detect the orifice of the canal.



Explorers

Explorer_Metal Handle

BEST

EXD5H

- Shepherd Hook Type + Orban Type
- Rigid type. For subgingival calculus and caries



BEST

EXD5W

- Shepherd Hook Type + Orban Type



EXD17-23

- Shepherd Hook Type + Orban Type



EXD54-17H

- Shepherd Hook Type + Orban Type



EXDG16

- Used to detect the orifice of the canal.



Video
Clip

EXD11-12

- For subgingival root examination



EXD3CH

- Cowhorn Type



EXD2

- Pigtail Type



Explorers

Explorer_Semi-Silicone Handle



EXDK

- Korean explorer
- It has two different types of tip at the ends for multiple uses. One is rigid and the other one is flexible.

1
RIGID TIP
To remove cement

2
SILICONE HANDLE
To lessen wrist fatigue

3
METAL HANDLE
Need delicate sense

4
FLEXIBLE TIP
For checking caries, calculus & margin

Explorer_Metal Handle

EXS96

Explorer, EXS96



EXS6
• Straight type
6.5mm

EXS6XL
• Straight type
12.5mm

EXS6L
• Straight type
9mm

EXS3
• Curved Type
• For detecting calculus of shallow periodontal pocket or gingival sulcus.

EXS23H
• Shepherd Hook Type

EXS54H
• Shepherd Hook Type

EXS23W
• Shepherd Hook Type

EXS3A
• Curved Type
Useful for detection of calculus on furcation and deep periodontal pocket.

Mirrors

Metal Handle

MHS (10pcs)

- 10 piece packet
- SS Type



MHC (10pcs)

- 10 piece packet
- CS Type




Fine materials & high-tech engineering process

Glass : High quality Japanese glass

Waterproof : Sealed by silicone to protect the reflection layer of glass

Durability : Double welded




What is front surface mirror ?

Front Surface Mirror

Front surface mirror


- Reflection layer exists on the front side of the mirror.
- As it has no secondary image, it is mainly used in endodontic treatment.



Rear Surface Mirror

Plain mirror (rear surface mirror)

- Reflection layer exists on the rear side of the mirror.
- It is durable and economic in price, but it has ghost image.




Choice of size of mirror

	#3	Dia 20mm, Suitable for children
	#4	Dia 22mm, Popular size
	#5	Dia 24mm, Used for wider view

Mirrors

Silicone Handle_Simple Stem (SS Type)

Autoclavable



2MHS1 (5pcs)

- 5 piece packet



2MHS2 (5pcs)

- 5 piece packet




2MHS3 (5pcs)

- 5 piece packet



Silicone Handle_Cone Socket (CS Type)

Autoclavable



2MHC1 (5pcs)

- 5 piece packet



2MHC2 (5pcs)

- 5 piece packet



2MHC3 (5pcs)

- 5 piece packet



Mirrors

Mouth Mirror(Rear Surface)

SS Type

DMSS3 (12pcs)

- Simple Stem No.3 (20mm)
- 12pcs

DMSS4 (12pcs)

- Simple Stem No.4 (22mm)
- 12pcs

DMSS5 (12pcs)

- Simple Stem No.5 (24mm)
- 12pcs



CS Type

DMCS3 (12pcs)

- Cone Socket No.3 (20mm)
- 12pcs

DMCS4 (12pcs)

- Cone Socket No.4 (22mm)
- 12pcs

DMCS5 (12pcs)

- Cone Socket No.5 (24mm)
- 12pcs



Mouth Mirror(Front Surface)

One-Sided

DMFSS5 (12pcs)

- Simple Stem No.5 (24mm)
- 12pcs

DMFSS4 (12pcs)

- Simple Stem No.4 (22mm)
- 12pcs

DMFCS5 (12pcs)

- Cone Socket No.5 (24mm)
- 12pcs

DMFCS4 (12pcs)

- Cone Socket No.4 (22mm)
- 12pcs

Double-Sided

DMDSS4 (5pcs)

- Simple Stem No.4 (22mm)
- 5pcs

DMDCS4 (5pcs)

- Cone Socket No.4 (22mm)
- 5pcs



Oversized Mirror

Autoclavable(Limited to 10 cycles max.)




DMS39

- Front surface
- Big size of Dia 39mm
- Convenient to use as the angle of the mirror can be adjustable.

Plastic Mouth Mirror

Autoclavable

DMAPA



Mirrors

Spoon Mirror

Rear Surface



DMPM

- Rear Surface mirror
- Use for multi-purpose



DMEM

- Rear Surface mirror
- Use for multi-purpose




Front Surface



DMFPM

- Front surface mirror
- Use for multi-purpose




Rear Surface



DMFEM

- Front surface mirror
- Use for multi-purpose



Front Surface

Para Mirror

Autoclavable

DMPRA

- Anterior

DMPRP

- Posterior



To check parallelism of an insert line of prosthetic appliance during preparation

Useful to check direction & implant site of a fixture during implant surgery



Posterior



Anterior

Photo Mirrors-Glass

Photo Mirrors-Glass


Ours has reflexivity of 94-97%

General Mirror : 85%
Other photo mirror for dentistry has reflexivity of 65-80%

Best quality and most competitive prices

Best quality and most competitive prices!

Type of Coating	Value of Reflectance
HR Coating	94~97%
German Ultra Bright Coating	Around 95% (Not Clearly Announced)
America and German Rhodium Coating and Titanium Coating	70~80%



The world-best reflectance having a range of 94-97% in the visible spectrum!


- Highly-reflective coating technology
- Coating technique that strengthens anti-scratch capability

Glass

• Has clear image but breakable

DME6G

- Occlusal
- X-Large (Adult Size)
- Unit(mm)




유리거울
Glass




DME1G

- Occlusal
- Large (Adult Size)
- Unit(mm)




DME3G

- Occlusal
- Medium (Adult Size)
- Unit(mm)




DME5G

- Occlusal
- Small (Pedo Size)
- Unit(mm)



DME2G

- Buccal
- Unit(mm)



DME4G

- Lingual
- Unit(mm)





Photo Mirrors-Metal

Metal

• Has less clear image than glass but not breakable.


DME1

- Occlusal
- Large (Adult Size)
- Unit(mm)




DME3

- Occlusal
- Medium (Adult Size)
- Unit(mm)




DME5

- Occlusal
- Small (Pedo Size)
- Unit(mm)




DME2

- Buccal
- Unit(mm)



DME4

- Lingual
- Unit(mm)



Handle Photo Mirror (Metal)

BDMHL

- Occlusal, Large
- Unit(mm)



BDMHM

- Occlusal, Medium
- Unit(mm)



BDMHS

- Occlusal, Small
- Unit(mm)



BDMHLT

- Lateral
- Unit(mm)




- Ordinary photo mirrors should be held in the middle area with fingertips as they have no handle. Handle photo mirrors can be held conveniently as they have a silicone handle.
- Sufficient mirror length covering the final tooth. As the handle keeps the surface of the mirror off the floor, the mirror becomes free from scratch.
- Safe metal type made of stainless steel and does not break.

Contraster

Contraster

Photo Contrastors are accessories for intraoral photos to form a background when taking macro shots of the anterior region. Photo Contrastors allow photos of the teeth and gingiva to be taken aesthetically. One side can be used for adults and the other side for children.

NEW
CTRM
Contraster
• CTRM



Guide for using Photo Contrastors



1. The lips cover the teeth, so open the lips outward using a Side Winder(RTSWH).
2. Position the Photo Contrastor at the rear of the tooth to be taken.
3. Fix it in the proper state and take shots.
4. Remove the Photo Contrastor after taking shots.

Before using the Photo Contrastor



After using the Photo Contrastor

Sterilization instruction for Photo Contrastors

- High-pressure sterilization is available. (Repeated sterilization may cause damage to the product.)
- Make sure to wrap in a sterilization pouch when disinfected with high-pressure sterilization.
- Disinfecting the Contrastors without being wrapped in a sterilization pouch may cause damage to the products.
- Please make sure not to touch the hot wire when disinfected with high-pressure sterilization. Exposure to excessive heat may cause damage to the products.
- High-pressure steam sterilization may cause slight water stains, but it can be wiped and used. The product itself is not defective, and the water stains can be removed by wiping with a soft cloth.




Photo Mirror_FF-photo

FF-Photo(Fog-Free Intraoral Photo Mirror)



Good image focus is difficult because mirrors tend to fog in mouth due to breathing and the light from a dental unit may be reflected to some extent by the mirror.

Even though these mirrors have already been specially coated to resist fogging, a combination of additional steps such as air syringes and hot water dipping are frequently taken in the clinic for defogging.

The air syringe or hot water dipping is only briefly effective and must be continuously applied, particularly if multiple photographs need to be taken.

Both methods require additional staff to help defog.

Conventional lighting may be lost by the mirror making a clear problem of reflection which needs to be overcome by an alternative light source.

We would like to introduce a new device which defogs and improves focus with fewer assistants and less photography time than the conventional method.

FF-photo used



No FF-photo used

FF-photo used

Photo Mirror, FF-Photo



Photo Mirrors_FF-Photo Slide

FF-Photo(Fog-Free Intraoral Photo Mirror)_Metal

- You must use a Dedicated Mirror made of metal.

FF-Photo was introduced to the Journal of Clinical Orthodontics (2008,2), an international journal of orthodontics as a patented invention.
 • It is commercially available products based on research and development data from the team of department of dentistry at the Catholic Medical College.

DMBF-220


FF-Photo Slide

- Components
- ① FF-Photo Body
- ② Micro 5-pin USB charging cable
- ※ Photo Mirror sold separately

Chargers are not provided. You can use your cell phone charger and more.




◀ Dedicated Mirror sold separately ▶



DMBFL

FF-Photo Mirror


- Occlusal
- Large
- 단위(mm)



DMBFM

FF-Photo Mirror


- Occlusal
- Medium
- 단위(mm)



DMBFS

FF-Photo Mirror

- Occlusal
- Small
- 단위(mm)



DMBFLT

FF-Photo Mirror

- Lateral
- 단위(mm)

Photo Mirrors_FF-Photo Lever

FF-Photo(Fog-Free Intraoral Photo Mirror)

- It is a product that complements the fastening part that can be compatible with one another for all metal and glass mirrors of our company and other companies.

DMBF1-220

FF-Photo Lever

- Components
- ① FF-Photo Body
- ② Micro 5-pin USB charging cable
- ③ Thickness control panel
- ※ Photo Mirror sold separately

Chargers are not provided. You can use your cell phone charger and more.



◀ Various types of glass mirrors and metal mirrors are available. ▶

Glass Mirrors and Metal Mirrors P.21-21



TIP

The fastening parts are basically designed to allow the use of glass mirrors, and metal mirrors can be used when necessary using the [thickness control panel].




Probes


Type of Probes



Novatec
Effective to approach molar area.
3-6-9-12



Goldman-Fox
Fiat tip
1-2-3-5-7-8-9-10



Williams
1-2-3-5-7-8-9-10



WHO
0,5 ball type's tip
3.5-5.5-8.5-11.5




Nabers
Useful for detecting root furcation involvement with long and bended tip
3-6-9-12




CP 2
Color graduation in 2-4, 6-8, 10-12mm
2-4-6-8-10-12



CP 11
Color graduation in 3-6, 8-11mm
3-6-8-11




CP 12
Color graduation in 3-6, 9-12mm
3-6-9-12




UNC 15
Color graduation in 1-15mm
1-2-3-4-5-6-7-8-9-10-11-12-13-14-15

Practice

Ball-End ?
Ball-End probe increases tactile sensitivity and gives a more accurate picture of gingival pocket depth. Also it doesn't damage tissue.



Probes



1. Determine the consistency of the gingival tissue
2. Survey sulcus and pocket depths
3. Measure clinical attachment levels
4. Measure the width of attached gingiva.
5. Verify whether bleeding and/or purulent exudate is in presence

Silicone Handle_Single-Ended

Autoclavable It measures the depth of periodontal pocket, height of attachment, anatomical shape and bleeding of gingiva.

2BPCP8



2BPUNC15



2BPCP12



Implant Probe

Autoclavable

3IBPCP12

3IBPCP12-5T Plastic Tip Only (5pcs)



It fits SS type mirror handle. The packet doesn't include any handle. 17 page

NEW

3IBPCP2

3IBPCP2-5T Plastic Tip Only (5pcs)



It fits SS type mirror handle. The packet doesn't include any handle. 17 page

NEW

3IBPCP2-12



Practice

- Excellent discrimination (yellow, black colors)
- More flexible than metal, which allows for better contour-following of the alveolar bone.
- Autoclaves available (less heat resistant than metals, but no problem with repeated autoclave use).
- You can measure the distance after bending to suit your needs(available multiple times).




Probes

Metal Handle_Single-Ended

BPWHO

Probe, WHO
• WHO



BPW
Probe, PW
• Williams



BPWB
Probe, PWB
• Williams



BPCP2
Probe, CP2



BPCP8
Probe, CP8



BPCP10
Probe, CP10



BPCP11
Probe, CP11




BPCP12
Probe, CP12



BPUNC15
Probe, UNC15

Probes



1. Determine the consistency of the gingival tissue
2. Survey sulcus and pocket depths
3. Measure clinical attachment levels
4. Measure the width of attached gingiva
5. Verify whether bleeding and/or purulent exudate is in presence

Metal Handle_Double-Ended

It measures the depth of periodontal pocket, height of attachment, anatomical shape and bleeding of gingiva.

BEST BPGF-W

Probe, GF-W
• Goldman-Fox + Williams




BPQ2N

• Nabers
• To measure horizontal and vertical pocket depth of multirooted teeth in furcation areas.



BPCP12N


Probe, 12N
• Novatec with CP12



Left tip has curved 90 is for anterior and helps to measure even little opened mouth. Right tip is adopted BPCP12 popular.

BPW15

Probe, PW15
• Probe X



The both of tips are bent 15 at 3mm point from tip end and this helps to measure depth following curved teeth in status of a little opened mouth.

EX-Probes

Ex-Probe is an excellent combination instrument for dental diagnostics. It has an explore tip on one end and a probe tip on the other.

Metal Handle_Double-Ended



XP23-WHO

Ex-Probe, XP23-WHO
• EXS23 with WHO



XP23-W

Ex-Probe, XP23-W
• EXS23 with Williams



XP23-8

Ex-Probe, XP23-8
• EXS23 with CP8



XP23-12

Ex-Probe, XP23-12
• EXS23 with CP12

Tweezers· Locking Pliers

stainless steel Made
Precise, steady serration on the tip gives a firm grip when taking material.

Tweezer


PC1

• Length : 164mm (± 5mm)



PCW150

• Length : 151mm (± 5mm)
• Wide grip for easy taking





PCU155

• Length : 155mm (± 5mm)
• Has shape of one bending angle.
• Useful for suturing for surgical knot or taking surgical sponge.
• Useful for placing transplant tissue.




Practice



Locking Plier


PCUL155

• Length : 155mm (± 5mm)
• Easy to take material as it has locking function.



EAPCUL155

• Length : 155mm (± 5mm)
• Easy to take material as it has locking function.
• Useful for moving material as it has furrow inside the tip.



Retractors · Lip Widens

Lip Retractor with Anterior Shield

1. Will prevent water from jumping up to the face and give the patient comfort feeling during the process of anterior maxillary tooth preparation.
2. Adult size

RTLAS (2pcs)

- 2pcs
- Autoclave is not available



Practice



1. Press retractor up & down by hand and retract cheek with mirror in order to put one side and then put the other side.
2. Wash and sterilize/disinfect with EO gas or antiseptic solution only after use.

Plastic Lip Wider

- 134°C It is made of plastic but autoclavable.
- Sterilize at autoclave (134°C)

- Place product on flat surface to prevent deformation during autoclave.
- Keep goods away from heat-source of autoclave during sterilization.

RTCPD1 (2pcs)

- Large Size, 2pcs
- Size A:130mm/B:92mm

RTCPD2 (2pcs)


- Medium Size, 2pcs
- Size A:118mm/B:85mm

RTCPD3 (2pcs)

- Small Size, 2pcs
- Size A:97mm/B:70mm



This Lip Wider is used for tooth whitening treatments



Retractors

Plastic Lip Wider

- 134°C It is made of plastic but autoclavable.
- Sterilize at autoclave (134°C)

- Place product on flat surface to prevent deformation during autoclave.
- Keep goods away from heat-source of autoclave during sterilization.

RTCPS1 (2pcs)

- Lip Wider
- Large Size, 2pcs
 - Size A:120mm/B:53.5mm

RTCPS2 (2pcs)

- Lip Wider
- Small Size, 2pcs
 - Size A:114mm/B:42mm



For having lip wide for intra-oral photography.

Side Wider


- 134°C Plastic products, but autoclave sterilizable
- 134 degrees autoclave sterilizable
- An assistant tool for helping to understand the oral health when capturing intraoral photos
- Do not place it near the internal hot line of the autoclave during sterilization
- Select the appropriate shape and size for the area you want to shoot.
- It's made transparent and solid, so there is no deformation.

NEW RTSWM (2pcs)

- Side Wider M
- Medium size, 2pcs
 - Size A : 139mm / B : 51mm

NEW RTSWSS (2pcs)

- Side Wider S
- Small size, 2pcs
 - Size A : 135mm / B : 41mm



NEW RTSWH (4pcs)

- Side Wider H
- Hook, 4pcs
 - Size A : 86mm / B : 15mm



Side Wider disinfection instructions


- High-pressure sterilization is available.
- Make sure to wrap in a sterilization pouch when disinfected with high-pressure sterilization.
- Disinfecting the Side Wider without being wrapped in a sterilization pouch may cause damage to the product.
- Please note that repeated disinfection may cause the product to be cloudy or crack.

Retractors

Lip Retractor


Parkman Design

- Designed to give you a better view of the oral cavity as a self-retaining lip retractor.




NEW
RTLW8050

Lip Wider
• Large Size




NEW
RTLW7045

Lip Wider
• Medium Size



NEW
RTLW6540


Lip Wider
• Small Size



Practice

- Less pain for the patient when opening the mouth with thicker material than other products.
- Designed for the pressure point of the lips for a better view of the oral cavity.
- Stainless steel for easy cleaning and sterilization

※ It is recommended to apply Vaseline, etc., lightly to the lips of patients when using a retractor.



Retractors

Lip and Cheek Retractor


RTCRL

Lip Retractor, CRL



RTCRC

Cheek Retractor, CRC
• Columbia



Lip and Tongue Retractor

- For earning easy impression / possible to use during implant or surgical surgery

RTLS

- Using Lip retractor, you don't need to use multiple mirrors during earning upper jaw impression.



RTTG


- Using both lip retractor and tongue retractor simultaneously, you can reduce amount of work during earning lower jaw impression. It is available during an implant operation or a surgical operation.



Mouth Props · Suction Tips



Mouth Prop

Autoclavable



Suction Tip (Stainless Steel)

• The body and the tip are made with stainless steel which would not cause any discoloration or peeling away of coating. This suction tip can be used semi-permanently.



Products for Dentistry

OSUNG Catalogue 2022/2023

Periodontal


Products for Dentistry

OSUNG Catalogue 2022*2023



PERIODONTAL

Ultrasonic Scaling	Ultrasonic Scaler Tip	036
Scaling	Sickle Scaler	037
Root Planing & Curettage	Gracey Curette	042
	Rigid Gracey Curette	050
	Mini Five Curette	052
	Universal Curette	054
	Special Curette	056
	Implant Curette	057
	Chisel Scaler	057
	Hoe Scaler	058
	File Scaler	059
Option	Sharpening Stone	060
	Perio Scaling Kit	060
Periodontal Treatment	Manual	061



Ultrasonic Scaler Tips

Ultrasonic scaler tip made by 100% Korean technique



- Improved quality thru structural analysis
- Lowered price thru process improvement
- No damage to handpiece.



Provide a Different Service



Recommend to change scaler tip when tip is worn by 50%.



Autoclavable



• Torque Wrench

USETW

- For EMS and SATELEC tip
- Free from infection as the tip do not touch hand during connecting to handpiece.
- Do not use Dry Heat

SATELEC EMS

Ultrasonic Scaler Tip

SATELEC



• Universal

USS1

- Compatible with SATELEC "No.1"



• Supragingival & Subgingival

USEP

- Compatible with EMS "Type P"



• General deposit removal

USEA

- Compatible with EMS "Type A"




• Interproximal & Subgingival

USEPS


- Compatible with EMS "Type PS"

Sickle Scalers



Curved Sickle Scaler

- Both cutting edges are focused on one point following the shape of rounded curved blade.



Straight Sickle Scaler = Jacquette Scaler

- Both cutting edges are focused on one point following the shape of straight blade.
- Jacquette scaler

Sickle Scaler_Silicone Handle

Autoclavable

BEST

2LSH5-33

Sickle Scaler, H5-33

- To remove calculus of interproximal & cervical in anterior.



시클 스케일러
사용동영상



2LSH6-H7

Sickle Scaler, H6-H7

- Anterior, Premolar
- To remove calculus of interproximal



2LSJAC30-33

- Anterior
- One side is standard size and the other side is small size



BEST

2LSJAC31-32

- Posterior
- Standard sized Jacquette tip



2LSJAC34-35

- Posterior
- Small sized Jacquette tip



Sickle Scalars

Sickle Scaler_Metal Handle

BEST

LSH5-33

- To remove calculus of interproximal & cervical in anterior.



Video
Clip



LSH6-H7

- Anterior, Premolar
- To remove calculus of interproximal



LSJAC30-33

- Anterior
- One side is standard size and the other side is small size



BEST

LSJAC31-32

- Posterior
- Standard sized Jacquette tip



LSJAC34-35

- Posterior
- Small sized Jacquette tip



Sickle Scalars

Sickle Scaler_Metal Handle

LS204

- To remove calculus on proximal of supragingiva in posterior.
- For removal of heavy calculus.



LS204S

- To remove supragingival calculus between the tooth in posterior.



LS204SD

- For posterior
- Standard sized Jacquette tip



LSSCM152

- Useful to remove stain with spoon shaped knife.




Science & Technology

The selection of proper steel and heat treatment is a very important thing for instruments. But how do you make a right product if there is no analytical technique for metal crystal structure and heat treatment?

Surprisingly, even the famous instrument manufacturers in developed countries are not able to secure these analytical skills. However, we have world-class technology and know-how in analytical engineering for metal as a result of many efforts for a long time.

Figure. SEM image for checking the crystal grain size, solid solubility of carbide and etc..



Sickle Scalers

Towner (U15)

- For removal of heavy calculus of interproximal, buccal and lingual

BEST

LSU15-30

- Anterior
- Towner-Jacquette
- For removal of heavy calculus



LSU15-33

- Anterior
- Towner-Jacquette
- For removal of heavy calculus




Micro Sickle Scaler

- Elongated terminal shank with sharp and slender blade.
- Very useful for tight proximal surfaces.
- Used for removal of supragingival calculus for all teeth surfaces.
- Used for removal of subgingival calculus near the edge of gum.

2LSMS1-2

- Silicone Handle

134°C 오토클레이브 사용가능



LSMS1-2

- Metal Handle




Mini Sickle Scaler

- Used for removal of supragingival calculus of all tooth surfaces, especially of proximal surfaces.
- Used for removal of subgingival calculus near the edge of gum.
- Used for removal of filling, adhesive and cement excesses.
- Used for removal of calculus and granulation tissue in flap operation.

2LSMS11-12

- Silicone Handle

134°C 오토클레이브 사용가능




LSMS11-12

- Metal Handle



Curettes


Four Types of Curettes



A curette is the most suitable instrument for removing subgingival calculus and for doing root planing. It is specially useful for deep periodontal pocket or furcation lesion.

A curette is designed to avoid tissue trauma & damage to the teeth.


- Two cutting edges meet to make round toe and cross section is round shape.
- Face & lateral surface meet to make cutting edge.



1. Gracey Curette

Used for fixed specific area according to each instruments

The lower cutting edge is used only and have 70° angle on the basis of terminal shank



2. Universal Curette


Used for all root conditioning

Both cutting edges are 90° angle on the basis of terminal shank

Gracey Curette and Universal Curette

	Gracey Curette	Universal Curette
Area for use	Designated to be used for the specific area.	Available to all area & root conditioning as one curette
Blade angle	70°	90°
Use cutting edge	One lower cutting edge	Both cutting edges
Curve of cutting	Curved toward the end & side of a tip	Curved toward the end of a tip
How to use	Use to each area according to curette number.	Narrow terminal angle & short length is for anterior and wide terminal angle & long length is for posterior.

Gracey Curettes-Standard



Gracey Curette
GR11-12 VS 15-16

It is difficult to place connecting part of GR11-12 on the mesial surface of the lower posterior parallelly. The angle of GR15-16 reaches the mesial surface of posterior when it is fixed in the mouth in front of a patient.

It has a blade that is laterally offset by 70 degrees relative to the shank and has a lower cutting edge and an upper non-cutting edge. Only one side of the blade can be used.

Standard Gracey Curette_Silicone Handle

Autoclavable

※OSUNG's silicone handle gives no-stress on wrist and provides an excellent grip.



BEST
2CGR1-2

• Anterior



Video Clip



BEST
2CGR11-12

• Mesial surface of all posterior teeth



Video Clip




BEST
2CGR13-14

• Distal surface of all posterior teeth



Video Clip

Gracey Curettes-Standard



Gracey Curette
GR13-14 VS 17-18

Designed to reach distal surface effectively and makes reach the mesial surface of posterior when it is fix in the mouth in front a patient.

Standard Gracey Curette_Silicone Handle

Autoclavable

2CGR3-4

• Anterior & premolar



2CGR5-6

• Anterior & premolar



2CGR7-8

• Premolar & molar
(facial and lingual surface)



2CGR9-10

• Molar (facial and lingual surface)



2CGR15-16

• Mesial surface of all posterior teeth
• Shank has the same angle with GR13-14 but useful for mesial surface of posterior



2CGR17-18

• Distal surface of all posterior teeth



Gracey Curettes-Standard



Gracey Curette
GR11-12 VS 15-16

It is difficult to place connecting part of GR11-12 on mesial surface of lower posterior parallelly. The angle of GR15-16 makes reach mesial surface of posterior when it is fix in the mouth in front of a patient.

Standard Gracey Curette_Plastic Handle

Autoclavable



BEST
3CGR1-2

• Anterior

Video
Clip



BEST
3CGR11-12

• Mesial surface of
all posterior teeth

Video
Clip



BEST
3CGR13-14

• Distal surface of
all posterior teeth

Video
Clip

Gracey Curettes-Standard



Gracey Curette
GR13-14 VS 17-18

Designed to reach distal surface effectively and makes reach mesial surface of posterior when it is fix in the mouth in front of a patient.

Standard Gracey Curette_Plastic Handle

Autoclavable

3CGR3-4

• Anterior & premolar



3CGR5-6

• Anterior & premolar



3CGR7-8

• Premolar & molar
(facial and lingual surface)



3CGR9-10

• Molar
(facial and lingual surface)




3CGR15-16

• Mesial surface of all posterior
teeth
• Shank has the same angle
with GR13-14 but useful for
mesial surface of posterior



3CGR17-18

• Distal surface of all posterior teeth



Gracey Curettes-Standard



Gracey Curette
GR11-12 VS 15-16

It is difficult to place connecting part of GR11-12 on mesial surface of lower posterior parallelly. The angle of GR15-16 makes reach mesial surface of posterior when it is fix in the mouth in front of a patient.

Standard Gracey Curette_Metal Handle



BEST
CGR1-2
• Anterior



Video Clip



BEST
CGR11-12
• Mesial surface of all posterior teeth



Video Clip




BEST
CGR13-14
• Distal surface of all posterior teeth



Video Clip

Gracey Curettes-Standard



Gracey Curette
GR13-14 VS 17-18

Designed to reach distal surface effectively and makes reach mesial surface of posterior when it is fix in the mouth in front of a patient.

Standard Gracey Curette_Metal Handle

CGR3-4

• Anterior & premolar



CGR5-6

• Anterior & premolar



CGR7-8

• Premolar & molar
(facial and lingual surface)



CGR9-10

• Molar
(facial and lingual surface)



CGR15-16

• Mesial surface of all posterior teeth
• Shank has the same angle with GR13-14 but useful for mesial surface of posterior



CGR17-18

• Distal surface of all posterior teeth



Gracey Curettes-Standard

Standard Gracey Curette_Silicone Handle

Autoclavable

※OSUNG's silicone handle gives no-stress on wrist and provides an excellent grip.



BEST
2CLGR1-2
• Anterior



BEST
2CLGR11-12
• Mesial surface of all posterior teeth



BEST
2CLGR13-14
• Distal surface of all posterior teeth

Gracey Curettes-Standard

Standard Gracey Curette_Silicone Handle

Autoclavable

2CLGR3-4

• Anterior & premolar



2CLGR5-6

• Anterior & premolar



2CLGR7-8

• Premolar & molar (facial and lingual surface)



2CLGR9-10


• Molar (facial and lingual surface)



We applied colour coding for curettes by ISO 13397-2:2005.

Type	Colour Coding	Area
GR 5/6	YELLOW	Anterior / Canine Teeth
GR 7/8	GREEN	Molar and Premolar, Buccal and Oral
GR 11/12	RED	Molar and Premolar, Mesial, Furcations
GR 13/14	BLUE	Molar and Premolar, Distal, Furcations

Rigid Gracey Curettes



The tip is strong and thick compared to gracey curette so it can remove much calculus without extra use of sickle scaler or hoe scaler.

Rigid Gracey Curette_Plastic Handle

Autoclavable

BEST
3CRGR1-2

• Anterior



3CRGR3-4

• Anterior & premolar



3CRGR5-6

• Anterior & premolar



3CRGR7-8

• Premolar & molar
(facial and lingual surface)



3CRGR9-10

• Molar (facial and lingual surface)



BEST
3CRGR11-12

• Mesial surface of all posterior teeth




BEST
3CRGR13-14

• Distal surface of all posterior teeth



Rigid Gracey Curettes



Rigid Gracey Curette_Metal Handle

BEST
CRGR1-2

• Anterior



CRGR3-4

• Anterior & premolar



CRGR5-6

• Anterior & premolar



CRGR7-8

• Premolar & molar
(facial and lingual surface)



CRGR9-10

• Molar
(facial and lingual surface)



BEST
CRGR11-12

• Mesial surface of all posterior teeth




BEST
CRGR13-14

• Distal surface of all posterior teeth



Mini Five Curettes



The terminal shank is 3mm longer than that of standard curette for access into deep periodontal pockets and root surfaces of 5mm or more. The blade length is reduced in half from the standard gracey curette, for a better adaptation in narrow pockets and furcations. The blade is a little thinner than that of the standard gracey curette to enable easy gingival insertion and reduce tissue damage.


Mini Five Curette_Silicone Handle

Autoclavable ※OSUNG's silicone handle gives no-stress on wrist and provides an excellent grip.

BEST

2CMGR1-2

• Anterior



2CMGR3-4

• Anterior & premolar



2CMGR5-6

• Anterior & premolar



2CMGR7-8

• Premolar & molar (facial and lingual surface)



2CMGR9-10

• Molar (facial and lingual surface)



BEST

2CMGR11-12

• Mesial surface of all posterior teeth



BEST

2CMGR13-14

• Distal surface of all posterior teeth



Mini Five Curettes

Mini Five Curette_Metal Handle

BEST

CMGR1-2

• Anterior



CMGR3-4

• Anterior & premolar



CMGR5-6

• Anterior & premolar



CMGR7-8

• Premolar & molar (facial and lingual surface)



CMGR9-10

• Molar (facial and lingual surface)



BEST

CMGR11-12

• Mesial surface of all posterior teeth



BEST

CMGR13-14

• Distal surface of all posterior teeth



Mini Five Rigid Curettes . After Five Rigid Curettes

Mini Five Rigid Curette_Plastic Handle

- The terminal shanks of the Rigid gracey curette are extended by 3 mm.
- The blade is 10% thinner than the Rigid gracey curette to ease gingival insertion and reduce tissue expansion.
- Mini Five Rigid Curettes feature a reduced blade that is to half the length of the After Five Rigid Curette, for better adaptation in narrow pockets and furcations.

NEW
3CMRGR11-12

Mini five rigid curette

- Easy access to use into deep periodontal pockets and root surfaces of 5 mm or more.
- Used for all of the mesial surfaces of posterior teeth



134°C
오토클레이브 사용가능

NEW
3CMRGR13-14

Mini five rigid curette

- Easy access to use into deep periodontal pockets and root surfaces of 5 mm or more.
- Used for all of the distal surfaces of posterior teeth



134°C
오토클레이브 사용가능

After Five Rigid Curette_Plastic Handle

- The terminal shanks of the Rigid gracey curette are elongated by 3 mm.
- Easy access to use into deep periodontal pockets and root surfaces of 5 mm or more.
- The blade is 10% thinner than the Rigid gracey curette to ease gingival insertion and reduce tissue expansion.

NEW
3CARGR11-12

After five rigid curette

- Easy access to use into deep periodontal pockets and root surfaces of 5 mm or more.
- Used for all of the mesial surfaces of posterior teeth



134°C
오토클레이브 사용가능

NEW
3CARGR13-14

After five rigid curette

- Easy access to use into deep periodontal pockets and root surfaces of 5 mm or more.
- Used for all of the distal surfaces of posterior teeth





134°C
오토클레이브 사용가능

Curettes

Detail of Curette


※ Design comparison by type



1. Rigid Gracey Curette

Used to remove heavy calculus due to thick and strong shank.


The terminal shank is thicker and stronger than the Gracey Curette, which is thicker and stronger than the Gracey type Curette and is used to remove heavy calculus without the use of additional sickle scalers or hoe scalers.



2. After Five Rigid Curette

Thick shank and long terminal shank make it easy to remove heavy calculus while improving access to root surfaces.

The terminal shank is 3 mm longer, and the blade width decreased by 10% compared to Rigid Gracey Curettes, which make it easy to access to root surfaces, and it is thicker than Standard After Five, which make it less bounced over the calculus removal.



3. Mini Five Rigid Curette

The long terminal shank and the short blade improve accessibility into periodontal pockets.

The blade is 1/2 shorter than the After Five Rigid Curettes, making it ideal for narrow pockets and furcations. (The tips are rigid rather than the standard Mini Five Curettes.)

Universal Curettes

Blades are sharpened on both sides. Blade curved at 90 degrees to the shank with a rounded toe. Designed so that the working ends can be adapted to all tooth surfaces of all regions of the mouth with one double-ended instrument.

Universal Curette_Silicone Handle

 Autoclavable

※OSUNG's silicone handle gives no-stress on wrist and provides an excellent grip.

BEST

2CUC13-14

- Used for removal of minor calculus of all teeth.
- For both mesial and distal surfaces.



2CU2L-2R

- Used for removal of minor calculus of incisors and premolars.
- For supra- and subgingival.



2CUGF3

- Used for removal of minor calculus of premolars and molars.
- Also for concave tooth surfaces and furcation lesions.



2CUMC13S-14S

- Rigid shank with strong blade is suitable for removal of heavy calculus deposits.
- For removal of supra- and subgingival calculus



2CUSYN15-16

- A combination of universal and finishing curette for the removal of supra- and subgingival calculus.
- For both concave and convex surfaces.



2CUSYN15-16M

- Excellent in tight, deep pocket



Universal Curettes

Universal Curette_Metal Handle

BEST

CUC13-14

- Used for removal of minor calculus of all teeth.
- For both mesial and distal surfaces.



CU2L-2R

- Used for removal of minor calculus of incisors and premolars.
- For supra- and subgingival.



CU4L-4R

- Used for removal of minor calculus of incisors and premolars.
- For supra- and subgingival.



CUGF3

- Used for removal of minor calculus of premolars and molars.
- Also for concave tooth surfaces and furcation lesions.
- Series #3 of Goldman fox



CUMC13S-14S

- Rigid shank with strong blade is suitable for removal of heavy calculus deposits.
- For removal of supra- and subgingival calculus



Universal Curettes · Special Curettes

Universal Curette_Metal Handle

Blades are sharpened on both sides. Blade curved at 90 degrees to the shank with a rounded toe. Designed so that the working ends can be adapted to all tooth surfaces of all regions of the mouth with one double-ended instrument.

CUYG7-8

- Younger-Good 7-8
- Premolar & molar



CUSYN15-16

- A combination of universal and finishing curette for the removal of supra- and subgingival calculus.
- For both concave and convex surfaces.



CUSYN15-16M

- SYN15-16 MINI
- Excellent in tight, deep pocket.



New Combination of Gracey Curette

CGR11-14

- Combine with gracey curettes 11 & 14
- Access to mesial & distal on posterior as one curette.



CGR12-13

- Combine with gracey curettes 12 & 13
- Access to mesial & distal on posterior as one curette.



SUB-ZERO Curette

CUSUB-0

- It is ideal for anterior flap surgery with a long shank that reaches furrow.
- The blade of short hook type removes calculus clearly.
- One sub-zero curette is enough for flap surgery.



Video Clip

Implant Curettes · Chisel Scaler

Implant Curette

Soft titanium is used as a material for scratch-free and contamination-free on the implant. And the implant curettes are more workable than plastic curettes. Also the curettes can be used for a second surgery. As the curettes have the same blade shapes and angles as those of standard curettes, they give a more comfortable feeling when we use them.

BEST

ICGR1-2

- Similar to Gracey 1-2
- For anterior

ICGR5-6

- Similar to Gracey 5-6
- For anterior

ICGR7-8

- Similar to Gracey 7-8
- For posterior

BEST


ICGR11-12

- Similar to Gracey 11-12
- For posterior

BEST

ICGR13-14

- Similar to Gracey 13-14
- For posterior



Chisel Scaler

CSZ

- Metal Handle/Single End
- Removes calculus on the mandibular anterior.
- It is push stroke type not like hoe scaler



Hoe Scalers

Used for removal of heavy supramarginal calculus.

Hoe Scaler

For anterior buccal and lingual surfaces.

HSA12-13
• Anterior Hoe Scaler



For the buccal and lingual surfaces of all teeth. Also can be used in furcation areas.

HSL34-35
• Lateral Hoe Scaler



For the mesial and distal surfaces of molar. Can be used in furcation areas.

HSP56-57
• Posterior Hoe Scaler



It is used for buccal and lingual of posterior.

HSO8-9
Hoe Scaler, O8-9
• Orban 8-9
• 폭1.8mm/폭1.8mm



File Scalers

File Scaler

Used for crushing large calculus deposits so that the deposit can be more easily removed by a curette. Can also be used for smoothing the margins of amalgam restorations.

FSH3-7
• Hirschfeld 3-7
• Buccal/Lingual



FSH5-11
• Hirschfeld 5-11
• Mesial/Distal



Periodontal File Scaler

For interproximal. To crush and remove heavy deposits from subgingival and supragingival interproximal areas.

File on both sides for using proximal. Push & pull.

PDS1-2S
• Mesial/Distal



Has file on one side. No harm on gingival during using at buccal & lingual.

PDS3-4S
• Buccal/Lingual



Used for crown lengthen procedure, implant surgery, removing torus and oarplasty

PDS9-10S
• Curved File
• Buccal/Lingual



Sharpening Stone · Perio Scaling Kit

Sharpening Stone

• Used for sharpening hand instruments.

SST-C3

- Ceramic Sharpening Stone #3C (Medium Grit)
- Brown
- 80 x 33 x 6.3H (mm)

 Autoclavable



Perio Scaling kit

 Autoclavable

3LSK01

- Perio Scaling Kit Part 1
- Diagnostic & Supragingival Scaling
 - Ex-Probe 1ea, Sickle Scaler 2ea

3LSK02

- Perio Scaling Kit Part 2
- Root Planning & Subgingival Curettage
 - Anterior Curette 1ea, Molar Curette 2ea



- 3XP23-WHO**
- Exploring calculus & measuring periodontal pocket

- 3LSU15-33**
- Used for removing of supragingival calculus in anterior.

- 3LSJAC34-35**
- Used for removing of supragingival calculus in posterior

- 3CGR1-2**
- Used on the anterior of teeth.

- 3CGR11-12**
- Used on the mesial portions of posterior teeth.

- 3CGR13-14**
- Used on the distal portions of posterior teeth.



Periodontal Treatment

Non-surgical treatment to maintain a healthy periodontal condition, to restore periodontally diseased tissue to a healthy state to prevent progression of periodontal disease.

Rigid curette

Rigid gracey curette has a thicker and stronger terminal shank than that of the standard gracey curette.

It is built for removal of heavier levels of calculus.


Light weight plastic handle design provides easy handling and reduces hand and wrist fatigue.

Plastic handle with embossed dot pattern gives a more positive grip and its rolling stopper on the handle allows minimizing tip damage from rolling or sliding when it is placed on the table.



Arrangement

01. Probe	BPWHO	P.026
02. Explorer	EXD11-12	P.014
03. Torque Wrench	USETW	P.036, 334
04. Ultrasonic Scaler Tip	USEA	P.036
05. Sickle Scaler	LSH5-33	P.038
06. Sickle Scaler	LSJAC31-32	P.038
07,08,09. Gracey Curette	CGR1-2	P.046
	CGR11-12	P.046
	CGR13-14	P.046



Process

BPWHO
EXD11-12



01. Measuring periodontal pocket depth
02. Detecting subgingival calculus

USFTW
USEA
LSH5-33



03. 04. Removal of supragingival calculus with ultrasonic scaler
05. Removal of supragingival calculus(anterior)

LSJAC31-32
CGR1-2



06. Removal of supragingival calculus(posterior)
07. Root planing(anterior)

CGR11-12
CGR13-14



08. Root planing(mesial surface of posterior teeth)
09. Root planing(buccal surface of posterior teeth)

Practice

01. Measuring periodontal pocket depth

Used
Periodontal probes are used to measure the depth and determine the configuration of a periodontal pocket, gingival bleeding response to the periodontal probing, gingival recession and clinical attachment loss. It can also be used for determining the extent of furcation involvement on multi rooted teeth and measuring the pathologic lesions and width of the attached gingiva.

Character
It has a ball end of diameter 0.5mm and a first colored band at 3.5-5.5mm. Blunt ball end makes the patient comfortable when inserting the periodontal probe into the gingival sulcus. Specially designed for detecting subgingival calculus and overhanging margin.

02. Detecting subgingival calculus

Used
Used for detecting the amount & distribution of subgingival calculus & plaque and examining the condition of tooth surface after receiving treatment of scaling and root planing. Detect the anatomic configuration of root and root anomalies.

Character
Angled like Gracey 11/12 Curette for improved calculus detection. Extra-long complex shank allows deeper insertion and better access into the periodontal pocket. Available for anterior or posterior application.

Probe_BPWHO

How to use

1. Correctly adapt the periodontal probe using a proper pen grasp.
2. While probing, the tip of the probe is kept vertically parallel to the long axis of the tooth and placed gently on the gingival margin until the junctional epithelium is contacted. The minimal force of around 20-25g should be used.
3. Proceed with walking stroke. The side of the probe tip should be kept in contact with the tooth surface.



The ball ended tip of 0.5mm in diameter minimizes patient discomfort due to probing.



Probing(posterior tooth): The side of the probe tip(1-2mm) is applied to the distofacial line the angle of the buccal/lingual surface.**Probing(anterior tooth):** The side of the probe tip is applied to the superjacent free gingival margin at the central region of the labial/lingual surface.



Adapt the tip to the tooth surface at the angle of 0 to 15 degrees vertically to the tooth and gently apply to a tooth until the junctional epithelium is contacted. Tilt the probe following the shape of Col.

Explorer _ EXD11-12

How to use

1. Use a modified pen grasp with finger rest on an adjacent tooth surface which makes it possible to provide stability and control.
2. Keep the shank parallel to the long axis of the tooth.
3. Insert a tip with a light pressure keeping in contact with the tooth surface.
4. Exploring with walking stroke. The side of the tip should be kept in contact with the tooth surface.



When choosing the correct working end, place the terminal shank parallel to the long axis of the tooth surface.



Apply 1-2mm of a tip(back) to the tooth surface at the superjacent free gingival margin to the tooth surface.



Adapt the tip to the tooth surface into a vertical position, gently insert the tip until the junctional epithelium is contacted and stroke in a vertical direction.

Wrong Position



The incorrect working end has been selected if the terminal shank is not parallel to the long axis of the tooth and it curves around the tooth surface when placing the point to the lingual surface from the buccal surface.

Wrong Position



If the point is directed toward the tooth surface, the wrong working angle will be set.



Proceed from the junctional epithelium to the gingival margin with overlapping short stroke.

Practice

03, 04. Removal of supragingival calculus with ultrasonic scaler

Used

Used for tightening the ultrasonic scaler tip and handpiece.

Character

Made of high strength steel for hexa head a connecting shaft of a scaler tip to ensure improved validity of the instrument.

Used

Ultrasonic scaler tips are used for removal of calculus, plaque and temporary sealing material rapidly from tooth surface during dental prophylaxis.

Character

Tips are designed for EMS scalers and allow various functions to be performed.

Torque Wrench _ USFTW

How to use

1. Insert the ultrasonic scaler tip into the wrench, screw it into the handpiece by turning clockwise to tighten.
2. For loosening, turn the torque wrench counter-clockwise direction.
3. After use, be sure to clean and sterilize the wrench completely with the scaler tip is attached.



Ultrasonic Scaler Tip _ USEA

How to use

1. The side of the tip should be applied 15 degree angle to the long axis of the tooth.
2. It is recommended that the tip be adapted to stroke parallel to the tooth surface in a sweeping-like motion.
3. After 50% of the tip is worn away, it should be discarded as the tip's efficiency is lost.



Tighten the ultrasonic scaler tip and handpiece.



Apply the side of the tip to the tooth surface of 15 degrees angle.



If you want to remove the calculus accumulated on the adjacent tooth surface, place the side of the tip to the tooth surface.

Sickle Scaler _ LSH5-33

How to use

1. Use a modified pen grasp with finger rest on an adjacent tooth surface which makes it possible to provide stability and control.
2. Maintain 45 to 90 degrees for working stroke.
3. Overlapping motions with tip 1/3 with a vertical or oblique stroke.



Tilt the shank slightly toward the tooth surface to establish correct angulation.



Adapt the tip 1/3 of the cutting edge to the center of the cervical line, directing the point toward the mesial surfaces.



Tilt the facial surface of the blade toward the tooth to achieve an approximate 70-80 angle between the tooth and blade. Apply lateral pressure against the tooth and pull the scaler firmly upward and diagonal with overlapping strokes.

05. Removal of supragingival calculus

Used

Designed for removal of moderate to a heavy accumulation of supragingival calculus on anterior teeth and subgingival calculus located just below free gingiva.

Character

Double-ended straight shank for use anyway and anywhere. Two cutting edges on a straight triangular-shaped blade. The internal angle of the blade is 70-80°.

06. Removal of supragingival calculus

Used

Designed for removal of moderate to a heavy accumulation of supragingival calculus on posterior teeth and subgingival calculus located just below free gingiva.

Character

The blade of a curette is correctly adapted when the terminal shank is parallel to the long axis of tooth surface when the blade tip directed toward the two adjacent teeth surface. Two cutting edges on a paired working end with a triangular cross-section and the internal angle of the curved blade is 70-80°.

Sickle Scaler _ LSJAC31-32

How to use

1. Use a modified pen grasp with finger rest on an adjacent tooth surface which makes it possible to provide stability and control.
2. Adapt the tip 1/3 to the distobuccal surface line angle. Oblique stroke on the buccal-lingual surface and vertical stroke on the proximal surface.
3. Use a short pull stroke.



The correct working end for scaling is evident when the terminal shank is parallel with the long axis of the tooth surface.



Adapt the tip 1/3 of the cutting edge to the distobuccal surface.



Wrong Position

It is a wrong making end if the terminal shank is not parallel to the long axis of the tooth when the pin is positioned buccolingual.

07. Root planing

Used

Curettes are fine instruments used for subgingival scaling, root planing and removal of soft tissue lining the pocket. It is used on anterior teeth surface.

Character

It can be adapted and provide good access to a deep pocket with minimal soft tissue trauma. The two cutting edges meet together and make a round toe. The cross section of a toe is a semicircle.

Gracey Curette _ CGR1-2

How to use

1. Use a modified pen grasp with finger rest on an adjacent tooth surface wherever possible to provide stability and control.
2. Position the face of the blade toward the tooth surface, and the angulation between the tooth should be as close to zero as possible. Adapt the tip 1/3 of the lower cutting edge against the tooth surface.
3. Use a short pull stroke.



The blade of a curette is correctly adapted when the terminal shank is parallel to the long axis of tooth surface when the blade tip directed toward the two adjacent teeth surface.



Position the face of the blade toward the tooth surface, and the angulation between the tooth should be as close to zero as possible. Adapt the tip 1/3 of the lower cutting edge against the tooth surface.



Proceed from the junctional epithelium to the gingival margin with overlapping short pull stroke.

Practice

08.09. Root planing

Used

Curettes are fine instruments used for subgingival scaling, root planing and removal of soft tissue lining the pocket. It is used on the mesial surface of posterior teeth.

Character

It can be adapted and provide good access to deep pocket, with minimal soft tissue trauma. The blade has a round toe and two cutting edges for scaling, which makes it an efficient design for a better adaption to the root surface unlike the straight design and pointed end, which can cause tissue laceration and trauma.

Used

Curettes are fine instruments used for subgingival scaling, root planing and removal of soft tissue lining the pocket. It is used on the mesial surface of posterior teeth.

Character

It can be adapted and provide good access to deep pocket, with minimal soft tissue trauma. The blade has a round toe and two cutting edges for scaling, which makes it an efficient design for a better adaption to the root surface unlike the straight design and pointed end, which can cause tissue laceration and trauma.

Gracey Curette _ CGR11-12

How to use

1. Select a suitable blade that can be used for a mesial application.
2. Correctly adapt the tip 1/3 of the working end to the tooth surface.
3. Gently insert the tip until the junctional epithelium is contacted with the angle of 0 degree maintained to the tooth surface.
4. Apply overlapping pull stroke towards root canal side with working angulation of between 70 and 80 degrees.



The blade of a curette is correctly adapted when the terminal shank is parallel to the long axis of tooth surface when the blade tip directed toward the two adjacent teeth surface.



Position the face of the blade toward the tooth surface, and the angulation between the tooth should be as close to zero as possible. Adapt the tip 1/3 of the lower cutting edge to the tooth surface.



Keep the handle slightly away from the tooth surface so there is an angle of 70 to 80 degrees for working stroke.



Wrong Position

It is a wrong making end if the terminal shank is not parallel to the long axis of the tooth when the pin is positioned buccolingual

Gracey Curette _ CGR13-14

How to use

1. Select a suitable blade that can be used for a buccal application.
2. Correctly adapt the tip 1/3 of the working end to the tooth surface.
3. Gently insert the tip until the junctional epithelium is contacted with the angle of 0 degree maintained to the tooth surface.
4. Apply overlapping pull stroke towards root canal side with working angulation of between 70 and 80 degrees.



The blade of a curette is correctly adapted when the terminal shank is parallel to the long axis of the tooth surface when the blade tip directed toward the two adjacent teeth surface.



Position the face of the blade toward the tooth surface, and the angulation between the tooth should be as close to zero as possible. Adapt the tip 1/3 of the lower cutting edge against the tooth surface.



Keep the handle slightly away from tooth surface so there is an angle of 70 to 80 degrees for working stroke.



Wrong Position

It is a wrong making end if the terminal shank is not parallel to the long axis of the tooth when the pin is positioned buccolingual

Products for Dentistry

OSUNG Catalogue 2022/2023

Surgery


Products for Dentistry

OSUNG Catalogue 2022°2023



SURGERY

Preparation for Surgery	Surgical Suction Tip	076
	Extension Hose	076
	Extension Hose Adapter	076
	Frazier Suction Tip	076
	Suction Tip	077
	Surgical Drape	078
	Wrapping Cloth	078
	Towel Clamp	078
	Anesthesia Syringe	079
	Scalpel Handle	080
Composite Scalpel Handle	081	
Oral Surgery	Periosteal Elevator	082
	NiTi-Periotome	085
	Periotome	086
	Root Picker	087
	Luxating Elevator	088
	Elevator	090
	Extraction Forceps(Adult)	092
	Extraction Forceps(Pedo)	096
	Surgical Curette	098
	Bone Rongeur	102
	Nipper	102
	Bone File	103
	Mallet	103
	Hemostat	104
	Needle Holder	105
Anatomic Dressing Forceps	106	
Tissue Plier	107	
Scissors	107	
Periodontal Surgery	Periodontal Knife	110
	Periodontal Chisel	110
	Periodontal Surgical Curette	111
Maxillofacial Surgery	Periosteal Elevators for Maxillofacial Surgery	112
	V-Notch Periosteal Elevator	112
	Channel Retractor	113
	Spatula Periosteal Chisel	113
	Retractor	114
	Tunneling Instrument	116
	Palatal Wedge	118
Simple Extraction Excision of Torus	Manual	119
	Manual	123




Surgical Suction Tips

Practice

Surgical Suction Tip Guide

Recommend to use extension hose to solve the problem of a short length of a built-in suction system of a unit chair.



Frazier Suction Tip

- As it has a function of controlling a suction force by closing & opening a hole, and it prevents soft tissue damage. (Especially mucous area under the tongue)
- For strong suction power, close the hole.

SNF20

- Cleaning Stylet (Suction cleaner) is not included.



BEST

SNF25

- Cleaning Stylet (Suction cleaner) is not included.



SNF30

- Cleaning Stylet (Suction cleaner) is not included.



SNKHS

- Extension Hose**
- Connected to the unit chair by 'Extension Hose Adapter'
 - Silicone made
 - Length : 1.5m
 - Autoclavable



SNKHS-1

- Extension Hose**
- Connected to the unit chair by 'Extension Hose Adapter'
 - Silicone made
 - Length : 1.5m
 - Autoclavable



SNKHA

Extension Hose Adapter



SNKGS

- Cleaning Stylet



Suction Tips


It gives a nice grip-feeling, and it is a highly available item for implant surgery.
*Recommended to use together with an extension hose.

Titanium Suction Tip

- Titanium tip is combined for implant surgery.

SN4TI

Titanium Suction Tip



BEST

SN3TI

Titanium Suction Tip



Stainless Steel Suction Tip

- Made of stainless steel
- Superior durability

SN4SUS

Stainless Steel Suction Tip



BEST

SN3SUS

Stainless Steel Suction Tip




SN3SUSL

- 15mm longer than SN3SUS



Surgical Drape · Wrapping Clothes · Towel Clamp

Surgical Drape



- A magnet is inserted into drape to stick hinge-typed dental instruments such as scissors, needle holders, and so on.
- Magnet can be removed.



WDMA

- Size : 900 x 900mm
- Hole diameter : 90mm

Wrapping Cloth



- No discolor during washing. Tenacious fabric. Made in Korea.

WR5050

- Wrapping cloth for sterilizing dental instruments.
- Size : 500 x 500(mm)

WR7575

- Wrapping cloth for sterilizing dental instruments.
- Size : 750 x 750(mm)




Video Clip

Towel Clamp

CPTC135

- Length 135mm(±5mm)



Anesthesia Syringes

Aspirating Syringe



Video Clip


- Very simple to use, and is made of high quality materials.
- Our unique hook design guarantees perfect aspiration.



BEST

SAA1

- Type A
- 1.8cc
- Hook Shape



- The hook shaped harpoon provides an excellent clamping force with the lidocaine ampoule, making it well secured



SAB1

- Type B
- 1.8cc
- Hook Shape



SAF1

- 1.8cc
- Arrow Shape

Cartridge Syringe



SAC1

- 1.8cc



Video Clip

- The end of the rod is wide enough to push a rubber plunger of a cartridge stably.


Scalpel Handles

Scalpel Handle

BEST

SHF

- With ruler for measuring length.



SHS

- Straight



SHC

- Curved
- For Posterior & Palatal Areas




- The rounded handle helps to make a curved incision.
- Useful for cutting the deepest posterior & palatal.


SH2S

- Straight

134°C
 오토클레이브 사용가능



- The silicon handle makes a comfortable feeling of grip (tactile sensitivity.)



SHTL

- Users can freely change the angle of the blade for one's needs.




Composite Scalpel Handle

Micro Surgical Scalpel Handle

NEW

SHMS

- Micro Blade Handle
- The microsurgical scalpel handle is a chucking type used by inserting various microsurgical blades (3mm in width).



Composite Scalpel Handle

SHCS


- Blade replaceable

Based on # 15 blade fastened



SHCS2

- Blade replaceable
- Blade length adjustable




It is designed to improve safety and convenience from the old way as inserting the blade into the lidocaine ampoule.

- Safety - Protect your hands from the blade
- Convenience - Easy blade fastening and save time
- Hygienic - Autoclave sterilizable
- Semi-permanent use is possible only by replacing the blade

When the adjacent tooth surface of the anterior or premolars is laminated with composite resin, no matter how well the metal matrix band or Mylar strip is sealed, the resin may overflow, or the overhanging margin may form on the gingival margin. At the end of the photopolymerization process, after finishing and polishing with a high speed bur, mostly adjust the lowermost part using the blade of #12 scalpel comes out of the composite scalpel handle and form the gingival margin and embrasure without touching the gingival as much as possible.

Practice



Periosteal Elevators

To separate tissue from tooth or bone. To hold tissue away from surgical site.

Periosteal Elevator

BEST

EP9

- Used the most for surgery.
- To elevate mucous periosteal with a wide tip.
- To elevate interdental papilla with a sharp tip.



EP9H

- Allen
- The hole in the spoon blade helps to hold tissue during suture.



EP9S

- To elevate mucous periosteal with a wide tip.
- To elevate interdental papilla with a sharp tip.



EP14

- Goldman-Fox



Periosteal Elevators

Periosteal Elevator

BEST

EP24G



EPKN1

- Small periosteal elevator with one rounded end and the other pointed end for delicate tissue retraction.



BEST

EPBUSER

- This is ideal for delicate tissue retraction. The spear shaped end is used for initiating the flap by retracting the interdental papilla. The paddle end is used for continuing flap retraction apically.



EP20

- To elevate mucous periosteal after cutting gingiva.
- One tip is straight shaped and the other one is contra-angle type with a thin and sharp edge



Periosteal Elevators

Used to lift the mucoperiosteal flaps after the incision of gingival tissues.

Periosteal Elevator



BEST
EPPR3

- Prichard
- To retract flap



EP15

- Freer



EP16

- Freer



EP23

- Selden
- To fix tissue by retracting during a flap surgery.



NiTi-Periotome

Used for teeth removal without damaging adjacent bone and peripheral tissues by cutting periodontal ligaments using thin blades on both sides. The tip of the NiTi material curves along adjacent bone tissue to minimize the discrepancy between the adjacent bone tissue and the surrounding tissue.

NiTi-Periotome_Straight

NEW
PRTMTS

- NiTi-Periotome**
- Special Alloy: Nickel + titanium alloy
 - Width 2.1 mm
 - Joint 3mm




- Nickel-titanium alloys provide excellent corrosion resistance, and the elasticity allows easy access to the roots.



- Sold separately P.81
- SHMS: the microsurgical scalpel handle (combined use) is inserted to use.

Practice

The tip of the NiTi material curves along adjacent bone tissue to minimize the discrepancy between the adjacent bone tissue and the surrounding tissue.



Periotomes

Used for atraumatic extraction of teeth.

Periotome

PRRS3

- Straight
- For anterior



BEST

PR2-2R

- Curved
- For posterior



PRR256

- Curved
- Serrated blade
- For anterior & posterior



PRR258

- Curved
- Serrated blade
- For anterior & posterior



PRM1

- Straight
- For malleting (Single End)




Root Pickers

Root Picker

- To remove root tips when the root is fractured during extracting.

Autoclavable



- ERP1 Root Picker, RP01 • Heidbrink
- 3ERP1
- 4ERP1



- ERP2 Root Picker, RP02 • Heidbrink
- 3ERP2
- 4ERP2



- ERP3 Root Picker, RP03 • Heidbrink
- 3ERP3
- 4ERP3



ERHB13-14

- Root Picker, HB13-14



Luxating Elevators




Plastic

Aluminium



Luxating Elevator

• Designed to easily extract the root of the tooth by cutting the periodontal ligament. The merit is less tissue damage and preservation of the alveolar bone. • Autoclavable



Regular

3ELLR301
4ELLR301

• For anterior

3ELLR303
4ELLR303

• For posterior

3ELLR302
4ELLR302

• For posterior

3ELLR501
4ELLR501

• For anterior

3ELLR503
4ELLR503

• For posterior

3ELLR502
4ELLR502

• For posterior

Ultra-Sharp

3ELL301
4ELL301

Ultra-Sharp
• Thinner and sharper blade than the regular type.

3ELL303
4ELL303

Ultra-Sharp
• Thinner and sharper blade than the regular type.

3ELL302
4ELL302

Ultra-Sharp
• Thinner and sharper blade than the regular type.

3ELL501
4ELL501

Ultra-Sharp
• Thinner and sharper blade than the regular type.


3ELL503
4ELL503

Ultra-Sharp
• Thinner and sharper blade than the regular type.

3ELL502
4ELL502

Ultra-Sharp
• Thinner and sharper blade than the regular type.

Luxating Elevators




Luxating Elevator Regular or Ultra-Sharp

Luxating Elevator Dual-Edge

The concept of the Luxating elevator is to luxate the dental root easily by incising the gingival ligament unlike an elevator that needs more power to elevate dental root. Because of the thinner and sharper blade than an elevator, it can get damaged easily when it is overused for operation but the luxating elevators have the merit of less damage of tissue and preservation of the alveolar bone. It is more effective using with other extracting instruments like forceps if needed. It is redesigned of the dentist's ideas based on their clinical experience. The grip is also specially designed to be controlled by minimum force.

Dual-Edge

• Tiny narrow blades can be used for deeply fractured teeth or dental caries.



Regular

3ELLD501
4ELLD501

• For anterior

3ELLD503
4ELLD503

• For posterior

3ELLD502
4ELLD502

• For posterior

3ELLD50K
4ELLD50K

• For posterior distal surface

Elevators which need strong and steady force to operate specified in ISO standard to test the tightness of the connector. However, no test equipment is available on the market to test these instruments. You must build your own torsion tester in accordance with what is presented in the ISO standard document. Prior to building the tester, finite-element analysis is required. Based on this analysis, we determine the materials to be used and the force to be applied to the specimen. In despite of its simple appearance, the mechanism has many complex implementations.



▲ Figure. Deformation analysis result of Torsion tester

Compound Curved



Regular

3ELLR30K
4ELLR30K

• For posterior distal surface

3ELLR50K
4ELLR50K

• For posterior distal surface

Ultra-Sharp

3ELL30K
4ELL30K

Ultra-Sharp
• Thinner and sharper blade than the regular type.

3ELL50K
4ELL50K



Ultra-Sharp
• Thinner and sharper blade than the regular type.

Elevators

Elevator

• Used for loosening the tooth from the periodontal ligament and making the extraction easy.

• Autoclavable



EL34 Elevator, E34
• Upper & Lower anterior roots

3EL34

4EL34

EL34S Elevator, E34S
• Upper & Lower anterior roots

3EL34S

4EL34S



EL31F Elevator, E31F
• Upper & Lower cuspids

3EL31F

4EL31F

EL32F Elevator, E32F
• Upper & Lower cuspids

3EL32F

4EL32F



EL41 Elevator, E41
• Upper & Lower anterior teeth and roots

3EL41

4EL41

EL81 Elevator, E81
• Upper anterior apices and small teeth

3EL81

4EL81

Elevators

Elevator



EL301 Elevator, E301
• Deeply seated roots

3EL301

4EL301

EL301A Elevator, E301A
• Deeply seated roots

3EL301A

4EL301A



EL304W Elevator, E304W
• Deeply seated roots


3EL304W

4EL304W

EL52 Elevator, E52
• Curved Backward
• 3rd molars

3EL52

4EL52



EL4 Elevator, E4
• Broken or deeply seated roots

3EL4

4EL4

EL5 Elevator, E5
• Broken or deeply seated roots
• Deeply seated roots

3EL5

4EL5



EL73 Elevator, E73
• Broken or deeply seated roots

3EL73

4EL73

EL74 Elevator, E74
• Broken or deeply seated roots

3EL74

4EL74

Extraction Forceps

Extraction Forceps_Adult

For upper jaw teeth

FXX2

• Upper jaw's left and right of anterior

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FXX7

• Upper jaw's left and right of premolar

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FXX17

• Upper jaw's right molar

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FXX18

• Upper jaw's left molar

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



Extraction Forceps

Extraction Forceps_Adult

For upper jaw teeth

BEST

FXX67A

• Upper jaw's left and right of 3rd molar

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FXX13

• Upper jaw's left and right of anterior and premolar

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FXX22

• Lower jaw's left and right of molar

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FXX79

• Lower jaw's left and right of 3rd molar

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



Extraction Forceps

Extraction Forceps_Adult

FX1

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FX150

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FX10S

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FX53R

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FX53L

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



Extraction Forceps

Extraction Forceps_Adult

FX151

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FX17

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FX222

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FX300

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



FX301

8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8
8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8



Pedo Extraction Forceps

It is specially designed for children patients. When a user grasps the forceps, the forceps are mostly covered by hand in visual. So children are not afraid of the instrument during extracting.



Asian Type (Pedo)

FXX29C

- Upper jaw's left and right of primary anterior




FXX7C

- Upper jaw's left and right of primary anterior



FXX51C

- Upper jaw's left and right of primary posterior




FXX33C

- Upper jaw's left and right of primary anterior



FXX13C

- Upper jaw's left and right of primary posterior




Pedo Extraction Forceps

It is specially designed for children patients. When a user grasps the forceps, the forceps are mostly covered by hand in visual. So children are not afraid of the instrument during extracting.

American Type (Pedo)

FX150S

- Upper jaw's primary anterior and posterior
- For root



FX151S

- Lower jaw's primary anterior and posterior
- For root



FX101

- Upper and lower jaw's left and right of permolar
- Upper and lower jaw's left and right of primary posterior



Surgical Curettes(Buccal-Lingual)

Surgical Curette_Plastic Handle

- Used for curettage and oral cyst removal, removal of alveolar bone necrosis.
- Same use as surgical curette but more precise removal with serrated tip.
- Autoclavable



3URCL85C

- Surgical Curette, CL85C
- Lucas • Curved
- Serrated blade



3URCL86C

- Surgical Curette, CL86C
- Lucas • Curved
- Serrated blade




3URCL87C

- Surgical Curette, CL87C
- Lucas • Curved
- Serrated blade



Practice

Features : Excellent soft tissue removal ability with serrated tip.
How to use : Remove remaining soft tissue with a light force.
Caution : The removal ability is excellent and there is a risk of bone loss when used with excessive force.
Clinical application : Removal of root apron granulomas and root adenoma cysts. When extracted due to periodontal disease, it effectively removes the remaining soft tissue to help later bone regeneration



Surgical Curettes(Mesial-Distal)

Surgical Curette_Metal Handle

URCL84MD

- Mesial-distal Curette
- Curved




URCL85MD

- Mesial-distal Curette
- Curved



Practice

The use of conventional buccolingual curette with mesiodistal curette makes it very easy to remove granulation tissue after extraction



Surgical Curettes

Surgical Curette_Metal Handle

URCM2-4

Surgical Curette, CM2-4
• Miller • Straight



URCM9

Surgical Curette, CM9
• Miller • Straight



URCM10

Surgical Curette, CM10
• Miller • Curved



URCM11

Surgical Curette, CM11
• Miller • Curved



URCL84

Surgical Curette, CL84
• Lucas • Curved



Surgical Curettes

Surgical Curette_Metal Handle

URCL85

Surgical Curette, CL85
• Lucas • Curved



URCL86

Surgical Curette, CL86
• Lucas • Curved



URCL87

Surgical Curette, CL87
• Lucas • Curved



URCL88

Surgical Curette, CL88
• Lucas • Curved



Bone Rongeurs · Nippers

Bone Rongeur

- To remove granulation tissue, fractured bone.
- It has sharp blades on both sides of the beak.


RNGF140

- Length : 142mm(±5mm)
- Small-sized



RNG165

- Length : 165mm(±5mm)
- Normal-sized



RNG178

- Length : 178mm(±5mm)
- Double action type for increasing the grip force.
- Less sliding and hand fatigue



Nipper

- To cut soft tissue and bone as well as tissue residues.

TNC100

- Length : 100mm(±5mm)



Bone Files · Mallets

To smoothen the surface of a bone. There are a string of blades on the file.

Bone File

BF22

- Miller
- Straight-cut blades
- Pull stroke



BF1X

- Miller
- Cross-cut blades
- Pull & Push stroke



BF45

- Miller
- Pull & Push stroke
- For the wide surgical area



Mallet

ML25



- Autoclavable
- 196g, Ø25, 180mm
- Replacement disk is optional
- Soft use with less shock than ML20



ML25D

(2pcs)



- Autoclavable
- set(2pcs)




ML20

- Made of stainless steel, 221g, Ø19.8, 165mm
- Stainless steel material gives full power even if it is small.



ML29

- 467g, Ø29mm, 185mm
- Used for bone crusher



Hemostats

Multiple use for taking or removing something or pressing blood vessel.

Hemostats_Mosquito

- Useful for taking small fibrous tissue.

BEST

HTM130


- Straight
- Length : 130mm(±5mm)



BEST


HTM130C

- Curved
- Length : 130mm(±5mm)



NEW

HTM120C



Hemostats_Kelly


HTK145

- Straight
- Length : 145mm(±5mm)



HTK145C

- Curved
- Length : 145mm(±5mm)



Needle Holders

Used for taking & guiding the needle.

Needle Holders

BEST

NH160

- Straight
- Length : 160mm(±5mm)



NH160TC

- Length : 160mm(±5mm)
- Tungsten Carbide on beak



BEST


NHC150TC

- Straight
- Length : 150mm(±5mm)
- Tungsten Carbide on beak



OPLF01

- Ligature Forcep



Needle Holders · Anatomic Dressing Forceps

Needle Holder

NHC130TC

- Needle Holder
- Length : 130mm(±5mm)



NH150TC-P

- Needle Holder
- Length : 150mm(±5mm)



NH160TC-P

- Needle Holder
- Length : 160mm(±5mm)



Anatomic Dressing Forceps

- Used for holding soft tissue.

BEST


PT41

- Tissue Plier, Tissue Forceps
- 길이 128mm (±5mm)



PTB19

- Tissue Plier, Bakey
- 길이 154mm (±5mm)




Scissors · Tissue Pliers

Tissue Plier

PTG1

- Length : 178mm (±5mm)




PT52C

- Length : 150mm (±5mm)



PTS22C

- Length : 150mm (±5mm)
- Hole : 2.2mm



PT42

- Length : 122mm (±5mm)



Metzemaum Scissors

SCMB130

- Metzemaum Scissors
- Length : 130mm(±5mm)



SCMB145

- Metzemaum Scissors
- Length : 150mm(±5mm)



Scissors

Scissors

BEST

SCD170

Dean Scissors

- Length : 170mm(±5mm)
- Serrated blade on one side
- It helps a suture not to slide.



BEST

SCC105

Crown Scissors

- Length : 105mm(±5mm)
- Cut or trim crown or gold metal



SCT115

Tissue Scissors

- Straight
- Length : 115mm(±5mm)
- Cut tissue



BEST

SCTC115

Tissue Scissors

- Curved
- Length : 114mm(±5mm)
- Cut tissue



Scissors


Scissors

SCLC115

Scissors, LaGrange

- Compound Curved
- Length : 115mm(±5mm)

Easy access to the mouth with curved twice.



SCGS130

Scissors, Goldman-Fox

- Straight
- Length : 130mm(±5mm)
- Remove granulation tissue from the interdental papilla and surgical flaps.



SCGC130

Scissors, Goldman-Fox

- Curved
- Length : 130mm(±5mm)
- Remove granulation tissue from the interdental papilla and surgical flaps.



SCLSS115

Scissors, Littauer, Suture

- Straight
- Length : 115mm(±5mm)
- For suture
- Useful for edema



SCLSA115

Scissors, Littauer, Suture

- Angled
- Length : 118mm(±5mm)
- For suture
- Useful for edema



Periodontal Knives · Periodontal Chisels

Periodontal Knives

BEST

KNK15-16

- Kirkland 15-16
- Used for initial bevel incision for gingivectomy or gingivoplasty procedures.
- Easy to access to the distal surface of the posterior.



KNO1-2

- Orban 1/2
- The blade & shank are properly angled for posterior use.
- Useful for making tunnels of the recipient site.



Periodontal Chisel

- Used for removing and shaping bone

BEST

CHCO1

- Ochsensbein & Fedi Curved
- Used for reshaping a bone.
- Can also be used to reflect flaps or remove secondary palatal flaps.
- Useful to get a small amount of autogenous bone during implant surgery.
- To remove soft and hard tissue strongly stick to the palate.
- To remove thin bone after supporting the backside of the knife against the root surface.



BEST

CHCO2

- Ochsensbein & Fedi Curved
- Has a knife of opposite direction against CHCO1



BEST

CHC36-37

- Ochsensbein & Fedi Curved
- Back-Action
- Pull-stroke, ideal for removing bone adjacent to tooth without injury
- Useful for molar's distal



Periodontal Chisels · Periodontal Surgical Curette

Periodontal Chisel

CHS13K-13KL

- Kirkland
- Curved
- For removing and reshaping bone 3,3mm



CHC13K-TG

- Curved
- Suitable for root forming



CHCP3-4

- Straight



CHBC1

- Used for splitting bone.



Periodontal Surgical Curette

Larger and heavier curette for the removal of granulation tissue and tenacious subgingival deposits.

BEST

URPR1-2

- Periodontal Surgical Curette, Prichard 1-2



Periosteal Elevators for Maxillofacial Surgery

For oral & maxillofacial surgery

Kang's Elevator

Designed by NARA KANG, DDS

MXP3S

• Straight



3.5mm


MXP3SC

• Slightly Curved



MXP3C

• Curved



MXP6S

• Straight



6.0mm


MXP6SC

• Slightly Curved



MXP6C

• Curved



MXP9S

• Straight



9.0mm


MXP9SC

• Slightly Curved




MXP9C

• Curved



V-NOTCH Periosteal Elevator



MXP75S

V-NOTCH Periosteal Elevator



Channel Retractors · Spatula Periosteal Chisels

MXS1

Kang's Elevator Kit
• Size 210 x 271 x 30(H) (mm)



- 1 MXP3S
- 2 MXP3SC
- 3 MXP3C
- 4 MXP6S
- 5 MXP6SC
- 6 MXP6C
- 7 MXP9S
- 8 MXP9SC
- 9 MXP9C
- 10 MXP75S
- 11 EFCL15

Channel Retractor

• Surgery for lower jaw

RTCN8



8mm

RTCN10



10mm

Spatula Periosteal Chisel

MXSP6

Spatula Periosteal Chisel
• 폭 6mm



6.0mm

Marked at 5, 10, 15, 20, 25, 30mm

MXSP9

Spatula Periosteal Chisel
• 폭 9mm



9.0mm

Marked at 5, 10, 15, 20, 25, 30mm

Retractors

Lip and Cheek Retractor


RTCRL

Lip Retractor, CRL



RTCRC

Cheek Retractor, CRC
• Columbia



Lip and Tongue Retractor

- New choice for earning a more comfortable impression.
- Can be used as a surgical retractor, too.

RTL5

Lip Retractor
• Using Lip retractor, you don't need to use multiple mirrors during upper jaw impression.



RTTG

Tongue Retractor
• Using both Lip retractor and Tongue retractor simultaneously, you can reduce the amount of work during a lower jaw impression.



Retractors

Cheek Retractor

Patent pending 2015-0173145

RTCRM

• Minnesota



RTAN20

Cheek Retractor,
ArmyNavy, RTAN-20
• Retract lip
• Parkman Design



Langenbeck Retractor

- Pull cheek or incised gum in order to secure a clear view during treatment.

RTP90-1

• Langenbeck Retractor



RTP30

• Retractor



Tunneling Instruments

- To separate coronal tissue while tunneling surgery. (periodontal plastic treatment or tunneling techniques for bone graft)

Tunneling Instrument

TITU1

- Anterior
- Tunneling for the flat area
- To take tissue or insert tissue into tunneling



TITU2

- Tunneling for the heavy curved area



TITU3

- Combination of TITU1 & TITU2



Tunneling Instruments

- To separate coronal tissue while tunneling surgery. (periodontal plastic treatment or tunneling techniques for bone graft)

Tunneling Instrument

TITU4

- 10° angle
- Similar to TITU1 but Knife is rounded.



TITU5

- 30° angle
- To keep expanding toward the curved area.



NEW TITU6



Palatal Wedge

An accessory to lift the palatal mucosa for incisional biopsy in the protruding connective tissues.

Palatal Wedge

NEW

PW16

• 10x7x1.6H (mm)
An accessory to lift the palatal mucosa for incisional biopsy in the protruding connective tissues



NEW

PW21

• 10x7x2.1H (mm)
An accessory to lift the palatal mucosa for incisional biopsy in the protruding connective tissues



Practice

Surgical method using Palatal Wedge




1. Cut the connective tissues to be collected and make space for the Palatal Wedge to be entered



2. Insert the Palatal Wedge into the space spread.



3. Make an incision in the connective tissue protruded by the Palatal Wedge and take it.



4. Graft the taken connective tissue where necessary.

Simple Extraction

Allowed to minimize damage to the surrounding bone during a tooth extraction for orthodontic treatment and implants.


Luxating Elevator

A luxating elevator has a thin and ultra-sharp blade than the standard elevator. It minimizes tissue trauma and preserves alveolar bone. Osung luxating elevator is redesigned based on field experiences and clinical practice to improve efficiency and safe clinical activities. An optimal handle design reduces the force required for gripping the instrument.



Arrangement

01. Anesthesia Syringe	SAF1	P.079
02. Periotome	PR2-2R, PRRS3	P.086
03. Luxating Elevator	3ELL303, 3ELL302	P.088
04. Forceps	FX151, FXX13	P.093,095



Process

SAF1
PR2-2R
PRRS3



01. Local anesthesia

02. Periodontal ligaments cutting

3ELL303
3ELL302



03. Luxation

FX151
FXX13



04. Extraction

Practice

01. Local anesthesia

Used

Local anesthesia syringe. Harpoon type syringe provides stable aspiration during nerve block anesthesia.


Character

Harpoon is designed to hold the rubber plunger of the cartridge and thumb ring is designed to make negative pressure for aspiration.

Anesthesia Syringe_SAF1

How to use

1. Choose a local anesthetic ampoule based on the patient's condition and the dentist's decision.
2. Check the validity period, whether the ampoule is cracked and the integrity of the rubber packing.
3. Make sure that the thumb ring of the syringe and the screw hub are tight.
4. Hold the syringe with one hand and pull the thumb ring back to insert a local anesthetic ampoule.
5. Remove the short side protection cap on the needle and secure it by screwing it onto the screw on the syringe hub. (At this time, the cap of the needle is not removed.)
6. Press the plunger that is thumb ring back so that the needle passes through the rubber septum. (Be careful not to bend the tip of the needle.)
7. Make sure that no air bubbles are generated.
8. Examine the treatment site.



Pull the rod back to insert the ampoule.



Push the harpoon firmly into the rubber membrane of the ampoule and attach the needle.



Pull the thumb ring to make negative pressure for aspirating.

Periotome_PR2-2R,PRRS3

How to use

Insert the blade into the periodontal sulcus along the root surface, severing the periodontal ligament directly below the alveolar crest.



Hold the PRRS3 with the pen grasp for the anterior teeth extraction.



Hold the PR2-2R with the pen grasp for the posterior teeth extraction.



Insert the blade into the periodontal sulcus along the root surface, severing the periodontal ligament directly below the alveolar crest.

02. Detecting subgingival calculus

Used

An extraction instrument used for cutting periodontal ligaments, preventing excess trauma to the interproximal papillae and marginal gingiva. It also can be used when considering the placement of a dental implant with minimal damage to the surrounding alveolar bone during the extraction process.

Character

Thin sharp blades to facilitate the removal of the teeth.

03. Luxation

❖ Used

Used for luxating the tooth from the periodontal socket reducing damage to a surrounding bone.


❖ Character

Luxating elevators have thin and sharp blades for cutting and separating the periodontal ligaments from the tooth. It reduces trauma during extraction as they used in a luxating motion, compared to a standard elevator that pries and lifts. Allows reduced pressure on adjacent teeth.

Luxating Elevator_3ELL303,3ELL302

How to use

Insert the tip between the root and alveolar bone and apply rotary motion to expand the socket. Cut the periodontal ligament and separate the root from the socket.



Hold it in the palm of the hand and support the shank with the index finger to control the forces applied to the elevator.

Curved shank provides easier access to the teeth in posterior region.

04. Extraction

❖ Used

Extraction forceps are used for extracting the teeth. Various sizes and dimensions are available as per the particular area of the mouth.

❖ Character

Plier type.

Forceps_FX151

How to use

1. Put the beak on the tooth surface with the handle opened wide.
2. Adapt the beak to the tooth with the handle closed.
3. Apply the force to allow the forcep to grasp the tooth and apply the constant rotatory force laterally to release the tooth from the alveolar bone.



Grip of the maxillary extraction forcep.



Grip of the mandibular extraction forcep.



Beaks should be adapted on the tooth and moved apically during extraction.

Forceps_FXX13

How to use

1. Put the beak on the tooth surface with the handle opened wide.
2. Adapt the beak to the tooth with the handle closed.
3. Apply the force to allow the forcep to grasp the tooth and apply the constant rotatory force to release the tooth from the alveolar bone.



Grip of the mandibular extraction forcep.



Beaks should be adapted on the tooth and moved apically during extraction.

Excision of Torus

Treatment that makes the alveolar bone shape gentle and soft when a patient feels uncomfortable due to benign osteophytosis on the buccal in the maxillary molar area, mandibular lingual and maxillary palate or exostosis on the basal bone, not the alveolar bone or in case of interfering with the denture stability and causing pain.


Needle Holder

Made of high-quality stainless steel which promotes corrosion resistance after a long period of use. Grasp more securely and open more smoothly. Tungsten carbide beak allows maximum grip and prevents needle rotation and slippage.



Arrangement

01. Anesthesia Syringe	SAF1	P.073
02. Scalpel Handle	SHF, SHS, SHC	P.074
03. Periosteal Elevator	EP9, EP9H	P.076
04. Bone Rongeur	RNG178	P.094
05. Periodontal Chisel&mallet	CHBC1, ML20	P.103, 095
06. Periodontal Chisel	CHCO1, CHCO2	P.102
07. Bone File	BF1X	P.095
08. Tissue Plier	PTS22C	P.099
09. Scissors	SCLSS115	P.101
10. Needle Holder	NHC150TC	P.097



Process

SAF1
SHF
SHS, SHC



01. Local anesthesia

02. Mucosal incision

EP9, EP9H
RNG178



03. Creating a mucoperiosteal flap

04. Removal of sharp bony edges

CHBC1, ML20
CHCO1, CHCO2
BF1X



05.06. Smoothing rough or sharp edges

07. Smoothing rough or sharp edges

PTS22C
SCLSS115
NHC150TC



08. 09. 10. Suture

Practice

01. Local anesthesia

Used

Local anesthesia syringe, Harpoon type syringe provides stable aspiration during nerve block anesthesia.


Character

Harpoon is designed to hold the rubber plunger of the cartridge, and thumb ring is designed to make negative pressure for aspirating.

Anesthesia Syringe_SAF1

How to use

1. Choose a local anesthetic ampoule based on the patient's condition and the dentist's decision.
2. Check the validity period, whether the ampoule is cracked and the integrity of the rubber packing.
3. Make sure that the thumb ring of the syringe and the screw hub are tight.
4. Hold the syringe with one hand and pull the thumb ring back to insert a local anesthetic ampoule.
5. Remove the short side protection cap on the needle and secure it by screwing it onto the screw on the syringe hub. (At this time, the cap of the needle is not removed.)
6. Press the thumb ring that is pulled back so that the needle passes through the rubber septum. (Be careful not to bend the tip of the needle.)
7. Make sure that no air bubbles are generated.
8. Examine the treatment site.



Buccal application on the maxillary posterior teeth



Palatal application on the maxillary posterior teeth



Labial application on the maxillary anterior teeth



Buccal application on the mandibular posterior teeth



Lingual application on the mandibular posterior teeth



Labial application on the mandibular anterior teeth

02. Mucosal incision

Used

Installing a blade, it can be used for the soft tissue incision or trimming the proximal restoration.


Character

It comes with a ruler making for measuring the size of your surgical sites.

Scalpel Handle_SHF

How to use

1. Insert a blade with Kelly holding the backside of a blade.
2. To remove the blade, lift the bottom of the blade carefully until it unlocks from the handle.



Use Kelly to hold the blade securely, align the blade into the grooves in the handle and gently slide it toward the scalpel handle until it locks in place. Always hold the back of the blade as it is extremely sharp.



To remove the blade, carefully lift the bottom of the blade until it unlocks from the handle.



Metric scale on the back side.

Used

Installing a blade, it can be used for the soft tissue incision or trimming the proximal restoration.

Character

Easy to mount the blade with a gentle push forward and backward in any direction. The cylindrical handle is advantageous for the curved incision. Easy to access to the maxillary posterior and palatal areas.

Scalpel Handle_SHS, SHC

How to use

1. Attach the blade to the handle by slipping the slit in the blade into the grooves on the handle using Kelly clamp.
2. To remove the blade, lift carefully bottom of the blade until it unlocks from the handle.



Cylindrical handle design for bidirectional blade mounting.



The cylindrical handle makes it possible to smoothly curve the incision with force applied to the handle with the pressure of finger movement without putting any pressure on your wrist.



SHC can easily access to the maxillary posterior and palatal areas.

03. Making a mucoperiosteal flap

Used

Used for detaching and lifting the flap.

Character

Combines a wide flat blade for raising or lifting muco-periosteal flap and a pointed sharp blade for detaching interdental papilla.

Periosteal Elevator_EP9, EP9H

How to use

1. Select the blade according to the incision size and shape.
2. Retract soft tissue with a not too strong force to minimize soft tissue injury.
3. The convex surface of the blade is toward soft tissue and detach the flap with a wide blade.



Detach interdental papilla using the sharp & narrow edge.



Retract soft tissue with not too strong force to minimize soft tissue injury.



A suture hole on EP9H anchors soft tissue flaps during suturing.

04. Removal of sharp bony edges

Used

Rongeur is used for the removal of bony fragments or soft tissue.

Character

It can deliver forces efficiently due to two hinges. Therefore, it lessens hand fatigue.

Bone Rongeur_RNG178

How to use

Place the bone to be cut between beaks and remove it by exerting the cutting force.



Remove the protuberant bony fragment.

05. The Removal of sharp bony edges

Used

The Removal of sharp bony edges

Character

Chisel with a beveled cutting blade on one side for reshaping a bone. Small stainless steel mallet provides effective striking.

Bone chisel&mallet_CHBC1, ML20

How to use

1. Fix the blade of a chisel outward.
2. Strike a flat striking surface of the chisel with a mallet to drive the blade.



Remove bony protrusions.



The beveled cutting blade on the one end.

06.07. Smoothing rough or sharp edges

Used

Used for trimming irregular bony surfaces. Ideal for the atraumatic removal of bone on adjacent teeth during the crown lengthening procedure. Allows easy removal of fractured tooth root and secondary inner flap.

Character

The semicircular blade is beveled on one side..


CHCO1-Up (The direction towards the occlusal surface)

CHCO2-Down (The direction towards the root)

Bone Chisel_CHCO1,CHCO2

How to use

Hold in the modified pen grasp to allow maximal control, use push stroke.



1. CHCO1-Up (The direction towards the occlusal surface)
2. CHCO2-Down (The direction towards the root)



Fix the flat surface of CHCO1 on the bony surface and remove the thin cortical bone.



Fix the flat surface of CHCO1 on the root surface and remove the thin cortical bone.

Used

Used for trimming the alveolar bone after tooth extraction. It also can be used for smoothing the sharp alveolar ridge away after the alveolar bone osteotomy or osteoplasty.

Character

There are a string of blades on the file. There are two kinds of blades, straight type and grid type.

Bone File_BFIX

How to use

Use a modified pen grasp, use push and pull motion to smooth the bone surface.



Straight cutting blades are used with a pull stroke, Grid type can easily applied to the area which is hard to access.



A rounded tip design of BFIX is useful for smoothing the small area.



BF45 is ideal for smoothing alveolar bone after the removal of many teeth at one time.

08, 09, 10. Suture

Used
Used for holding soft tissue to guide a suture needle.

Character
Useful for suturing soft tissue which is relatively firm. Wide contact area minimizes pressure on the tissue and a hole on the tips allows suturing needle to pass through soft tissue without slipping.

Tissue Plier_PTS22C

How to use

1. Place the plier between the thumb and index finger, while index finger helps guide.
2. Gently pick up the tissue in place.
3. Hold the suturing needle and insert to the hole of a tip.



PTS22C allows easy insertion of suturing needle into suture space. PT41 is ideal for holding thin soft tissue.



Used for holding soft tissue which is relatively firm.



Used for holding soft tissue to guide a suture needle.

Used
Scissors are used for cutting silk during suturing. Used for eliminating excessive soft tissue & granulation tissue.

Character
Similar to Iris.


Scissors_SCLSS115

How to use

Hold the scissors by placing the thumb and ring finger through the rings in the handles with the index finger supporting the shank, hook the silk on looked edge and cut it at an appropriate length.



Hook the silk into the U-shaped area.



Used
Used for holding a suture needle.

Character
A tungsten carbide beak with a grid-shaped blade makes the suture needle non-slippery.

Needle Holder_NHC150TC

How to use

1. Proper length of a needle holder for easy handling is 7~8 in. Fix the suturing needle in the jaws.
2. Hold the scissor by placing the thumb and ring finger through the rings in the handles with the index finger supporting the shank.



Proper length of needle holder for easy handling is 18~20cm. FIX the suturing needle in the jaws.



A tungsten carbide beak for the easy control of a needle



Products for Dentistry

OSUNG Catalogue 2020/2021

Implant

Products for Dentistry

OSUNG Catalogue 2022°2023



IMPLANT / 임플란트


Fixture Implantation	Lindemann Drill	132
	Implant Depth Gauge	134
	Caliper	134
Bone Graft	Bone Spreader	135
	Bone Expander Hand Kit	136
	Bone Expander Engine Kit	137
	Micro Saw Shield	140
	Micro Saw	142
	Trephine Bur	144
	Narrow Bone Removing Bur	145
	Ridge Split Bur	145
	Lateral Approach Bur	145
	Surgi-Drill Stand	145
	Convex Osteotome	146
	Concave Osteotome	147
	Bone Scraper	148
	Block Bone Clamp	148
	Bone Collect Chisel	149
	Bone Collector	149
	Hexa Wrench	149
	Bone Mill	150
	Bone Crusher	150
	Bone Crusher Mallet	150
Bone Syringe	151	
Bone Well	152	
Bone Carrier	152	
Bone Packer	152	
Membrane Forceps	153	
Sinus Rongeur	153	
Sinus Lift	Crestal Approach Kit	154
	Lateral Approach Kit	156
	Sinus Lift	158
	Bone Screw	160
	Bone Tack	161
Bone Tack offset Holder	161	
Implant 2nd Surgery	Hand Tissue Punch	162
	Tissue Punch	163
Implant Crown Setting & Maintenance	Screw Removal Kit	164
	Implant Curette	166
Option	PRF & GRF Box	167
	New Product	168




Lindemann Drills

Lindemann Drill

LLS21




LLS26




LLL21

As it is long enough, it is very useful when a drill cannot access to the aiming position easily, due to the proximal teeth.




LLL26

As it is long enough, it is very useful when a drill cannot access to the aiming position efficiently due to the proximal teeth.



Practice

- Specially-designed blades make excellent cutting power in verticality and horizontality.
- Useful for multiple implants procedures. Specially used for relocation and redirection.
- Effective for thick cortical bone cutting.
- Effective for the site preparation of socket for an immediate implant.



Side cutting at the ridge of a socket

Change the path and cut sidewall of a socket

Rotating Instrument List

Rotating Mechanism List

Components

No	Product	Shape		Page
1	Lindemann Drill	Ø2.1 LLS21	The drill capable of deletion in the lateral direction as well as in the vertical direction.	132
2	Bone Expander		Used for the ridge expansion	137
3	Micro Saw Shield	NEW	Saw shield for safe bone cutting	140-141
4	Micro Saw		For bone cutting	142-143
5	Trephine Bur	Ø4.0/Ø3.0	For bone harvesting	144
6	Narrow Bone Bur	NEW	Used to flatten a narrow or irregular bone width	145
7	Ridge Split Bur	NEW	Used to remove the remaining bones on both ends during Ridge split	145
8	Lateral Approach Bur	NEW	Used for window formation during Sinus craft	145
9	A.I. Dill		A functional drill that stops itself when it encounters the maxillary sinus.	154-157
10	Tissue Punch		Used to remove wipes	163
11	Crew Removal	Ø1.5	Tool for removing the broken screw	164-165
12	Diamond Bur		Bur for tooth preparation	224-286
13	Zirconia Removing Bur	NEW	Bur for removing Zicornia	295

Implant Depth Gauge · Caliper

Implant Depth Gauge

DG1

- Used for checking the depth of the socket.
- Can also be used for checking if the sinus membrane is perforated.
- Gradation scale by 24mm with 3mm interval.



Caliper

LPC90


- Caliper, Castroviejo**
- Castroviejo 30 Angled
 - Length 90mm the scope of measurement is 0-40mm



COMING SOON

LPC90S

- Caliper, Castroviejo**
- Castroviejo 30 Angled/Castro Viejo measuring instrument with the front part be to a curved shape for about 30 degrees.
 - The total length 95mm, the measurement range 0-40mm
 - Used to measure dimensions such as bone or tooth size, spacing between teeth, and equal spacing of teeth




Bone Spreaders

Bone Spreader

BOC22

- First, make alveolar bone fractured with a mess or a chisel. Then use the spreader for expanding.



BOC28



BOC35




BOC35R




Character

1. Easy to operate with supporting bar on the handle
2. Offset design for easy access in posterior
3. Gradation mark to measure the depth
4. Use a mallet if necessary

CODE	A	B	C
BOC22	0.5mm	1.6mm	2.2mm
BOC28	0.5mm	1.6mm	2.8mm
BOC35	0.5mm	1.8mm	3.5mm
BOC35R	0.8mm	2.55mm	3.5mm



Practice



Cut alveolar bone with a disk (saw)

1. Malleting the spreader at the incision site, decide how deeply you need to put the spreader according to the size of the implant fixture.
2. Sway the spreader back and forth to expand the incision bone.
3. Twist the spreader to expand the incision bone wider.

Bone Expander Hand Kit

BEPD

• Size 160 x 85 x 65H(mm)

- Used for expanding bone width in case of narrow alveolar bone.
- Used for getting stable osseous tissue by compressing bone toward lateral when it is difficult to get Initial fixation due to low alveolar bone density.



Ratchet Wrench
Used for tightening or loosening the expander.

Engine Adapter
For using the expander with an engine.

Length Extension
Used for extending the length of the expander.

Point Drill Ø1.3
Used for making a position where a fixture is to be placed.

Hand Driver
Used for tightening or loosening the expander by finger instead of a ratchet wrench.

Expander

(mm)

	<ul style="list-style-type: none"> • After the pilot drill used, work first with a 2.3mm expander checking the path & depth. • Then gradually use a larger expander. • The final size of the expander has to be decided according to the size of a fixture. 						
		Apical	Ø1.0	Ø1.4	Ø1.8	Ø2.2	Ø3.0
	Diameter	Ø2.3	Ø2.7	Ø3.1	Ø3.5	Ø4.3	Ø5.1

Bone Expander Engine Kit

BEPD2

• Size 160 x 85 x 65H(mm)

- How to choose the final size of an expander?
1) Normal bone density : choose the same size as the fixture size.
2) Poor bone density : choose two step smaller size than the size of the fixture.



Length Extension
Used for extending the length of the expander

Hand Driver
Used for tightening or loosening the expander by finger instead of a handpiece.

Point Drill Ø1.8
Used for making a position where a fixture is to be placed.

Ratchet Adapter
Connect an expander with a ratchet wrench

Expander

(mm)


	<p>First, use a point drill to make a position. Then use the smallest size, 2.3mm expander with in & out motion at 50rpm. After that gradually use a larger expander.</p>						
		Apical	Ø1.0	Ø1.4	Ø1.8	Ø2.2	Ø3.0
	Diameter	Ø2.3	Ø2.7	Ø3.1	Ø3.5	Ø4.3	Ø5.1

Bone Expander Hand Kit · Bone Expander Engine Kit

Practice

Used for expanding bone width in case of narrow alveolar bone

1. Drill to the proper depth using a point drill at the implant site. (800 ~ 1,200rpm)
2. Expand the bone to the desired size using gradually from a small expander to a large expander.
 - After drilling, cut the bone with a disk(saw). Then use the bone expander.
 - Use a ratchet wrench for expanding and a hand wrench when taking out an expander.



Implantation is impossible due to narrow bone width

Initial drilling on the implantation site


Expand bone width using from a small expander to a large expander gradually

Fixture placement

Practice

Used for stable osseous tissue with pressing bone toward lateral when it is difficult to get Initial fixation due to low bone density

1. Drill to the proper depth using a point drill at the implant site (Engine speed : 800-1,200rpm)
2. Insert a fixture after expanding the bone to the desired size using from small a screw to a large screw gradually
3. Please check whether the bone is pressed too much. (need to cut the bony bone surface using a reamer if bleeding doesn't occur.)



Worry about initial fixation due to low bone density

Decide the path and implant position with a point drill


Press bone toward lateral using an expander

Getting better initial fixation with higher bone density

Bone Expander Hand Kit · Bone Expander Engine Kit

Practice

Diagram of ridge split






Micro Saw Shield

Micro Saw Shield

Design Application: 30-2020-0015736
Patent Application: 10-2020-0042753

Safe & excellent surgical procedure
Patent Application

- The shield cover with a 360 - degree revolution can approach in the oral cavity freely
- The shield cover rotates to fit on the surface of the bone. It enables you to do the safe surgical procedure.



Micro Saw Shield



It can be combined with various handpieces such as KaVo, W&H, NSK, etc.

Shield




Micro Saw - Contra Angle Type



Micro Saw

Micro Saw



Contra Angle




Straight Angle



Micro Saw


Micro saw is a rotary cutting instrument for removing unnecessary bone torus or extremely thin ridge crest. Also it can be used for ridge split or block bone surgery.

Micro Saw



Use a contra angle handpiece for cutting bone in widthwise direction.

In case of expanding a little, the lengthwise cutting is not necessary.




Use a straight handpiece for cutting bone in lengthwise direction in supplement for wider expanding.

Do make lengthwise cutting for wide expansion.

Advantages compared to competitors


- Holes for lessening bone heating.
- Saw thickness : 0.35mm including diamond coating layer.
- No shaking and superior cutting power.
- The vibration is minimized during bone cutting because of the diamond coating surface.
- Possible to change direction gradually during bone cutting.

Contra angle type



Order No.	MICSA08	MICSA10	MICSA13
Diameter	ø8	ø10	ø13

Straight angle type




Order No.	MICSS08	MICSS10	MICSS13
Diameter	ø8	ø10	ø13

Trephine Burs


Trephine Bur

- 3.0/2.0mm trephine bur (THB30) can be used as an initial drill harvesting bone.
- Make a path with reverse rotation and use a normal rotation at 800-1,200rpm
- Used for the harvesting of mandibular bone.
- Choose a proper size of trephine bur for failed fixture removing.


Trephine bur can be initially used for harvesting bone without using a drill. In this case, it is the required size of the Trephine bur.




THB30
• For bone harvesting




THB40
• For bone harvesting




THB50
• For bone harvesting




THB60
• For bone harvesting



THB70
• For bone harvesting




THB80
• For bone harvesting




THB90
• For bone harvesting

To remove implant



THR50 ₩ 70,000
• To remove implant




THR60 ₩ 70,000
• To remove implant

Implant Surgical Bur · Surgi-Drill Stand

Implant Surgical Bur


NEW
LABEC62

- Lateral approach bur**
- For Angle Handpiece
 - Used for window formation during Sinus graft
 - Less splash than cutting burs due to the electrodeposited diamonds.




NEW
NBRBEC50

- Narrow bone removing bur**
- For Angle Handpiece
 - Used to flatten the narrow bone width after pulling teeth.
 - Less splash than cutting burs due to the electrodeposited diamonds.



NEW
RSBTRM

- Ridge split bur**
- For Angle Handpiece
 - Used for horizontal osteotomy to remove the remaining bone at both ends after Ridge split using a saw.
 - More suitable for vertical osteotomy due to its accessibility.




Surgi-Drill Stand

- Surgi-Drill Stand is a perfect solution for managing surgical drills and burs.
- The stand has 16 multi-silicone-holders and the multi-holder can hold any kinds of drills, low speed burs and high speed burs.
- Use the cover as a mini-tray as a stand is opened.

DSTA16

- SIZE 99 x 40 x 45H(mm)




P.358

Convex Osteotomes

BOVXSET


Convex Osteotome 5ea + Cassette



It has a stopper for safe osteotome technique.

A convex osteotome is used for the stable initial fixation of an implant through condensing alveolar bone laterally.

Convex Osteotome Parkman Design




Concave Osteotomes

BOCVSET

Concave Osteotome 6ea + Cassette

EFS8B

Instrument Cassette



A concave osteotome is used for the sinus elevation surgery. It can be used when the alveolar bone is very thin up to 1-2mm. The floor of the sinus is then lifted by tapping the sinus floor with the use of osteotomes.

Concave Osteotome Parkman Design



Bone Scraper · Block Bone Clamps

Bone Scraper


BSF

- Sterilized product. Single-use.
- Size 148 x 17 x 16H(mm)
- If the blade becomes dull, rotate the blade with any 1.2mm hex driver to expose a new sharp edge.
- 5 pcs



The bone scraper comes packed as sterile products. It is very convenient to use.

Practice




Rotate the blade of a bone scraper with a 1.2mm hex driver

Rotate the blade by 90 degree to have a new sharp edge.

Harvesting bone

Harvested bone particle in the chamber.

Clinical case



Gingival flap opened

Taking bone

Flap gingival

Fill bone

Fix with bonetack after covering membrane


Fixed with bonetack

Sutured

Block Bone Clamp

RCA197

- Anterior
- Used for holding block bone
- The length is 197mm



RCP200

- Posterior
- Used for holding block bone
- The length is 200mm
- Drilling available through the slot of the beak

Bone Collect Chisel · Bone Collector · Hexa Wrenches

Bone Collect Chisel

STSB-1

- Used for collecting osseous coagulum, Back-action



Bone Collector

ST1

- Bone Collector
- Used for collecting bone particles cut off while drilling.
- To take osseous coagulum inhaled by suction using a filter.
- Length 235mm
- The filter is single-use.



ST1-F

- Bone Collector Filter
- Disposable



Practice

Advantages of using bone collector

- Collect usable bones during implant drilling.
- Reduce the extra surgery for getting bone graft from a patient.
- Can save bone graft material.
- Can save surgery time.

Notice

- Avoid inhaling saliva to the collector through suction.
- Use two suctions for the bone collector and inhaling saliva.
- Keep collected osseous coagulum in gauze wet with saline solution and place it to the bone defect area as soon as possible.

Hexa Wrench

IDH5-15N

- 15N Torque Wrench
- Alert not to give excess power as it is bent if power is over 15N



IDH7-15N

- 15N Torque Wrench
- Alert not to give excess power as it is bent if power is over 15N




Bone Mill · Bone Crusher · Bone Crusher Mallet

Bone Mill

BMH


- Hinge type
- Serrated disks are made of titanium.
- Place a small bock bone between two serrated disks for grinding so the particles are minutely broken.
- Easy & simple to use.



Bone Crusher

BCR3


- The punching head has a round shape to effectively crush bone with a relatively small force. And the punching tool doesn't get jammed to the bowl due to its rounded shape.



Bone Crusher Mallet

ML29

- Durable and strong stainless steel mallet. Suitable for bone crushing procedure.



Bone Mill · Bone Syringes

Bone Syringe

- Placing graft material into recipient site / stainless steel made
- To carry grained bone into the recipient site


BSY35

- For the socket approach
- Useful for fine & high-flowable bone particles.
- Inner Diameter 2.5mm / Outer Diameter 3.5mm



BSY47

- Most popular size.
- Inner Diameter 3.7mm / Outer Diameter 4.7mm



BSY70

- Used for lateral approach. It is able to carry large volume of material.
- Inner Diameter 6.0mm / Outer Diameter 7.0mm



Practice



1. Firstly pull the bone syringe knob and move up and down the bone syringe to fill up.



2. Lastly, push the knob to extract the bone from the bone syringe.

Implant

Bone Well · Bone Carrier · Bone Packers

Bone Well

BWSUS1

- Used for mixing bone graft material
- Made of stainless steel
- Size : Diameter 38 x 27.5H(mm)



Bone Carrier



BSC3539

- To carry bone graft material
- Titanium made

Bone Packer

GP2530



GP3340



Implant

Membrane Forceps · Sinus Rongeur

Membrane Forceps

MF01

- Used for holding & carrying membrane effectively
- Length 121mm / Width 5.0mm



Sinus Rongeur

RNGSK100

- Kerrison
- To open sinus window
- Length of beak 10cm
- Width of bite 4mm



Crestal Approach Kit




OCAKIT

• Size 260 x 165 x 75H(mm)

- Everyone gets the same result in operation as the drill stops by itself when it touches the sinus membrane.
- Riskless in the matter of membrane-tear as the membrane is lifted by using hydraulic pressure



Video Clip



Aquainjector

- Use together with a disposable plastic syringe.
- Designed to generate strong hydraulic pressure enough for a sinus lift.
- Also, this device is intended for injecting liquid slowly to give membrane recovery -time from the strain formed inside.
- Possible to know the volume of space for bone grafting.



1	Pilot Drill $\phi 1.8$		• To mark a position of an implant after confirmation of X-ray			
2	Stopper (2mm ~ 10mm)		• Control depth			
3	A.I. Drill		• Drill stops automatically when the drill touches sinus membrane (by 1,200rpm)			
4	Cannon Drill $\phi 2.4$		• Initial drill to access the sinus before the use of A.I. drill.			
	Cannon Drill $\phi 3.2$		• Used for expanding a hole size. • Can be used with a stopper.			
5	Aqua tap		• Used for injecting a saline solution or contrast medium into the perforated maxillary sinus			
			Diameter	3.15S	3.15	3.5
			Length	28mm	33mm	33mm
6	Hand Adapter					
7	Engine Adapter		• Adapter for connecting Aquatap to the implant handpiece.			
8	Hand Wrench					
9	Sensor Gauge		• To check perforation of the sinus membrane by using compression of spring			
10	Bone Condensor		• Condensing bone material into the maxillary sinus			
11	Bone Carrier		• Carry bone material to the maxillary sinus			

Crestal Approach Kit

Practice

01 Pilot Drilling

- $\phi 1.8$ 1,200rpm
- Mark a drilling position on the cortical bone with a pilot drill.

02 Attaching a stopper

- Take a stopper which is 1mm longer in length than the length measured from alveolar bone to sinus membrane.
- Then attach the stopper to the Canon drill.

03 Cannon Drilling

- Run the Canon drill 2.4mm with a stopper at 1,200rpm.
- Use 3.2mm drill if necessary for a large hole.
- Drill by the spot which is 1mm away to the sinus membrane.

04 Cannon A.I. Drilling

- $\phi 2.8$ 1,200rpm
- Select the stopper which is 2mm shorter than the stopper used with Canon drill.
- Attach the stopper to the A.I. drill, and run it at 1,200rpm.
- The A.I. drill stops automatically when the end of the drill reaches the membrane.

05 Checking the perforation

- Attach the stopper which is used with A.I. drill.
- Check the perforation with the sensor gauge.

06 Aqua Tapping

- $\phi 3.15 \sim \phi 3.5$ 30rpm
- Attach the same stopper which is used with A.I. drill.
- Fix the Aqua tap into the bone by using a handpiece.

07 Adjusting the height of Aquatap – Hand wrench

- To adjust the height of an Aquatap by finger minutely.
- Use with a tube.

08 Connecting with Injector

- Attach a 10cc plastic syringe to the Aqua injector and connect the Aqua tap to the plastic syringe by a flexible hose.

09 Lifting sinus membrane

- Inject liquid slowly with a speed of 1 click/10sec to give membrane recovery time from the formed inside.
- To check penetration and the strain amount of sinus membrane lifting, use contrast medium.
- Decide the amount of Bone powder by calculating the difference of liquid volume in a syringe from beginning to end.

10 Condensing bone graft material

- Carry and push graft material into the socket

Lateral Approach Kit




OLAKIT

• Size 260 x 165 x 75H(mm)

- Everyone gets the same result in operation as the drill stops by itself when it touches the sinus membrane.
- Riskless in the matter of membrane-tear as the membrane is lifted by using liquid.



Video Clip



- Aqua Lifter**
- Use together with a disposable syringe.
 - Inject 1.0~1.5cc of saline solution or contrast medium.



- A.I. Drill_Straight**
- For straight angle.
 - It stops automatically when the end point of the drill touches the membrane. (Running speed : 6,000rpm)



- A.I. Drill_Contra Angle**
- Extra shank for the contra-angle handpiece.
 - Combine the shank with the upper structure of Lateral A.I.Drill

- Sinus Lift**
- Useful for separating the sinus membrane after making a window.



- Burin Drill**
- A smooth ball-ended tip minimizes the damage of the membrane.
 - Cut the alveolar bone using side blades after detaching the membrane.




- A.I. Drill Dismantler**
- Tool for changing Lateral A.I.Drill to Contra Angle.




- Wide Burin Drill**
- A smooth ball-ended tip minimizes the damage of the membrane.
 - Cut the alveolar bone using side blades after detaching the membrane.

Lateral Approach Kit

Practice




- Lateral AI Drilling**
- Position a hole in the front and lower area if possible as it is the best position to lift membrane safely and effectively.




- Lateral Aqua Lifter**
- Fill 1.0~1.5cc saline solution in the chamber of a plastic syringe.
 - Attach the Aqua lifter to the plastic syringe and then engage the end part of the Aqua lifter into the hole.
 - Inject the solution to elevate the sinus membrane.



- Burin Drilling**
- After Detaching the sinus membrane, use the Burin drill to open the window.



- Enlarging the size of window**
- Use wide Burin drill to enlarge



- Sinus Membrane And Elevation**
- Lift sinus membrane from the side and the lower area through the window.



- Bone Graft**
- Fill graft material.

ISSL1 Left Tip



1. Detach the sinus membrane and tidy up the area around the window.

ISSL2 Left & Right Tips



2. Detach the sinus membrane of the lower area of a hole.

ISSL3 Left & Right Tips



3. Separate sinus membrane of right & left of a hole.

ISSL4 Left & Right Tip



4. Separate sinus membrane by long distance in the left and right of a hole.

ISSL1 Right Tip



5. Carrying bone graft material

ISSL3 Left & Right Tips



6. Put bone graft material into the window.

Sinus Lifts

Sinus Lift

• It is for detaching the sinus membrane from the bone.

IS6577SC5

• Acute angled



IS65785

• Obtuse angled



ISPKN152

• Used for beginning the delicate separation of the sinus lining.



ISSC1



ISSC2

• Long blade type of ISSC1 to access to the deep inner space.




Sinus Lifts

Sinus Lift

• It is for separating the sinus membrane from the bone.

ISSLKIT

Sinus Lift Kit
• ISSL1/ISSL2/ISSL3/ISSL4
• 72 x 192 x 22H(mm)



ISSL1

Sinus Lift
• The left tip is used for detaching the sinus membrane.
• The right tip is used for carrying bone graft material.



ISSL2

Sinus Lift
• Used for detaching the sinus membrane on the mesial & distal area.



ISSL3

Sinus Lift
• Used for detaching membrane around a hole.



ISSL4

Sinus Lift
• Used for detaching membrane by long distance of left & right surrounding a hole.



Bone Screw

Bone Screw Kit

This has a double-thread structure that ensures quick insertion and good settlement in a bone. Further, ultra-precision machining that cuts threads at an accuracy of 3/1000mm will warrant accurate engagement with the driver.


BSKIT

- This product is made of titanium GR5 ELI.
- You can buy these bone screws as a full-set kit or several pieces of your desired sizes separately.

- 1 Bone Screw ϕ 1.5mm
- 2 Bone Screw ϕ 1.7mm
- 3 Bone Screw ϕ 2.0mm
- 4 Drill
- 5 Driver
- 6 Driver Holer
- 7 Case



Products Constitution



Consist of product

Name	Size	Code	(ea)
Bone Screw ϕ 1.5mm	4mm	BSW15-004	9
	5mm	BSW15-005	6
	6mm	BSW15-006	3
	8mm	BSW15-008	3
	10mm	BSW15-010	3
	12mm	BSW15-012	3
Bone Screw ϕ 1.7mm	4mm	BSW17-004	6
	5mm	BSW17-005	4
	6mm	BSW17-006	2
	8mm	BSW17-008	2
	10mm	BSW17-010	2
	12mm	BSW17-012	2
Bone Screw ϕ 2.0mm	4mm	BSW20-004	6
	5mm	BSW20-005	4
	6mm	BSW20-006	2
	8mm	BSW20-008	2
	10mm	BSW20-010	2
	12mm	BSW20-012	2
Drill	1.0mm	BSWDR1.0	1
	1.3mm	BSWDR1.3	1
	1.6mm	BSWDR1.6	1
Driver	For hand driver	BSWDTL	1
	For contra-angle	BSWDTS	1
Driver Holer	Hand Driver Holder	BSWDH	1
Case	Aluminium Case	BSWDC	1

Bone Tack

Bone Tack Kit

This Bone tack is designed to be driven in using a mallet at the time of insertion while being screwed out using a driver at the time of removal, with a view to fixing membranes or foils.


BTKIT

- This Bone Tack has a screw structure directly under its head, which allows for easy removal using the driver supplied with the product.
- It is made of titanium GR5 ELI and is available in 3mm and 5mm in size. 3mm Tacks are for normal bone while 5mm Tacks for soft bone.

- 1 Bone Tack Holder
- 2 Case
- 3 Bone Tack M0.85 x 5mm
Bone Tack M0.75 x 3mm
- 4 Hexa Driver 0.9mm



Products Constitution



Consist of product

Name	Code	(ea)
Bone Tack Holder	BTSHC(Straight)	1
Case	BTSC	1
Bone Tack M0.85x5mm	BTS85-50	5
Bone Tack M0.75x3mm	BTS75-30	9
Hexa Driver 0.9mm	BTIDH09	1

Bone Tack Offset Holder

BTSHCO

- TOffset
- For posterior



Tissue Punches

Hand Tissue Punch

It is useful to operate as a surgeon's intention.

TPH35S
• Hand Tissue Punch, Straight



TPH35C
• Hand Tissue Punch, 90 Angled



TPH40S
• Hand Tissue Punch, Straight



TPH40C
• Hand Tissue Punch, 90 Angled



TPH45S
• Hand Tissue Punch, Straight



TPH45C
• Hand Tissue Punch, 90 Angled



TPH50S
• Hand Tissue Punch, Straight



TPH50C
• Hand Tissue Punch, 90 Angled



Tissue Punches

Center Guide Inclined Tissue Punch

TPI40G
• For Handpiece (15rpm)
• Inner dia 4.0mm



TPI50G
• For Handpiece (15rpm)
• Inner dia 5.0mm



TPI40
• For Finger
• Inner dia 4.0mm



TPI50
• For Finger
• Inner dia 5.0mm



Center Guide Tissue Punch

• This tool is mainly used in the primary surgery and especially useful in flapless implant surgery.
• When a path is created by a 2mm guide drill, use the center guide to set the tissue punch in the right position at the speed of 15rpm.

TP40CJ
• Inner dia 4.0mm



TP50CJ
• Inner dia 5.0mm



Tissue Punch

Used for removing tissue without unnecessary trauma.

TP35
• Inner dia 3.5mm



TP40
• Inner dia 4.0mm



TP45
• Inner dia 4.5mm



TP50
• Inner dia 5.0mm




Practice

Inclined Tissue Punch

• Inclined blade can cut the tissue clearly which a normal tissue punch is unable to do. It is generally hard for a normal tissue punch to cut clearly because the bone is not flat structurally.

▼ The picture shows when the tissue has been clearly removed using a center guide tissue punch and an inclined tissue punch.




Screw Removal Kit

Used for removing a broken screw in the fixture safely & speedily.

Screw Removal Kit

OSRKIT

• Size 160 x 85 x 65H(mm)



Tap Repair
• Consist of M1.6, M1.8, and M2.0

Ratchet Wrench
• For abutment hex removal and tap repair

Internal Guide Handle
• Used for holding an internal guide firmly.

Claw Removal
• Consist of Ø1.1S, Ø1.1L, Ø1.5S, and Ø1.5L for removing a broken screw inside a fixture

	Ø1.1S	Screw for dia Ø1.6
	Ø1.1L	
	Ø1.5S	Screw for dia Ø2.0
	Ø1.5L	



Hand Driver

• For using a slot driver and claw removal by hand



Slot Driver

• For damaged Hex of healing abutment, cover screw and abutment screw




Abutment Hex Removal

• used for removing octa or broken hex of abutment



Internal Guide

• Consist of IH2.4, IH2.5, and IO3.1
• It is used with claw removal Ø1.1L or Ø1.5L as a guide



Screw Removal Kit





Practice

1. Remove fractured one-piece Abutment



01. Make a linear slot using over 0.8mm bur.



02. Fit a slot driver into the preformed linear slot and turn counterclockwise to remove the broken abutment.


2. Remove fractured Abutment Screw



01. Fractured screw




02. Hold guide handle not to move internal guide.



03. Put claw removal to reach the screw through an internal guide.

3. Remove fractured Abutment Hex



01. Fractured abutment hex



02. Turn Abutment Hex removal clockwise until it hangs to fractured Abutment hex (Fixture can be moved if you put excessive power)




03. Remove Hex.

4. Renew of interior Tap for damaged fracture



01. Do not harm by heat injecting saline continuously.



02. Refine internal tap with a force of 10-30N.



03. Remove residue in the tap. (2-3 times repeat)



Before



After

Implant Currettes

Implant Curette

• Soft titanium is used as a material for scratch-free and contamination-free on the implant surface.

BEST

ICGR1-2

• For Anterior



ICGR5-6

• For Anterior



ICGR7-8

• Used on the buccal and lingual portions of posterior teeth



BEST

ICGR11-12

• Used on the mesial portions of posterior teeth



BEST

ICGR13-14

• Used on the distal portions of posterior teeth



PRF & GRF Box

It is joint work with dentists and designed to handle growth factors that have been extracted by a centrifuge.

PRF & GRF Box

Parkman Design

GRF

Platelet-Rich-Fibrin & Growth Factor-Rich-Fibrin Box, PRF & GRF Box

- Used for handling PRF and GRF
- Size 155 x 108 x 60H(mm)



Practice



Gained PRF by a centrifuge.



Flattened like membrane by the press.



Also can make a ball shape for socket treatment.




A ball shaped PRF.



The liquid plasma is gathered in the bath.

Implant Basic Kit



No	Product	Name	Code	Page
1	DIAGNOSIS	Mirror(Rear Surface)	DMCS4	18p
2		EX-Probe	XP23-8	28p
3	DEBRIDEMENT	Surgical Curette	URCM10	92p
4		Periodontal Curette	URPR1-2	103p
5	INCISION & FLAP	Silicone Scalpel Handle	SH2S	74p
6		Periosteal Elevator	EP24G	77p
7		Surgical Curette	URCM2-4	92p
8	OTHERS	Periosteal Elevator	EPPR3	76p
9		Titanium Suction Tip	SN4TI, SN3TI	71p
10	OTHERS	Bone Well	BWSUS1	140p
11	ANESTHESIA	Anesthesia Syringe	SAF1	74p
12	DIAGNOSIS	Wide Tweezer	PCW150	29p
13	SUTURE	Dean Scissor	SCD170	100p
14	OTHERS	Towel Clamp	CPTC135	72p
15	SUTURE	Needle Holder	NHC150TC	97p
16	DEBRIDEMENT	Hemostat	HTM130C	96p
17	OTHERS	Caliper	LPC90	122p
18	SUTURE	Tissue Plier	PT42	99p
19	OTHERS	Saline Bowl	SALB-10	
20		Instrument Cassette	EFCL1	329p
21		Wrapping Cloth	WR7575	332p
22		Surgical Drape	WDMA	332p

Products for Dentistry

OSUNG Catalogue 2022/2023

Restorative

Products for Dentistry

OSUNG Catalogue 2022°2023



RESTORATIVE

/ 수복

Cavity Preparation	Excavator	172
	Gingival Retractor	175
	Margin Trimmer	176
Amalgam	Amalgam Carrier	177
	Amalgam Well	177
	Amalgam Plugger	178
	Carver	179
	Amalgam Burnisher	180
Composite Resin	Measuring Instrument	181
	Placement	181
	Composite Instrument	182
	Composite Instrument Kit	190
Amalgam Filling	Manual	191
	Resin Filling	197



Excavators

Used for removing carious dentin. Also used for carving amalgam and direct wax pattern.

Excavator_Plastic Handle

Autoclavable

3EXC38-39



BEST
3EXC17



3EXC18



3EXC63-64



3EXC65-66



Excavators

Excavator_Meatal Handle

EXC38-39



BEST
EXC17



BEST
EXCE1



EXC18



EXCE2



EXC63-64



EXC65-66



Excavators

Excavator_Plastic Handle

Autoclavable

3EXCL61-62



BEST

3EXCL63-64



3EXCL65-66



Excavator_Meatal Handle

EXCL61-62



BEST

EXCL63-64



EXCL65-66



Gingival Retractors

Useful for separating and protecting gingiva during cavity preparation or resin filling.

Gingival Retractor_Plastic Handle

Autoclavable

3GRM1

• For lower incisors



3GRM2

• For canine & premolar in upper & lower



3GRM3

• For upper central incisor and molar teeth in upper & lower. It can be used for a wide canine.



Gingival Retractor_Meatal Handle

GRM1

• For lower incisors



GRM2

• For canine & premolar in upper & lower



GRM3

• For upper central incisor and molar teeth in upper & lower. It can be used for a wide canine.



Practice



• Used for protecting gingival tissue while cavity preparation or resin filling.

Margin Trimmers

Margin Trimmer

- Used for making proper bevel on enamel margins

MT26

Margin Trimmer, MT26
• Distal



BEST

MT27

Margin Trimmer, MT27
• Mesial



BEST

MT28

Margin Trimmer, MT28
• Distal



MT29

Margin Trimmer, MT29
• Mesial



Amalgam Carriers • Amalgam Well

Amalgam Carrier

PM1520

- Mini/Regular



PM2025

- Regular/Large



PM2030

- Regular/Jumbo




Amalgam Well

- To store amalgam before it is placed in preparation.
- Designed for easy amalgam handling.

PLGWL1

- Can be available as a bone well.



Amalgam Pluggers

Amalgam Plugger

• Pluggers are used to condensing filling materials into cavity preparations.

BEST

PLG0-1

• Serrated tip-end



BEST

PLG1-2

• Serrated tip-end



PLGOR1

• Oregon 1
• Easy to access inside wall of the cavity
• Plain tip-end



PLGOR3

• Oregon 3
• Easy to access inside wall of the cavity



Carvers

Carvers

• To carve occlusal shape or remove excessive condensed material.

BEST

CVCD89-92

• Cleoid Discoid



CVCD3-6

• Cleoid Discoid



BEST

CV3S

• Hollenback 3S



CV3

• Hollenback 3



CV74-75

• Can be used to cut soft tissue during surgery



CV76-77



Amalgam Burnishers

Amalgam Burnishers

• Used to condense, smooth, and polish amalgam.

BB26-27S

• Ball Burnisher



BEST

BB27-29

• Ball & Egg (Football) Burnisher



BB31-32

Burnisher, BB31-32
• Ball Burnisher
• 직경 1.8mm/직경 2.4mm



BBL3

• Ladmore 3



Measuring Instruments · Placements

Calibration Instrument

NEW

MIVD-01

Calibration instrument
• Instrument for measuring depth and width of cavity.
• For the dental hygiene students
• 0.4 mm diameter x 1.5 mm length
• 0.75 mm diameter x 1.5 mm length



NEW

MIVD-02

Calibration instrument
• Instrument for measuring depth and width of cavity.
• For the dental hygiene students
• 1.0 mm diameter x 2.5 mm length
• 1.5 mm diameter x 4.0 mm length



Practice



Depth, thickness measurements



Clearance measurement



Depth measurement

Placement

• Used to apply calcium hydroxide or liner in the cavity.

• Also useful as a small burnisher

PIS

• Metal Handle/
Single-End



BEST

PICH

• Calcium Hydroxide Placer
• Dycal Applicator
• Metal Handle / Double-Ended



Practice




Mixing dycal



Apply base and liner like calcium hydroxide or glass ionomer to cavity

Composite Instruments

Science & Technology



The technical core of the Composite Instrument is the surface roughness of the working end. Plainly said, the smooth surface makes the resin materials do not stick on it. This is the operating principle of the Composite Instrument.

All Composite Instruments on the market today are made with this principle. They are highly polished and coated for better performance. In fact, a more important factor is the polishing. Since the coating is too thin that it is about 1/400 of the hair thickness, the role of the coating is only an additional part.

In reality, the quality of the material depends on the surface condition. The quality is not a part that can be confirmed with the naked eye, because it is extremely fine to use an electron microscope. The technique of the metal polishing for composite instruments not simple.

OSUNG's products are differentiated from the products of other companies by implementing the highest level of surface polishing technology. We, OSUNG are continuously making efforts to improve quality by using quality analysis system with state-of-art testing equipment.

◀ Figure. RA value(Arithmetical mean deviation of the profile) of Composite Instrument

Composite Instrument

2CSAT6

- Silicone Handle/Double-Ended



134°C Autoclavable

CSAT6

- Metal Handle/Double-Ended



Flowable

- Useful to operate flexible composite resin and glass ionomer / Used for precise reappearance of occlusal groove / Used to manage minute connection area

2CSCT15

- Silicone Handle/Double-Ended



134°C Autoclavable

CSCT15

- Metal Handle/Double-Ended



Composite Instruments

As the tool for the composite resin, the polished tip surface is smooth and lubrication-coated, so the composition resin is not sticky. The thin profile allows access to narrow interdental areas, especially for the aesthetic resin treatment.

Metal Handle

NEW

CSF1W

- Composite instrument
- 폭 1.5mm/폭 1.5mm



NEW

CSF2

- Composite instrument
- 폭 1.5mm/폭 1.5mm



NEW

CSF3

- Composite instrument
- 폭 1.5mm/폭 1.5mm



Practice



Composite Instruments

Plastic Handle

Autoclavable

BEST

3PFWDS2

• Woodson 2



3PFWDS3

• Woodson 3



BEST

3PFIW3

• Combination of middle sized paddle blade and small sized condenser tip



3PF43-47

• For anterior



3PF48-702

• For posterior



Composite Instruments

Composite Instrument_Metal Handle

BEST

PFWDS2

• Woodson 2



PFWDS3

• Woodson 3



BEST

PFIW3

• Combination of middle sized paddle blade and small sized condenser tip



PF43-47

• For anterior



PF48-702

• For posterior



Composite Instruments

Non-stick-instruments for fast and efficient placement and shaping of composite and easy cleanup. The tips are highly polished and titanium-coated.

Composite Instruments_Silicone Handle

Autoclavable

BEST

2CSCT1

- Combination of paddle and rounded condenser tip. For placement and contouring.



2CSCT6

- For proximal contouring



2CSCT7

- For placement and contouring



2CSCT8

- For condensing and contouring



2CSCT10

- For condensing and contouring



BEST

2CSCOM11

- Blade type for universal use
- Straight type



2CSCOM13

- Corn type for occlusal use



Composite Instruments

Non-stick-instruments for a fast and efficient placement and shaping of composite and an easy clean up. The tips are highly polished and titanium-coated.

Composite Instruments_Metal Handle

BEST

C SCT1

- Combination of paddle and rounded condenser tip. For placement and contouring.



C SCT6

- For proximal contouring



C SCT7

- For placement and contouring



C SCT8

- For condensing and contouring



C SCT10

- For condensing and contouring



BEST

C SCOM11

- Blade type for universal use
- Straight type



C SCOM13

- Corn type for occlusal use



Composite Instruments

Composite Instruments_Silicone Handle

Autoclavable



2CSCOMKIT

Composite Instruments Set
• Including sterilization case
• Size 180 x 99 x 18H(mm)

2CSCOM1

Composite Instrument
• For placement




2CSCOM2

Composite Instrument
• For shaping and carving



2CSCOM3

Composite Instrument
• For margin trimming



2CSCOM4

Composite Instrument
• For margin arrangement



2CSCOM5

Composite Instrument
• For margin trimming



2CSCOM6

Composite Instrument
• For occlusal shaping



Composite Instruments

Composite Instruments_Metal Handle



CSCOMKIT

Composite Instruments Set
• Including sterilization case
• Size 180 x 99 x 18H(mm)


CSCOM1

Composite Instrument
• For placement




CSCOM2

Composite Instrument
• For shaping and carving



CSCOM3

Composite Instrument
• For margin trimming



CSCOM4

Composite Instrument
• For margin arrangement



CSCOM5

Composite Instrument
• For margin trimming



CSCOM6

Composite Instrument
• For occlusal shaping



Composite Instrument Kit

Composite Instrument Kit

Autoclavable

3CSK01

Composite Instrument Kit Part 1
• Composite Resin Placement

- 3CSCOM11**
- Ideal for placement.
 - Thin tip provides maximum comfort with accurate control.

- 3CSCT8**
- Excellent for condensing restorative prostheses.

- 3CSCT1**
- Useful for placement and condensing.

3CSK02

Composite Instrument Kit Part 2
• Contouring Instrument

- 3CSCT6**
- Thin tip allows the easier restorations in narrow proximal surface.

- 3CSCOM13**
- Make it easy to create the ideal occlusal anatomy.

- 3CSCT15**
- Optimized for flowable resin handling. Useful for reproducing the proper anatomy of the finer points of the tooth such as fit and fissure.




Amalgam Filling

Treatment to filling the mixed amalgam in the cavity after removing the carious dentin.


Amalgam Carrier

To place the prepared amalgam to the cavity preparation and properly condense it.



Arrangement

01. Anesthesia Syringe	SAF1	P.073
02. Excavator	EXC18	P.161
03. Placement	PICH	P.164
04. Amalgamwell	PLGWL1	P.165
05. Carrier	PM1520	P.165
06. Plugger(=Condenser)	PLG1-2	P.166
07. Burnisher	BB27-29	P.168
08. Carver(Discoid-Cleoid)	CVCD89-92	P.167
09. Carver(Hollenback)	CV3	P.167
10. Burnisher	BB27-29	P.168



Process



Practice

01. Local anesthesia

- ❖ **Used**
Local anesthesia syringe, Harpoon type syringe provides stable aspiration during nerve block anesthesia.
- ❖ **Character**
Harpoon is designed to hold the rubber plunger of the cartridge and thumb ring is designed to make negative pressure for aspirating.

Anesthesia Syringe_SAF1

How to use


1. Choose a local anesthetic ampoule based on the patient's condition and the dentist's decision.
2. Check the validity period, whether the ampoule is cracked and the integrity of the rubber packing.
3. Make sure that the thumb ring of the syringe and the screw hub are tight.
4. Hold the syringe with one hand and pull the plunger back to insert a local anesthetic ampoule.
5. Remove the short side protection cap on the needle and secure it by screwing it onto the screw on the syringe hub. (At this time, the cap of the needle is not removed.)
6. Press the plunger that is pulled back so that the needle passes through the rubber septum. (Be careful not to bend the tip of the needle.)
7. Make sure that no air bubbles are generated.
8. Examine the treatment site.



Excavator_EXC18

How to use

Remove the carious lesion along the outer wall of the cavity using spoon shaped working end.



02. Cavity preparation

- ❖ **Used**
Suitable for removal of small cavities, carious dentin and temporary sealing materials after cavity preparation.
- ❖ **Character**
Have a spoon-shaped cutting surface. Various size of Small, Medium and etc.

03. Pulp protection

❖ Used

For mixing and applying a base and liner (calcium hydroxide, glass ionomer and etc) at the cavity.

❖ Character

Double ended. Ball-ended tips with different diameters. Each length of the two shanks is different.

Placement_PICH

How to use

1. Hold it with a pen grasp and mix the ingredients.
2. Coat the ball end with a small amount and apply it in the cavity.



Weight the same amount of base and catalyst on the mixing paper and start mixing.



Apply the base and liner such as calcium hydroxide, glass ionomer and etc. at the cavity.



It can also be used to fill materials in narrow grooves such as premolars.

04. Trituration and mulling of amalgam

❖ Used

Container that holds triturated amalgam prior to its being delivered to the cavity preparation.

❖ Character

Concave bowl-shaped well. It is more convenient to use than rubber sheet and is safe from the danger of mercury.

Well_PLGWL1

How to use

1. Hold the well with a left hand and hold the carrier with the other hand.
2. Fill the Amalgam carrier with the amalgam along the inside wall of the well.



For placing of triturated amalgam before it transferred to the carrier.



It can also be used to mix bone during implant surgery.

05. Amalgam placement

❖ Used

To place the prepared amalgam to the cavity preparation and properly condense it.

❖ Character

Consists of cylindrical pellets of different sizes.

Carrier_PM1520

How to use

1. Select the appropriate carrier according to the size of the cavity
2. Fill the freshly mixed amalgam in the cylindrical pellet and apply it while pressing the lever when placing into the prepared cavity to fill it.
3. Keep in mind that once amalgam has been triturated it immediately begins to harden. Use amalgam in the pellet immediately to prevent hardening of amalgam.



Hold the carrier with a Palm Grasp and place your index finger between the levers to fill the carrier with the mixed amalgam along the inside wall of well.



Hold the carrier with a Palm Grasp and place your index finger on the lever and press it.

06. Amalgam condensation

❖ Used

Also known as amalgam condenser, amalgam plugger compacts and condenses amalgam into the cavity preparation.

❖ Character

The shape of the tip is variable. Round, flat and diamond working end that can be serrated or plain. Diamond shaped working end for packing amalgam into larger portions of preparations.

Plugger(=Condenser)_PLG1-2

How to use

1. Hold the plugger with its tip held 90 degrees to the interface of the tooth and start the compaction at the center.
2. Use the smaller face plugger end and apply lateral pressure for condensing the corner of the cavity. (sidewall)
3. Amalgam should be condensed into the cavity with a minimal amount at several times.



In case of packing amalgam into wide portions of preparation, start condensing from the center to the sidewall pushing with large faced plugger end.



Smaller face plugger end is designed to compact amalgam with greater force.

07. Pre-carve burnishing

❖ Used

To smooth amalgam after condensing, used to create occlusal anatomy.

❖ Character

Double ended. Egg-ball and ball type.

Burnisher_BB27-29

How to use

Start stroking from the center to the cusp pushing with constant pressure.



Using Ball type, gently stroke with minimal force.



Using egg-ball type. Start burnishing with forming central groove.

08, 09. Carving

❖ Used

To carve occlusal anatomy (pits and fissures) into amalgam restorations.

❖ Character

Discoid is disk shaped, Cleoid is pointed, sharp. Designed for removing excess amalgam from the occlusal surface.

Carver_CVCD89-92
(Discoid-Cleoid)

How to use

To reduce the removal of large amounts of amalgam, place the blade on the adjacent tooth and pull it from the distal side to the mesial side.



Use discoid-cleoid type to remove excess amalgam from the occlusal surface.

❖ Used

To contour and carve occlusal and interproximal anatomy in amalgam restorations.

❖ Character

Sharp stiff metal blade. The angles of working ends are different.

Carver_CV3
(Hollenback)

How to use

Start carving by place the lateral side of blade in the inclination of cusp. And continue to carve the surface along the margin.



Use a hollenback type for carving central groove.

10. Post-carve burnishing

❖ Used

To smooth amalgam after carving, and burnish to obtain adequate adaptation.

❖ Character

Double ended. Egg-ball and ball type.

Burnisher_BB27-29

How to use

Start stroking from the center to the cusp pushing with constant pressure.




Using the burnisher, produce polished flat surface by stroking the amalgam surface.

Resin Filling

Treatment of filling the mixed amalgam in the cavity after removing the carious dentin.


Composite Instrument

Used to remove excess composite or resin cement from the interproximal surface of a tooth. A polished and lubricate coated surface tip prevents composite resin from sticking to the instrument.



Arrangement

01. Excavator	EXC18	P.173
02. (echant, 3way syringe)		
03. (bonding, 3way syringe, Light curing unit)		
04. Carrier Placement	CSCT7	P.187
05. Condenser	CSCT8	P.187
06. Carver	CSCOM2	P.189
07. Occlusal shaping	CSCOM13	P.187
08. (Light curing unit)		
09. Paper holder	PHNS	P.287



Process

EXC18 ▶
PICH ▶



01. Cavity preparation

02. Acid etching



03. Applying the composite resin primer

04, 05. Composite resin filling

CSCT7 ▶
CSCT8 ▶
CSCOM2 ▶
CSCOM13 ▶



06, 07. Composite resin filling

08. Light curing

PHNS ▶



09. Shaping and occlusal adjustment

Practice

01. Cavity preparation

✦ **Used**
Used to remove soft carious decay and temporary sealing materials after cavity preparation.

✦ **Character**
Spoon-shaped cutting surface. Various size of Small, Medium and etc.

02. Cavity preparation

03. Applying the composite resin primer

04. Composite resin filling

✦ **Used**
Carry composite resin to the cavity preparation.

✦ **Character**
Paddle shaped working end.

05. Composite resin filling

✦ **Used**
For compacting the dental restorative materials into a prepared tooth cavities (as temporary sealing materials).

✦ **Character**
Rounded working end with different diameter.

Excavator_EXC18

How to use

Remove the carious lesion along the outer wall of the cavity using spoon shaped working end.



Use small size excavator to remove the carious dentin.



Using spoon-shaped excavator, remove the carious decay.



EXC65-66 blade type is suitable for removing carious decay.

(Etching, 3Way Syringe)

(Bonding, 3Way Syringe, Light curing unit)

Carrier_CSCT7

How to use

Apply the appropriate composite resin into the cavity.



Take the composite resin out using the paddle end as much as desired.



Place the composite resin on the occlusal surface.



Used for contouring buccal, lingual surface.

Condenser_CSCT8

How to use

1. Hold the plugger with its tip held 90 degrees to the interface of the tooth and start the compaction at the center.
2. Use the smaller face plugger end and apply lateral pressure for condensing the corner of the cavity (sidewall).
3. Composite resin should be condensed into the cavity with a minimal amount at several times.



Compact the composite resin into the narrow cavities.



Compact the composite resin into the wide cavities.

06. Composite resin filling**Used**

For carving or contouring amalgam restorations to obtain optimal occlusion.

Character

Curved paddle shaped and Acorn shaped working end.

07. Composite resin filling**Used**

Used for shaping occlusal surface during the posterior composite resin restoration.

Character

Triangular ridge shaped.

08. Light curing**09. Shaping and occlusal adjustment****Used**

It is used for handling articulating paper during the occlusal adjustment

Character

Serrated Jaw, Tweezer type.

Carver_CSCOM2**How to use**

1. Shaping the buccal surface using paddle shaped tip.
2. Shaping the occlusal surface using acorn shaped tip.



Carve the composite resin by pushing it on the buccal surface.



Carve the composite resin by pushing it on the occlusal surface.

Occlusal Shaping_CSCOM13**How to use**

Reproduce the occlusal surface of the posterior.



Reproduce the occlusal surface of the posterior.

(Ligth curing unit)**Paper Holder_PHNS****How to use**

1. Fix the articulating paper to the end of the paper holder about 5mm longer from the tip of the paper holder.
2. Place the paper holder on the buccal and check the occlusion.



The serrated shape is designed to fix the articulating paper.



Fix the articulating paper to the end of the paper holder about 5mm longer from the tip of the paper holder.



Place the paper holder on the buccal and check the occlusion.

Products for Dentistry

OSUNG Catalogue 2022/2023

Endodontic