



WARNING: Follow all instructions carefully and use only as directed. Always consult with your healthcare professional before treating any serious injury or if pain persists beyond treatment with the Kanjo Seat Cushion.

WARNING: Risks are associated with any exercise program. Please consult with your healthcare professional before performing this or any physical activity.

Please visit
www.gokanjo.com/warranties for
information on our Limited Express
Product Warranty.

Questions? Concerns?

Please contact us at support@gokanjo.com
www.gokanjo.com

Kanjo Memory Foam Acupressure **SEAT CUSHION**


Welcome to Kanjo, and thank you for your purchase. Kanjo reminds you to make time for you. Relieve pain. Destress. Relax and rejuvenate.

Kanjo Acupressure products have stimulators that target acupressure points to reduce pain throughout the body, much like ancient Eastern methodologies used in acupuncture, but non-invasive and easy to use at home and on the go.

You are on your way to experiencing the benefits of your **Kanjo Memory Foam Acupressure Seat Cushion**, including relief of spinal, tailbone, seat and hip pain caused by long hours of sitting.

KANJO

Relief for Real Life™



A few instructions:

- Place the Kanjo Acupressure Pain Relief Seat Cushion on a chair with the opening toward the back of the chair.
- When you sit on the Cushion, your tailbone will be in the open space, preventing compression of the tailbone and spine.

A few notes:

If it's your first time using the Kanjo Seat Cushion, we recommend a 15 minute session to start. First-time users may feel strange or even uncomfortable on the plastic stimulators. Any discomfort will subside as your body responds to the stimulators' effects and pain relief begins, usually after a few minutes of contact..

How do I clean this?:

To clean fabric cover, remove it from the memory foam core, spot clean with water or a mild detergent, and let dry completely before use.

Cushion on chair in position



Correct seated position

