



WARNING: This product is not intended for sleeping purposes or for use while driving.

WARNING: Follow all instructions carefully and use only as directed. Always consult with your healthcare professional before treating any serious injury or if pain persists beyond treatment with the Kanjo Neck Pain Relief Support Cradle.

WARNING: Risks are associated with any exercise program. Please consult with your healthcare professional before performing this or any physical activity.

Please visit www.gokanjo.com/warranties for information on our Limited Express Product Warranty.

Questions? Concerns?

Please contact us at support@gokanjo.com
www.gokanjo.com

Kanjo Neck Pain Relief SUPPORT CRADLE

Welcome to Kanjo, and thank you for your purchase. Kanjo reminds you to make time for you. Relieve pain. Destress. Relax and rejuvenate.

Kanjo wellness products reduce pain throughout the body and are easy to use at home and on the go.

You are on your way to experiencing the benefits of your **Kanjo Neck Pain Relief Support Cradle**, including neck and shoulder pain relief, and neck traction and spinal alignment.

KANJO

Relief for Real Life™



A few instructions:

- Lie down on the floor or on a firm bed or couch, and place the Kanjo Neck Pain Relief Support Cradle beneath your neck and head. Position it so that the flexible curve cradles your neck.
- Adjust the position if necessary to maximize comfort.
- Can also be used in an upright, seated position.

A few notes:

If it's your first time using the Neck Cradle, we recommend a 15 minute session to start, followed by longer sessions as you become accustomed to the Cradle.

How do I clean this?:

Wipe plastic clean with water and/or a mild detergent.

