



AIR P03M Pulse Oximeter

1 Device Installation

1. Download the free iHealth MyVitals app on the App Store or Google Play depending on the model of your smartphone then create your iHealth account. Download the free iHealth MyVitals app from the App Store.
2. Charge the pulse oximeter battery. While your device is charging, the battery icon flashes. The battery is fully charged when the icon is steady.
3. Start the iHealth MyVitals app and sign into your account.
4. Go to MENU > Set up my devices > Select a new device > Pulse oximeter and follow the instructions on the screen.
 - Make sure that Bluetooth is activated on your smartphone.
 - Press the pulse oximeter button to start the oximeter and initialise the connection. The Bluetooth icon on the pulse oximeter appears steady once the connection is established.
 - The pulse oximeter's date and time are automatically set when connected to your smartphone.

2 Connecting to SNUG

1. Download Snug from the Google Play or Apple App Store and create your Snug account
2. Access "My Devices" through the app menu in the top left-hand corner
3. Select "Add a data source", choose "Human API", then press "Okay"
4. Choose "iHealth", add your iHealth account email and password (created in Step 1), then "Sign In"
5. Once Human API has processed, you can close this window manually to return to Snug
6. "iHealth" will appear in your list of Data Sources
7. Allow up to 30 minutes for data to flow into Snug from iHealth
8. For any connectivity issues contact support@snughealth.com.au

