



SMART BG5 Gluco-Monitoring System

1 Device Installation

1. Charge your glucometer's battery with the supplied USB cable (about 2-4 hrs).
2. Install the iHealth Gluco-Smart app on your smartphone.
3. Create your iHealth account and follow the instructions on the screen.
4. Connect your glucometer to your smartphone:
5. Select the Gluco + glucometer

2 Connecting to SNUG

1. Download Snug from the Google Play or Apple App Store and create your Snug account
2. Access "My Devices" through the app menu in the top left-hand corner
3. Select "Add a data source", choose "Human API", then press "Okay"
4. Choose "iHealth", add your iHealth account email and password (created in Step 1), then "Sign In"
5. Once Human API has processed, you can close this window manually to return to Snug
6. "iHealth" will appear in your list of Data Sources
7. Allow up to 30 minutes for data to flow into Snug from iHealth
8. For any connectivity issues contact support@snughealth.com.au

