
no touch forehead thermometer

INSTRUCTIONS FOR USE

(please read and keep for future reference)

tommeetippee.com

WARNINGS AND PRECAUTIONS

- This thermometer is not intended as a substitute for consultation with your doctor.
- Keep out of reach of children aged under 12 years.
- Do not allow children to take their own temperature unattended.
- To avoid risk of choking from swallowing small parts or batteries please keep the device out of the reach of children and pets.
- The device should not be submerged into water or other liquids (not waterproof).
- The scanner lens is the most delicate part of the thermometer. Use with care when cleaning the probe lens to avoid damage. See care and cleaning section.
- Do not use the thermometer if there are signs of damage to the scanner or the thermometer itself. If damaged do not attempt to repair the product.
- No modification of this equipment is permitted.
- Intended for household use only.
- Forehead temperature is displayed in oral mode. This mode converts the forehead temperature to display its "oral-equivalent" value.
- The "Clinical Bias" is -1,4 ~ -1,7°C
- The "Limits of Agreement" is 0,98
- The "Repeatability" is 0,20°C
- Holding the thermometer too long may cause a higher ambient temperature reading of the probe. This could make the body temperature measurement lower than usual.
- Intended operator: At least 11 years old (5 years intensive reading experience), no maximum.
- Intended use: The infrared forehead thermometer can measure the infrared heat generated by central forehead and surrounding tissues to reflect patient's body temperature accurately.
- Accuracy: Forehead mode: $\pm 0.2\text{ }^\circ\text{C}$ (0,4 °F) within 35–42°C (95–107,6°F) (Ambient Temp: 15–35°C), $\pm 0.3\text{ }^\circ\text{C}$ (0,5 °F) for other range.

FUNCTIONS

Forehead temperature

See "how to use" section to learn how to measure body temperature.

Surface temperature

See "how to use" section to learn how to measure the temperature of objects (for example baby milk).

Fever indication

If the thermometer detects a temperature equal to or greater than 37,5°C (or 99,5°F) under forehead mode, three short beeps will sound followed by one long beep to warn for potential fever.

Memory locations

The thermometer memory holds the latest 25 readings. Once you have powered on, press the **STAND-BY (ON/MEM)** button to see temperature records with the icon.

FUNCTIONS

How to switch between °C / °F

In **power off** mode, press and hold the **START** button then press the **STAND-BY (ON/MEM)** button for 3 seconds. °C will switch to °F. You can use the same process to switch back from °F to °C. **NOTE:** when you switch between °C and °F the memory will be cleared.

Mute mode

The device is automatically set with buzzer on but you can set it **ON** or **OFF** using the **MUTE** mode. In **power on** mode, press and hold the **STAND-BY (ON/MEM)** button for 3 seconds. This icon will flash on the screen. Release the **STAND-BY (ON/MEM)** button to set **MUTE**.

You can use the same process to turn the **MUTE** function off.

NOTE: If you hold the **STAND-BY (ON/MEM)** button down for 5 seconds after this mute icon starts to flash the device will **power off** WITHOUT setting Mute.

HOW TO USE

Power on

Press and release the **STAND-BY (ON/MEM)** button once. The backlight will come on and the start-up sequence begins.

Measuring body temperature on the forehead (forehead mode is the default mode). When you **power on** the unit you will see this icon on the screen and you will hear 2 beeps.

Position the thermometer on or up to 4cm away from the centre of the forehead just between the eyebrows and hold steady.

Press the **START** button and hold the thermometer in place until you hear a final beep. Your temperature is ready to read. After each forehead measurement, wait for this icon to stop flashing before taking the next measurement.

Measuring surface temperature

After **power on**, press and hold the **STAND-BY (ON/MEM)** button, and press the **START** button at the same time for "Infrared thermometer" mode. You will see icon on your LCD display. In this mode, you can get the target surface temperature.

When you press the **START** button, you will get the real time temperature immediately. If you press and hold the **START** button, the reading of measurement will be continuously updated.

Applications include temperature measurements for water, milk, cloth, skin or other objects.

NOTE: This mode shows the actual and unadjusted surface temperature which is different from the body temperature.

Power Off

The device will automatically power off after 1 minute of being idle to extend the battery life. To manually power off the device, press the **STAND-BY (ON/MEM)** button.

TIPS FOR TAKING AN ACCURATE READING

- Patient must be inside for 30 minutes before taking a measurement.
- The patient and device must be in the same stable ambient (room) temperature for 15 minutes before operating.
- Always make sure the scanner lens is clean and undamaged.
- If the eyebrow area is covered with hair, sweat, make-up or dirt, please clean the area beforehand to improve reading accuracy.
- Patients should not drink, eat or be physically active before / while taking the measurement. Remove hats and wait 10 minutes before taking a measurement.
- Always hold the thermometer and forehead steady when taking a reading. Do not move the thermometer until you hear the final beep.
- Don't take a measurement during or immediately after nursing a baby.

HOW TO INTERPRET YOUR READING

Understanding the meaning of your child's temperature.

READING	MEANING
>35,7° - 37,4° (>96,3° - 99,4°F)	Acceptable
>37,4° - 38,5°C (>99,4° - 101,3°F)	Elevated
>38,5° - 42,2°C (>101,3° - 108,0°F)	Possible High Fever

REMEMBER if you are concerned about your child's health you should see a doctor. Even in the absence of fever those who exhibit a normal temperature may still need to receive medical attention.

CARE AND CLEANING

- Make sure the probe is clean to ensure an accurate reading.
- Use a 70% alcohol swab or cotton wool moistened with 70% alcohol to gently wipe the lens clean.
- Allow the lens to fully dry for at least 1 minute.
- Never insert a sharp object into the scanner area or any other open surface on the thermometer.
- Keep the thermometer dry and away from any liquids and direct sunlight.
- It should be stored at room temperature between -20 ~ +50°C, RH < 85%.

TROUBLESHOOTING

	Error 5-9, the system is not functioning properly.
	Measurement before device stabilization.
	The ambient temperature is not within the range between 10°C and 40°C (50°F-104°F).
Solution: Allow the thermometer to rest in a room for at least 15 minutes at room temperature: 10°C and 40°C (50°F-104°F).	
	(1) In forehead mode: Temperature taken is higher than +42.2°C (108°F). (2) In surface mode: Temperature taken is higher than +80°C (176°F).
Solution: Please select the target within specifications. If a malfunction still exists, please contact our UK Careline.	
	(1) In forehead mode: Temperature taken is lower than +34°C (93.2°F). (2) In surface mode: Temperature taken is lower than -22°C (-7.6°F).
Solution: Please select the target within specifications. If a malfunction still exists, please contact our UK Careline.	
	Device cannot be powered on to the ready stage.
Solution: Change with new batteries.	

REPLACING THE BATTERIES

- When the "low battery" icon indicates the batteries are low, the batteries should be replaced immediately with 2x AAA batteries.
- Open the battery cover: use the thumbs to push battery cover out (see fig. 1).
- Insert the new batteries (see fig. 2).
- Replace the battery cover (see fig. 3).

fig. 1

fig. 2

fig. 3

BATTERY WARNINGS

Batteries are to be inserted with the correct polarity. Ensure that the + and - ends of the batteries are positioned correctly. Requires 2 AAA batteries. **DO NOT USE RECHARGABLE BATTERIES.** The use of rechargeable batteries could cause your Thermometer to malfunction and will void the warranty. Non rechargeable batteries are not to be recharged. Different types of batteries or new and used batteries are not to be mixed. Only batteries of the same or equivalent type as recommended are to be used. Exhausted batteries are to be removed from the Thermometer. Batteries should always be removed or replaced by an adult. Keep batteries away from children. Remove batteries if the Thermometer will not be used for some time.

PRODUCT SPECIFICATIONS

- Temperature measurement range:
Forehead mode: 34–42.2°C (93.2–108°F).
Surface mode: -22–80°C (-7.6–176°F).
- Operating temperature range: 10–40°C (50–104°F), 15%–85% RH.
- Storage temperature range: It should be stored at room temperature between -20~+50°C, RH equal to or less than 85%.
Transportation temperature shall be less than 70°C, RH equal to or less than 95%.
- Atmospheric pressure: 800–1013 hPa.
- Comply with ASTM E1965-98, EN ISO 80601-2-56, IEC/EN60601-1-2(EMC), IEC/EN60601-1(Safety) standards, ISO10993, RoHS.
- Accuracy:
Forehead mode: $\pm 0.2\text{ }^\circ\text{C}$ (0,4 °F) within 35–42°C (95–107,6°F) (Ambient Temp: 15–35°C), $\pm 0.3\text{ }^\circ\text{C}$ (0,5 °F) for other range.
Surface mode: $\pm 0.3\text{ }^\circ\text{C}$ (0,5°F) within 22–42,2°C (71.6–108°F), others $\pm 4\%$ or $\pm 2^\circ\text{C}$ (4°F) whichever is greater.
- Fever indicator, memory and °C/°F switch function.

- Battery: AAA x 2 pcs.
- Battery life: around 3,000 continuous readings.
- This thermometer converts the forehead temperature to display its "oral equivalent." (according to the result of the clinical evaluation to get the offset value).
- Blue LED Backlight.
Blue LED Backlight will be automatically turned on after measurement, and automatically turned off after 5 seconds.
- Enclosure Rating: **IP22**.
- Dimensions: 158.0 x 48.0 x 40.2 mm.
- Weight: 100 grams including batteries.

Manufacture Date: as the serial number (please open the battery cover, it is shown on the inside of the device).

Ex.SN:E912A000001, the first "E" is External, the second number "9" is the manufacture year 2009, the third and the fourth number "12" is the manufacture month, the others is the serial number.

Note: The thermometer is calibrated at the time of manufacture. If at any time you question the accuracy of temperature measurements, please contact the dealers or nearest service address.

SYMBOL DESCRIPTIONS

- Stand-by
- The CE mark and Notified Body Registration Numbers, the requirement of Annex II from Medical Device Directive 93/42/EEC are met.
- CAUTION.
- BF type applied part.
- IP22** Classification for water ingress and particulate matter.
Indicates this device is subject to the Waste Electrical and Electronic Equipment Directive in the European Union. To protect the environment, dispose of useless device at appropriate collection sites according to national or local regulations.
- Please read the instructions for use.
- Battery recycling.

EC REP Authorized representative in the European community
Medical Technology PromedT Consulting GmbH
Add: Altenhofstrasse 80, D-66386 St. Ingbert, Germany.

Paper recycling.

Manufactured by:
Radiant Innovation Inc.
Http://www.radiantek.com.tw
Add: 1F, No.3, Industrial East 9th Road,
Science-Based Industrial Park, HsinChu, Taiwan 300.

Distributed by:
Jackel International Limited trading as Mayborn Group
Northumberland Business Park West
Cramlington
Northumberland
NE23 7RH, UK

Any questions? Visit tommeetippee.com

Guidance and manufacturer's declaration – electromagnetic emissions

The THD2FE is intended for use in the electromagnetic environment specified below. The customer or the user of the THD2FE should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The THD2FE uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The THD2FE is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration – electromagnetic immunity

The THD2FE is intended for use in the electromagnetic environment specified below. The customer or the user of the THD2FE should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the THD2FE, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2 \cdot \sqrt{P}$ $d = 2,3 \cdot \sqrt{P}$ 80 MHz to 800 MHz $d = 2,3 \cdot \sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range b. Interference may occur in the vicinity of equipment marked with the symbol:
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the THD2FE is used exceeds the applicable RF compliance level above, the THD2FE should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the THD2FE.
b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer's declaration – electromagnetic immunity

The THD2FE is intended for use in the electromagnetic environment specified below. The customer or the user of the THD2FE should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV contact 8 kV air	6 kV contact 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV line(s) to line(s) 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.

Guidance and manufacturer's declaration – electromagnetic immunity

The THD2FE is intended for use in the electromagnetic environment specified below. The customer or the user of the THD2FE should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the THD2FE requires continued operation during power mains interruptions, it is recommended that the THD2FE be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM.

The THD2FE is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the THD2FE can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the THD2FE as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1,2 \cdot \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \cdot \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \cdot \sqrt{P}$
0,01	N/A	0,12	0,23
0,1	N/A	0,38	0,73
1	N/A	1,2	2,3
10	N/A	3,8	7,3
100	N/A	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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