**drive**

**DeVilbiss**  
HEALTHCARE

1025 Series  
Serie de 1025  
Série 1025

CE 0044

## **EN** DeVilbiss® 10-Liter Oxygen Concentrator Instruction Guide

**WARNING**—Read instruction guide before operating this equipment.

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 **NE PAS FUMER**

[Click here to view detailed product descriptions and prices.](#) Call 800-397-5899 to order today!

**TABLE OF CONTENTS**

Symbol Definitions . . . . .	EN - 3
Important Safeguards . . . . .	EN - 3
Introduction . . . . .	EN - 4
Intended Use . . . . .	EN - 4
Indications For Use . . . . .	EN - 4
Why Your Physician Prescribed Supplemental Oxygen . . . . .	EN - 4
How Your Concentrator Works . . . . .	EN - 4
Important Parts of Your Concentrator . . . . .	EN - 5
Setting Up Your Concentrator . . . . .	EN - 6
Operating Your Concentrator . . . . .	EN - 6
DeVilbiss OSD® Operation . . . . .	EN - 7
Reserve Oxygen System . . . . .	EN - 7
Caring for Your Concentrator . . . . .	EN - 7
Troubleshooting . . . . .	EN - 8
Overview of Alarms . . . . .	EN - 9
Warranty . . . . .	EN - 9
Return and Disposal . . . . .	EN - 9
Provider's Notes . . . . .	EN - 9
Specifications . . . . .	EN - 10
Guidance and Manufacturer's Declaration . . . . .	EN - 11

**WARNING**

Under certain circumstances, oxygen therapy can be hazardous. Seek medical advice before using an oxygen concentrator.

**Physician Information**

Physician Name: \_\_\_\_\_  
 Telephone: \_\_\_\_\_  
 Address: \_\_\_\_\_

**Prescription Information**

Name: \_\_\_\_\_  
 Oxygen liters per minute  
 at rest: \_\_\_\_\_ during activity: \_\_\_\_\_ other: \_\_\_\_\_  
 Oxygen use per day  
 Hours: \_\_\_\_\_ Minutes: \_\_\_\_\_  
 Comments: \_\_\_\_\_


**DeVilbiss 10-Liter Oxygen Concentrator w/OSD Serial Number:** \_\_\_\_\_

**DeVilbiss Equipment Provider Information**

Set-Up Person: \_\_\_\_\_



This instruction guide was reviewed with me, and I have been instructed on the safe use and care of the DeVilbiss Oxygen Concentrator.

Signature: \_\_\_\_\_ Date: \_\_\_\_\_

DeVilbiss 10-Liter Series

## SYMBOL DEFINITIONS

	It is mandatory to read and understand the operating instructions prior to use.	<b>REF</b>	Catalog Number		Danger - No smoking near patient or device.
	Electric Shock Hazard. Cabinet to be removed by authorized personnel only.	<b>SN</b>	Serial Number		Do NOT use Oil, Grease or Lubricants
	Off On	$\checkmark O_2$	Normal Oxygen		Do not use near heat or open flames
	Reset	$\downarrow O_2$	Low Oxygen		Manufacturer
	Alternating Current		Service Required	<b>EC REP</b>	European Representative
	Type B applied part		TUV Rheinland C-US approval mark	<b>CE</b> 0044	European Rep CE mark
	Double Insulated		TUV Rheinland Certified approval mark	 <p>≤ 10 LPM</p> <p>Maximum recommended flow rate: 10 LPM</p>	
	Atmospheric Pressure Range 840 to 1010 hPa (Approximate sea level to 5000 ft)	<b>R<sub>X</sub> ONLY</b>	Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.		
	Operating Temperature Range +5 to +35°C (+41 to +95°F)	<b>LOT</b>	LOT Number		
	Ingress Protection - Protected against finger access to hazardous parts; protected against vertically falling water drops.				
	This device contains electrical and/or electronic equipment that must be recycled per EU Directive 2012/19/EU- Waste Electrical and Electronic Equipment (WEEE)				

## IMPORTANT SAFEGUARDS

Read this entire guide before using your DeVilbiss concentrator. Important safeguards are indicated throughout this guide. Pay special attention to all safety information.

### READ ALL INSTRUCTIONS BEFORE USING.

#### WARNING:

*A warning indicates the possibility of injury to the user or the operator.*

- Oxygen causes rapid burning. Do not smoke while your oxygen concentrator is operating, or when you are near a person utilizing oxygen therapy.
- Oxygen makes it easier for a fire to start and spread. Do not leave the nasal cannula or mask on bed coverings or chair cushions if the oxygen concentrator is turned on but not in use. The oxygen will make the materials flammable. Turn the oxygen concentrator off when not in use to prevent oxygen enrichment.
- Keep the oxygen concentrator and cannula at least 2 m (6.5 feet) from hot, sparking objects or naked sources of flame.
- Open flames during oxygen therapy are dangerous and are likely to result in fire or death. Do not allow open flames within 2 m (6.5 feet) of the oxygen concentrator or any oxygen carrying accessories.
- Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- Smoking during oxygen therapy is dangerous and is likely to result in facial burns or death. Do not allow smoking within the same room where the oxygen concentrator or any oxygen carrying accessories are located.
  - If you intend to smoke, you must always turn the oxygen concentrator off, remove the cannula and leave the room where either the cannula or mask or the oxygen concentrator is located. If unable to leave the room, you must wait 10 minutes after you have turned off the oxygen concentrator before smoking.
- In order to prevent a fire propagating from the patient through the cannula towards the unit, a means of protection should be located as close to the patient as practicable. Please contact your provider for this means of protection.
- DeVilbiss 10-liter oxygen concentrators are equipped with a fire mitigating outlet fitting that prevents propagation of fire into the unit.
- To avoid electric shock, do not remove the concentrator cabinet. The cabinet should only be removed by a qualified DeVilbiss technician. Do not apply liquid directly to the cabinet or utilize any petroleum-based solvents or cleaning agents.
- Improper use of the power cord and plugs can cause a burn, fire or other electric shock hazards. Do not use the unit if the power cord is damaged.
- The surface temperature of the exhaust vents on the bottom of the unit may exceed 105.8°F (41 °C) under certain conditions.
- When device is used under extreme operating conditions, the temperature near the exhaust vents on the bottom of the unit may reach 138.2°F (59°C). Keep body parts a minimum of 30" (76.2 cm) away from this area.
- Use of this device at an altitude above 5000 feet (1524 meters) or above a temperature of 95°F (35°C) or greater than 93% relative humidity may affect the flow rate and the percentage of oxygen and consequently the quality of the therapy. Refer to specifications for details regarding parameters tested.
- It is recommended that the homecare provider lock the flow control knob to prevent inadvertent adjustment. A flow setting other than prescribed may affect the patient therapy.
- To ensure you receive the therapeutic amount of oxygen delivery according to your medical condition, the Oxygen Concentrator must:
  - be used only after one or more settings have been individually determined or prescribed for you at your specific activity levels.
  - be used with the specific combination of parts and accessories that are in line with the specification of the concentrator manufacturer and that were used while your settings were determined.
- If you feel discomfort or are experiencing a medical emergency while undergoing oxygen therapy, seek medical assistance immediately to avoid harm.
- Geriatric, pediatric or any other patient unable to communicate discomfort can require additional monitoring and/or a distributed alarm system to convey the information about the discomfort and/or the medical urgency to the responsible caregiver to avoid harm.
- For your safety, the oxygen concentrator must be used according to the prescription determined by your physician.
- Your delivery settings of the oxygen concentrator should be periodically reassessed for the effectiveness of therapy.

**NO  
SMOKING**

- Under certain circumstances, oxygen therapy can be hazardous. Seek medical advice before using an oxygen concentrator.
- Before attempting any cleaning procedures, turn the unit “Off.”
- Do not service or clean this device while in use with a Patient.
- Use of harsh chemicals (including alcohol) is not recommended. If bactericidal cleaning is required, a non-alcohol based product should be used to avoid inadvertent damage.
- Use only water-based lotions or salves that are oxygen-compatible before and during oxygen therapy. Never use petroleum or oil-based lotions or salves to avoid the risk of fire and burns.
- Do not lubricate fittings, connections, tubing or other accessories of the oxygen concentrator to avoid the risk of fire and burns.
- Use only spare parts recommended by the manufacturer to ensure proper function and to avoid the risk of fire and burns.
- Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the Electromagnetic Compatibility [EMC] information provided in the accompanying documents.
- Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.
- Do not use a low-output flow meter with this concentrator.
- The Device is classified as IP21 which means it is protected against finger access to hazardous parts and protected against vertically falling water drops.
- The equipment or system should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.
- This device contains electrical and/or electronic equipment. Follow local governing ordinances and recycling plans regarding disposal of device components.

**CAUTION:**

*A Caution indicates the possibility of damage to the device.*

- It is very important to follow your oxygen prescription. Do not increase or decrease the flow of oxygen – consult your physician.
- To prevent product damage, do not attempt to operate the unit without the air filter, or while the filter is still damp.

# SAVE THESE INSTRUCTIONS.

## INTRODUCTION

---

This instruction guide will acquaint you with your DeVilbiss oxygen concentrator. Make sure that you read and understand this guide before operating your unit. Important safeguards are indicated throughout this guide. Pay special attention to all safety information. Contact your DeVilbiss equipment provider should you have any questions.

### Intended Use

The DeVilbiss 10 Liter Oxygen Concentrator intended use is to provide supplemental low flow oxygen therapy for patients suffering from COPD, cardiovascular disease, and lung disorders. The oxygen concentrator is used in home type environments, homes, nursing homes, patient care facilities, etc.

### Indications For Use

The DeVilbiss Oxygen Concentrator is intended for use as an oxygen concentrator to provide supplemental low flow oxygen therapy in the home, nursing homes, patient care facilities, etc.

**NOTE**– *The device is not intended to be life supporting or life sustaining.*

### Why Your Physician Prescribed Supplemental Oxygen

Today, many people suffer from heart, lung and other respiratory diseases. Many of these people can benefit from supplemental oxygen therapy. Your body requires a steady supply of oxygen to function properly. Your physician prescribed supplemental oxygen for you, because you are not getting enough oxygen from room air alone. Supplemental oxygen will increase the amount of oxygen that your body receives.

Supplemental oxygen is not addictive. Your physician prescribed a specific oxygen flow to improve symptoms such as headaches, drowsiness, confusion, fatigue or increased irritability. If these symptoms persist after you begin your supplemental oxygen program, consult your physician.

The oxygen delivery setting has to be determined for each patient individually with the configuration of the equipment to be used, including accessories.

The proper placement and positioning of the prongs of the nasal cannula in the nose is critical to the amount of oxygen delivered to the respiratory system of the patient.

Your Delivery settings of the oxygen concentrator should be periodically reassessed for the effectiveness of therapy.

### How Your DeVilbiss Oxygen Concentrator Works

Oxygen concentrators are the most reliable, efficient and convenient source of supplemental oxygen available today. The oxygen concentrator is electrically operated. The unit separates oxygen from room air which allows high-purity supplemental oxygen to be delivered to you through the oxygen outlet. Although the concentrator filters the oxygen in a room, it will not affect the normal amount of oxygen in your room.

## IMPORTANT PARTS OF YOUR CONCENTRATOR

Please take time to familiarize yourself with your DeVilbiss oxygen concentrator before operating.

**Front View (Figure A)**

1. Operating instructions
2. Power Switch  
| = ON  
O = OFF
3. Flow meter knob
4. Flow meter
5. Circuit breaker – resets the unit after electrical overload shutdown
6. Oxygen outlet – oxygen is dispersed through this port
7. Normal Oxygen (green) light (see page 7)
8. Low Oxygen (yellow) light (see page 7)
9. Red Service Required (red) light – when illuminated contact your DeVilbiss provider
10. Hour meter

**Back View (Figure B)**

11. Handgrip
12. Exhaust  
**WARNING –** When device is used under extreme operating conditions, the temperature near the exhaust vents on the bottom of the unit may reach 138.2°F (59°C). Keep body parts a minimum of 30” (76.2 cm) away from this area.
13. Power cord and/or IEC power connector.
14. Line cord strap
15. Air filter – prevents dirt, dust and lint from entering your unit.
16. Oxygen Outlet Port: Your concentrator is equipped with an auxiliary oxygen port that can be used to fill oxygen cylinders with an FDA-cleared cylinder filling device that is designed to use oxygen from a concentrator to fill a cylinder. The port is only for use with FDA-cleared filling devices with compatible oxygen input specifications. Refer to the cylinder filling device instruction guide for the oxygen input/output specifications, connection and operating instructions.

### Accessories

Transfiller Caddy .....	DeVilbiss 525DD-650
High Flow (6-15 LPM) Bubble Humidifier.....	Salter Labs 7900 or equivalent
Low Flow (up to 6 LPM) Bubble Humidifier .....	Salter Labs 7600 or equivalent
High Flow Nasal Cannula.....	Salter Labs 1600HF or equivalent

There are many types of humidifiers, oxygen tubing and cannulas/masks that can be used with this device. Certain humidifiers and accessories may impair the device's performance. A mask or any nasal cannula can be used with continuous flow delivery and may be sized according to your prescription as recommended by your homecare provider who should also give you advice on the proper usage, maintenance and cleaning.

### WARNING:

**When using the Transfiller Caddy with a Transfill device, always keep the system on a flat surface. Disassemble the system prior to moving.**

**NOTE–** The bubble humidifier should be supplied with a permanent fire stop device. If a bubble humidifier needs to be used without a permanent fire stop device, a secondary fire stop device must be used and placed as close to the humidifier as possible. Failing to do so could increase the risk of fire.

**NOTE–** The device should only be used with bubble humidifiers that are designed for use with flows up to 10 liters per minute and 20 psi pressure.

**NOTE–** A maximum of 50 feet (15 meters) of crush-proof oxygen tubing, plus 7 feet (2.1 meters) of cannula, plus a bubble humidifier is allowed between the concentrator and the patient.

**NOTE– The oxygen supply accessory (patient tubing) shall be equipped with a means that, in case of fire, stops the delivery of oxygen to the patient. This means of protection should be located as close to the patient as practicable.**

**NOTE–** Your healthcare provider should verify the compatibility of the oxygen concentrator and all of the parts used to connect to the patient before use.

## SETTING UP YOUR OXYGEN CONCENTRATOR

1. Position your unit near an electrical outlet in the room where you spend most of your time.

**NOTE**– Do not connect to an electrical outlet controlled by a wall switch.

### WARNING

Oxygen causes rapid burning. Do not smoke while your oxygen concentrator is operating, or when you are near a person utilizing oxygen therapy. Keep the oxygen concentrator and cannula at least 2 m (6.5 feet) from hot, sparking objects or naked sources of flame.

2. Position your unit on a flat surface at least 6 inches (16 cm) from walls, draperies or any other objects that might prevent the proper flow of air in and out of your oxygen concentrator. The oxygen concentrator should be located in a well-ventilated area to avoid pollutants or fumes.  
**NOTE**– To move the unit, firmly grasp the handle located on the top of the unit, rolling and/or lifting the unit over pathway obstacles.
3. Before operating your unit, always check to be sure the air filter (located on the back of your unit) is clean. Proper cleaning of this filter is discussed in the Caring For Your Concentrator section on page 7.
4. Attach the appropriate oxygen accessories to the oxygen outlet.

#### Oxygen Tubing Connection:

- a. Thread the supplied oxygen outlet connector onto the oxygen outlet.
- b. Attach the oxygen tubing directly to the connector (Figure 1).

#### Oxygen Tubing Connection With Humidification:

If your physician has prescribed an oxygen humidifier as part of your therapy, follow these steps (If using a prefill, go to step b.):

- a. Fill the humidifier bottle as per manufacturer's instructions.
- b. Thread the wing nut located on the top of the humidifier bottle to the oxygen outlet so that it is suspended (Figure 2). Make sure it is securely tightened.
- c. Attach the oxygen tubing directly to the humidifier bottle outlet fitting (Figure 3).

**NOTE**– Your physician has prescribed either a nasal cannula or face mask. In most cases, they are already attached to the oxygen tubing. If not, follow the manufacturer's instructions for attachment.

5. Remove the power cord completely from the line cord strap. Make sure the power switch is in the "Off" position and insert the plug into the wall outlet. The unit is double insulated to guard against electric shock.

**NOTE**– (only 120 volt units) The plug on the DeVilbiss oxygen concentrator has one blade wider than the other. To reduce the risk of electric shock, this plug is intended to fit in a wall outlet only one way. Do not attempt to defeat this safety feature. Your healthcare provider should verify the compatibility of the oxygen concentrator and all of the parts used to connect to the patient before use.

**NOTE**– To check your oxygen concentrator and accessories for proper operation; 1. Check the output flow by placing the end of the nasal cannula under the surface of a half-full cup of water and look for the bubbles. 2. Check the system for leaks by bending the nasal prongs over and squeeze tight to stop the flow of oxygen. Look at the flow meter to see that the indicator ball on the flow meter drops to zero. If the indicator ball does not drop to zero, check all connections for possible leaks. Parts to check for leaks are: tubing connections, humidifier bottle and other accessories like firebreaks. Repeat these steps until the flow meter ball drops to zero. Contact your provider or service supplier immediately if you encounter any problems.

### WARNING

Improper use of the power cord and plugs can cause a burn, fire or other electric shock hazards. Do not use the unit if the power cord is damaged.

## OPERATING YOUR DEVILBISS OXYGEN CONCENTRATOR

### WARNING

Oxygen causes rapid burning. Do not smoke while your oxygen concentrator is operating, or when you are near a person utilizing oxygen therapy. Keep the oxygen concentrator and cannula at least 2 m (6.5 feet) from hot, sparking objects or naked sources of flame.

The DeVilbiss 10-liter oxygen concentrators are equipped with a fire mitigating outlet fitting that prevents propagation of fire into the unit.

In order to prevent a fire propagating from the patient through the cannula towards the unit, a means of protection should be located as close to the patient as practicable. Please contact your provider for this means of protection.

1. Press the power switch to the "On" position. When the unit is turned "On," all three lights (Service Required, Low Oxygen and Normal Oxygen) on the front panel will illuminate briefly and an audible signal will briefly alarm confirming that the LEDs and audible signal are functioning properly. The unit will then operate in "start up" mode with the Low Oxygen light lit until a normal oxygen level is achieved, at which time the Normal Oxygen light will remain lit. The "start up" may take up to 15 minutes.

**NOTE**– DeVilbiss recommends for optimal service life that the DeVilbiss Oxygen Concentrator be operated for at least 30 minutes after it is powered on. Shorter periods of operation, operating in extreme temperature/humidity conditions or in the presence of contaminants, and/or handling and storage conditions outside those specified, may affect the long term reliable operation of the product.

### WARNING

Oxygen makes it easier for a fire to start and spread. Do not leave the nasal cannula or mask on bed coverings or chair cushions if the oxygen concentrator is turned on but not in use, the oxygen will make the materials flammable. Turn the oxygen concentrator off when not in use to prevent oxygen enrichment.

**NOTE**– If the audible signal alarms but the unit is not operating, there is no power to the unit. Refer to the Minor Troubleshooting chart on page 8 and contact your DeVilbiss provider if necessary.

**NOTE**– If an audible low-frequency vibration sound is detected, the unit is not operating properly. Refer to the Minor Troubleshooting chart on page 8, and contact your DeVilbiss provider if necessary.

2. Check the flow meter to make sure that the flow meter ball is centered on the line next to the prescribed number of your flow rate.

**NO  
SMOKING**

**CAUTION**– It is very important to follow your oxygen prescription. Do not increase or decrease the flow of oxygen – consult your physician.

**NOTE**– Your DeVilbiss provider may have preset the flow meter so that it cannot be adjusted.

**NOTE**– If the flow meter knob is turned clockwise, the flow increases (and eventually will shut off the oxygen flow). If the knob is turned counterclockwise, the flow decreases.

**NOTE**– For prescriptions of 10 LPM, be sure the ball is centered on the 10 liter line. The ball should not touch the red line. Setting the flow higher than 10 may cause the oxygen purity level to drop.

**NOTE**– The low-flow alarm may activate if the flow meter ball is set below 2 lpm. The unit will continue to run; however, the Service Required light will come on accompanied by an audible alarm. Adjust the flow meter to your prescribed flow. Do not use a low-output flow meter with this concentrator.

3. Your DeVilbiss concentrator is now ready for use, properly position the cannula with the nasal prongs facing upward, insert the prongs into nose. Wrap the cannula tubing over the ears and position in front of body (Figure 4). Allow 15 minutes for the oxygen concentrator to reach stated performance.

### DeVilbiss OSD® Operation

The OSD (Oxygen Sensing Device) is a device within your concentrator that monitors the oxygen produced by your unit.

The OSD lights on the top panel are defined as follows:

- Green Normal Oxygen light–acceptable oxygen level.
- Yellow Low Oxygen light– below an acceptable oxygen level.

If the oxygen purity falls below the acceptable level, the green Normal Oxygen light will shut off, and the yellow Low Oxygen light will illuminate. Switch to your reserve oxygen system. Refer to the Minor Troubleshooting section in this guide on page 8 and contact your DeVilbiss provider.

As an added safety feature, should the oxygen purity continue to drop, an intermittent audible signal will sound. Contact your DeVilbiss provider immediately. Do not attempt any other maintenance.

## RESERVE OXYGEN SYSTEM

As a precaution, your DeVilbiss provider may supply you with a reserve oxygen system. If your unit loses electrical power or fails to operate correctly, the Patient Alert System will sound to signal you to switch to your reserve oxygen system (if provided) and contact your DeVilbiss provider.

## CARING FOR YOUR DEVILBISS OXYGEN CONCENTRATOR

DeVilbiss recommends using only original DeVilbiss parts and filters in order to guarantee reliable operation of the product.

### WARNING

Do NOT use lubricants, oils or grease.

Before attempting any cleaning procedures, turn the unit “Off.”

### Cannula/Mask, Tubing and Humidifier Bottle

Clean and replace the cannula/mask, tubing and humidifier bottle according to the manufacturer’s instructions.

### Air Filter and Oxygen Outlet Connector

The air filter and connector should be cleaned at least once a week. To clean, follow these steps:

1. Remove the air filter, located on the back of the unit. Remove the oxygen outlet connector (if used).
2. Wash in a solution of warm water and dishwashing detergent (Figure 5).
3. Rinse thoroughly with warm tap water and towel dry. The filter should be completely dry before reinstalling.

**CAUTION**– To prevent product damage, do not attempt to operate the unit without the air filter or while the filter is still damp.

### Exterior Cabinet

Clean the concentrator exterior cabinet weekly by using a damp cloth and wiping dry.

### Cleaning

	Recommended cleaning interval	Number of cleaning cycles *	Compatible cleaning method
Outer Cabinet	7 days	260	Water, use only a damp cloth
Air Filter	7 days	104	Mild dish soap and warm water
Oxygen Outlet Connector	7 days	104	Mild dish soap and warm water

\* number of cleaning cycles determined by recommended cleaning interval and expected service life

### WARNING

To avoid electric shock, do not remove the concentrator cabinet. The cabinet should only be removed by a qualified DeVilbiss technician. Do not apply liquid directly to the cabinet or utilize any petroleum-based solvents or cleaning agents.

Use of harsh chemicals (including alcohol) is not recommended. If bactericidal cleaning is required, a non-alcohol based product should be used to avoid inadvertent damage.

## TROUBLESHOOTING

The following troubleshooting chart will help you analyze and correct minor oxygen concentrator malfunctions. If the suggested procedures do not help, switch to your reserve oxygen system and call your DeVilbiss homecare provider. Do not attempt any other maintenance.

### WARNING

To avoid electric shock hazard, do not remove the cabinet. The cabinet should only be removed by a qualified DeVilbiss homecare technician.

### Minor Troubleshooting Chart

SYMPTOM	POSSIBLE CAUSE	REMEDY
A. Unit does not operate. All lights are off when the power switch is "On." Audible alert is pulsing.	1. Power cord not properly inserted into wall outlet.	1. Check power cord connection at the wall outlet. On 230 volt units, also check the mains connection on the back of the unit.
	2. No power at wall outlet.	2. Check your home circuit breaker and reset if necessary. Use a different wall outlet if the situation occurs again.
	3. Oxygen concentrator circuit breaker activated.	3. Press the concentrator circuit breaker reset button located below the power switch. Use a different wall outlet if the situation occurs again. If the above remedies do not work, contact your DeVilbiss provider.
B. Unit operates. Red Service Required light is illuminated. Audible alert may be sounding.	1. Air filter is blocked.	1. Check the air filter. If the filter is dirty, wash it following the cleaning instructions on page 7.
	2. Exhaust is blocked.	2. Check the exhaust area and make sure there is nothing restricting the unit exhaust.
	3. Blocked or defective cannula, face mask, or oxygen tubing.	3. Detach cannula or face mask. If proper flow is restored, clean or replace if necessary. Disconnect the oxygen tubing at the oxygen outlet. If proper flow is restored, check oxygen tubing for obstructions or kinks. Replace if necessary.
	4. Blocked or defective humidifier bottle.	4. Detach the humidifier from the oxygen outlet. If proper flow is obtained, clean or replace humidifier.
	5. Flow meter set too low.	5. Set flow meter to prescribed flow rate. If the above remedies do not work, contact your DeVilbiss provider.
C. Unit operates. Audible low-frequency vibration sound is detected.	1. Electronic Assembly Malfunction.	1. Turn your unit "Off." Switch to your reserve oxygen system and contact your DeVilbiss provider immediately.
D. Yellow Low Oxygen light is on.	1. Unit in "start up" mode.	1. Allow unit up to 15 minutes to complete start up period.
E. The yellow Low Oxygen light is on and the intermittent audible signal is sounding.	1. Flow meter is not properly set.	1. Ensure the flow meter is properly set to the prescribed number. (The maximum flow meter setting is 6 LPM when an oxygen bottle is being filled with oxygen from the auxiliary port.)
	2. Air filter is blocked.	2. Check the air filter. If the filter is dirty, wash it following the cleaning instructions on page 7.
	3. Exhaust is blocked.	3. Check the exhaust area and make sure there is nothing restricting the unit exhaust. If the above remedies do not work, contact your DeVilbiss provider.
F. Red Service Required light is on and an intermittent audible signal is sounding.	1. Flow meter is not properly set.	1. Ensure the flow meter is properly set to the prescribed number. (The maximum flow meter setting is 6 LPM when an oxygen bottle is being filled with oxygen from the auxiliary port.)
	2. Air filter is blocked.	2. Check the air filter. If the filter is dirty, wash it following the cleaning instructions on page 7.
	3. Exhaust is blocked.	3. Check the exhaust area and make sure there is nothing restricting the unit exhaust. If the above remedies do not work, contact your DeVilbiss provider.
	4. Electronic Assembly Malfunction.	4. Contact your DeVilbiss provider.
G. If any other problems occur with your oxygen concentrator.		1. Turn your unit "Off." Switch to your reserve oxygen system and contact your DeVilbiss provider immediately.

## OVERVIEW OF ALARMS

This device contains an alarm system which monitors the state of the device and alerts of abnormal operation, loss of essential performance or failures. Alarm conditions are shown on the LED display. The alarm system functions are tested at power up by lighting all visual alarm indicators and sounding the audible alarm (beep).

All alarms are Low Priority Technical Alarms.

Alarm Condition	LED Indicator	Visual Alarm Signal Meaning	Audible Alarm Signal	Visual Alarm Signal Cleared by	Action to take
Start-up Period		YELLOW Low O <sub>2</sub> LED ON	No	After startup period, O <sub>2</sub> reaches at least 87%	Wait for unit to finish startup period, up to 15 minutes
Low Oxygen Concentration		YELLOW Low O <sub>2</sub> LED ON when O <sub>2</sub> is <87%	Yes Before O <sub>2</sub> drops below 82%	Turn unit off	See Troubleshooting table
Malfunction		RED Service Required LED ON	Yes	Turn unit off	Return unit to provider for service

## WARRANTY

DeVilbiss provides warranty for this device. Use only DeVilbiss approved original spare parts for maintenance and repair to retain warranty. Using unapproved or non-original spare parts would void warranty.

## RETURN AND DISPOSAL

This device may not be disposed of with household waste. After use of the device, please return the device to the provider for disposal. This device contains electrical and/or electronic components that must be recycled per EU Directive 2012/19/EU-Waste Electrical and Electronic Equipment (WEEE). Non-infectious used accessories (e.g. nasal cannula) can be disposed of as residential waste. The disposal of infectious accessories (e.g. nasal cannula from an infected user) must be made via an approved waste disposal company. Names and addresses can be obtained from the local municipality.

## PROVIDER'S NOTES - Cleaning and Disinfection When There is a Patient Change

When medical devices have already been used with a patient, contamination with human pathogenic germs should be assumed (unless there is evidence to the contrary). The next patient, user or third party should be protected by appropriate handling and preparation.

Therefore, when there is a patient change, people must be protected during the transport and handling of the device, and the device must be fully processed, i.e., cleaned and disinfected, by suitably trained personnel before reuse to protect the next patient. The complete processing may only be done by the manufacturer or by a qualified DeVilbiss provider/service technician.

**NOTE**– If the following described complete processing of the concentrator by a qualified DeVilbiss provider/technician is not possible, the device must not be used by another patient!

DeVilbiss Healthcare recommends that at least the following procedures be carried out by the manufacturer or a qualified third party between uses by different patients.

**NOTE**– If preventive maintenance is due at this time, these procedures should be carried out in addition to the servicing procedures.

1. Dispose of all accessory components that are not suitable for reuse, i.e., particularly the oxygen tubing, the nasal cannula/mask, oxygen outlet connector and humidifier bottle.
2. **The concentrator must be disconnected from the power supply for this step:** Open the concentrator and remove all dust deposits inside the cabinet with an appropriate vacuum cleaner.
3. Clean and disinfect all parts of the cabinet inside and outside and the power cord with a suitable disinfecting agent, e.g., Microbac Forte or Terralin®.
4. Check the cord, the plug on the back of the device, the power switch, the fuse holder and the indicator light for possible damage.
5. Replace all damaged or worn components.
6. Replace the cabinet air filter on the back of the device.
7. Check the oxygen concentration. If the device is within specifications, the extended life intake bacteria filter does not need to be replaced between patients. If the concentration is not within specifications, the provider should refer to the service manual section on Troubleshooting.

**NOTE**– There is no portion of the gas pathways through the concentrator that could be contaminated with body fluids under normal conditions.

The device patient connection may unintentionally become contaminated with expired gases for a single fault condition i.e., a hose internal to the device becomes disconnected. This condition will cause no flow out of the device and/or an alarm condition. Should this occur refer to the Service Manual for additional instructions.

## Disinfection

### WARNING

**Electric shock hazard. Do not attempt to open or remove cabinet. There are no user-serviceable internal components.**

**NOTE**– The disinfection process can only be completed by the manufacturer or by a qualified DeVilbiss provider/service technician.

	Recommended disinfection interval	Number of disinfection cycles *	Compatible disinfection method
Cabinet parts, inside and out, power cord	Between patients	20	Microbac Forte, Terralin®, Aldahol, Cidex OPA or Peract, 1:10 chlorine bleach (5.25%) and water solution
Oxygen Tubing, nasal cannula, Oxygen outlet connector, humidifier bottle	Do not clean, replace between patients	N/A	N/A
Inside cabinet	Between patients	N/A	Remove dust with a vacuum cleaner

# SPECIFICATIONS

DEVILBISS 10-LITER SERIES			
Catalog Number	1025DS	1025KS	1025UK
Delivery Rate	2 to 10 LPM	2 to 10 LPM	2 to 10 LPM
Maximum Recommended Flow (@ nominal outlet pressures of zero and 7 kPa)**	10 LPM	10 LPM	10 LPM
Outlet Pressure	20.0 ± 1.0 psi (138 kPa +/- 7 kPa)	20.0 ± 1.0 psi (138 kPa +/- 7 kPa)	20.0 ± 1.0 psi (138 kPa +/- 7 kPa)
Auxiliary Oxygen Port **	Outlet Pressure: 8.5±0.5 psig (58.6±3.5 kPa) Outlet Flow: 2.5±0.5 LPM Continuous	Outlet Pressure: 8.5±0.5 psig (58.6±3.5 kPa) Outlet Flow: 2.5±0.5 LPM Continuous	Outlet Pressure: 8.5±0.5 psig (58.6±3.5 kPa) Outlet Flow: 2.5±0.5 LPM Continuous
Electrical Rating	120 V, 60 Hz, 6.1 Amp	230 V~, 50 Hz, 3.2 Amp	240 V~, 50 Hz, 3.2 Amp
Operating Voltage Range	102-132 V~, 60 Hz	195-253 V~, 50 Hz	204-264 V~, 50 Hz
Oxygen Percentage	2-10 LPM=87%-96%	2-10 LPM=93% +3%/-6%	2-10 LPM=93% +3%/-6%
Operating Altitude	(tested at 70°F {21°C} only) 0-1500 M (0-5000 ft)		
	Across the voltage range: No degradation of performance	Across the voltage range: No degradation of performance	Across the voltage range: No degradation of performance
Operating Environment Range*	41°F (5°C) to 95°F (35°C), humidity range of 15% to 93% non-condensing		
	No degradation in performance across the operating voltage range.	No degradation in performance across the operating voltage range.	No degradation in performance across the operating voltage range.
Power Consumption	120 vac, 60HZ: 639 watts average	230 vac, 50Hz: 664 watts average	240 vac, 50Hz: 670 watts average
Weight	42 lbs. (19 Kilograms)	42 lbs. (19 Kilograms)	42 lbs. (19 Kilograms)
Safe Working Load	53 lbs. (24 Kilograms)	53 lbs. (24 Kilograms)	53 lbs. (24 Kilograms)
Sound Pressure Level at 3 and 10LPM	<59 dBA	<57 dBA	<57 dBA
Sound Power Level at 3 and 10 LPM	<69 dBA	<67 dBA	<67 dBA
Dimensions	24.5"H x 13.5"W x 12"D (62.2 x 34.2 x 30.4 cm)	24.5"H x 13.5"W x 12"D (62.2 x 34.2 x 30.4 cm)	24.5"H x 13.5"W x 12"D (62.2 x 34.2 x 30.4 cm)
Maximum Limited Pressure	35 PSIG (241 kPa)	35 PSIG (241 kPa)	35 PSIG (241 kPa)
Operating System	Time Cycle / Pressure Swing	Time Cycle / Pressure Swing	Time Cycle / Pressure Swing
The visible "low oxygen" indicator will activate at the following level	86% ± 3% (The audible alarm will alert at approximately 85%. At less than 60%, the red "service required" light will activate.)	86% ± 3% (The audible alarm will alert at approximately 85%. At less than 60%, the red "service required" light will activate.)	86% ± 3% (The audible alarm will alert at approximately 85%. At less than 60%, the red "service required" light will activate.)
Storage Conditions	-13°F (-25°C) to 158°F (70°C), humidity range of 15% to 93% non-condensing, including condensation	-13°F (-25°C) to 158°F (70°C), humidity range of 15% to 93% non-condensing, including condensation	-13°F (-25°C) to 158°F (70°C), humidity range of 15% to 93% non-condensing, including condensation
Equipment Class and Type	<input type="checkbox"/> Class II Equipment Double Insulated; ⚠ Type B Applied Part Ordinary Equipment, IP21	<input type="checkbox"/> Class II Equipment Double Insulated; ⚠ Type B Applied Part IP21	<input type="checkbox"/> Class II Equipment Double Insulated; ⚠ Type B Applied Part IP21
Approval Body and Safety Standard	TUV IEC 60601-1 3rd Edition	TUV IEC 60601-1 3rd Edition IEC 60601-1 IEC 60601-1-11 ISO 80601-2-69	TUV IEC 60601-1 3rd Edition IEC 60601-1 IEC 60601-1-11 ISO 80601-2-69
CE mark	No	Yes	Yes
EMC Compliance To	EN60601-1-2	EN60601-1-2	EN60601-1-2

\*NOTE- The OSD performance at 41°F (5°C) to 95°F (35°C), 93% R.H. through voltage range on the 1025 verified at 670m.

\*\*The maximum recommended flow is 6 LPM when an oxygen bottle is being filled with oxygen from the auxiliary oxygen port.

Specifications subject to change without notice.

## Oxygen Concentration vs Flow Rate

Flow L/m	%O <sub>2</sub>
10	87% - 92%
9	87% - 93%
8	87% - 95%
7	87% - 96%
6	87% - 96%
5	87% - 96%
4	87% - 95%
3	87% - 95%
2	87% - 94%

## DEVILBISS GUIDANCE AND MANUFACTURER'S DECLARATION

### WARNING

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the Electromagnetic Compatibility [EMC] information provided in the accompanying documents.

Portable and Mobile RF Communications Equipment can affect Medical Electrical Equipment.

The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

**NOTE**– The EMC tables and other guidelines provide information to the customer or user that is essential in determining the suitability of the Equipment or System for the Electromagnetic Environment of use, and in managing the Electromagnetic Environment of use to permit the Equipment or System to perform its intended use without disturbing other Equipment and Systems or non-medical electrical equipment.

<b>Guidance and Manufacturer's Declaration – Electromagnetic Emissions</b>			
The 1025DS, 1025KS and 1025UK Oxygen Concentrators are intended for use in the electromagnetic environment specified below. The customer or the user of the Concentrator should assure that it is used in such an environment.			
<b>Emissions Test</b>	<b>Compliance</b>	<b>Electromagnetic Environment – Guidance</b>	
RF Emissions CISPR 11	Group 1	The 1025 (series) Oxygen Concentrator uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.	
RF Emissions CISPR 11	Class B	The 1025DS, 1025KS and 1025UK Oxygen Concentrators are suitable for use in all establishments including domestic and those directly connected to the public low-voltage power supply network power supply that supplies buildings used for domestic purposes.	
Harmonic emissions IEC/EN 61000-3-2	Class B		
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies		
<b>Guidance and Manufacturer's Declaration – Electromagnetic Immunity</b>			
The 1025DS, 1025KS and 1025UK Oxygen Concentrators are intended for use in the electromagnetic environment specified below. The customer or the user of the Concentrator should assure that it is used in such an environment.			
<b>Immunity Test</b>	<b>IEC 60601 Test Level</b>	<b>Compliance Level</b>	<b>Electromagnetic Environment - Guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±15kV air	Complies	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Radiated RF IEC 61000-4-3	10V/m 80MHz to 2.7GHz & Proximity fields from RF wireless communications equipment as per section 8.10	Complies	Field strengths outside the shielded location from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than 3 V/m. Interference may occur in the vicinity of equipment marked with the following symbol: 
Conducted RF IEC 61000-4-6	3 Vrms 150kHz to 80MHz 6Vrms ISM Band frequencies between 150kHz-80MHz	Complies	
Electrical fast transient IEC 61000-4-4	±2kV power line ±1kV I/O lines	Complies	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV differential ±2kV common	Complies	
Power frequency magnetic field IEC 61000-4-8	30 A/m	Complies	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interrupts and voltage variations on power supply input lines IEC 61000-4-11	100% dip 0.5 cycle 100% dip 1 cycles 30% dip 25 cycles 100% dip 5 secs.	Complies	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Concentrator requires continued operation during power mains interruptions, it is recommended that the Concentrator be powered from an uninterruptible power supply or battery.

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DeVilbiss® and OSD® are registered trademarks of DeVilbiss Healthcare.
© 2017 DeVilbiss Healthcare LLC. 08.17 All Rights Reserved.
**DeVilbiss Healthcare Ltd**
**DeVilbiss Healthcare Pty. Limited**
**DeVilbiss Healthcare SAS**
**DeVilbiss Healthcare GmbH**