



Vinyl Cleaning
& Disinfecting
Guide



Table
Maintenance

Contemporary III

Operating Instructions & Manual For Elevation and Stationary Models

MANUAL CONTENTS:

Safe Operation	Page 2
Warranty / Return Policy / Maintenance	Page 3
Headpiece Operation	Page 4
Table Body Operation	Page 5
Trouble Shooting	Page 6

CONTACT INFORMATION:

Website: sunchirotable.com

Phone: (863) 414-2893

Email: admin@sunchirotable.com

Our Facility:

Sun Chiro Table
4226 Commercial Dr
Sebring, FL 33870

Facility Hours:

8am to 4pm, Monday—Friday

THESE INSTRUCTIONS ARE CURRENT AS OF JANUARY 2023.

Please contact us if you have any questions – our online customer support is always available!



Thank you for your purchase. With proper care your new table should provide many years of dependable service. If you ever have any issues or questions please contact us at 863-414-2893. We are always here to improve your experience with our equipment.

Your table has been shipped to you fully assembled and ready to use. Please note that during shipment the table may have been jarred and require re-adjusting. If anything does not seem quite right do not hesitate to contact us immediately at the number above. We will attempt to rectify any issues over the phone in a timely manner. Please read this manual in its entirety to ensure your safe use and operation of the table.

SAFE OPERATION



CRUSH HAZARD:

Improper use of elevation equipped tables can result in personal injury. It is your responsibility as the operator to ensure that tables are disabled when not in active use treating patients. Tables should only be operated by instructed personnel that have read the contents of this manual.



TRIP / SHOCK HAZARD

This table is powered by 120 V AC current. It must be plugged into a grounded 120V receptacle of at least 15 amps. If an extension cord is to be used it should be of industrial quality, 14 gauge, and no longer than 25 feet. Any cord running across the floor should be covered so as to prevent a trip hazard and protect exposed wires from damage.



DISABLE TABLE BY REMOVING KEY

Table should be disabled by removing key when operator not in the room.



WARRANTY:

- All Standard Line Tables carry a one-year warranty on materials and workmanship (upholstery is warranted for workmanship only).
- Prestige Line and Contemporary Tables carry a two-year warranty on materials and workmanship (upholstery is warranted for workmanship only).
- Elevation columns installed on Contemporary tables are warranted for one year by the manufacturer of the provided columns.
- As a courtesy service we will ship replacement parts to you at no expense, however you will be responsible for the installation of the provided parts. The exception to this service is tables which have been shipped outside the Continental US. In such situations the customer will be responsible to cover applicable shipping costs less an equivalent Domestic shipping allowance.

At Sun Chiropractic Tables, LLC sole discretion "hard parts" (defined as those mechanical table parts subject to wear and breakage) may be provided to the original customer at no charge except for the cost of shipping, at any time both before and after the termination of the warranty periods stated above.

10 DAY RETURN POLICY:

We offer a 10 day, no-questions asked return policy for all the standard tables we offer. ***Please note that this policy does NOT include tables with custom colors or other customized options. Tables with a custom color or other custom options may not be returned for any reason.***

If not satisfied with your table for any reason you may return it at your expense for a complete refund excluding original shipping paid. You must contact us within 10 days of arrival to let us know that you will be returning the table. This allows you 10 full days to decide if you are happy with your purchase.

Returns will not be accepted if you fail to contact us within 10 days of the table's delivery. Tables must be returned in "new" condition. Tables not returned in "new" condition will have their refund decreased by the stated cost of noted damage.

For the most updated warranty and return policy information, please see our website at sunchirotable.com.

TABLE MAINTAINANCE:

Maintenance:

Owners are encouraged to inspect their tables at least once every six months. Tables should be examined for loose nuts, bolts, frayed wires, etc. A log of such inspections should be kept. Drops should be lubricated at this time, if not more frequently, with a spray of WD40 on the drop plungers.

Cleaning Table Naugahyde:

The genuine Naugahyde Vinyl should be cleaned using a solution of common bleach at 2% to 98% water. This meets all disinfection requirements in the United States at the time of this writing. For heavier soil or contamination with bodily fluids please use a solution of 10% household bleach to 90% water and allow the solution to sit on the stain for a few minutes. **COVID, until the pandemic is over tables should be disinfected with the 10% bleach to 90% water solution between patients. To avoid staining patient's clothes keep a bottle of plain water available to rinse bleach residue from the table after cleaning.**

Face Paper:

Tables are designed to be used with 8.5 inch wide 150 foot long rolls commonly available from most chiropractic suppliers.

**Vinyl Cleaning
& Disinfecting
Guide**



**Table
Maintenance
Log**



General Operation

Below please find the controls and general operation instructions. For more in depth assistance with use please contact us at the number provided on the first page of this guide.

Toggle Arm, lift to unlatch and switch to toggle. There is an arm on both sides. Latch to restore forward cervical drop

Headpiece



Cocking handle

Cervical lateral bending release handle.




Cervical Tension Adjustment

Rear headpiece Up/Down Adjust release

Flexion/Extension Release Handle

Main Body Operation Elevation and Stationary*



* Stationary Tables use kick drops to activate the thoracic and lumbo/pelvic drops. All other operations are the same as the Elevation version.

Troubleshooting

Problem	Solution
Bolt spins or won't tighten	Thread has been stripped. Call support .
Drop won't cock.	Some surface rust might have formed in shipment. Easy solution is to rapidly and gently move the cocking arm up and down until the drop does cock. We then suggest spraying a very small amount of WD40 on the plunger to lubricate it.
With time the drop seems stiff	Spray the drop with a small amount of WD40.
Feel a sharp object through a cushion	Every effort has been made to ensure that all of the cushion attachment screws were not over driven into the cushion. Should you think that you feel a sharp object through a cushion, look under the table at the point where you feel the object. If you see an attachment screw, use a Phillips head screw driver to back the screw out several turns. This should solve the problem.
Drop won't cock after prolonged use.	Spray between the lifting arm and the white block some WD 40. Then using a quick tapping method, move the cocking rod up and down until the drop seats itself and returns to normal function.
Elevation Does not work	Check electrical connection to table control box, check control box to column plug. Call support.

We are always available to help! Contact us via phone, email or chat on our website (sunchirotable.com) at any time and we will be happy to assist you!