UltraLite Patient Chair with Scissor Base Model 4012

Maximum Weight: 300 Pounds

Setup



Place chair on floor.



Unfold chair.



Insert red leg into red bracket. Insert green footrest support into green bracket



Place bracket notch on crossbar.



Place chair on its feet.

Adjustment



Raise backrest to appropriate height by pushing up while stepping on leg.



To lower backrest: Hold backrest with one hand.



Slowly lift up lever on clamp block.



Adjust optional headrest.



Insert optional armrest rods.