## Olive Union



# Olive SmartEar Plus

#### Olive SmartEar Plus Instruction Manual

## Olive SmartEar Plus Instructional

Scan the QR Code to watch our instructional "How to..." Olive SmartEar Plus Videos!





	Battery	tery Olive: Lithium-ion battery 65mAh Charging Case: Lithium-ion battery 420mAh	
	Power input	Olive: DC 5V, 130mA Charging Case: DC 5V, 500mA	
=	Usage time	- Usage time (when battery is 100%) Music 6h / Hearing 9h 30m / Call 5h  * Battery life may vary depending on environmental and usage factors such as Bluetooth device distance and volume level. (Operating temperature: 0 ~ 40°C) Charge time (when battery is 0%)  - Olive within 1 hour / Charging Case within 1 hour 40 minutes	
	Size and weight	Olive EA: 0.78in(20mm) × 0.78in(20mm) × 1.22in(31mm) / 0.22oz(6.2g) Charging Case: 2.60in(66mm) × 0.90in(23mm) × 1.77in(45mm) / 1.25oz(35.5g	
	Main unit and accessories	Olive Left Earbud x 1 Olive Right Earbud x 1 Charging Case x 1 - USB Type-C cable - Urethane Eartips: XS, S, M, L, XL (2 for each size)	
	Compatible models	Compatible with iOS 12 or later. Compatible with Android 5.1 or later. Usage at a distance of 10 meters or further is not guaranteed to work. Compatibility with older or newer versions of the App is not guaranteed. Compatibility with	

## **Check Your Olive & Accessories**

Please make sure to set up indoors

old Bluetooth standards is not guaranteed. Bluetooth radio waves may cause interference with other radio waves.


Charging Case


9

LED

USB Type-C

#### Olive (Left ear & Right ear)









\* Please make sure to use the included cable.




XS/S/L/XL size 2 each

\* Size M is attached to the left and right sides

at the time of purchase. \*\* XS-size Yellow / S-size Red / M-size Green / L-size Blue / XL-size White

## **Wearing Your Olive**


To experience the best hearing with your Olive SmartEar Plus, make sure you wear them correctly!



Remove your Olive from the case and grasp it by the sides. Avoid touching the touchpad when inserting them into your

## How to distinguish left and right

The [R] at the base of one Olive indicates the right earbud, and the [L] on the other indicates the left earbud



## Pinch the eartip with your fingers, compressing it.

\* The eartip will expand once inserted into your ear canal, creating a tight seal

Insert the Olive into your ear, gently twisting it slightly as you do to ensure a deep and secure fit.


## What if my Olive doesn't fit my ears?

Finding the right fit for your Olive is very important. Make sure to try all 5 sizes (XS, S, M, L, and XL) to find the best comfort and fit.



## 3 Charging Your Olive

- (1) Put your Olive into the Charging Case and close the lid.
- (2) Connect the Charging Case to a USB AC adapter with the included USB Type-C Cable and charge for about an hour. (\* The charging time may vary slightly depending on your environment.)



## How do I charge my Olive on the go?

Your Olive can be charged by simply returning them to the Charging Case.

The Charging Case can hold around 2 full charges you can check the Charging Case's remaining battery life by opening the lid and looking at the LED lights.

## Remaining Battery Life in Charging Case



## 4 Install The Dedicated Olive Application

While your new Olive SmarEar Plus may perform well on its own, use it with the "My Olive" App to experience its full potential! Follow the instructions below to download the "My Olive" App.

\* Application is free to use and download, but user is responsible for data or packet communication fees incurred during download or use.



(1)

Open your smartphone's camera app.



Aim the camera at the QR code above. Once read, the "My Olive" app download screen should appear.



Tap the pop-up or link that appears to transition to the "My Olive" app installation screen.

#### What if I'm having trouble reading the QR Code?

Search for "My Olive" in the AppStore or Google Play and download the corresponding app shown to the right. \* Alternatively, please search for "My Olive App" in the search engine and install the app below.








## 5 Registering Your Olive With The "My Olive" App

Start the "My Olive" app and follow the instructions on the screen to proceed. When Bluetooth Connection appears on your screen, tap the "Go to Settings" text. When the Bluetooth settings screen appears, follow the steps below to connect to Bluetooth.

## If you are using an iPhone




Turn on Bluetooth and Tap "Olive" from "OTHER DEVICES"



When the "Olive" of "MY DEVICES" is connected, tap "My Olive" to return to the app



Move the Olive you want to pair close to the phone, and select the "Olive App" that has the strongest signal




Tap "Pair" from "Bluetooth Pairing Request" gu-gog

## 5 Registering Your Olive With The "My Olive" App

## If you are using an Android smartphone

\* The display may vary depending on the model




## What if the connection failed and I can't complete registration?

Please refer the "Bluetooth Connection Troubleshooting" section of this manual.

## 6 Customizing Your Olive's Touch Controls


You can customize various touch control settings via "My Olive" App. Change them for your convenience to maximize the Olive SmartEar Plus experience.





The Olive's default Touch Control Setting is a "Tap & Hold" Touch Control Action assigned to the

To assign a function to a Touch Control Action, first tap the menu bar next to the Touch Control Action (Double Tap in the example above)



Choose the function you would like to assign to the Touch Control Action (Volume Control function in the

Confirm that the Volume Control function has been assigned to the "Double Tap" Touch Control Action

Single Tap

Double Tap

R No

Q Volu

Tap & Hold Hear-Thro

On/Off
 On/Off


Additional customizable touch controls include Song Control, Volume Control, Mode Control, and Voice Assistant. Please note that the Call Control actions (answer, hang up, cancel) can only be assigned to the "Tap & Hold" Touch Control Action.

## **Bluetooth Connection Troubleshooting**

If you are having issues with pairing, please try Manual Pairing and/or resetting your smartphone's Bluetooth connection.


#### Before resetting Bluetooth connection.

- 1. Return Olives to the Charging Case.
- 2. Press & hold the Charging Case button for 4 seconds to activate manual pairing.
- 3. The Charging Case's LED will flash 3 times and blink left to right to indicate manual pairing has activated.




## For iPhone

- "Settings" > Touch "Bluetooth"
- 2 Touch the (i) icon in "MY DEVICES" for "Olive"
- Touch "Forget This Device"
- 4 Touch the (i) icon in "MY DEVICES" for "Olive App"
- Touch "Forget This Device"
- 6 Uninstall the "My Olive" app.
- Restart your smartphone.
- Please install the app again and try connecting.



- "Settings" > "Connected Devices" > "Bluetooth"
- Touch the icon in "Paired devices"
- 3 Touch "Unpair"
- 4 Uninstall the "My Olive" app.
- 6 Restart your smartphone.
- 6 Please install the app again and try connecting.



#### **Customer Support Services**

Have guestions or need further assistance? Call or email us with the contacts below! Tel. (+1) 702-848-4837

E-mail. support\_us@oliveunion.com

### **Important Safety Precautions** Please read and keep all safety, security, and use instructions.

The following section contains important information regarding safe use. Please read the following carefully before use. In addition, this instruction manual contains excerpts of particularly important items. Please check our website for the full text. If you have trouble finding the full text,

Please contact customer support at (+1) 702-848-4837

#### Safety Precautions

Marning (	Symbol indicates information that, if ignored, could possibly result in serious illness or even death.
<u></u> Caution	Symbol indicates information that, if ignored, could possibly result in personal or property damage.

### Description of Symbol

Symbol indicates usage or action that should be avoided.

■ Handling of the main unit, charger, and other accessories

## Warning

- Never attempt to disassemble, remodel, or repair the product yourself. There is a possibility of fire and/or injury.
- Do not use at unnecessarily loud volume. exacerbate existing hearing loss.
- Do not expose this product to water, liquids, or moisture. Do not use wet hands to handle the product. Contact with liquids may result in electrical shock, failure, and injury.
- Do not attempt to use with outlets or outlet strips that would exceed the safe, rated value. Overloading an electrical circuit can result in fires



Do not hit, crush, or use excessive force with the product Also, do not handle the parts roughly, such as forcibly pulling them.

Rough and improper use can cause damage, electric shock, fire, and accidents.

## Storage in the following places may cause electrical shock and failure and should be avoided.

Areas with extremely high or low temperatures/Close to heat sources and fire such as heating appliances/Hot or humid places such as bathrooms/Very dusty areas/Outdoors or in direct sunlight/Near harmful substances/Near slippery substances (grease, oil, etc.)

If you experience or notice any of the following, stop using the product immediately. Ignoring the following and continuing to use the product may cause unexpected accidents The product is broken or degraded. The product is damaged by chemicals. The product is

Use of the product in the following ways can cause unexpected accidents, injuries, or physical damage. Do not allow children or people who require supervision to use. Furthermore the product contains small pieces which should be kept out of reach of children due to choking hazard.

Do not lend or borrow this product. Do not continue to use if broken or damaged. Do not use

### Check before use

If you fall under any of these eight categories, please refrain from use and ask an otolaryngologist for a diagnosis.

- 1) Has damage or deformity due to surgery, injury, etc
- 2) Otorrhea has occurred within the last 90 days due to otitis media
  3) Has experienced sudden or progressive hearing loss within the last 90 days
  4) Has experienced hearing loss in one ear or the other within the last 90 days

- Has experienced acute or chronic dizziness
   Has bots/unusual amounts of earwax
   Has eczema, pain, or discomfort in the external auditory meatus. 8) Has abnormal air bone conduction in hearing at 500,1000 and 2000 Hz

### Care Procedures

Before and after use, clean the main body with a dry and clean cloth or brush.

#### About Product Disposal

Please dispose of the product in an appropriate way according to the instructions of the local



#### ECC ID : 2AOI HOSE300M

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause

## Warranty

Product	Olive SmarEar Plus					
Warranty Period	Purchase Date	Month	Day	Year		
Retailer	Name of Retailer :					
	Address :					
	Phone Number :					

#### **Warranty Policy**

■ In the unlikely event that a problem or failure occurs within the warranty period,

If the product breaks down under normal usage during the warranty period, we will replace it for free. However, consumable parts and accessories are not covered by the warranty.

Consumable parts: Eartips

Accessories: USB cable

- For the following cases, even if within the warranty period, replacement and repairs will require payment. If warranty information is not properly filled out
   If this document has been lost, or the wording of this document has been modified
- 3. Failure or damage caused by misuse or improper repair or modification
- Natural disasters such as fire, earthquake, flood, lightning, etc. Furthermore, pollution, abnormal voltage, or damage from the usage of an undesignated power supply
   Damage caused by carelessness of the customer or failure caused by improper repair,
- decomposition, or modification Failure or deterioration due to use that is not in accordance with this manua
- 7. Exchange of consumables and accessories
- This warranty is only valid in the United States of America. This manual & warranty will not be reissued, so please keep it safe and secure.

#### [Contact Us] Olive Customer Support

Tel. (+1) 702-848-4837 E-mail. support\_us@oliveunion.com