

OPERA® SOLO END
WALL PROTECTOR

Installation and User Guide





Before you begin! We recommend that you conduct a risk assessment to ensure that the Opera Solo End Wall Protector meets your requirements and does not become a trip hazard.

1. General Information and Safety

1.1 Intended & Non Intended Uses

The Opera Solo End Wall Protector is designed to reduce the risk of an Opera Solo Profiling Bed being placed too close to a wall, plug socket, skirting board or fixture, which may be damaged by the bed. The wall protector gives a 100mm clearance from the bed surrounds to the wall or fixture.

The Opera Solo Wall Protector is designed to be used with the Opera Solo Profiling beds. Opera will not be liable for any hazards, injury, or damage arising from using the wall protector incorrectly or with a different bed.

2. User Instructions



Always read the instruction manual before using the wall protectors, pay particular attention to the instructions below.



Always carry out a risk assessment to ensure that the wall protector is suitable for the required application. Ensure that the installation of wall protectors does not create risks for both staff and users, in particular, trip hazards.



The wall protectors should only be fitted to the bed at points where the bed is against a wall, socket, fixture or skirting board, non compliance with this may result in creating trip hazards. This can result in serious injury.

If the bed is moved away from the wall without the wall protectors being removed, this can result in creating trip hazards, and may result in serious injury. Always remove the wall protectors before moving and relocating the bed.



Check the fixings monthly to ensure that there are no loose parts.

It is recommended that wall protectors are attached to the chassis of the bed's head or foot end that face a wall, socket or skirting board. To attach the wall protectors, please follow the instructions below:

1. Place the end tube of the wall protector into the hole on the castor surround.
2. Place the opposite end of the tube into the opposite castor surround.
3. Adjust the wall protector if required and ensure that the bed can still lower to its floor position.

When in doubt, consult Opera at support@operabeds.com.

1 | Align the end of the wall protector with the hole.



2 | Place the wall protector into the hole.



3 | Place the opposite end of the wall protector into the castor hole and ensure the bed can still lower to its lowest position.





Azure House, Connaught Road, Kingswood, Hull, HU7 3AP

0333 222 8584 | support@operabeds.com | operabeds.com

Version 1 | Published in August 2022