i-FIX®

EN	Glass ionomer luting cement
LT	Stiklo jonomerinis cementas cementavimui
DE	Glasionomer-Befestigungszement
BG	Стъклено -иономерен цемент за цементиране
CS	Skloionomerní tmelovací cement
DA	Glasionomer -cementering af cement
ET	Klaasionomeerist liimtsement
FR	Ciment de scellement verre ionomère
EL	Τσιμέντο ιοντομερούς γυαλιού
HU	Üvegionomer tömítő cement
IT	Cemento cementizio vetroionomerico
LV	Stikla jonomēru blīvējošais cements
NO	Glasionomer luting sement
PL	Cement glasionomerowy do osadzania
PT	Cimento de cimentação de ionómero de vidro
RO	Ciment ionizant de sticlă ionomer
SK	Skloionomerný cement na tmelenie
SL	Stekloionomerni cement za cementiranje
ES	Cemento de fijación de ionómero de vidrio
SV	Glasjonomer som luter cement
TR	Cam iyonomer yapıştırma simanı
RU	Стеклоиономерный цемент для фиксации
UA	Склоіономерний цемент для цементування

EN ISO 9917-1 4.2a Glass Polyalkenoate Cement

INSTRUCTION FOR USE EN

DESCRIPTION

i-FIX is a chemical curing, bioinert, glass ionomer fixing cement with high adhesion to dentin, enamel and non-precious alloys. Has low film thickness, good translucency and excellent biocompatibility.

COMPOSITION

Fluoroaluminosilicate glass 75-95%, water 20-25%, polyacrylic acid 5-25%, excipient 0-2%, pigment 0-1%.

i-FIX does not contain medicinal substance, including human blood or plasma derivative; tissues or cells, or their derivatives, of human origin; tissues or cells of animal origin, or their derivatives, as referred to Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

compressive strength	126.1 MPa
film thickness	15.4 µm

INTENDED PURPOSE AND CLINICAL BENEFITS

i-FIX restores/improves aesthetic appearance of restorable tooth; restores/maintains dental function of restorable tooth; protects biological structures of restorable tooth and adjacent tissues.

CLINICAL INDICATIONS

- For cementation of crowns, bridges, inlays, onlays, orthodontic brackets and bands.

CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS

None known.

UNDESIRABLE SIDE EFFECTS

In susceptible individuals, i-FIX may cause allergic or irritation reactions (eye).

RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OR BODY FLUIDS

Part of the body – mouth. Tissues or body fluids contacted by the device – tooth, oral mucosa, saliva.

INTENDED USER

i-FIX is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental cements. There is no need for specific training.

STERILITY

i-FIX is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT

i-FIX is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of cement is suitable for single use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

CONSUMABLE COMPONENTS AND ACCESSORIES

No accessories are supplied with the device. Consumables, such as measuring scoop, mixing pads, are supplied with the device.

INSTRUCTION FOR USE

1. Prepare surface as always. For deep cavities use calcium hydroxide liner.

2. Mixing ratio: 1 level spoon of powder / 2 drops of liquid. Lightly tap the bottle of powder against the hand before dispensing. Do not shake or invert. Hold liquid bottle vertically upright and squeeze gently when dispense a liquid.
3. For mixing use metal or plastic spatula and glass slab or paper mixing pad. Divide powder into two parts. Mix the first half of powder with liquid for 10 seconds, and then add remainder. Mix until exact consistency is reached. Total mixing time is 30 seconds. Do not add powder in small portions.
4. The working area must be kept dry. Apply the mixed cement with suitable instrument to both restoration and prepared tooth immediately after mixing. Total working time is 2 minutes including mixing. Excess cement may be removed immediately or after cement has set. Do not attempt to remove excess cement in the phase of setting (1.5-6 minutes from start of mix) as adhesion may be impaired.
5. After 8 minutes from start of mixing restoration is ready for finishing.

WARNINGS

As i-FIX highly adhesive, clean all instruments before cement set on. Higher amount of powder or higher temperature shorten the working/setting time and lower temperature, or higher amount of liquid prolong the working/setting time. Do not use i-FIX for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. i-FIX does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS

It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

Wash hands thoroughly after handling. Avoid contact with eyes. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

SHELF-LIFE

Shelf-life of i-FIX is 3 years of powder / 5 years of liquid from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

STORAGE

Keep product tightly closed in dry well-ventilated place at 5-25°C. Protect from moisture, direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

i-FIX is safe and performs as intended if it is used in accordance to manufacturer's instruction for use. Summary of safety and clinical performance will be introduced in EUDAMED as soon as it will start work.

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

PACKAGING

REF IIXTP 20g powder, 13g liquid, measuring scoop, mixing pad

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