Instructions to User

Dear Customer,

Thank you for purchasing this quality product. Please read the manual very carefully before using this device. Failure to follow these instructions can cause measuring abnormality or damage to the Oximeter.

No part of this manual may be photocopied, reproduced or translated into another language without the prior written consent. We reserve the right to revise and amend it at any time without prior notice.

Version of the Manual: Ver 1.6

Issued Date: November 6, 2018

Copyright CMI Health Inc.

All rights reserved.

Notes:

- The contents contained in this manual are subject to change without prior notice.
- Information furnished by **CMI Health** is believed to be accurate and reliable. However, no responsibility is assumed by **CMI Health** for its use, or any infringements of patents or other rights of third parties that may result from its use.

3502-2490003

Instructions for Safe Operation

- Check the device to make sure that there is no visible damage that may affect user's safety and measurement performance. It is recommended that the device should be inspected minimally before each use. If there is obvious damage, stop using the device.
- Necessary service must be performed only by qualified technicians. Users are not permitted to service this device.
- The oximeter must not be used with the devices and accessories not specified in User Manual.

Cautions

- Explosive hazard—DO NOT use the oximeter in environment with inflammable gas such as some ignitable anesthetic agents.
- **DO NOT** use the oximeter while the Patient is under MRI or CT scanning. This device is NOT MRI Compatible.

Warnings

- Discomfortable or pain may occur if using the sensor of this device continuously on the same location for a long time, especially for the patients with poor microcirculation. It is recommended that the Oximeter should not be applied to the same location for longer than 2 hours or less if any abnormal condition is found. Frequently check and re-position the Oximeter sensor.

- For the individual patients, there should be a more prudent inspecting in the placing process. The sensor can not be placed on the edema and tender tissue.
- The local law should be followed when disposing of the expired device or its accessories.

Attentions

- Keep the Oximeter away from dust, vibration, corrosive substances, explosive materials, high temperature and moisture.
- If the Oximeter gets wet, please stop operating it and do not resume operation until it is dry and checked for correct operation. When it is carried from a cold environment to a warm and humid environment, please do not use it immediately. Allow at least 15 minutes for the Oximeter to reach ambient temperature.
- **DO NOT** operate the button on the front panel with sharp materials or sharp point.
- **DO NOT** use high temperature or high pressure steam disinfection on the oximeter and probes. Refer to related chapter for instructions regarding cleaning and disinfection.
- The intended use of this device is not for therapy purpose.

Declaration of Conformity

The manufacturer hereby declares that this device complies with the following standards:

IEC 60601-1, ISO 9919:2009, ISO 80601-2-61:2011

and follows the provisions of the council directive MDD93/42/EEC.

Table of Contents

1 Overview.....	1
1.1 Appearance.....	1
1.2 Product Name and Model.....	2
1.3 Structure.....	3
1.4 Features.....	3
1.5 Intended Use.....	3
1.6 Working Environment.....	3
2 Battery and Stand.....	4
3 SpO₂ Probe Connection.....	6
4 Operation.....	8
4.1 Power on/off the Oximeter.....	8
4.2 Default Display Screen.....	8
4.3 Display Screen with PI Value.....	10
4.4 Menu Setup.....	10
4.5 Data Recall.....	17
5 Technical Specifications.....	21
6 Over-limit Indication.....	23
6.1 Limit settings.....	23
6.2 Over-limit indication sound mute setting.....	23
7 Packing List.....	24
8 Repair and Maintenance.....	25
8.1 Maintenance.....	25
8.2 Cleaning and Disinfecting Instruction.....	25
9 Troubleshooting.....	26
Appendix.....	27

1 Overview

1.1 Appearance


Figure 1-1

1. Display screen: display SpO₂ plethysmogram and parameter values.

2. Navigation keys:

▲ : Up/Left/Increase

Press this key, the default screen can be shifted to display pulse rate (PR) or perfusion index (PI). If on the system setup screen, press it to move the cursor upwards or to the left and adjust parameter values.

▼ : Down/Right/Decrease

Its function is similar with the key “▲ : Up/Left/Increase”.

3.  **(Mode/OK):** press this key, the screen can be shifted between default screen display and alternative screen display; longtime press it, the menu screen will be displayed; when you finish parameter setting, press this key to confirm.
4.  **(Data interface):** used for uploading data (Optional function).
5.  **(DC power input):** used for connecting external DC power input for recharging the built-in rechargeable battery.
6.  **(Power/Back):** Power on/off the device by longtime pressing; short time press it to back to the previous level of menu while setting menu operation.
7.  **(Recall/Mute):** Longtime press it to enter SpO₂ trend data recall screen; when the device is beeping, short time pressing will mute the indication sound, the mute state will persist for about 90s. After this mute period (90s), then the indication sound will resume.

8. **Icon: “SpO₂”:** SpO₂ Probe Connector.

1.2 Product Name and Model

Name: Pulse Oximeter

Model: PC-66L

1.3 Structure

It consists of the main unit and an adult SpO₂ probe.

1.4 Features

- ❖ It is lightweight, small in size and easy to carry
- ❖ Color LCD to display plethysmogram and parameters
- ❖ Monitor SpO₂ and Pulse Rate simultaneously
- ❖ PI (Perfusion Index) display is available
- ❖ Up to 384 hours storage and recall of SpO₂ and PR data.
- ❖ Audible and visual alert function is available
- ❖ Data transmission to PC for view and analysis
- ❖ Power saving mode is available

1.5 Intended Use

This Pulse Oximeter is intended for checking and recording the pulse rate and functional oxygen saturation (SpO₂) of adult and children for wellness use. This device is not intended to diagnose and/or treat any medical condition or disease.

1.6 Working Environment

Operating temperature: 5~40°C

Operating humidity: 30~80%

Atmospheric pressure: 70kPa~106kPa

2 Battery and Stand

- 1) The battery compartment is located at the back of the device and can be open with a coin or a flat-head screwdriver, as shown in Figure 2-1. **DO NOT open the battery compartment unless told by CMI tech support.**

Figure 2-1

- 2) Insert the metal stand as shown in Figure 2-2.

Figure 2-2 Metal Stand

Note: When the device shows low battery, the user should recharge the battery in time. To recharge the battery, connect one end of the charging cable to the device's top DC power input, and connect the other end to the AC power adapter. When the screen shows the rolling icon "", it means the battery is charging; When the screen shows icon "", it means the battery is fully charged. Unplug the charging cable from the device, then the battery icon will go back to normal.

Note: When changing the rechargeable lithium battery, unplug the old lithium battery and plug in the new battery into the corresponding connector in the battery compartment. Securely place the new battery in the battery compartment.
Do Not attempt unless told by CMI tech support.

Safety instruction for operation:

- Do not throw the device or battery into water, liquid or fire.
- Keep the battery out of the reach of the child.
- Do not disassemble the battery.
- The local law should be followed when disposing of the expired device or its accessories in order to protect environment from being polluted.

- ⌚ If the device does not power properly or the battery is damaged, please contact CMI customer service.
- ⌚ In order to prolong the lithium battery's using life, please pay attention to the battery maintenance.

3 SpO₂ Probe Connection

Connect the SpO₂ probe to the connector labeled “SpO₂” at the right side of the Oximeter. After starting up the Oximeter, insert one finger (index finger, middle finger or ring finger with proper nail length) into the probe according to the demonstration shown in the following figure.


Figure 3-1 Illustration of using SpO₂ probe (adult finger clip)

Instructions of Operation

1. The finger should be put in properly and correctly.
2. Do not shake the finger and keep at ease during measurement.
3. Do not put wet finger directly into sensor.

4. Avoid placing the sensor on the same limb which is wrapped with a cuff for blood pressure measurement or during venous infusion.
5. Do not let anything block the emitting light from the sensor.
6. Vigorous exercise and electrosurgical device interference may affect the measuring accuracy.
7. Using enamel or other makeup on the nail may affect the measuring accuracy.
8. If the first reading appears with poor waveform (irregular or not smooth), then the reading is unlikely true, the more stable value is expected by waiting for a while, or a re-inserting finger is needed when necessary.

4 Operation

4.1 Power on/off the Oximeter

- Press and hold "Ⓐ" power/back key for 2 seconds to power the device ON or OFF.
- During measurement, if the "Probe off" status is showing for longer than one minute, and if there is no keypad operation for one minute, then the device will power off automatically.

4.2 Default Display Screen

Press and hold “Ⓐ” power key for 2 seconds to power on the Oximeter. The display screen will show: "**Please use appropriate sensor according to the patient type and change the related settings**". Then the display will show the default screen, as shown in Figure 4-1.


Figure 4-1 Default Display Screen

Description:

Sound indicator "🔇" means the device is on mute status, the user can press recall/mute "🔇" key to turn on the indication sound. During the measurement, over-limit event or probe off event can trigger the alarm.

"꽉": Memory full icon; if the memory is full, the memory full icon will appear on the screen. No display of this icon means the current storing space is not full. If the memory is full, the data storing will continue in such way the new record will overwrite the oldest record, so that it's recommended to upload the stored data into the computer before the memory gets full.

During the measurement, press recall/mute "🔇" key to turn off (or resume) the device sound (including pulse beep sound, audible alert and key click), while the pulse symbol "❤️" still blinks. For alert sound mute function. Refer to Section 6.2 for detail of the sound mute function.

Note: The pitch tone of pulse beep is modulated by the SpO₂ value, that means the pitch tone changes when the measured SpO₂ changes. The higher the SpO₂ value is, the higher the tone frequency of pulse beep (sound becomes sharper); The lower the SpO₂ value is, the lower the tone frequency of pulse beep (sound becomes flatter).

4.3 Display Screen with PI Value

On the default display screen, press “ Δ/∇ ” Navigation key to shift screens between default screen and display screen with PI value. The display screen with PI value is shown below.



Figure 4-2 Display Screen with PI Value

4.4 Menu Setup

On the above mentioned screens, longtime press “ \square/\leftarrow ” key for entering into setup menu screen (as shown in Figure 4-3).


Figure 4-3 Setup Menu Screen

Screen Description

- “**Setting**”: set parameter values, refer to Chapter 4.4.1 for details.
- “**Upload data**”: enter into data uploading state, refer to Chapter 4.4.2 for details.
- “**Default**”: enter into the factory default setting, refer to Chapter 4.4.3 for details.
- “**espanol**”: this Oximeter provides the display with two languages: English and espanol(Spanish).
- “**Version**”: for viewing version number of the software, refer to Chapter 4.4.5 for details.

4.4.1 Setting

On the menu screen, select “Setting” and then press “” key for entering into system setup screen. The setup screen is as shown in Figure 4-4.

Setting	
Patient	ADU 
SpO2 Low-limit	99%
PR Hi-limit	100
PR Low-limit	30
Date	2013-10-22
Time	13:20:22
Recording	Interval 1s
Power saving	ON
Volume	2

Figure 4-4 System Setup Screen

Operation Instructions:

- **Patient:** two options:"ADU " for adult, and "PED " for pediatric.

- **SpO₂ Lo-Limit:** SpO₂ low limit setting; range: 50%~99%, the step is 1%. The factory default value for adult is 90% and 95% for pediatric.
- **PR Hi-Limit:** High limit setting of pulse rate; range: 100~240bpm. From 100 to 150, the step is 1bpm, and from 150 to 240, the step is 5bpm. The factory default value for adult is 120bpm and 160bpm for pediatric.
- **PR Lo-Limit:** Low limit setting of pulse rate; range: 30~99bpm, and the step is 1bpm. The factory default value for adult is 50bpm and 60bpm for pediatric.

Note: When the SpO₂ reading is equal to or below the preset alert setting or the PR reading is equal to or above the preset alert setting, then the over-limit alert sound will trigger, and the corresponding reading(s) will blink. When measured on pediatric patients, the alert sound and blinking display will trigger if the SpO₂ reading is off the limit for 10 seconds.

➤ **Date:** Date setting

- 1) When cursor stays on the Year of the date, press “” (Mode/OK) key to active Year option, the cursor flashes on the Year of the date;
- 2) Press / (Navigation key) to adjust year.
- 3) Press “” (Power/back) key or “” (Mode/OK) key to

confirm and exit from date setting.

- 4) The procedures of adjusting Month value and Day value are the same with Year adjustment.

Date Format: yyyy-mm-dd

Note: The setting operations of other parameters (such as TIME, PATIENT, RECORDING INTERVAL, POWER SAVING etc.) are the same with date setting.

- **Time:** Time setting
- **Recording:** Time interval for recording data (SpO₂ & PR), there are five options: "1s, 2s, 4s, 8s" and "OFF".
 - 1) "1s": the least length of data record is set to 30 seconds, and the maximal length for one record is limited to 1 hour. The total storage time is up to 48 hours.
 - 2) "2s": the least length of data record is set to 60 seconds, and the maximal length for one record is limited to 2 hours. The total storage time is up to 96 hours.
 - 3) "4s": the least length of data record is set to 120 seconds, and the maximal length for one record is limited to 4 hours. The total storage time is up to 192 hours.
 - 4) "8s": the least length of data record is set to 240 seconds, and the maximal length for one record is limited to 8 hours. The total storage time is up to 384 hours.
 - 5) When the option is set to "OFF", the device will not store the measuring data.

- **Power saving:** power saving setting; two options: "on" and "off". The factory default setting is "on". When the power saving mode is set to "on" option, during the measurement, if there is no key operation for 2 minutes, the screen display will be dim for power saving. The display brightness will resume to normal condition by pressing any key.
- **Volume:** 3 levels setting: "1", "2" and "3", corresponding to the speaker volume icons of "", "" and "" respectively. The speaker volume icon indicates the volume of any sound generated by the device for key press, over-limit alarm, pulse beep and audible probe-off indication etc..

4.4.2 Upload Data

On the menu screen, select “UPLOAD DATA” and then press key for entering into connecting status (as shown in Figure 4-5).

Note: Make sure the provided USB data cable is well connected between the device and PC before uploading data.

When you transmit data (SpO_2 and PR values) to your computer, please let the oximeter stay in connecting status, then data uploading will be activated. Refer to the instruction in “Oximeter Data Manager User Manual”for detailed operation.


Figure 4-5 Connecting Status Screen

4.4.3 Default

On the menu screen, select "Default" and press " " key to enter the default setting screen (as shown in Figure 4-6). Press the Navigation " / " key to choose "Yes" or "No", and press " " mode/OK key to confirm or exit. Press " " power/back key to return to the previous menu screen.


Figure 4-6 Default Setting Screen

4.4.4 Language Selection

- In Spanish version: on main menu screen, select "English" and then press " " key to change display

language changes to English.

- In English version: on main menu screen, select "espanol" and then press "≡" key to change display language to Spanish.

4.4.5 Version

On the menu screen, select “VERSION” and then press “≡” key for entering into version screen (as shown in Figure 4-7).


Figure 4-7


4.5 Data Recall

On the default display screen, press and hold “⊗” (Recall/mute) key to enter the record list screen.

2013-01-09	12:09:35
2013-01-09	15:07:35
2013-01-09	10:03:35
2013-01-09	12:50:35

Figure 4-8 Record List

4.5.1 Data Recall

Choose record in the record list, then press “

The figure shows a dual-axis trend graph. The top axis represents SpO₂ levels, ranging from 60% to 100%. The bottom axis represents PR (Pulse Rate) in beats per minute (bpm), ranging from 60 to 140. A yellow line represents the SpO₂ trend graph, which remains relatively stable around 98-100%. A green line represents the PR trend graph, showing a fluctuating pattern between approximately 75 and 100 bpm. The graph is set against a background grid. At the top left, there is a small silhouette icon. To the right of the icon are the text "12:50:30" and "2s". To the right of the PR axis, the text "PR" is displayed. Labels on the left side identify the "SpO₂ trend graph (yellow line)" and the "SpO₂ ordinate". Labels on the right side identify the "Recording interval", "PR ordinate", and "PR trend graph (green line)".

Figure 4-9A Trend Graph Display Screen

Operation Instructions:

- Press “18


Figure 4-9B Trend Graph Display Screen


Figure 4-9C Trend Graph Display Screen

4.5.2 Data Deletion (all device data)

On the menu screen, longtime press “” key and the records list will appear. At this time, longtime pressing “” key again, an message “Are you sure to delete all?” prompts on the screen, as shown in Figure 4-10.


Figure 4-10

At this time, press / navigation key to select "Yes" or "No", and press mode/OK key to confirm or exit. Or short time press power/back key to return to record list screen.

5 Technical Specifications

A. Display Mode: Color dot-matrix LCD;

B. Power Supply:

Supply voltage: d.c. 3.6V (Lithium rechargeable battery)

Operating current: $\leq 180\text{mA}$

C. SpO₂ Measurement

Transducer: dual-wavelength LED sensor with wavelength:

Red light: 663 nm, Infrared light: 890 nm.

Maximal average optical output power: $\leq 2\text{mW}$

Display range: 35~99%

Measuring range: 35~100%

Measuring accuracy:

A_{rms} is not greater than 3% for SpO₂ range from 70% to 100%

*NOTE: A_{rms} is defined as root-mean-square value of deviation according to ISO 9919 / ISO 80601-2-61.

D. Pulse Rate Measurement

Display and measuring range: 30bpm~240bpm

Accuracy: $\pm 2\text{bpm}$ or $\pm 2\%$ (whichever is greater)

Alarm setting range: 25bpm~250bpm

Default Alarm Setting: High -- 120bpm, Low -- 50bpm

E. Perfusion Index Display

Range: 0.2%~20%

F. Operating Environment

Operating Temperature: 41°F ~104°F

Operating Humidity: 30%~80%

Atmospheric pressure: 70kPa~106kPa

Note: portable and mobile RF communications equipment may affect the performance of the Oximeter.

G. Data Update

8 beats moving average for both SpO₂ and pulse rate readings

H. Data Storage

Recording SpO₂ and pulse rate data every 1/2/4/8 second(s), up to 384-hour records can be stored.

I. Low Perfusion Performance

The accuracy of SpO₂ and PR measurement will meet the precision described above when the modulation amplitude is at least 0.5%.

J. Resistance to interference of surrounding light:

The difference between the SpO₂ value measured in the condition of indoor natural light and that of darkroom is less than $\pm 1\%$.

K. Dimensions: 145 mm (L) × 74 mm (W) × 29 mm (H)

Net Weight: 210g (including batteries)

L. Classification

Type of protection against electric shock:

Internally powered equipment

Degree of protection:

Type BF applied parts.

Degree of protection against harmful ingress of liquids:

IP22 .

The device is not intended for use in the environment with rich oxygen or ignitable gas.

Mode of operation: Continuous operation.

Electro-Magnetic Compatibility: Group I, Class B

6 Over-limit Indication

6.1 Limit settings

- SpO₂ low limit setting range: 50% ~ 99%.
- Pulse Rate limits setting range:
High: 100bpm--240bpm Low: 30bpm--99bpm

During the measurement, if the measured value exceeds the preset value, the alert beeping sound will be activated, the value that is over-limit will blink at the same time.

6.2 Over-limit alarm mute setting

- During the measurement, if the over-limit alarm is on, press "⊗" recall/mute key, then the over-limit alarm will mute for 90 seconds, but the over-limited value will still keep blinking. At this moment, the speaker volume icon becomes "🔈". If this alert event persists over 90 seconds, then the over-limit alarm will be activated again.
- During the measurement, if the probe is off or disconnected, the message "Check Probe" shows on the

display screen. The alert sound starts (interval is 5 seconds) and lasts for about 1 minute. If the probe is still off, the Oximeter will power off automatically.

7 Packing List

1. The PC-66L oximeter
2. An adult SpO₂ probe
3. A SpO₂ extension cord
4. A metal stand
5. Charging cable and adapter
6. User manual
7. Quality inspection certificate
8. A data cable

Note: The accessories are subject to change. See the package in your hand for detailed items and quantity.

8 Repair and Maintenance

8.1 Maintenance

The expected service life(not a warranty) of this device is 5 years. In order to ensure its long service life, please pay attention to the maintenance;

- The recommended storage environment of the device:
Ambient temperature: -4°F ~140°F
Relative humidity 10%~95%
Atmospheric pressure: 50kPa~107.4kPa
- The oximeter is calibrated in the factory before sale, there is no need to calibrate it during its life cycle. However, if it is necessary to verify its accuracy routinely, the user can do the verification by means of SpO2 simulator, or it can be done by the local third party test house.

8.2 Cleaning and Disinfecting Instruction

- Surface-clean sensor with a soft cloth by wetting with a solution such as 75% isopropyl alcohol, if low-level disinfection is required, use a 1:10 bleach solution.
- Then surface-clean by a dampened cloth and let it air dry or wipe it with a cloth.

⚠ High-pressure disinfection cannot be used on the device.

⚠ Do not immerse the device in liquid.

9 Troubleshooting

Trouble	Possible Reason	Solution
Device kept showing “Check probe!”	The cord is disconnected or broken circuit.	Check the probe to see if the light source is blinking red light. If not, either the connection cord or the sensor is malfunctioning. Please replace or repair the cord or/and the sensor.
	The device could not detect signals.	Misplaced finger or foot in the sensor for a period could lead to undetected SpO2 or pulse rate signals. Turn off the devices, place the finger or the foot securely in the sensor, and turn on the device again.
	The device is multifunctioning.	Please contact customer service.
Infant foot warp sensor kept coming off	The Velcro loosens.	Replace the Velcro every 2-3 weeks depending on user scenario.
Suspect inaccurate reading	Finger or foot misplaced in the sensor	Replace finger or secured foot wrap.
	Light disruption	Secure finger or foot in the sensor. Light from other sources may affect the sensor and the reading signals.
	Body movement	Try to stabilize the measuring site when measuring SpO2 and pulse rate.

	Improper sensor/improper body type setting	Please use the suggested sensor according to user's body type and usage. Improper sensor or body type setting may cause inaccurate readings.
	The device is multifunctioning.	Please contact customer service.
No signal	Device shows two dashes	If the monitor still shows pulse plethysmogram, wait for few more seconds for the device to show numeric readings. If no pulse plethysmogram either, check sensor placement and sensor integrity.
	The cord is disconnected or broken circuit.	Check the probe to see if the light source is blinking red light. If not, either the connection cord or the sensor is broken circuit. Please replace or repair the cord or/and the sensor.
	The device could not detect signals.	Sensor misplacement could lead to undetected SpO2 or pulse rate signals. Turn off the devices, place the finger or the foot securely in the sensor, and turn on the device again.
	Finger or foot misplaced	Replace finger or secured foot wrap.
	The device is multifunctioning.	Please contact customer services.

FAQ

	Question	Answer
Battery	How long does the battery last with full charge?	About 10 hours.
	Can I replace the rechargeable batteries?	No. Do not attempt to replace the rechargeable battery. Doing so will void the device warranty.
Data storage	How is the data stored?	There's a build-in memory inside the device that could store up to 384 hours of data.
	How do I manage the data?	CMI Health's offers a free PC software, Oximeter Data Manager, which could easily download from our website at: www.cmihealth.com/download
	What is the data file format?	The data could be export as .csv format.
PC Software	I have a Mac computer; can I use the software?	No. The software is only compatible with PC computer. CMI Health currently does not offer computer software program for Mac users.
	Is the software free? Does the software come with the purchase?	Yes. The software is free. Simply download the software from CMI Health's website: www.cmihealth.com/download

	What computer system does the PC software support?	CMI Health's Oximeter Data Manager software support Windows 7/8/10.
	How can I transfer data from one computer to another?	CMI Health's Oximeter Data Manager software provides with data import and export feature. Simply export data from one computer and import data to another computer.
	Where can I find the user manual of the software?	The user manual could be found after installing the software. Click on "Help" and "user manual"

Appendix

Key of Symbols

	Symbol	Description
Symbols on the screen	%SpO ₂	The oxygen saturation
	PI%	Perfusion Index
	bpm	Pulse rate (Unit: beats per minute)
		Low battery voltage
		Battery fully charged
		Speaker mute icon
		Speaker volume icon
		Memory full
		(Pediatric/Adult) Patient type

Symbol	Description
Symbols on the panels	SpO₂ SpO ₂ probe connector
	 Power/Back Key
	 Mode/OK Key
	 Recall/Mute Key
	 Navigation Key
	 Data Interface
	 CE mark
	 Serial number
	 Date of manufacture
	 Manufacturer (including address)
	 With Type BF applied part
	 See User Manual
	 Disposal of this device according to WEEE regulations
	 No alarm
	 Do not litter at will

Pulse Oximetry Knowledge

1 Meaning of SpO₂

SpO₂ is the saturation percentage of oxygen in the blood, so called O₂ concentration in the blood; it is defined by the percentage of oxyhemoglobin (HbO₂) in the total hemoglobin of the arterial blood. SpO₂ is an important physiological parameter to reflect the respiration function; it is calculated by the following method:

$$\text{SpO}_2 = \text{HbO}_2 / (\text{HbO}_2 + \text{Hb}) \times 100\%$$

HbO₂ are the oxyhemoglobins (oxygenized hemoglobin), Hb are those hemoglobins which release oxygen.

2 Principle of Measurement

Based on Lamber-Beer law, the light absorbance of a given substance is directly proportional with its density or concentration. When the light with certain wavelength emits on human tissue, the measured intensity of light after absorption, reflecting and attenuation in tissue can reflect the structure character of the tissue by which the light passes. Due to that oxygenated hemoglobin (HbO₂) and deoxygenated hemoglobin (Hb) have different absorption character in the spectrum range from red to infrared light (600nm~1000nm wavelength), by using these characteristics, SpO₂ can be determined. SpO₂ measured by this oximeter is the functional oxygen saturation -- a percentage of the hemoglobin that can transport oxygen. In

contrast, hemoximeters report fractional oxygen saturation – a percentage of all measured hemoglobin, including dysfunctional hemoglobin, such as carboxyhemoglobin or metahemoglobin.

Clinical application of pulse oximeters: SpO₂ is an important physiological parameter to reflect the respiration and ventilation function, so SpO₂ monitoring used in clinical becomes more popularly, such as monitoring the patient with serious respiratory disease, the patient under anesthesia during operation, premature and neonate. The status of SpO₂ can be determined in time by measurement and find the hypoxemia patient earlier, thereby preventing or reducing accidental death caused by hypoxia effectively.

3 Normal SpO₂ Range and Default Low Limit

In campagna area, healthy people's SpO₂ value is greater than 94%, so the values below 94% are determined as hypoxia. SpO₂<90% is considered as the default threshold for determining anoxia by most researchers, so SpO₂ low limit of the oximeter is set as 90% generally.

4 Factors affecting SpO₂ accuracy (interference reason)

- ✧ Intravascular dyes such as indocyanine green or methylene blue
- ✧ Exposure to excessive illumination, such as surgical lamps, bilirubin lamps, fluorescent lights, infrared heating lamps, or direct sunlight.
- ✧ Vascular dyes or external used color-up product

such as nail enamel or color skin care

- ✧ Excessive patient movement
- ✧ Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line
- ✧ Exposure to the chamber with High pressure oxygen
- ✧ There is an arterial occlusion proximal to the sensor
- ✧ Blood vessel contraction caused by peripheral vessel hyperkinesias or body temperature decreasing

5 Factors causing low SpO₂ value (pathology reason)

- ✧ Hypoxemia disease, functional lack of HbO₂
- ✧ Pigmentation or abnormal oxyhemoglobin level
- ✧ Abnormal oxyhemoglobin variation
- ✧ Methemoglobin disease
- ✧ Sulfhemoglobinemia or arterial occlusion exists near sensor
- ✧ Obvious venous pulsations
- ✧ Peripheral arterial pulsation becomes weak
- ✧ Peripheral blood supply is not enough

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