


SPER SCIENTIFIC

INSTRUCTION MANUAL



Operation

Count up setting:

Power on, press START button to count up.

Count down setting:

Power on, press the number buttons to set the desire time, once finished press START button to start the count down. When the countdown time has finished, there will be Beeping sound for 60 second, with a RED LED indicator light flashing. Press any key to cancel the alarm and the timer will return back to the set time.

Note:

1. While timer is running, press STOP/RESET button, the timer will pause, press START button to count time again.
2. While timer is running, press STOP/RESET button, timer will pause, press STOP button again, timer will be reset into 00:00.
3. After 10 minutes of no operation, the timer will enter a sleep mode automatically. Press any key to return to wake up the timer.
4. When the battery is low, LCD display may go dark or not display any digits. If this happens, replace the two AAA batteries.