



cimilre S7 Dual Pro

Quick Start Guide

Assembly



Touch Controls



Tap the Cimilre logo to turn on the LED night light.
Tap again to adjust brightness and power off.

Press and hold both the up and down level arrows of either side for 3 seconds to mute or unmute sound.

Night Light



Pump Operation


To begin pumping, press and hold the power control for Pump 1 or Pump 2 for 3 seconds. Once the pump is powered on, tap the power control for the other side once to begin double pumping.



For single pumping, only use one motor. Tap the power control to pause the motor for the side you will not be pumping.

There are two modes. Levels 1-5 are Massage Mode. Levels 6-12 are Expression Mode.

Para esta información en español visite :
cimilrebreastpumps.com/manuals

Pump Operation



Use the   controls on each side to adjust suction levels independently. Tap the SYNC control to use the same suction level on both sides. The SYNC light will turn on.



In SYNC Mode, you have 3 cycle speeds in Massage and Expression Modes. To adjust the cycle speed, tap the CYCLE control. When you are not in SYNC Mode, the pump will default to the highest cycle speed (refer to Operations Manual for chart).



To adjust the suction levels for one side while in SYNC Mode, tap the SYNC control to turn SYNC off. SYNC light will turn off.

To end session in SYNC Mode, tap either power control.
To end session when pumping each side independently, tap the power control for each side.

Thanks for choosing Cimilre for your journey.



Register your pump!

for educational resources and product updates