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Filter rate: 82 % by weight = 65 % by volume. The variation of inorganic filler particle is between 0.02-3 µm.

Indications:

- Caries
- Old restorations
- Lengthening of the incisal edge
- Misaligned teeth
- Tooth fractures
- Tooth discoloration
- Veneering of old crowns
- Incorrect chromaticity
- Anatomical deformities
- Diastemas

Technical data:

Filler rate: 82 % by weight = 65 % by volume. The variation of inorganic filler particle is between 0.02-3 µm.

Instructions for use:

- Attrition
- Abrasion
- Erosion
- Cosmetic corrections

Recommended application:**1. Cleaning the teeth**

Clean the tooth to be treated as well as the mesial and distal adjacent teeth using fluoride-free cleaning paste. Roughen/clean interdental areas using abrasive strips.

2. Selection of the tooth shape

Selection of the tooth shape is performed using the sizing template provided (e.g. edelweiss VENEER Sizing Guide). The transparent sizing template is placed over the teeth to be restored in order to provide a visible profile when selecting the best fitting edelweiss VENEER. Necessary corrections can also be marked on the sizing template. Adjustments to the shape of the edelweiss VENEER can be best performed using a rough surfaced disc at a slow speed without water.

3. Selection of the enamel and dentin shades

Shade selection is performed before drying using the edelweiss Dentine & Enamel Shade Guide in broad, natural daylight or by using a light that represents natural daylight conditions.

4. Preparation of the tooth

Shade selection is performed before drying using the edelweiss Dentine & Enamel Shade Guide in broad, natural daylight or by using a light that represents natural daylight conditions.

5. Selection of shade tone

Determine shade tone in the cervical area with the dentine color rod, since color tone is most pronounced here.

6. Insert the dentine rod into the skin bleached enamel shell.

(Before doing this, apply the glycerin for the refraction in the shell)

7. Selection of shade intensity and opacity

Cervical area of the dentine shade / dentin composite

8. Incisal area of the enamel shade / enamel composite

Incisal area of the enamel shade / enamel composite

9. Selection of shade brightness

General preparation

10. Isolation (dental dam)

Dry work area provides the basis for optimal results.

11. Preparation of the tooth

Depending on the indication, minimally-invasive as possible.

12. Matrix / interdental wedges

To protect and separate the adjacent teeth, a thin, transparent matrix band should be seated. If necessary, fix the matrix interproximally using interdental wedges.

13. Pretreatment of edelweiss VENEERs

1. Adjust the fit of the edelweiss VENEER.

2. To improve adhesion, the edelweiss VENEER should be sandblasted on the backside using 25 µm or 50 µm aluminum oxide; or roughened using a diamond composite instrument.

3. Coat the prepared surface of the edelweiss VENEERs using edelweiss VENEER Bond. Apply one drop of edelweiss VENEER Bond onto the edelweiss VENEER using a brush, and massage it in for at least 30 s.

4. Gently dry edelweiss VENEER Bond using an oil- and water-free stream of air, so that a slightly moist layer remains on the inner surface of the edelweiss VENEER.

5. Light cure edelweiss VENEER Bond for at least 20 s.

Note:

The use of a dental dam is recommended. Before use, disinfect the edelweiss VENEER using alcohol.

6. Pretreatment of prepared teeth:

1. Clean tooth or the preparation of the edelweiss VENEER.

2. Etch using 37% phosphoric acid for at least 30 s to 60 s on the enamel and 15 s on the dentin.

3. Thoroughly rinse for approx. 20 s and dry using a gentle stream of oil- and water-free compressed air.

7. Application of adhesive:

1. Dispense one drop of Dentin/Enamel Bond (Peak Univ., Optibond, etc.) onto the microbrush, and massage it onto the tooth / preparation for 20 s.

2. Briefly dry and pre-cure for 30 s using a light intensity of less than 1000 mW/cm².

Note:

Avoid excessive residue, since this can impact the proper fit of the edelweiss VENEER.

8. Cementation of the edelweiss VENEER

In order to guarantee esthetic success, it is highly recommended to use edelweiss Nano-Hybrid Composite.

9. Select the composite shade. For substrates that are severely broken down or discolored a dentin shade is recommended to provide a homogeneous color. For shape corrections, an enamel composite can be used (Enamel). Too much enamel can make the restoration gray and too transparent.

10. Take roughly 1 cm of the composite strip in the fingertips and roll it into a ball. Use the finger tips to press the ball into the inside of the VENEER.

11. Use the spatula instrument to evenly distribute the composite into the edelweiss VENEER.

12. Press the central VENEERs in desired positions.

<p