Hinweis:

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Instruction for use**English****Definition**

Effect Shade is a light-curing, thin-flowing composite which has been specifically developed for the characterization of composites. Due to the special pigmentation, Effect Shades allow color individualization of restorations.

Indication

Light-curing, thin-flowing composite for the characterization of composite restorations

Recommended use

1. Select colors needed for the desired characterizations.
2. Place a small amount of the selected color on a mixing palette.
3. Apply the material using a thin brush or any other appropriate instrument in a layer thickness of ≤ 1 mm.
4. Light-cure Effect Shades using commercially available light-curing units (e.g. LED polymerization device) for at least 20 s.
5. As a general rule, Effect Shades should be covered with a layer of restorative composite, flowable composite or a restoration such as a veneer. If this is not done, the effect shade might be worn away when the restoration is polished or as well when submitted to brushing and function.
6. Effect Shades may only be mixed with one another (to obtain new color shades); they are not recommended to be mixed with other products.

Please note

- Light-curing units should have a light intensity of at least 500 mW/cm², emit at 450 nm and their output should be regularly checked. The light should be placed as close as possible to the material.
- To prevent premature polymerization due to ambient light, the syringes should be closed immediately after use and the prepared materials on the mixing palette should be covered with a black lid.
- It is important that the surface of the restoration is not contaminated before applying Effect Shades.
- Ensure that the covers of the color stains are not mixed up with one another.

Contraindications

Do not use the material in the case of a known allergy to one of the components.

Interactions

Avoid moisture, materials containing eugenol and oleaginous air; they can inhibit polymerization at the contact point. Discolorations can occur if cationic mouthwashes, plaque revealers or chlorhexidine are used.

Side effects

There are no known side effects to date. Isolated cases of hypersensitivity to components of the material cannot be excluded. Avoid skin contact with the paste. If there is accidental contact, wash with soap and water. In the event of eye contact, immediately rinse with plenty of water and seek medical attention.

Composition

Inorganic fillers in a matrix of dental resins. Filler content: 64% by weight = 38% by volume. The variation width of the inorganic filler particles is between 0.02-3 µm.

Storage

Store at room temperature (15-25°C / 59-77°F). Immediately close syringe after use so that the material does not cure due to ambient light. Do not use after the expiration date.

Note

Only to be supplied to dentists and dental laboratories or upon their request. Keep out of the reach of children! For dental use only!

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